



Covidien Adult Briefs

Recommended Application Techniques

Standing Application



Fold brief in half with soft inner lining exposed and leg ruffles folded out. With fold at top, place between legs.




Open back of brief over buttocks. Position tab closure at waist and hold brief in position.



Bring front of brief across abdomen. Smooth out for a comfortable fit.



Fasten tabs to front of brief.



Gently pull brief upwards until the crotch area fits snugly to the groin and inner leg. Re-adjust fastening tabs as required, applying pressure to the entire length and width of the tab.

Seated Application



Place brief on chair with back open and slightly lower than waist level: front of brief closed with leg ruffles folded out.




Have patient/resident sit on brief and bring front of brief through legs ensuring leg ruffles remain folded out.



Open across abdomen. Smooth out for a comfortable fit.




Fasten tabs to front of brief.




Gently pull brief upwards until the crotch area fits snugly to the groin and inner leg. Re-adjust fastening tabs as required, applying pressure to the entire length and width of the tab.


Bed/Contractured Application



Roll individual onto side. Prepare brief by tri-folding edge closest to individual. Slide tri-folded edge under individual, ensuring proper positioning. Roll individual back onto the opened brief.




Roll front of brief with backing together into a "log roll". Bring the "log roll" of the brief through the legs.



Open and smooth across abdomen for a comfortable fit.



Fasten tabs to front of brief.



Gently pull brief upwards until the crotch area fits snugly to the groin and inner leg. Re-adjust fastening tabs as required, applying pressure to the entire length and width of the tab.

Removal



Hold the front of the brief near the tab with one hand. Slowly peel the tab back with the other hand.

Trouble Shooting Guide

Problem	Cause	Correction
Redness in groin area	<ul style="list-style-type: none"> Poor fit Leg ruffles not turned out away from skin 	<ul style="list-style-type: none"> Proper fit, with brief pulled up snugly in groin area Leg ruffles turned out away from skin
Redness in inner thigh area	<ul style="list-style-type: none"> Poor cleansing techniques Baggy, too large brief 	<ul style="list-style-type: none"> Keep skin clean and dry Use correct size brief
Redness in buttocks area	<ul style="list-style-type: none"> Infrequent changes Poor cleansing techniques 	<ul style="list-style-type: none"> More frequent changes Keep skin clean and dry
Blisters	<ul style="list-style-type: none"> Tape tab touching skin Plastic backing touching skin 	<ul style="list-style-type: none"> Secure tape tabs to backing Ensure proper application of brief—no plastic touching skin
Wet linen / clothing	<ul style="list-style-type: none"> Brief is saturated Brief not properly fitted / applied 	<ul style="list-style-type: none"> More frequent changes Snug fit, proper size

Sizing

Youth (Green)
20" – 31" (23"L)

Small (Peach)
20" – 31" (26.5"L)

Medium (White)
32" – 44"

Large (Blue)
45" – 58"

X-Large (Beige)
59" – 64"

XX-Large (Green)
65" – 69"

Bariatric (White)
70" – 95"

Size requirements may vary based on the individual.

Wetness Indicator



Change the brief when the wetness indicator turns green.