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
IFU No. KM6454340

COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.

© 2011 Covidien.

Covidien Inc, 15 Hampshire Street, Mansfield, MA 02048 USA.

www.covidien.com

REV 02/2011

Catheterization Un-Complicated

Tips and Techniques in Preventing Complications

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Always Follow Your Facility's Standard Protocol

1. Set up clean work area

1a Perform hand hygiene. Put on non-sterile gloves.

1b Pick up under-buttocks drape by edge, without contaminating contents, and place under the area to be prepped with shiny side down

2. Prepare patient

2a Place the fenestrated drape shiny side down over the genital area with the urethral area visible.

2b Remove non-sterile gloves. Put on sterile gloves.

2c Squeeze bellows syringe to apply lubricating jelly. Fully lubricate catheter prior to insertion.

Catheterization Tips and Techniques in Preventing Complications

3. Insert catheter

3a Take forceps and separate 5 cotton balls. Open prepping solution and pour over cotton balls. Prep patient using each soaked cotton ball only once.

3b Proceed with catheterization according to facility guidelines.

3c After urine return, insert 1 inch. Inflate balloon with water-filled syringe to anchor catheter. Gently pull to feel for resistance.

3d Remove backing and apply securement device to inner thigh.

4. Secure collection device

4a Secure drainage bag to bed frame below the level of the bladder to ensure unimpeded flow. Do not allow bag to touch floor.

4b Adjust the sheet clip to ensure the extra tubing is looped on the bed, providing constant downward flow.

4c Remove gloves, dispose of supplies, perform hand hygiene.

4d Fill in name and date of insertion on tag and wrap around drain tubing above the area of the drainage bag.

Tips To Reduce Catheter-Associated Urinary Tract Infection

- Cleanse hands before and after any manipulation of the catheter or site.
- **Do not touch anything** which is non-sterile once you put on sterile gloves.
- Make sure the tip of the catheter is well lubricated for easy insertion and to help prevent damage to the urethra.
- **Do not** reinsert catheter if first insertion was unsuccessful.
If the catheter is inserted into the female patient's vagina by mistake, leave it there as a marker until a new catheter is properly placed in the urethra.
- Whenever possible, maintain a closed sterile drainage system after insertion.
- Make sure the catheter drains. Verify that tubing is not kinked or twisted.

Patient Education

1. How does a catheter drain urine from the body?

After the catheter has been inserted, a small balloon at the tip of the catheter is inflated by a syringe filled with sterile water. This maintains the proper position of the catheter inside the bladder. There are drainage holes in the tip of the catheter so urine will drain out of the bladder through the drainage tubing and down into the drainage bag. This system of catheter, tubing, and drainage bag is called a "closed system". It is specifically designed to help prevent urinary tract infections.

2. What care must I take?

The nursing staff knows what they must do to care for you, but there are a few things you can do to help; if the catheter and drainage tubing accidentally disconnects; or the bag falls or touches the floor, **notify your nurse at once** so the necessary steps may be taken to correct the situation.

Since the urine flow from the bladder to the drainage bag depends on gravity, a "downhill" position must be maintained at all times. If you are standing or walking you **must be sure** that the drainage bag stays below your waist.

Notify your nurse if the drainage bag becomes full. **Do not empty the bag yourself**, as a special technique must be followed. The staff will empty the bag and clean the catheter, tubing, and skin around the catheter daily.

3. What points should I remember?

- Don't disconnect, pull on, or twist the catheter and drainage tubing.
- **Notify your nurse immediately** if you notice any leaks or kinks in the system.
- Keep the drainage bag below the level of your waist, whether you are lying, standing or sitting.
- Notify your nurse if the drainage bag becomes full. **Do not empty the bag yourself.**
- Remember to ask your nurse daily if a Foley catheter is necessary.

4. What is a catheter-associated urinary tract infection (UTI)?

When a catheter introduces outside germs inside the body, it can cause an infection. To prevent infections hands **must** be washed before handling a catheter and catheters must be removed as soon as possible.

5. What are the symptoms of a urinary tract infection?

- Burning or painful urination, lower stomach pain
- Bloody urine
- Sudden fever
- Frequent or more urgent urination after catheter is removed