




Step-by-step guide

ProGrip™ self-gripping polyester mesh for incisional hernia prevention



1 Ensure closure of the muscular aponeurotic layers according to usual recommended practices.¹




2 European Hernia Society (EHS) suggestion: Use a small bite suture technique to wound length (SL/WL) ratio. Note: 4/1 is the current suggested method of fascial closure.¹



3 Dissect the subcutaneous space, creating a lodge for implant placement.



4 Place the implant, with self-gripping side on the muscular aponeurotic layer, to reinforce the suture line. Deploy the mesh without folding and ensure a minimum 3 cm mesh placement overlap.²



5 Fixate the mesh.[†]



6 Close the superficial layers.

Hernia care

Mesh. Fixation.
Biologics. Dissection.

Our comprehensive product portfolio can enhance your hernia repair procedures.

† Additional fixation is at the surgeon's discretion.

1. Muysoms FE, Antoniou SA, Bury K, et al; European Hernia Society. European Hernia Society guidelines on the closure of abdominal wall incisions. *Hernia*. 2015;19(1):1-24. doi: 10.1007/s10029-014-1342-5.
2. Jairam AP, Timmermans L, Eker HH, et al; PRIMA Trialist Group. Prevention of incisional hernia with prophylactic onlay and sublay mesh reinforcement versus primary suture only in midline laparotomies (PRIMA): 2-year follow-up of a multicentre, double-blind, randomised controlled trial. *Lancet*. 2017;390(10094):567-576. doi: 10.1016/S0140-6736(17)31332-6.

Important: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2022 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic.

All other brands are trademarks of a Medtronic company. EMEA-HR-2200110-progrip-prophylaxis-in-service-guide-en-we-7584091