

DAR™
Breathing System

en	Instructions For Use
fr	Mode d'emploi
it	Istruzioni per l'uso
de	Gebrauchsanweisung
nl	Gebruiksaanwijzing
es	Instrucciones de uso
sv	Bruksanvisning
da	Brugsanvisning
no	Bruksanvisning
fi	Käyttöohjeet
pt	Instruções de uso
ru	Инструкция по применению
pl	Instrukcja użytkowania
cs	Návod k použití
sk	Návod na použitie
hu	Használati útmutató
el	Οδηγίες χρήσης
tr	Kullanma Talimatı


Fig. A


Fig. B

	Adult-Pediatric Mechanical Filter HME Compact		Adult-Pediatric Mechanical Filter HME Large	
Recommended Vt (tidal volume) indications	200 - 1500 ml		300 - 1500 ml	
Resistance to flow (ISO 9360-1) (ISO 23328-2)				
before use	0.6 cm H ₂ O 1.1 cm H ₂ O 2.4 cm H ₂ O 4.2 cm H ₂ O	at 15 l/min at 30 l/min at 60 l/min at 90 l/min	1.1 cm H ₂ O 2.5 cm H ₂ O 4.2 cm H ₂ O	at 30 l/min at 60 l/min at 90 l/min
after use (24 h)	0.6 cm H ₂ O 1.2 cm H ₂ O 2.7 cm H ₂ O 4.6 cm H ₂ O	at 15 l/min at 30 l/min at 60 l/min at 90 l/min	1.3 cm H ₂ O 2.8 cm H ₂ O 4.7 cm H ₂ O	at 30 l/min at 60 l/min at 90 l/min
Moisture loss (ISO 9360-1)	5 mg H ₂ O / l 5 mg H ₂ O / l 6 mg H ₂ O / l 7 mg H ₂ O / l 7 mg H ₂ O / l 10 mg H ₂ O / l	at Vt 200 ml at Vt 250 ml at Vt 500 ml at Vt 750 ml at Vt 1000 ml at Vt 1500 ml	5 mg H ₂ O / l 5 mg H ₂ O / l 6 mg H ₂ O / l 7 mg H ₂ O / l 7 mg H ₂ O / l 9 mg H ₂ O / l	at Vt 300 ml at Vt 500 ml at Vt 750 ml at Vt 1000 ml at Vt 1500 ml
Gas leakage (ISO 9360-1)				
before use	0 ml / min	at 70 cm H ₂ O	0 ml / min	at 70 cm H ₂ O
after use (24 h)	0 ml / min	at 70 cm H ₂ O	0 ml / min	at 70 cm H ₂ O
Internal Volume (ISO 9360-1) (ISO 23328-2)	66 ml		96 ml	
Compliance (ISO 9360-1)	0,66 ml / KPa		1,14 ml / KPa	
Weight (approx.)	36 g		49 g	
Connectors	22M/15F -22F/15M Luer lock gas sampling port ISO	(ISO 5356-1) (ISO 594-1) (ISO 594-2)	22M/15F-22F/15M Luer lock gas sampling port ISO	(ISO 5356-1) (ISO 594-1) (ISO 594-2)

CLEAN/STERILE (CHECK LABEL) AND READY FOR USE: unless package has been opened or damaged. Single use device to be used on a single patient. Discard after use.

Identification of a substance that is contained or present within the product or packaging.

Identification of a substance that is not contained or present within the product or packaging.

This product contains DEHP. When used as indicated, very limited exposure to trace amounts of DEHP may occur. There is no clear clinical evidence that this degree of exposure increases clinical risk. However, in order to minimize risk of DEHP exposure in children and nursing or pregnant women, this product should only be used as directed. Refer to primary labeling.

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

DESCRIPTION:

DAR™ Adult-Pediatric Mechanical Filter HME Large and Adult-Pediatric Mechanical Filter HME Compact are mechanical bacterial Filters with Heat and Moisture Exchanger (Filter/HMEs) for breathing circuits. These Filter/HMEs help to protect patients, equipment, and hospital staff from cross-contamination while minimizing heat and moisture loss from patient airways. These devices also include a Luer lock fitting, luer lock cap, and a non-threaded cap mount for using during CO₂ monitoring.

INDICATIONS:

These Filter/HMEs are indicated for single use on patients under anesthesia or in intensive care as follows:

- Adult Patients: Adult-Pediatric Mechanical Filter HME Large, Adult-Pediatric Mechanical Filter HME Compact
- Pediatric Patients: Adult-Pediatric Mechanical Filter HME Compact

The Filter/HMEs features a passive heat and moisture exchanger changing the humidifying performance in relation to the patient's heat and moisture output. Alternative active humidification systems/devices can be used at the clinician's discretion. The Filter/HMEs must be placed at the patient's side, between the tracheal tube and breathing circuit in the positions indicated in figures A and B. Consult the attending physician for tidal volumes (Vt).

CONTRAINDICATIONS:

- Do NOT use these devices in conjunction with heated humidifiers or nebulizers.
- Do NOT use these devices in a position other than described.

DIRECTIONS FOR USE:

1. Remove the Filter/HME from its protective package.
2. Firmly connect the device to the breathing circuit as close to the patient as possible. Check all connections for an airtight seal.
3. Adjust tidal volume to compensate for dead space.
4. For CO₂ monitoring, attach the monitoring line to the Luer lock port.

Periodically inspect any increasing in the resistance to flow and replace the device if necessary.

Maximum Filter/HME use: 24 hours. Do not reuse. Discard after use.

Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

NOTE: the device is latex-free and non conductive.

WARNINGS:

- **Prescription device to be used only under medical supervision.**
- **The product is guaranteed only if the package has not been damaged. Check integrity of the filter package; any damage or opening may compromise sterility and/or performance of the device. In this case, do not use the filter concerned.**
- **Excessive secretions/bleedings can block the Filter/HME, making breathing difficult; in this case, replace the device immediately.**
- **Ensure that the connections to the breathing circuit are secure to reduce risk of accidental disconnection.**
- **Ensure the patency of the device.**
- **Do not use on patients having a tidal volume (Vt) lower than those indicated in the table.**
- **If the device is used with spontaneously breathing patients on assist-control or intermittent mandatory ventilation (IMV), keep the ventilation parameters monitored.**
- **If the device is used during anesthesia with circle breathing systems with soda lime canister, it is necessary to periodically monitor the resistance to flow. It is not recommended to use the device in this manner longer than 6 hours.**
- **Before use, check that the Luer lock cap on the CO₂ monitoring port is in place and secure, if the port is not in use.**
- **Do not connect any tubes to the Luer lock port other than CO₂ sampling lines. Connections to other devices, e. g. gas administration tubes, could be dangerous for the patient.**
- **The hydrophobic membrane does not allow the passage of liquid in either direction under normal pressure conditions.**
- **To reduce the risk of damage and obstruction by small objects, keep wrapped until required for actual use.**

PRECAUTIONS:

- Do not try to refurbish. Do not reuse. Do not soak, rinse, wash, sterilize or treat with any disinfectants (in particular phenolic and alcoholic-based solutions are to be avoided).
- Do not attach CO₂ monitoring line to the non-threaded Luer lock cap mount.
- Adjust ventilator to account for added dead space, especially for those patients with small tidal volumes (Vt). The effect of mechanical dead space should be evaluated individually for each patient.
- Use only with devices having ISO standard connectors.

POTENTIAL ADVERSE REACTIONS:

The following adverse reactions may result from the use of Filter/HMEs: mucus plugging and/or complications due to device occlusion including difficult breathing, dyspnea, hypercapnia, and hypoxia. The order of listing is alphabetical and does not reflect anticipated frequency or severity.

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© 2011 Covidien.
COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.
Covidien Ireland Limited,