

# DAR™ Breathing System

en	Instructions For Use
fr	Mode d'emploi
it	Istruzioni per l'uso
de	Gebrauchsanweisung
nl	Gebruiksaanwijzing
es	Instrucciones de uso
sv	Bruksanvisning
da	Brugsanvisning
no	Bruksanvisning
fi	Käyttöohjeet
pt	Instruções de uso
ru	Инструкция по применению
pl	Instrukcja użytkowania
cs	Návod k použití
hu	Használati útmutató
el	Οδηγίες χρήσης
tr	Kullanma Talimatı

	Mechanical Filter Large		Mechanical Filter Compact		Mechanical Filter Small	
<b>Recommended Vt (tidal volume) indications</b>	300- 1500 ml		200- 1500 ml		150- 1200 ml	
<b>Resistance to flow (ISO 9360-1 EN 23328-2) before use</b>	0.8 cm H <sub>2</sub> O at 30 l/min		0.8 cm H <sub>2</sub> O at 30 l/min		0.5 cm H <sub>2</sub> O at 15 l/min	
<b>after use (24 h)</b>	2 cm H <sub>2</sub> O at 60 l/min		1.9 cm H <sub>2</sub> O at 60 l/min		1.2 cm H <sub>2</sub> O at 30 l/min	
	3.6 cm H <sub>2</sub> O at 90 l/min		3.2 cm H <sub>2</sub> O at 90 l/min		2.7 cm H <sub>2</sub> O at 60 l/min	
	0.8 cm H <sub>2</sub> O at 30 l/min		0.8 cm H <sub>2</sub> O at 30 l/min		0.6 cm H <sub>2</sub> O at 15 l/min	
	2 cm H <sub>2</sub> O at 60 l/min		2 cm H <sub>2</sub> O at 60 l/min		1.4 cm H <sub>2</sub> O at 30 l/min	
<b>Gas leakage (ISO 9360-1 ISO 23328-2)</b>	3.6 cm H <sub>2</sub> O at 90 l/min		3.4 cm H <sub>2</sub> O at 90 l/min		2.8 cm H <sub>2</sub> O at 60 l/min	
	0.8 cm H <sub>2</sub> O at 30 l/min		0.8 cm H <sub>2</sub> O at 30 l/min		4.8 cm H <sub>2</sub> O at 90 l/min	
	2 cm H <sub>2</sub> O at 60 l/min		2 cm H <sub>2</sub> O at 60 l/min		2.8 cm H <sub>2</sub> O at 60 l/min	
<b>before use</b>	0 ml / min at 70 cm H <sub>2</sub> O		0 ml / min at 70 cm H <sub>2</sub> O		0 ml / min at 70 cm H <sub>2</sub> O	
<b>after use (24 h)</b>	0 ml / min at 70 cm H <sub>2</sub> O		0 ml / min at 70 cm H <sub>2</sub> O		0 ml / min at 70 cm H <sub>2</sub> O	
<b>Internal Volume (ISO 23328-2)</b>	92 ml		66 ml		42 ml	
<b>Weight (approx.)</b>	47 g		39 g		24 g	
<b>Connectors</b>	22M/15F - 22F/15M 22M-22F	(ISO 5356-1)	22M/15F - 22F/15M 22M/15F-22F	(ISO 5356-1)	22M/15F - 22F/15M	(ISO 5356-1)
	CO <sub>2</sub> Luer monitoring port	(ISO 594-1) ISO 594-2)	CO <sub>2</sub> Luer monitoring port	(ISO 594-1) ISO 594-2)	CO <sub>2</sub> Luer monitoring port	(ISO 594-1) ISO 594-2)

**CLEAN/STERILE (CHECK LABEL) AND READY FOR USE:** unless package has been opened or damaged. Single-use device to be used on a single patient. Discard after use.

Identification of a substance that is not contained or present within the product or packaging.

Identification of a substance that is contained or present within the product or packaging.

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

This product contains DEHP. When used as indicated, very limited exposure to trace amounts of DEHP may occur. There is no clear clinical evidence that this degree of exposure increases clinical risk. However, in order to minimize risk of DEHP exposure in children and nursing or pregnant women, this product should only be used as directed. Refer to primary labeling.

#### DESCRIPTION:

The large, compact and small mechanical filters are mechanical bacterial/viral filters for breathing circuits, that help to protect patients, equipment, and hospital staff from cross-contamination. These devices also include a Luer lock fitting, Luer lock cap, and a non-threaded cap mount for use during CO<sub>2</sub> monitoring.

#### INDICATIONS:

These filters are indicated for single use on adult patients (the small mechanical filter is indicated on pediatric patients) under anaesthesia or in intensive care. They must be connected on patient's side, between the tracheal tube and breathing circuit in position A or for ventilator protection on the machine's side, as indicated in B. For tidal volumes (Vt), physician judgement prevails. These filters are not equipped with heat and moisture exchanger (HME) for this reason they have limited humidification functionality. For proper humidification and for long-term application, it is advisable to evaluate whether active humidification systems/devices or Filter/HME's are needed.

#### CONTRAINDICATIONS:

- Do NOT use these devices placed between the patient and the heated humidifiers or nebulizers.
- Do NOT use this device in a position different from the intended one.

#### DIRECTIONS FOR USE:

1. Remove the Filter from its protective package.
2. Firmly connect the device to the breathing circuit and the tracheal tube. Check air seal of all connections.
3. Adjust tidal volume to account for dead space.
4. For CO<sub>2</sub> monitoring, attach the monitoring line to the Luer lock port.

Periodically inspect for any increase in the resistance to flow, and replace the device if necessary.

Maximum Filter use: 24 hours. Do not reuse. Discard after use.

Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

**NOTE:** These devices are latex-free and non conductive.

#### WARNINGS:

- Prescription device to be used only under medical supervision.
- The product is guaranteed only if the package has not been damaged. Check integrity of the filter package; any damage or opening may compromise sterility and/or performance of the device. In this case, do not use the filter concerned.
- Excessive secretions/bleedings can block the Filter, making breathing difficult; in this case, replace the device immediately.
- Ensure that the connections to the breathing circuit are secure to reduce risk of accidental disconnection.
- Ensure the patency of the device.
- If these Filters are used on the patient's side, do not use them on patients having a tidal volume lower than those indicated in the table.
- If the device is used during anaesthesia with circle breathing systems with soda lime canister, it is necessary to periodically monitor the resistance to flow. It is recommended not to use the device in this condition for more than six hours.
- Before use, check that the Luer lock port on the CO<sub>2</sub> monitoring port is in place and secure.
- Do not connect to the Luer lock port tubes different from monitoring CO<sub>2</sub> sampling line. Connections to other devices, i.e. gas administration tubes, could be dangerous for the patient.
- The hydrophobic membrane does not allow the passage of liquid in either direction under normal pressure conditions.
- To reduce the risk of damage and obstruction by small objects, keep wrapped until required for actual use.

#### PRECAUTIONS:

- Do not try to refurbish. Do not reuse. Do not soak, rinse, wash, sterilize or treat with any disinfectants (in particular phenolic and alcoholic-based solutions are to be avoided).
- Do not attach CO<sub>2</sub> monitoring line to the non-threaded Luer lock cap mount.
- Adjust ventilator to account for added dead space, especially for those patients with small tidal volumes (Vt). The effect of mechanical dead space should be evaluated individually for each patient.
- Use only with devices having ISO standard connectors.

#### POTENTIAL ADVERSE REACTIONS:

The following adverse reactions may result from the use of Filters: mucus plugging and/or complications due to device occlusion including difficult breathing, dyspnea, hypercapnia, and hypoxia. The order of listing is alphabetical and does not reflect anticipated frequency or severity.

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© 2011 Covidien.
COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.
Covidien Inc, 15 Hampshire Street,
Covidien Ireland Limited,