



Products for Non-Intubated Patients Directions for Use

Omnistream™ EtCO₂ Sampling for Non-Intubated Patients

Single patient use, not to be re-processed

Non-sterile / non-pyrogenic

Does not contain natural rubber latex

Warning: Loose or damaged connections may compromise ventilation or cause an inaccurate measurement of respiratory gases. Securely connect all components and check connections for leaks according to standard clinical procedures.

Warning: Do not cut or remove any part of the sample line. Cutting the sample line could lead to erroneous readings.

Warning: Carefully route the FilterLine to reduce the possibility of patient entanglement or strangulation.

Warning: Check CO₂ and O₂ tubing regularly during use to ensure that no kinks are present. Kinked tubing may cause inaccurate CO₂ sampling or affect O₂ delivery to patient.

Caution: Ensure that tubing is not stretched during use.

Caution: Do not attempt to clean, disinfect, sterilize or flush any part of the sampling line as this can cause damage to the monitor.

Caution: Dispose of sampling lines according to standard operating procedures or local regulations for the disposal of contaminated medical waste.

Caution: Use of the Smart OmniLine H Plus O₂ CO₂ sampling line (the H in its name indicates that it is for use in humidified environments) during MRI scanning may cause interference. The use of non H sampling lines is advised.

Note: Replace the sampling line according to hospital protocol or when a blockage is indicated by the device. Excessive patient secretions or a build-up of liquids in the airway tubing may occlude the sampling line, requiring more frequent replacement.

Note: When connecting a sampling line to the monitor, insert the sampling line connector clockwise into the monitor CO₂ port and turn the connector until it can no longer be turned, to ensure that it is connected securely to the monitor. This will assure that there is no leak of gases during measurement at the connection point and that measurement accuracy is not compromised.

Note: Following connection of the CO₂ sampling line, check that CO₂ values appear on the monitor display.

Note: Sampling lines with H in their names include a moisture reduction component (Nafion®* or its equivalent) for use in higher humidity environments where long duration use of CO₂ sampling is required.

Note: Sampling lines are recommended for use with oxygen provided at up to 5l/min. At higher levels of oxygen provision, dilution of CO₂ readings may occur, leading to lower CO₂ values.

Federal law (U.S.) restricts this device to sale by or on the order of a physician.

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 Single patient use	For use with a single patient
Rx ONLY	Rx Only
 Consult instructions for use	Consult instructions for use
 Caution, consult accompanying documents	Caution, consult accompanying documents
 -4°F -20°C 131°F 55°C Temperature limitation	Storage temperature
 PHT DEHP	DEHP
 This product does not contain natural rubber latex	Does not contain natural rubber latex
 MR Unsafe	Not safe for use in MR environments
 Fire hazard	Fire Hazard
 CE 0482	CE Mark

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Oridion®, Microstream®, FilterLine®, Smart CapnoLine®, CapnoLine®, Smart CapnoLine Guardian™, and Omnistream™ are trademarks or registered trademarks of Oridion Medical 1987 Ltd. ™* and ®* are trademarks of their respective owners.

Oridion Medical 1987 Ltd.
7 Hamarpe St.
P.O. Box 45025
Jerusalem 9777407, Israel
US Toll Free: 1-888-674-3466
1-888-ORIDION

Obelis S.A.
Bd Général Wahis 53
1030 Brussels
Belgium

www.oridion.com