Nellcor™

Oxygen Saturation Accuracy Specification Grid

en	Oxygen Saturation Accuracy Specification Grid	REF N-65
fr	Grille de spécifications de précision de la saturation en oxygène	REF N-85
ar	جدول مواصفات دقة التشبع بالأكسجين	REF N-560
bg	Таблица със спецификации за точност на кислородна сатурация	REF N-600
zh	氧饱和和精确度规格坐标网格	REF N-600X
hr	Specifikacijska mreža preciznosti zasićenja krvi kisikom	REF 10005941
cs	Předpokládaná přesnost naměřené saturace kyslíkem	REF GR101704
da	Nøjagtighedsspecifikationskema for iltmætning	REF CS08651
nl	Specificatierooster voor zuurstofverzadigingsnauwkeurigheid	REF CS08652
et	Hapniku saturatsiooni täpsuse tehnilised andmed	REF CS08657
fi	Happikyllästeisyyden tarkkuustaulukko	REF CS08658
ka	ჟანგბადის გაჯერების სიზუსტის მახასიათებლების სქემა	
de	Raster zur Genauigkeitsspezifikation der Sauerstoffsättigung	
el	Πίνακας προδιαγραφών ακριβείας κορεσμού οξυγόνου	
hu	Az oxigén-szaturáció túréshatárainak táblázata	
it	Griglia dei dati tecnici sull'accuratezza della saturazione d'ossigeno	
ja	酸素飽和精度仕様グリッド	
ko	산소 포화 정확도 사양표	
lv	Skābekļa piesātinājuma mērījumu precizitātes tabula	
lt	Deguonies saturacijos tikslumo specifikacijų lentelė	
no	Tabell med nøyaktighetsspesifikasjoner for oksygenmetning	
pl	Tabela dokładności określenia nasycenia krwi tlenem (saturacji)	
pt	Grelha de especificação da precisão do nível de saturação do oxigénio	
ro	Grila specificațiilor de precizie pentru saturația oxigenului	
ru	Таблица точности показаний насыщения кислородом	
sr	Tablica za specifikaciju preciznosti zasićenja kiseonikom	
sk	Tabuľka s údajmi o presnosti merania saturácie kyslíkom	
sl	Tabela s podatki o natančnosti saturacije krvi s kisikom	
es	Tabla de especificaciones de precisión de la saturación de oxígeno	
sv	Specifikationstabell över noggrannhet avseende syremättnad	
tr	Oksijen Doymunluğ Doğruluk Özellikleri Tablosu	

Nellcor™

Oxygen Saturation Accuracy Specification Grid

Accuracy Specifications

Accuracy specifications are based on controlled hypoxia studies with healthy non-smoking adult volunteers over the specified saturation SpO₂ range(s).¹ Pulse oximeter SpO₂ readings were compared to SaO₂ values of drawn blood samples measured by hemoximetry. All accuracies are expressed as ± "X" digits. Pulse oximeter equipment measurements are statistically distributed about two-thirds of pulse oximeter measurements can be expected to fall in this accuracy (ARMS) range. Because scatter and bias of pulse oximeter SpO₂ and blood SaO₂ comparisons commonly increase as the saturation decreases, and accuracy specifications are calculated from data spanning the stated range, different accuracy values may result when describing partially overlapping ranges.

Oxygen saturation accuracy can be affected by certain environmental, equipment, and patient physiologic conditions (as discussed in the operator's manual for the monitor) that influence readings of SpO₂, SaO₂, or both. Accordingly, observations of clinical accuracy may not achieve the same levels as those obtained under controlled laboratory conditions.

Functional testers or patient simulators cannot be used to assess or validate the accuracy of pulse oximeter equipment. SpO₂ measurement accuracy can only be evaluated in vivo by comparing oximeter readings with SaO₂ measurements made using a laboratory CO-oximeter obtained from simultaneously sampled arterial blood.

Compatibility Information

Use Nellcor sensors only with Nellcor patient monitors and monitors containing Nellcor oximetry, or with instruments licensed to use Nellcor sensors (Nellcor-compatible instruments). Consult individual manufacturers for accuracy specifications and compatibility information of particular instruments and Nellcor sensor models.

Range of Peak Wavelength and Maximum Output Power

Use Nellcor sensors only with Nellcor patient monitors and monitors containing LEDs that emit red light at a wavelength of approximately 660 nm and infrared light at a wavelength of approximately 900 nm. The total optical output power of the sensor LEDs is less than 15 mw.

	SpO ₂ Range ¹	
	LoSat™ 60% - 80%	70% - 100%
MAXA, MAXAL	± 3	± 2
MAXN ² (Adult & Neonate)	± 3	± 2
MAXP	± 3	± 2
MAXI	± 3	± 2
MAXFAST	± 3	± 2
SC-A (Adult)		± 2
SC-PR ³ (Neonate)		± 2
SC-NEO ³ (Neonate)		± 2
MAXR ⁴		± 3.5

1. Subjects used to validate SpO₂ measurement accuracies were healthy and recruited from local population. Comprised of both men and women, subjects spanned a range of skin pigmentations and ranged in age from 18-50 years old.

2. MAXN

Clinical functionality has been demonstrated on a population of hospitalized neonate patients. The observed SpO₂ accuracy was 2.5% in a study of 42 patients with ages of 1 to 23 days, weight from 750 to 4,100 grams, and 63 observations made spanning a range of 85 to 99% SaO₂.

3. Softcare® SC-PR, SC-NEO

Clinical functionality has been demonstrated on a population of hospitalized neonate and infant patients. The observed SpO₂ accuracy was 3.0% in a study of 57 patients with ages of 24 to 40 weeks, weight from 710 to 5,000 grams, and 185 observations made spanning a range of 63 to 100% SaO₂.

4. The accuracy specification has been determined between saturations of 80%-100%.

Nellcor™

Oxygen Saturation Accuracy Specification Grid

Accuracy Specifications

Accuracy specifications are based on controlled hypoxia studies with healthy non-smoking adult volunteers over the specified saturation SpO₂ range(s).¹ Pulse oximeter SpO₂ readings were compared to SaO₂ values of drawn blood samples measured by hemoximetry. All accuracies are expressed as ± "X" digits. Pulse oximeter equipment measurements are statistically distributed; about two-thirds of pulse oximeter measurements can be expected to fall in this accuracy (ARMS) range.

Neonatal Accuracy

When sensors are used on neonatal subjects as recommended, the **specified accuracy range is increased by ± 1 digit**, as compared to adult usage, to account for the theoretical effect on oximeter measurements of fetal hemoglobin in neonatal blood. For example, N accuracy on neonates is ± 3.5 digits, rather than ± 2.5.

Oxygen saturation accuracy can be affected by certain environmental, equipment, and patient physiologic conditions (as discussed in the operator's manual for the monitor) that influence readings of SpO₂, SaO₂, or both. Accordingly, observations of clinical accuracy may not achieve the same levels as those obtained under controlled laboratory conditions.

Functional testers or patient simulators cannot be used to assess or validate the accuracy of pulse oximeter equipment. SpO₂ measurement accuracy can only be evaluated in vivo by comparing oximeter readings with SaO₂ measurements made using a laboratory CO-oximeter obtained from simultaneously sampled arterial blood.

Compatibility Information

Use Nellcor pulse oximeter sensors only with Nellcor patient monitors and monitors containing Nellcor oximetry, or with instruments licensed to use Nellcor sensors (Nellcor-compatible instruments). Consult individual manufacturers for accuracy specifications and compatibility information of particular instruments and Nellcor sensor models.

Range of Peak Wavelength and Maximum Output Power

Use Nellcor pulse oximeter sensors only with Nellcor patient monitors and monitors containing LEDs that emit red light at a wavelength of approximately 660 nm and infrared light at a wavelength of approximately 900 nm. The total optical output power of the sensor LEDs is less than 15 mw.

	SpO ₂ Range
Nellcor™ (Single Patient Use)	70% - 100%
A	± 2.5
P	± 2.5
N (Adult)	± 2.5
N (Neonate)	± 3.5
I	± 2.5
Nellcor™ Reusable Sensors	70% - 100%
D-YS (Infant to Adult)	± 3
D-YS (Neonate)	± 4
D-YS with D-YSE Ear Clip	± 3.5
D-YS with D-YSPD SpotClip	± 3.5
DS-100A	± 3
OXI-A/N (Adult)	± 3
OXI-A/N (Neonate)	± 4
OXI-P/I	± 3

1. Subjects used to validate SpO₂ measurement accuracies were healthy and recruited from local population. Comprised of both men and women, subjects spanned a range of skin pigmentations and ranged in age from 18-50 years old.

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© 2012 Covidien.
COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.
Covidien Ireland Limited,