# Nellcor™ Neonatal SpO<sub>2</sub> Sensor Non-Adhesive

REF SCNEO

REF SCNEO-I

- en** Instructions For Use
- fr** Mode d'emploi
- de** Gebrauchsanweisung
- nl** Gebruiksaanwijzing
- es** Instrucciones de uso
- it** Istruzioni per l'uso
- sv** Bruksanvisning
- da** Brugsanvisning
- pt** Instruções de uso
- ru** Инструкция по применению
- zh** 使用说明

Identification of a substance that is contained or present within the product or packaging.

Do not use if package is opened or damaged

Identification of a substance that is not contained or present within the product or packaging.

**IP22** Protected against access to hazardous parts with a finger and protected against vertically falling water drops when sensor tilted up to 15°.

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

This device is not made with natural rubber latex or DEHP.

## Neonatal SpO<sub>2</sub> Sensor Non-Adhesive Directions for Use

### Indications/Contraindications

The Nellcor™ Neonatal SpO<sub>2</sub> Sensor Non-Adhesive, model SCNEO, is indicated for single patient use when continuous noninvasive arterial oxygen saturation and pulse rate monitoring are required for neonates weighing 1.5 kg to 5 kg.

Use this sensor only with Nellcor instruments and instruments containing Nellcor OxiMax™ technology, or with instruments licensed to use Nellcor OxiMax sensors (Nellcor-compatible instruments). When connected to an OxiMax - compatible instrument, this sensor uses OxiMax technology to provide additional advanced sensor performance features. Consult individual manufacturers for features and compatibility of particular instruments and sensor models.

Each manufacturer of Nellcor-compatible instruments is responsible for determining whether and under what conditions its instruments are compatible for safe and effective use with each Nellcor sensor model. This may include different specifications and/or warnings, cautions, or contraindications. Refer to the instrument operator's manual or consult manufacturer for complete instructions for use of this sensor with their Nellcor-compatible instrument.

**Note:** Non-adhesive sensor consists of sensor, cable, wrap and loop fastener.

### Instructions for Use

- 1) Locate transparent windows on the sensor. Windows cover optical components. This is the side that contacts the patient. Note corresponding alignment marks (a) on the outside (b) between the marks **A**.
- 2) Orient the SCNEO so the graphics are on the lateral edge of the site:  
  
The preferred site is a foot. Alternatively, use a hand. The window next to the cable goes on the sole of the foot as shown **B**.  
  
Note: When selecting a sensor site, priority should be given to an extremity free of an arterial catheter, blood pressure cuff, or intravascular infusion line.
- 3) Wrap the SCNEO firmly, but not too tightly around the foot. Windows must oppose each other. Secure the hook and loop fastener.
- 4) Use the cable wrap—attached to the sensor cable—to secure the sensor cable to the patient's ankle. **C** Secure with the hook and loop fastener.
- 5) Plug the SCNEO into the oximeter and verify proper operation as described in the oximeter operator's manual.

Note: If the sensor does not track the pulse reliably, it may be incorrectly positioned—or the sensor site may be too thick, thin, or deeply pigmented, or otherwise deeply colored (for example, as a result of externally applied coloring such as nail polish, dye, or pigmented cream) to permit appropriate light transmission. If any of these situations occurs, reposition the sensor or choose an alternate Nellcor sensor for use on a different site.

### Reapplication

The SCNEO can be reused *on the same patient* as long as the sensor attaches without slippage.

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Zusätzliche Exemplare dieser Gebrauchsanweisung sind kostenlos von Covidien bzw. von autorisierten Vertretungen von Covidien erhältlich. Unter dem Copyright von Covidien dürfen Käufer von Produkten, die von Covidien oder zugelassenen Händlern gekauft wurden, zur eigenen Nutzung Kopien dieser Anleitung erstellen.
COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.
Covidien Ireland Limited, IDA Business & Technology Park, Tullamore.