




Adult Flexible Tracheostomy Tube Cuffless
Disposable Inner Cannula

en	Instructions For Use	REF 4UN65A
fr	Mode d'emploi	REF 5UN70A
de	Gebrauchsanweisung	REF 6UN75A
nl	Gebruiksaanwijzing	REF 7UN80A
it	Istruzioni per l'uso	REF 8UN85A
es	Instrucciones de uso	REF 9UN90A
sv	Bruksanvisning	REF 10UN10A
da	Brugsanvisning	
no	Bruksanvisning	
fi	Käyttöohjeet	
pt	Instruções de uso	
ru	Инструкция по применению	
zh	使用说明	
pl	Instrukcja użytkowania	
cs	Návod k použití	
sk	Návod na použitie	
sl	Navodila za uporabo	
hu	Használati utasítás	
el	Οδηγίες χρήσης	
tr	Kullanma Talimatları	
ar	تعليمات الاستخدام	
bg	Инструкции за употреба	
ro	Instrucțiuni de utilizare	
ko	사용 설명서	
lt	Naudojimo instrukcija	
th	คำแนะนำการใช้งาน	
vi	Hướng dẫn sử dụng	
id	Petunjuk Penggunaan	

								
REF	(mm)	(mm)	θ (°)	C (mm)	(mm)	(mm)	(mm)	(fr)
4UN65A	6.5	9.4	73	62	5.5	62	21	✓
5UN70A	7.0	10.1	79	68	6.0	68	21	✓
6UN75A	7.5	10.8	85	74	6.5	74	24	✓
7UN80A	8.0	11.4	88	77	7.0	77	26	✓
8UN85A	8.5	12.2	90	79	7.5	79	27	✓
9UN90A	9.0	12.7	90	79	8.0	79	n/a	✓
10UN10A	10.0	13.8	90	79	9.0	79	n/a	✓

REF		Part Number, Numéro de pièce, Teilenummer, Onderdeelnummer, Codice, Número de pieza, Artikelnummer, Varennummer, Delenummer, Osanumero, Número da peça, Καταложный номер, 零件号, Numer katalogowy, Číslo dílu, Číslo dielu, Številka dela, Termékszám, Кодіток, εξεργίματος, Paçça Numarası, رقم القطعة, Καταложен номер, Код пієса, 부품 번호, Dalies numeris, 零件編號, Sô hiệu chi tiết, Normor Komponenten	ID – Inside Diameter (mm), DI – Diamètre interne (mm), ID – Innendurchmesser (mm), ID – Binnendiameter (mm), ID – Diametro interno (mm), DI – Diámetro interior (mm), ID – Innerdiameter (mm), ID – Indvendig diameter (mm), ID – Indre diameter (mm), ID – Sisäläpimitta (mm), DI – Diámetro interno (mm), BD – Внутренний диаметр (мм), ID – 內徑 (毫米), ID – Srednica wewnętrzna (mm), ID – Vnitřní průměr (mm), Vnut. pr. – Vnutorný priemer (mm), ID – Notranji premer (mm), ID – Belső átmérő (mm), ID – Εσωτερική διάμετρος (mm), ID – İç Çap (mm), ID – (قطر الداخلي (مم), ID – Вършен диаметър (mm), ID – Diametru interior (mm), ID – 내부 직경 (mm), ID – Vidinis skersmuo (mm), ID – (قطر الداخلي (مم), ID – Đường kính trong (mm), ID – Diameter Dalam (mm)
------------	---	--	---

	<p>OD – Outside Diameter (mm), DE – Diamètre externe (mm), AD – Außendurchmesser (mm), OD – Buitendiameter (mm), OD – Diametro esterno (mm), DE – Diámetro exterior (mm), YD - Ytterdiameter (mm), OD – Udvendig diameter (mm), OD – Ytre diameter (mm), OD – Ulkoläpmitta (mm), DE – Diámetro Externo (mm), HD – Наружный диаметр (мм), OD – 外径 (毫米), OD – Średnica zewnętrzna (mm), OD – Vnější průměr (mm), Von. pr. – Vonkajási priemer (mm), OD – Zunanji premer (mm), OD – Külső átmérő (mm), OD – Εξωτερική διάμετρος (mm), OD – Diş Çap (mm), OD – قطر الخارجي (mm), OD – Вышний диаметр (mm), OD – Diametru exterior (mm), OD – 외부 직경 (mm), OD – isorinis skersmuo (mm), OD – یشورینیس کشرموو (mm), OD – Đường kính ngoài (mm), OD – Diameter Luar (mm)</p>
	<p>C – Tube Length (mm), C – Longueur de la canule (mm), C – Kanülenlänge (mm), C – Tubelengte (mm), C – Lunghezza cannula (mm), C – Longitud do tubo (mm), C – Tubens længde (mm), C – Tubens længde (mm), C – Tubelengde (mm), C – Putken pituus (mm), C – Comprimento do tubo (mm), C – Длина трубки (мм), C – 套管长度 (毫米), C – Długość rurki (mm), C – Délka trubice (mm), C – Džika kanyly (mm), C – Dožina cevke (mm), C – Tubus hossza (mm), C – Мјхкос σωλήνα (mm), C – Tűp Uzunluğu (mm), C – طول الأنبوب (مم)، C – Дължина на тръбата (mm), C – Lungime tub (mm), C – 관 길이 (mm), C – Vamzdello ilgis (mm), C – یشورینیس کشرموو (mm), C – Panjang tabung (mm)</p> <p>θ – Tube Angle (°), θ – Angle de la canule (°), θ – Krümmung der Kanüle (°), θ – Tubehoek (°), θ – Angulo cannula (°), θ – Angulo del tubo (°), θ – Tubwinkel (°), θ – Tubens vinkel (°), θ – Ubel trubice (°), θ – Uhol kanyly (°), θ – Kot cevke (°), θ – Tubus szöge (°), θ – Γωνία σωλήνα (°), θ – Tűp Açısı (°), θ – زاوية الأنبوب (°), θ – брыл на тръбата (°), θ – Ughi tub (°), θ – 관 각도 (°), θ – Vamzdello kampas (°), θ – یشورینیس کشرموو (mm), θ – Góc của ống (°), θ – Sudut tabung (°)</p>
	<p>Disposable Inner Cannula Length (mm), Longueur de la chemise interne jetable (mm), Länge Innenkanüle zum Einmalgebrauch (mm), Lengte wegwerpbare binnencanule (mm), Lunghezza contro cannula monouso (mm), Longitud de la cánula interna desechable (mm), Längd på innerkanyl för engångsbruk (mm), Længden af inderkanylen til engangsbrug (mm), Lengde på engangs innerkanyle (mm), Kertakäyttöisen sisäkanylin pituus (mm), Comprimento da cânula interna descartável (mm), Длина одноразовой внутренией канюли (mm), 一次性內套管长度 (毫米), Długość wewnątrznej kaniiuli jednorazowego użyciu (mm), Délka jednorázové vnitřní kanyly (mm), Džika jednorazovej vnútornej kanyly (mm), Dožina notranje kanile za enkratno uporabo (mm), Eldobható belső kanül hossza (mm), Μηκος εσωτερικής κάνουλας μίας χρήσης (mm), Tek Kullanımlık İç Kanül Uzunluğu (mm), طول الكانويلا الداخلية للتحفظ مرة واحدة (مم)، 1 회용 내부 캐놀라 길이 (mm), Vienkartines vidines kanulės ilgis (mm), Vienkartinis kanūlas Sekali Pakai (mm), Đò dài nòng trong dùng một lần (mm), Panjang Kanula Bagian Sekali Pakai (mm)</p>
	<p>Loading Dilator (fr); not included, Dilatateur de charge (fr); non fourni, Lade-dilatator (fr); nicht inbegriffen, Laaddilatator (fr), niet inbegrepen, Introduttore di cateteri a palloncino (fr); non incluso, introduttore a palloncino (fr); no incluido, Dilator (fr), ingår inte, Indføringsdilator (fr), medfølger ikke, Innføringsdilator (fr); ikke inkludert, Asetusdilattori (fr), ei kuulu mukaan, Dilatator de carregamento (fr); não incluído, Dilatator (fr); ne входит в комплект, 負載擴張器 (fr); 不含, Rozszerzadlo (fr), brak w zestawie, Zaváděcí dilatátor (fr); není součástí balení, Dilatátor zavádzania (fr); nie je súčasťou dodávky, Polinili dilatator (fr); ni priložen, Behelyezési dilatátor (fr); nem képezi a csomagolás részét, Διεστολέας τοποθέτησης (fr), δεν περιλαμβάνεται, Yükleme Dilatörü (fr); dahil değildir, غرض تحميل (fr), موسع التجميل، Dilatator za zapredanje (fr); ne e vključen, Dilatator de încărcare (fr); neinclus, 부하 확장기 (fr), 불포함, Plečiamasis įvediklis (fr); nepridedamas, 負載擴張器 (fr), Que nong dân dùng (fr); không bao gồm, Dilator Pemuatan (fr); tidak termasuk</p>

	<p>Obturator, Introducteur, Obturator, Otturatore, Obturador, Obturátor, Obturátor, Obturador, Obturaattori, Obturador, Obturátor, 插管芯, Obturátor, Obturátor, Tesnilo, Obturátor, Електроутрікó, Obtüratör, السداد, Обруратор, Obturator, 달개, [kişikis, умкһтһа, Nüt biç, Obturator</p>
	<p>Neck Strap, Attache cervicale, Halteband, Nekband, Fettuccia di fessaggio, Correa para el cuello, Halsrem, Halsrem, Halsstropp, Kiinnitysnauha, Faixa para pescoço, Лента для крепления на шею, 固定頸帶, Pasek szyjny, Kісні popruh, Popruh na krk, Trak za pritrđitev na vrat, Nyakszj, Лопіба ұлооттіріңіс ауһево, Boyun kayışı, حریط الرقیة, Лента за шия, Suport pentru gât, 脖带, Kaklo raištis, мунд-тэма, Dây đeo cổ, Pengikat Leher</p>

Shiley™

Adult Flexible Tracheostomy Tube Cuffless Disposable Inner Cannula

Identification of a substance that is contained or present within the product or packaging.

Identification of a substance that is not contained or present within the product or packaging.

REF	4UN65A	6.5 mm
REF	5UN70A	7.0 mm
REF	6UN75A	7.5 mm
REF	7UN80A	8.0 mm
REF	8UN85A	8.5 mm
REF	9UN90A	9.0 mm
REF	10UN10A	10.0 mm

Instructions for Use

POST FOR USE BY ALL TRACHEOSTOMY CARE PERSONNEL

The primary user of this instructional material is the caregiver in a clinical environment such as the adult or critical care units of a hospital, non-critical care units of a hospital, long term care facilities and home care.

Description

The Shiley™ Adult Flexible Tracheostomy Tube Cuffless (“tracheostomy tube”) is a dual cannula tracheostomy tube with a disposable inner cannula. It contains an integrated 15 mm connector for use with standard ventilation and anesthesia equipment.

The tracheostomy tube has a radiopaque, biocompatible outer cannula constructed of polyvinyl chloride. A pliable neck plate allows conformity to individual neck anatomies and contains two suture holes. The tapered tracheostomy tube tip is compatible with percutaneous insertion. The smooth, rounded tip obturator facilitates insertion. The disposable inner cannula is translucent for easy inspection and allows for airway maintenance.

This tracheostomy tube can be used with or without the disposable inner cannula included in the package (and also available separately; reference 10006780 for part numbers). However, use of the inner cannula is recommended.

Indications for Use

The Shiley™ Adult Flexible Tracheostomy Tube Cuffless, Disposable Inner Cannula is intended for use in providing tracheal access for airway management.

The Shiley™ Adult Flexible Tracheostomy Tube Cuffless, Disposable Inner Cannula is also intended for use with Cook® Percutaneous Dilatational Tracheotomy (PDT) procedures.

WARNINGS

- When using a laser or electrosurgical device in conjunction with this or other tracheal tubes, care must be taken to avoid contact with the tube. Such contact, especially in the presence of oxygen-enriched mixtures or nitrous oxide containing mixtures, could result in rapid combustion of the tube with harmful thermal effects and with emission of corrosive and toxic combustion products, including hydrochloric acid (HCl).
- The tracheostomy tube, neck strap, and disposable inner cannula cannot be adequately resterilized by the user in order to facilitate safe reuse and are therefore intended for single patient use. Attempts to resterilize this device may result in product failure and increased risks to the patient.
- Sterile only if the protective tray and cover are not open, damaged, or broken. Do not resterilize. Do not store product at temperatures above 120°F (49°C).

- Prior to use, handle with due care to prevent user contamination.
- The disposable inner cannula is designed for single use and should not be cleaned or reused.
- When using a disposable inner cannula, it can only be replaced by a Shiley disposable inner cannula of the same size.
- It is recommended that a spare disposable inner cannula be with the patient at all times.
- When using a water-soluble lubricant, verify that the lubricant does not enter and occlude the tube lumen, which may restrict ventilation.
- Patients in the home care environment should be carefully instructed by a home health care provider in the proper use and handling of the tracheostomy tube and accessory products.

CAUTIONS

- Federal law (USA) restricts this device to sale by or on the order of a physician.
- The tracheostomy tube and obturator are single patient use medical devices. Duration of single patient use should not exceed twenty-nine (29) days. Covidien does not recommend and has not substantiated use of these devices beyond the 29-day timeframe. Decisions about tracheostomy tube changes should be made by the responsible physician or designate using accepted medical techniques and judgment.
- Clinicians placing this device percutaneously should be competent in the procedure prior to using this tube.
- Prior to use, check that the obturator can be inserted and removed without difficulty.
- To ease placement of the tracheostomy tube, generously lubricate the surface of the loading dilator.
- Do not use solutions or chemical agents other than sterile saline, hydrogen peroxide, or water and mild detergent, as this may result in tube damage.
- Precautions should be taken when discarding the tracheostomy tube and accessories. Disposal should be made in accordance with current institutional standards and applicable national regulations for biologically hazardous waste.
- Duration of use should not exceed twenty-nine (29) days.

NOTES

- The tracheostomy tube and accessories are safe to use for X-ray, ultrasound, Positron Emission Tomography (PET), and Radiation Therapy.
- The tracheostomy tube and accessories are composed of medical grade plastics which do not contain and have not been made with natural rubber latex material or DEHP.
- Sterilized with ethylene oxide.
- The product can be thermally cleaned in water at a temperature of 65 degrees Celsius without the addition of chemical additives.

Contents of Package: 1 tracheostomy tube with pliable neck plate; 1 disposable inner cannula; 1 obturator; 1 neck strap. Reference page 2 for sizes and dimensions.

Accessories Available from Covidien: CAP, Disposable Inner Cannula (reference 10006780 for part numbers).

Accessories Available Separately: Some of the tube sizes (reference page 2) are compatible with the Percutaneous Dilatational Loading Dilator from the Cook® Critical Care Ciaglia Percutaneous Introducer Sets and Trays.

NOTE: The sizing dimensions listed on page 2 refer to the ID (the inside diameter of the tube or inner cannula at its narrowest point), the OD (the outside diameter of the tube), and the length (the distance from the neck plate to the distal tip on the tube's or inner cannula's centerline).

Tube Preparation and Insertion Using the Obturator or a Percutaneous Loading Dilator

The tube sizes that have a compatible loading dilator available (reference page 2) can be used in percutaneous dilatational tracheostomy (PDT) procedures.

NOTE: As this tracheostomy tube can be placed percutaneously, sterile technique should be followed for the handling and placement of the tracheostomy tube.

1. Insert the **obturator or compatible loading dilator** into the tracheostomy tube.
 - If using the **obturator**, ensure it is fully seated. If indicated to ease insertion, apply a thin film of water soluble lubricant to the outside of the tracheostomy tube and protruding portion of the obturator.
 - If using a **loading dilator for a PDT procedure**, generously lubricate the surface of the appropriately sized loading dilator and load the tracheostomy tube onto the dilator. Ensure that the tapered section of the loading dilator clears the distal tip of the tracheostomy tube by approximately 2 cm.

NOTE: The loading dilator is available separately.

2. If using a guide wire with the obturator, first insert the guide wire in the patient, then thread it all the way through the obturator, starting at the distal tip of the obturator.

WARNING

If the guide wire is damaged during the procedure, it may be difficult to continue insertion without trauma to the trachea and misplacement of the tracheostomy tube. In these cases, it may be possible to continue only if:

- a. the damaged section can be advanced into the trachea, and
- b. there is a sufficient length of undamaged guide wire remaining.

If the guide wire cannot be salvaged, use a new guide wire.

3. Insert the tube into the patient's trachea.
4. After the tube is properly in place, remove the **loading dilator or obturator** and guide wire (if used) immediately.
5. Insert the disposable inner cannula into the tracheostomy tube and press to lock into position.

Tube Removal

1. Remove the neck strap from around the patient's neck.
2. While holding on to the tracheostomy tube's 15 mm connector, gently twist and pull the breathing circuit connector to disconnect the circuit from the tracheostomy tube.
3. Slowly remove the tracheostomy tube from the patient.

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© 2013 Covidien.
COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG.
Other brands are trademarks of their respective owners.
Covidien Inc, 15 Hampshire Street, Mansfield, MA 02048 USA.
EC REP Covidien Ireland Limited, IDA Business & Technology Park, Tullamore.