Shiley™

Oral/Nasal RAE Tracheal Tube with TaperGuard™ Cuff Murphy Eye

en	Instructions For Use	REF 115-50OR	REF 115-85OR
fr	Mode d'emploi	REF 115-55OR	REF 115-90OR
it	Istruzioni per l'uso	REF 115-60OR	REF 119-60NR
nl	Gebruiksaanwijzing	REF 115-65OR	REF 119-65NR
da	Brugsanvisning	REF 115-70OR	REF 119-70NR
no	Bruksanvisning	REF 115-75OR	REF 119-75NR
es	Instrucciones de uso	REF 115-80OR	REF 119-80NR
sv	Bruksanvisning		
de	Gebrauchsanweisung		
fi	Käyttöohjeet		
pt	Instruções de uso		
ru	Инструкция по применению		
el	Οδηγίες χρήσης		
pl	Instrukcja użytkowania		
cs	Návod k použití		
hu	Használati útmutató		
tr	Kullanma Talimatı		

Shiley™

Oral/Nasal RAE Tracheal Tube with TaperGuard™ Cuff Murphy Eye

Identification of a substance that is not contained or present within the product or packaging.

Identification of a substance that is contained or present within the product or packaging.

Cuff Resting Diameter

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

This product contains DEHP. When used as indicated, very limited exposure to trace amounts of DEHP may occur. There is no clear clinical evidence that this degree of exposure increases clinical risk. However, in order to minimize risk of DEHP exposure in children and nursing or pregnant women, this product should only be used as directed.

DESCRIPTION

The Shiley RAE Tracheal Tube with TaperGuard Cuff and Murphy Eye is a sterile, single-use device. The clear tube, available in both oral and nasal styles, is supplied with a standard 15 mm connector. The tube design incorporates a preformed curve at the point where the tube emerges from the patient's mouth or nose, a self-sealing valve with attached pilot balloon and a low pressure cuff. The beveled tip is of hooded Murphy type (with eye). The RAE tube has a radiopaque line to assist in radiographic visualization. The tube material meets the requirements of the implantation test as set forth in the United States Pharmacopeia. Sterilization is by ethylene oxide.

INDICATIONS

The RAE Tracheal Tube with TaperGuard cuff is indicated for use in airway management, during those surgical procedures involving the head, neck, or mouth, where it would be to the surgeon's advantage to remove all connections from the operative field. The RAE Tracheal Tube with TaperGuard cuff functions in the same manner as a standard tracheal tube with the Nasal indicated for nasal intubation and the Oral indicated for oral intubation.

CONTRAINDICATIONS

Use of RAE Tracheal Tube with TaperGuard cuff in procedures which will involve the use of a LASER or electrosurgical active electrode in the immediate area of the device is contraindicated. Contact of the beam or electrode with the tracheal tube, especially in the presence of oxygen-enriched or nitrous oxide containing mixtures could result in the rapid combustion of the tube with harmful thermal effects and with emission of corrosive and toxic combustion products including hydrochloric acid (HCl). It has been reported by Hirshman and Smith that mixtures of nitrous oxide and oxygen support combustion about the same as pure oxygen and that in addition to ignition by direct contact with the beam, the interior of the tube can also be ignited by contact with flaming tissue in close proximity to the tip of the tracheal tube (Hirschman, C.A. and Smith, J.: Indirect Ignition of the Endotracheal Tube During Carbon Dioxide Laser Surgery. Arch Otolaryngol Vol. 106:639-641, 1980).

WARNINGS/PRECAUTIONS (Cuff-Related)

- If the duration of intubation is expected to be more than 24 hours or is not predictable, consider using the TaperGuard Oral/Nasal Tracheal Tube to reduce the potential for tissue trauma.
- As these devices may have been subjected to handling, storage conditions or preparation which compromised functional integrity, each tube's cuff, pilot balloon and valve should be tested by inflation prior to use. If dysfunction is detected in any part of the inflation system, the tube should not be used. Initiating treatment using a tube already shown to have a dysfunction in the inflation system could unnecessarily subject the patient to the untoward effects of extubation, re-intubation, or loss of respiratory support. Furthermore, the integrity of the inflation system should be monitored both initially and periodically during the intubation period. Uncorrected failure of the inflation system could result in death.

- The use of Lidocaine Topical Aerosol has been associated with the formation of pinholes in PVC cuffs (Jayasuriya, K.D., and Watson, W.F.: P.V.C. Cuffs and Lignocaine-based Aerosol. Brit.J. Ann. 53:1368, 1981). Expert clinical judgment must be used when prescribing treatment involving use of this substance to help prevent situations of cuff leaks due to pinholes. The same authors report that lidocaine hydrochloride solution does not have this effect.
- Various bony anatomical structures (e.g., teeth, turbinates) within the intubation routes or any intubation tools with sharp surfaces present a threat to maintaining cuff integrity. Care must be taken to avoid damaging the thin-walled cuffs during insertion which would create the need to subject the patient to the trauma of extubation and re-intubation. If cuff is damaged, the tube should not be used.
- Diffusion of nitrous oxide mixture, oxygen or air may either increase or decrease cuff volume and pressure. Inflating the cuff with the gas mixture which will contact its external surface is recommended as a means to reduce the extent of such diffusion.
- Inflation of the cuff by “feel” alone or by using a measured amount of air is not recommended since resistance is an unreliable guide during inflation. Intracuff pressure should be closely monitored with a pressure measuring device. The pilot balloon is only intended to indicate the presence of pressure or vacuum in the cuff and is not intended to provide an indication of pressure level.
- Do not overinflate cuff. Ordinarily, the cuff pressure should not exceed 25 cm H₂O. Carroll and Grenvik recommend maintaining a seal pressure at or below 25 cm H₂O (Carroll, R.G., and Grenvik, A: Proper use of large diameter, large residual cuffs. Critical Care Medicine Vol. 1, No. 3: 153-154, 1973). Overinflation can result in tracheal damage, rupture of the cuff with subsequent deflation, or in cuff distortion which may lead to airway blockage.
- Minimal Occluding Volume or Minimum Leak techniques should be used in conjunction with an intracuff pressure measuring device in selecting the sealing pressure. Cuff pressure should continue to be monitored thereafter, and any deviation from the selected seal pressure should be investigated and corrected immediately.
- Deflate cuff prior to repositioning the tube. Movement of the tube with cuff inflated could result in patient injury, requiring possible medical intervention or damage to the cuff, requiring a tube change. When complete evacuation of the air from the cuff is accomplished, a definite vacuum will be noted in the syringe and the tracheal tube pilot balloon is collapsed. Verify correct placement of the tube after each repositioning.
- Syringes, three-way stopcocks or other devices should not be left inserted in the inflation valve for extended periods of time. The resulting stress could crack the valve housing and allow the cuff to deflate.

WARNINGS/PRECAUTIONS (General)

- RAE Tracheal Tube with TaperGuard cuff are designed to fit the average patient, and accordingly, the distance from the distal tip to bend may be either too long or too short for a given patient. Expert clinical judgment should be exercised in selection of the appropriate size tracheal tube for each given patient with the awareness that in a patient whose airway is shorter than average, endobronchial intubation may result, or in a patient whose airway is longer than average, impingement of the cuff on the cords may result. The RAE Tracheal Tube with TaperGuard cuff design utilizes an increase in the tip-to-bend length as the I.D. size increases, to correspond to the observed correlation between patient size and length of the airway. If the tip-to-bend length appears inappropriate for a given patient, consideration should be given to use of a RAE tube with a smaller or larger I.D. to change the length correspondingly.
- It is essential to verify that the tube position remains correct after intubation, especially when a patient's position or the tube placement is altered. Any tube mal-position should be corrected immediately.
- Exposure to elevated temperatures and ultraviolet light should be avoided during storage.
- Should extreme flexing (chin-to-chest) of the head or movement of the patient (e.g., to a lateral or prone position) be anticipated after intubation, use of a reinforced tracheal tube should be considered.
- The 15 mm tracheal tube connector is loosely seated. **Always assure the connector is firmly seated in both the tracheal tube and the breathing circuit to prevent disconnection during use.**
- Non-standard dimensioning of some connectors on ventilator or anesthesia equipment may make secure mating with the tracheal tube 15 mm connector difficult. Use only with equipment having standard 15 mm connectors.
- Minimize the traction exerted by the breathing circuit on the 15 mm connector.
- Intubation and extubation should be performed following currently accepted medical techniques.
- If lubricating jellies are used in conjunction with the tracheal tube, follow manufacturer's application instructions. Excessive amounts of jelly can dry on the inner surface of the tracheal tube resulting in either a lubricant plug or a clear film that partially or totally blocks the airway.
- During an MRI scan the pilot balloon should be secured near the Y connector of the ventilator circuit at least 3 cm from the area of interest to prevent movement and image distortion.

Note: The performance information shown on page 71 was collected using a Rigid Cylinder Trachea Model bench test that is intended to provide a comparison of the sealing characteristics of tracheal tube cuffs only in a laboratory setting and is not configured or intended to predict performance in the clinical setting.

ADVERSE REACTIONS

The following adverse reactions have been reported to be associated with the use of cuffed tracheal tubes during the intubation procedure, during the intubation period, or subsequent to extubation. The order of listing is alphabetical and does not indicate frequency or severity.

Reported adverse reactions* include: abrasion of the arytenoid cartilage vocal process; cartilage necrosis; cicatrix formation; consequences of failure to ventilate including death; damage to the perichondrium; development of dense or diffuse fibrosis invading the entire glottic area; emphysema; endobronchial aspiration; endobronchial intubation (hypoxemia); endotracheobronchial aspiration; epistaxis; esophageal intubation (stomach distention); excoriated membranes of the pharynx; eye trauma; fibrin deposition; formation of subglottic web; fracture-luxation of cervical column (spinal injury); fragmentation of cartilage; glottic edema (supraglottic, subglottic, retroarytenoidal); granuloma of the inner arytenoid area; infections (laryngitis, sinusitis, abscess, respiratory tract infection); inflammation; intermittent aphonia and recurrent sore throat; laryngeal fibrosis; laryngeal granulomas and polyps; laryngeal obstruction; laryngeal stenosis; laryngeal ulcers; laryngotracheal membranes and webs; membranous glottic congestion; membranous tracheobronchitis; mild edema of the epiglottis; mucosal changes; mucosal sloughing; paresis of the hypoglossal and/or lingual nerves; perforation of esophagus; perforation of the trachea; pneumothorax; replacement of the tracheal wall with scar tissue; respiratory obstruction; retrobulbar hemorrhage; retropharyngeal abscess; retropharyngeal dissection, rupture of the trachea; sore throat, dysphagia; stricture of nostril; stridor; subglottic annular cicatricial stenosis; submucosal hemorrhage, submucous puncture of the larynx; superficial epithelial abrasion; swallowed tube; synechia of the vocal cords; teeth trauma; tissue burns; tracheal bleeding; tracheal stenosis; trauma to lips, tongue, pharynx, nose, trachea, glottis, palate, tonsil, etc.; traumatic lesions of the larynx and trachea; ulcerations exposing cartilaginous rings and minor erosions at cuff site; ulceration of the lips, mouth, pharynx, ulcers of the arytenoid; vocal cord congestion; vocal cord paralysis, and vocal cord ulcerations.

*Helpful references for more detailed discussions of tracheal tube adverse reactions include the following:

1. Baron, S.H., and Kohlmoos, H.W.: Laryngeal sequelae of endotracheal anesthesia. *Ann. Otol.* 60:767-792, 1951.
2. Blanc, V.F., and Tremblay, N.A.G.: The complications of trachea following intubation. *Anesth. Analg.* 53:202-213, 1974.
3. Campbell, D.: Trauma to larynx and trachea following intubation and tracheostomy. *J. Laryng.* 82-981-986, 1968.
4. Cooper, J.D., and Grillo, H.C.: The evolution of tracheal injury due to ventilatory assistance through cuffed tubes. *Ann. Surg.* 169:334-348, 1969.
5. Gross, C.W., and Gros, J.C.: Rare complications after prolonged translaryngotracheal intubation. *Ann. Otol.* 80:582-584, 1971.
6. McGovern, F.H., Fitz-Hugh, G.S., and Edgemon, L.J.: The hazards of endotracheal intubation. *Ann. Otol.* 80-556-564, 1971.

SUGGESTED DIRECTIONS FOR USE

1. Remove the sterile RAE Tracheal Tube with TaperGuard cuff from its protective package and perform the following inspections:
 - A. Test the cuff, pilot balloon and valve of each tube by inflation prior to use. Insert a luer tip syringe into the cuff inflation valve housing and inject enough air to fully inflate cuff. After test inflation of the cuff, completely evacuate the air.
 - B. Carefully inspect the tube and 15 mm connector prior to use. If the 15 mm connector is not seated at all or is loose in the package, do not use.
2. After orienting the 15 mm connector position, firmly seat it as completely as possible.
3. Verify the security of this connection before intubation. If the 15 mm connector can be pulled out easily, use a new tracheal tube and repeat the procedure.
4. Intubate the patient following currently accepted medical techniques with consideration given to the specific cuff-related WARNINGS and CAUTIONS stated in this product insert. Follow currently accepted medical techniques to ensure that the endotracheal tube has been placed into the trachea, and not inadvertently placed into the esophagus or a mainstem bronchus.

5. Once the patient is intubated, inflate cuff only with enough gas mixture to provide an effective seal at the desired lung inflation pressure. The use of Minimal Occluding Volume, Minimum Leak techniques and monitoring (measuring) of cuff pressure can help reduce occurrence of many of the adverse reactions associated with the use of cuffed tracheal tubes.
6. Remove syringe from the valve housing after cuff inflation. Leaving the syringe attached will keep the valve open, permitting the cuff to deflate.
7. Always assure the 15 mm connector is firmly seated in both the tracheal tube and the breathing circuit to prevent disconnection during use.
8. Check to verify inflation system is not leaking. Integrity of the system should be verified periodically during the intubation period. Uncorrected failure of the inflation system could result in death. Cuff pressure should be closely monitored and any deviation from the selected seal pressure should be investigated and corrected immediately.
9. The device should be secured in the patient's airway following currently accepted medical techniques.
10. Minimize the traction exerted by the breathing circuit on the 15 mm connector.
11. If patient position is altered while intubated, verify that the tracheal tube placement remains correct.
12. Verify the integrity of the circuit-connector-RAE Tracheal Tube connections throughout the case.
13. Prior to extubation, deflate cuff by inserting syringe into valve housing and removing gas mixture until a definite vacuum is noted in the syringe and the pilot balloon is collapsed.
14. Extubate patient following currently accepted medical techniques.
15. Discard tracheal tube.

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1. Baron, S.H., and Kohlmoos, H.W.: Laryngeal sequelae of endotracheal anesthesia. *Ann. Otol.* 60:767-792, 1951.
2. Blanc, V.F., and Tremblay, N.A.G.: The complications of trachea following intubation. *Anesth. Analg.* 53:202-213, 1974.
3. Campbell, D.: Trauma to larynx and trachea following intubation and tracheostomy. *J. Laryng.* 82-981-986, 1968.

4. Cooper, J.D., and Grillo, H.C.: The evolution of tracheal injury due to ventilatory assistance through cuffed tubes. *Ann. Surg.* 169:334-348, 1969.
5. Gross, C.W., and Gros, J.C.: Rare complications after prolonged translaryngotracheal intubation. *Ann. Otol.* 80:582-584, 1971.
6. McGovern, F.H., Fitz-Hugh, G.S., and Edgemon, L.J.: The hazards of endotracheal intubation. *Ann. Otol.* 80-556-564, 1971.

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1. Baron, S.H., and Kohlmoos, H.W.: Laryngeal sequelae of endotracheal anesthesia. *Ann. Otol.* 60:767-792, 1951.
2. Blanc, V.F., and Tremblay, N.A.G.: The complications of trachea following intubation. *Anesth. Analg.* 53:202-213, 1974.
3. Campbell, D.: Trauma to larynx and trachea following intubation and tracheostomy. *J. Laryng.* 82:981-986, 1968.
4. Cooper, J.D., and Grillo, H.C.: The evolution of tracheal injury due to ventilatory air tubossistance through cuffed tubes. *Ann. Surg.* 169:334-348, 1969.
5. Gross, C.W., and Gros, J.C.: Rare complications after prolonged translaryngotracheal intubation. *Ann. Otol.* 80:582-584, 1971.

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1. Baron, S.H., and Kohlmoos, H.W.: Laryngeal sequelae of endotracheal anesthesia. *Ann. Otol.* 60:767-792, 1951.
2. Blanc, V.F., and Tremblay, N.A.G.: The complications of trachea following intubation. *Anesth. Analg.* 53:202-213, 1974.
3. Campbell, D.: Trauma to larynx and trachea following intubation and tracheostomy. *J. Laryng.* 82-981-986, 1968.
4. Cooper, J.D., and Grillo, H.C.: The evolution of tracheal injury due to ventilatory assistance through cuffed tubes. *Ann. Surg.* 169:334-348, 1969.
5. Gross, C.W., and Gros, J.C.: Rare complications after prolonged translaryngotracheal intubation. *Ann. Otol.* 80:582-584, 1971.

6. McGovern, F.H., Fitz-Hugh, G.S., and Edgemon, L.J.: The hazards of endotracheal intubation. Ann. Otol. 80-556-564, 1971.

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1. Baron, S.H., and Kohlmoos, H.W.: Laryngeal sequelae of endotracheal anesthesia. Ann. Otol. 60:767-792, 1951.
2. Blanc, V.F., and Tremblay, N.A.G.: The complications of trachea following intubation. Anesth. Analg. 53:202-213, 1974.
3. Campbell, D.: Trauma to larynx and trachea following intubation and tracheostomy. J. Laryng. 82-981-986, 1968.
4. Cooper, J.D., and Grillo, H.C.: The evolution of tracheal injury due to ventilatory assistance through cuffed tubes. Ann. Surg. 169:334-348, 1969.
5. Gross, C.W., and Gros, J.C.: Rare complications after prolonged translaryngotracheal intubation. Ann. Otol. 80:582-584, 1971.
6. McGovern, F.H., Fitz-Hugh, G.S., and Edgemon, L.J.: The hazards of endotracheal intubation. Ann. Otol. 80-556-564, 1971.

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- Campbell, D.: Trauma to larynx and trachea following intubation and tracheostomy. J. Laryng. 82-981-986, 1968.
- Cooper, J.D., and Grillo, H.C.: The evolution of tracheal injury due to ventilatory assistance through cuffed tubes. Ann. Surg. 169:334-348, 1969.
- Gross, C.W., and Gros, J.C.: Rare complications after prolonged translaryngotracheal intubation. Ann. Otol. 80:582-584, 1971.
- McGovern, F.H., Fitz-Hugh, G.S., and Edgemon, L.J.: The hazards of endotracheal intubation. Ann. Otol. 80-556-564, 1971.

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1. Baron, S.H., and Kohlmoos, H.W.: Laryngeal sequelae of endotracheal anesthesia. *Ann. Otol.* 60:767-792, 1951.
2. Blanc, V.F., and Tremblay, N.A.G.: The complications of trachea following intubation. *Anesth. Analg.* 53:202-213, 1974.
3. Campbell, D. Trauma to larynx and trachea following intubation and tracheostomy. *J. Laryng.* 82-981-986, 1968.
4. Cooper, J.D., and Grillo, H.C.: The evolution of tracheal injury due to ventilatory assistance through cuffed tubes. *Ann. Surg.* 169:334-348, 1969.
5. Gross, C.W., and Gros, J.C.: Rare complications after prolonged translaryngotracheal intubation. *Ann. Otol.* 80:582-584, 1971.
6. McGovern, F.H., Fitz-Hugh, G.S., and Edgemon, L.J.: The hazards of endotracheal intubation. *Ann. Otol.* 80-556-564, 1971.

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(A) Tracheal Tube Cuff Performance by Tracheal Tube Size, Cuff Pressure 27 hPa (cmH₂O) [Per ISO 5361 method]						
(D) Tracheal Tube Size (mm)	(B) Minimum Trachea Model Diameter (mm)			(C) Maximum Trachea Model Diameter (mm)		
	(E) Trachea Diameter (mm)	(F) Leak Rate Range (ml/h)		(E) Trachea Diameter (mm)	(F) Leak Rate Range (ml/h)	
		(G) 50th Percentile	(H) 90th Percentile		(G) 50th Percentile	(H) 90th Percentile
5.0	15.0	0.0	1.5	17.0	0.0	1.1
5.5	15.0	0.0	1.5	17.0	0.0	1.0
6.0	16.0	0.0	1.2	19.0	0.0	0.9
6.5	16.0	0.0	1.4	19.0	0.0	0.7
7.0	19.0	0.0	1.5	25.0	0.0	0.6
7.5	19.0	0.5	1.9	25.0	0.0	0.3
8.0	19.0	0.0	0.7	25.0	0.0	0.0
8.5	21.0	0.0	0.5	26.0	0.0	0.0
9.0	21.0	0.0	0.7	26.0	0.0	0.0

A	<p>Tracheal Tube Cuff Performance by Tracheal Tube Size, Cuff Pressure 27 hPa (cmH₂O) [Per ISO 5361 method], Performances du ballonnet de tube endotrachéal par taille de tube, avec une pression de ballonnet de 27 hPa (cmH₂O) [conformément à la méthode ISO 5361], Prestazioni della cuffia del tubo endotracheale in base alla dimensione del tubo endotracheale, pressione della cuffia 27 hPa (cmH₂O) [Per metodo ISO 5361], Prestatie van de tracheacanulecuff per tracheacanule-afmeting, cuffdruk 27 hPa (cmH₂O) [met de ISO 5361-methode], Trakealtubemanchetteres ydeevne efter trakealtubestørrelse, Manchetrykk 27 hPa (cmH₂O) [Per ISO 5361-metode], Ytelsen til trakealtubemansjetter etter trakealtubestørrelse, manjettrykk 27 hPa (cmH₂O) [Etter ISO 5361-metode], Rendimiento del manguito del tubo traqueal por el tamaño del tubo traqueal, presión del manguito de 27 hPa (cmH₂O) [Para el método ISO 5361], Trakealtubskuffens egenskaper efter trakealtubsstorlek, kufftryck 27 hPa (cmH₂O) [per ISO 5361-metod], Leistung eines Trachealtubus mit Manschette gemäß Größe des Trachealtubus, Manschettendruck 27 hPa (cmH₂O) [gemäß ISO 5361-Methode], henkitorviputken kuffin suorituskyky henkitorviputken koon mukaan, kuffipaine 27 hPa (cmH₂O) [ISO 5361 -menetelmän mukaisesti], Desempenho do manguito do tubo traqueal, de acordo com o Tamanho do tubo traqueal, Pressão do manguito 27 hPa (cmH₂O) [método Per ISO 5361], Рабочие характеристики манжеты трахеальной трубки, определяемые размером трубки, при давлении в манжете 27 гПа (см H₂O) [по методу ISO 5361], Απόδοση αεροθαλάμου τραχειοσωλήνα βάσει μεγέθους τραχειοσωλήνα, πίεσης αεροθαλάμου 27 hPa (cmH₂O) [Σύμφωνα με τη μέθοδο ISO 5361], Sprawność mankietów rurek intubacyjnych wg rozmiarów rurek intubacyjnych, przy ciśnieniu w mankiecie wynoszącym 27 hPa (cm H₂O) [zgodnie z metodą normy ISO 5361], Výkon manžety tracheální kanyly podle velikosti tracheální kanyly, tlak manžety 27 hPa (cm H₂O) [podle ISO 5361], Tracheális tubus mandzsetta teljesítmény tracheális tubusméret szerint, mandzsettanyomás 27 hPa (cmH₂O) [ISO 5361 módszer szerint], Trakeal Tüp Büyüklüğüne Göre Trakeal Tüp Kafının Performansı, Kaf Basıncı 27 hPa (cm H₂O) [ISO 5361 yöntemi uyarınca]</p>
B	<p>Minimum Trachea Model Diameter (mm), Diamètre minimal du modèle de trachée (mm), Diametro minimo del modello per trachea (mm), Minimum diameter tracheal model (mm), Minimum diameter for tracheamodel (mm), Minimum trakemodell diameter (mm), Diámetro mínimo del modelo de tráquea (mm), Minsta trakeamodellsdiameter (mm), Mindstdurchmesser des Tracheamodells (mm), Pienin trakeamallin läpimitta (mm), Diámetro mínimo do modelo da traqueia (mm), Минимальный диаметр модели трахеи (мм), Ελάχιστη διάμετρος μοντέλου τραχείας (mm), Minimalna średnica modelu tchawicy (mm), Minimální průměr modelu trachey (mm), Minimális trachea modell átmérő (mm), Minimum Trakea Model Çapı (mm)</p>
C	<p>Maximum Trachea Model Diameter (mm), Diamètre maximal du modèle de trachée (mm), Diametro massimo del modello per trachea (mm), Maximum diameter tracheal model (mm), Maksimum diameter for tracheamodel (mm), Maksimum trakeamodell diameter (mm), Diámetro máximo del modelo de tráquea (mm), Största trakeamodellsdiameter (mm), Höchstdurchmesser des Tracheamodells (mm), Suurin trakeamallin läpimitta (mm), Diámetro máximo do modelo da traqueia (mm), Максимальный диаметр модели трахеи (мм), Μέγιστη διάμετρος μοντέλου τραχείας (mm), Maksymalna średnica modelu tchawicy (mm), Maximální průměr modelu trachey (mm), Maximális trachea modell átmérő (mm), Maksimum Trakea Model Çapı (mm)</p>
D	<p>Tracheal Tube Size (mm), Taille du tube endotrachéal (mm), Dimensione del tubo endotracheale (mm), Afmeting tracheacanule (mm), Trakealtubestørrelse (mm), Trakealtubestørrelse (mm), Tamaño del tubo traqueal (mm), Trakealtubsstorlek (mm), Größe des Trachealtubus (mm), Henkitorviputken koko (mm), Tamanho do tubo traqueal (mm), Размер трахеальной трубки (мм), Μέγεθος τραχειοσωλήνα (mm), Rozmiar rurki intubacyjnej (mm), Velikost tracheální kanyly (mm), Tracheális tubus mérete (mm), Trakeal Tüp Büyüklüğü (mm)</p>
E	<p>Trachea Diameter (mm), Diamètre de la trachée (mm), Diametro trachea (mm), Tracheadiameter (mm), Tracheadiameter (mm), Trakeadiameter (mm), Diámetro de la tráquea (mm), Trakeadiameter (mm), Tracheadurchmesser (mm), Trakean läpimitta (mm), Diámetro da traqueia (mm), Диаметр трахеи (мм), Διάμετρος τραχείας (mm), Średnica tchawicy (mm), Průměr trachey (mm), Trachea átmérő (mm), Trakea Çapı (mm)</p>

F	Leak Rate Range (ml/h), Plage de taux de fuite (ml/h), Range di frequenza delle perdite (ml/h), Lekratiobereik (ml/u), Lækhastighed (ml/t), Lekkasjeområde (ml/t), Amplitud de fugas (ml/h), Läckageintervall (ml/h), Leckagewert (ml/h), Vuotonopeusalue (ml/h), Taxa de vazamento (ml/h), Диапазон скоростей утечки (мл/ч), Εύρος ρυθμού διαρροής (ml/h), Zakres przecieku (ml/h), Rozsah míry úniku (ml/h), Szivárgási sebesség tartomány (ml/óra), Sızıntı Hız Aralığı (ml/s)
G	50th Percentile, 50e centile, 50° percentile, 50e percentiel, 50-percentilen, 50. persentil, Percentil 50, 50:e percentilen, 50 %-Perzentil, 50. prosenttipiste, 50° Percentil, 50-й процентиль, 50ο εκατοστημόριο, 50. percentyl, 50. percentil, 50-edik percentilis, 50. Yüzdelik Kısım
H	90th Percentile, 90e centile, 90° percentile, 90e percentiel, 90-percentilen, 90. persentil, Percentil 90, 90:e percentilen, 90 %-Perzentil, 90. prosenttipiste, 90° Percentil, 90-й процентиль, 90ο εκατοστημόριο, 90. percentyl, 90. percentil, 90-edik percentilis, 90. Yüzdelik Kısım

Page Left Intentionally Blank

Page Left Intentionally Blank

STERILE EO

Single use

Do not use if
package is opened
or damaged

Do not expose to UV

Not made
with natural
rubber latex

DEHP

Caution, consult
accompanying
documents

0123

Part No. PT00045062 Rev A 2016-05

COVIDIEN, COVIDIEN with logo and
Covidien logo are U.S. and internationally
registered trademarks of Covidien AG.
Other brands are trademarks of a
Covidien company.

© 2012 Covidien.

 Covidien llc, 15 Hampshire Street,
Mansfield, MA 02048 USA.

 Covidien Ireland Limited,
IDA Business & Technology Park, Tullamore.

www.covidien.com