
Puritan Bennett™

Re/X800 and D/X800 Expiratory Filter/ Drain System Instructions 800 Series Ventilators

EN US	Instructions For Use
FR	Mode d'emploi
DE	Gebrauchsanweisung
NL	Gebruiksaanwijzing
IT	Istruzioni per l'uso
ES	Instrucciones de uso
SV	Bruksanvisning
DA	Brugsanvisning
NO	Bruksanvisning
FI	Käyttöohjeet
PT	Instruções de uso
RU	Инструкция по применению
EL	Οδηγίες χρήσης
PL	Instrukcja użytkowania
CS	Návod k použití
SK	Návod na použitie
HU	Használati útmutató
TR	Kullanma Talimatı

Re/X800 and D/X800

Expiratory Filter/Drain System Instructions 800 Series Ventilators

Identification of a substance that is not contained or present within the product or packaging.

Identification of a substance that is contained or present within the product or packaging.

 "Use by" mark; anticipated filter retirement date.

 Indicates flow direction; filter is designed for use only in this orientation.

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

This product contains DEHP. When used as indicated, very limited exposure to trace amounts of DEHP may occur. There is no clear clinical evidence that this degree of exposure increases clinical risk. However, in order to minimize risk of DEHP exposure in children and nursing or pregnant women, this product should only be used as directed.

Directions for Use

Description

The expiratory filter reduces the particles and bacteria in the patient's exhaled gas, and protects the ventilator's exhalation and spirometry components. Part of a comprehensive condensate drainage system that includes a collector vial and optional drain bag, the expiratory filter is available in reusable (*Re/X800*) filter and single-patient use (*D/X800*) filter versions.

Installation (Illustrations on page 2)

1. *Reusable filter and vial only:* Filters are shipped nonsterile. Steam-sterilize before the first use, and after each patient use. (Do not attempt to sterilize single-patient use filters.)
2. *Collector vial to filter:* Insert the collector vial (A) into the bottom of the expiratory filter (B) as shown. Be sure the lock ring is in the open position, and be sure to line up the expiratory limb connector (C) and the collector vial drain port (D) by turning; then turn the lock ring (E) to secure the collector vial into place. To unlock, turn the lock ring counterclockwise.

3. *Filter and collector vial to ventilator:* With the exhalation filter latch (F) in up position, slide filter rim onto the tracks (G) in the housing area (H) with drain port (D) facing toward you. Push the exhalation latch down: it will position the filter properly. Attach the expiratory limb of the ventilator breathing circuit to the filter's conical connector (C).
4. *If you are not using the drain bag:* The collector vial drain port (D) must be capped (J). If you are using the drain bag: Install clamp (K) on tubing (L). Uncap collector vial drain port and install tubing to collector vial drain port. Connect other end of tubing to drain bag (M). If the ventilator is mounted on the cart, place the drain bag in the ventilator drawer.

Note: The drain bag is designed to lie flat, and is not designed to be suspended.

WARNING

1. **To ensure that all ventilator breathing circuit connections are leak-tight, perform a circuit leak test every time you install the filter on a ventilator.**

Use

When the collector vial fills: Empty the collector vial before liquid reaches the maximum fill line (N). If you are using a drain bag, squeeze the clamp to drain liquid from the collector vial to the drain bag. Close clamp as soon as liquid level drops below drain port. Do not allow air to fill the drain bag. When the drain bag is full, disconnect and discard.

When the drain bag is full: Disconnect drain bag from tubing. Install bag fitting (P) onto tab to seal the bag. Discard bag.

WARNING

1. **Empty the collector vial before fluid reaches the maximum fill line. Collector vial overflow can allow fluid to enter the filter or the ventilator breathing circuit, and can cause increased flow resistance.**
2. **The use of nebulized medication can cause a buildup of exhalation flow resistance and may even block the filter. Inspect and test expiratory filters at patient setup and frequently during use.**
3. **Removing the collector vial while the patient is connected to the ventilator can cause loss of circuit pressure, ventilator autocycling, or direct contact with liquid.**
4. **Always follow your institution's infection control guidelines.**

Cleaning

Tubing, drain bag, clamp: If necessary, wipe down with a damp cloth.

Filter housing, collector vial: If necessary, wipe down. Use only cleaning solutions that are compatible with polysulfone plastic (Re/X800) filter. Do not attempt to clean filter medium.

Sterilization

CAUTION

1. Do not attempt to sterilize or reuse single-patient use items.

Reusable filter, reusable collector vial: Separate the filter from collector vial and wrap each in muslin or equivalent paper before autoclaving. Effective sterilization of the filter occurs by steam autoclaving at 132 °C (270 °F) for 20 minutes for gravity displacement cycles, or 132 °C (270 °F) for 4 minutes for prevacuum cycles (requires minimum 20 minute drying cycle after sterilization). Follow the steam sterilizer manufacturer's instructions.

Clamp: Heat-pasteurize according to manufacturer's instructions. Do not autoclave.

WARNING

1. **To avoid increased flow resistance, do not immerse filter in liquid. To avoid reducing filtration efficiency, do not attempt to scrub or touch the filter medium. To avoid retaining harmful residues, do not immerse in liquid or sterilize using ETO.**

Inspection

Filter: Inspect for any visible damage to the housing or filter medium before each use. Check filter flow resistance before every reuse. Discontinue use if the flow resistance is more than 4 cmH₂O at 100 L/min.

Collector vial, tubing, drain bag, clamp: Inspect for any visible damage.

Replacement

Single-patient use (D/X800) filter, single-patient use collector vial, tubing, drain bag: These are single-patient use items. Replace at each circuit change. Do not attempt to sterilize or reuse.

Reusable (Re/X800) filter: The filter is designed to withstand up to 100 autoclave sterilization cycles at temperatures not to exceed 135 °C, 275 °F. Replace the filter after one year of service or 100 autoclave cycles, whichever occurs first. For your convenience a space is provided on the filter to note the anticipated retirement date. Calculate an expected filter retirement date based on your typical use rates.

Replace filter seal if cracked or damaged.

Reusable collector vial: Replace if collector vial is cracked or damaged. Replace drain cap if cracked or damaged.

Clamp: Replace if visibly damaged.

Note: Dispose of all items according to your institution's protocol. Sterilize before nondestructive disposal. Follow local governing ordinances and recycling plans regarding disposal or recycling of device components.

Specifications

Filter

- Viral/bacterial filtration efficiency: Greater than 99.999%.
- Particle filtration efficiency: Minimum of 99.97% of 0.3- μ m nominal particle size at 100 L/min flow.
- Inlet connector: 22-mm ISO connector
- Resistance to flow: Less than 2.5 cmH₂O at 100 L/min when new.
- Flow leakage: Less than 0.01 L/min at 140 cmH₂O internal pressure.
- Size: 11.7 cm diameter x 12.7 cm long (4.6 in. diameter x 5 in. long).
- Internal volume: Approximately 400 mL.

Collector vial

- Capacity: Approximately 200 mL usable volume (to maximum fill line).

Tubing

- Length: 17.8 cm (7 in.)
- Internal diameter: 0.48 cm (3/16 in.)

Drain bag

- Capacity: Approximately 1 L.

Ordering information

4-070305-00	Expiratory filter, Re/X800, reusable
4-074647-00	Collector vial, reusable
4-070315-00	Expiratory filter and collector vial, D/X800, single-patient use (carton of 12)
4-048491-00	Drain bag, single-patient use (package of 25)
4-048493-00	Tubing, drain bag, single-patient use (package of 10)
4-048492-00	Clamp, reusable (package of 5)
4-074613-00	Drain cap
4-070311-00	Seal, expiratory filter

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© 2011 Covidien.
COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.
Covidien Inc, 15 Hampshire Street,
Covidien Ireland Limited,