

Clinical Evidence Guide

MONITORING CEREBRAL PERFUSION IN CAROTID ENDARTERECTOMY

With the INVOS™ cerebral/somatic oximeter



The common methods of monitoring bilateral cerebral perfusion adequacy and identifying the need for shunting during carotid endarterectomy (CEA) under general anesthesia are electroencephalogram (EEG), stump pressure (SP), or transcranial Doppler (TCD) monitoring.

This clinical evidence guide explores how INVOS™ monitoring technology can be a practical and continuous method to monitor perfusion and identify patients who may benefit from shunting, either alone or in combination with other monitoring modalities.

Medtronic
Further, Together

INVOS™ monitoring technology can be used to identify patients at risk for cerebral ischemia during CEA.

Near infrared spectroscopy monitoring during carotid endarterectomy: which threshold value is critical?

Mille T, Tachimiri ME, Klersy C, et al.

Study design	Single-center, retrospective study (Italy)
Arms	Single arm: INVOS™ monitoring technology
Objective	Identify the relative change from baseline rSO ₂ values 2 minutes after clamping that is predictive of impending cerebral ischemia
N	594
Population	CEA patients under general anesthesia
Threshold for intervention	No interventions were described and no shunts were placed
Results	<ul style="list-style-type: none">▪ 20/594 patients (3.4%) presented with neurological complications▪ 16/594 patients (2.7%) showed an early decrease in rSO₂ >20% from baseline<ul style="list-style-type: none">– 6/16 patients (37.5%) showing an early decrease in rSO₂ >20% from baseline also presented with neurological complications– 14/578 of the remaining patients (2%) showing a decrease in rSO₂ <20% from baseline presented with neurological complications (p = 0.0001)▪ An rSO₂ decrease of 11.7% from baseline at 2 minutes after clamping was identified as the optimal threshold to identify patients with a neurological complication:<ul style="list-style-type: none">– Positive predictive value (PPV) = 10%– Negative predictive value (NPV) = 99%– Sensitivity = 75%– Specificity = 77%▪ A threshold of 20% was associated with lower sensitivity (30%), higher specificity (98%), with a PPV and NPV of 37% and 98%, respectively, for identifying patients with neurological complications
Conclusions	Near-infrared spectroscopy (NIRS) monitoring technology is a useful method to detect cerebral ischemia during CEA.

INVOS™ technology may be an effective monitoring method compared to EEG to inform shunt selection during CEA.

Near-infrared spectroscopy to indicate selective shunt use during carotid endarterectomy

Pennekamp CWA, Immink RV, den Ruijter HM, et al.

Study design	Single-center, prospective observational study (The Netherlands)
Arms	Single arm: concomitant monitoring with INVOS™ technology, EEG, and TCD
Objective	Compare the ability of INVOS™ monitoring technology and TCD to EEG in detecting cerebral ischemia and determining the need for shunting
N	151
Population	CEA patients under general anesthesia with a sufficient TCD window
Definition of desaturation	rSO ₂ ≥16% from baseline value
Threshold for intervention	Shunt was placed when EEG detected the occurrence of new delta or theta activity
Results	<ul style="list-style-type: none"> ▪ 17/151 (11%) of patients were shunted based on EEG changes <ul style="list-style-type: none"> – In 16/17 cases, rSO₂ also dropped to ≥16% from baseline – In 5/134 of the remaining patients not shunted, rSO₂ dropped to ≥16% from baseline value without EEG changes – 129 patients showed no change in EEG or rSO₂ drop of ≥16% from baseline value ▪ Compared to EEG, rSO₂ demonstrated the following regarding indications for shunting: <ul style="list-style-type: none"> – PPV = 76% – NPV = 99% – Sensitivity = 94% – Specificity = 96% ▪ Similar results were seen with TCD, although the PPV and specificity were lower at 53% and 90%, respectively ▪ A trend toward higher incidence of stroke in shunted patients was seen (11.8% vs. 1.5%; p = 0.06)
Conclusions	NIRS monitoring technology may be effective for shunt selection. The optimal threshold for shunt selection requires further study.

Compared to transcranial Doppler, NIRS monitoring proved safe and reliable during CEA and reduced unnecessary shunt use.

Near-infrared spectroscopy versus transcranial Doppler-based monitoring in carotid endarterectomy

Cho JW, Jang JS.

Study design	Retrospective review of medical records of patients who underwent carotid endarterectomy (CEA) at a single center (Korea)																										
Arms	TCD group: transcranial Doppler (TCD) used for cerebral monitoring during CEA NIRS group: near-infrared spectroscopy (NIRS) used for cerebral monitoring during CEA (INVOS™ 5100C system)																										
Objective	To compare the safety and reliability of NIRS to that of TCD for cerebral monitoring during CEA																										
N	74 (45 TCD, 29 NIRS)																										
Population	Patients who underwent CEA using selective shunt-based TCD or NIRS																										
Threshold for intervention	TCD group: selective shunt inserted when middle cerebral artery blood flow decreased to <30% of baseline, or in patients with poor temporal view NIRS group: shunt inserted when regional cerebral oxygen saturation decreased to <80% of baseline																										
Results	<table border="1"> <thead> <tr> <th>Outcome</th> <th>TCD group</th> <th>NIRS group</th> <th>p value</th> </tr> </thead> <tbody> <tr> <td>30-day mortality, n (%)</td> <td>1 (2.2)</td> <td>0</td> <td>1</td> </tr> <tr> <td>30-day stroke, n (%)</td> <td>4 (8.9)</td> <td>0</td> <td>0.15</td> </tr> <tr> <td>Hypoglossal nerve injury, n (%)</td> <td>1 (2.2)</td> <td>0</td> <td>1</td> </tr> <tr> <td>Other neurologic symptoms, n (%)</td> <td>10 (22.2)</td> <td>0</td> <td>0.005</td> </tr> <tr> <td>Shunt used, %</td> <td>44.4</td> <td>10.3</td> <td>0.002</td> </tr> </tbody> </table> <p>Shunt use had no significant effect on 30-day mortality or 30-day stroke but researchers reported an effect of shunt use on other neurologic symptoms ($p = 0.079$)</p>			Outcome	TCD group	NIRS group	p value	30-day mortality, n (%)	1 (2.2)	0	1	30-day stroke, n (%)	4 (8.9)	0	0.15	Hypoglossal nerve injury, n (%)	1 (2.2)	0	1	Other neurologic symptoms, n (%)	10 (22.2)	0	0.005	Shunt used, %	44.4	10.3	0.002
Outcome	TCD group	NIRS group	p value																								
30-day mortality, n (%)	1 (2.2)	0	1																								
30-day stroke, n (%)	4 (8.9)	0	0.15																								
Hypoglossal nerve injury, n (%)	1 (2.2)	0	1																								
Other neurologic symptoms, n (%)	10 (22.2)	0	0.005																								
Shunt used, %	44.4	10.3	0.002																								
Conclusions	<ul style="list-style-type: none"> Regarding postoperative stroke and neurologic symptoms, NIRS appears to be safe and reliable for monitoring cerebral perfusion during CEA. NIRS monitoring reduces unnecessary shunt usage compared to TCD monitoring. 																										

Preoperative cerebral oxygen saturation values and changes in oxygen saturation during carotid clamping predict the risk of postoperative neurologic complications after CEA.

Brain oxygen supply parameters in the risk assessment of cerebral complications during carotid endarterectomy

Kamenskaya OV, Yu I, Lomivorotov VV.

Study design	Prospective, nonrandomized cohort study
Arms	One arm: NIRS monitoring (INVOS™ 5100 system)
Objective	To determine if the risk of neurologic complications from carotid endarterectomy (CEA) can be predicted by the preoperative cerebral oxygen saturation (rSO ₂) and the decrease in rSO ₂ during carotid clamping
N	466
Population	Adult patients who underwent CEA with bilateral monitoring of rSO ₂
Threshold for intervention	No intervention
Results	<ul style="list-style-type: none"> ▪ Patients with neurological complications had: <ul style="list-style-type: none"> – Lower preoperative rSO₂ on the surgical side (55% vs. 62%, p = 0.018) – Greater ΔrSO₂ during clamping on the surgical side (11% vs. 22%, p = 0.029) ▪ Multivariate regression analysis showed the following: <ul style="list-style-type: none"> – Preoperative rSO₂ and ΔrSO₂ are independent predictors of stroke – ΔrSO₂ is an independent predictor of cognitive disorders ▪ Preoperative rSO₂ and ΔrSO₂ cutoffs were significantly correlated with outcomes: <ul style="list-style-type: none"> – Preoperative rSO₂ <50% increased the risk of stroke by 6-fold (OR 6.0; 95% CI 1.30 – 27.68) – ΔrSO₂ >20% during clamping increased: <ul style="list-style-type: none"> ▪ Risk of stroke by 10-fold (OR 10.8; 95% CI 2.05 – 56.77) ▪ Risk of cognitive disorders by 8-fold (OR 8.7; 95% CI 2.55 – 29.59)
Conclusions	<ul style="list-style-type: none"> ▪ During temporary carotid artery clamping, ΔrSO₂ >20% below baseline significantly increases the risks of ischemic stroke (10-fold) and cognitive disorders (8-fold) after CEA. ▪ Preoperative rSO₂ <50% significantly increased the risk of ischemic stroke (6-fold) after CEA.

INVOS™ monitoring technology paired with an interventional algorithm can inform the need for shunt placement during CEA.

Evaluation of an intraoperative algorithm based on near-infrared refracted spectroscopy monitoring, in the intraoperative decision for shunt placement, in patients undergoing carotid endarterectomy

Zogogiannis ID, Iatrou CA, Lazarides MK, et al.

Study design	Multicenter, prospective, randomized controlled trial (Greece)
Arms	GROUP A: INVOS™ monitoring technology and use of an interventional protocol to restore rSO ₂ GROUP B: INVOS™ monitoring technology without the use of an interventional protocol GROUP C: Control, no INVOS™ monitoring technology used
Objective	Examine the effect of using NIRS monitoring with an interventional protocol on the decision to place a shunt
N	253
Population	CEA patients under general anesthesia
Threshold for intervention	rSO ₂ <80% of baseline value prompted use of interventional protocol, including consideration to place a shunt
Results	<ul style="list-style-type: none"> ▪ The incidences of shunt placement were significantly different between groups (p < 0.001): <ul style="list-style-type: none"> – Group A = 27.7% – Group B = 59.5% – Group C = 100% ▪ Compared to group A, patients in groups B and C had a 3.7 and 70.6 times greater likelihood of receiving a shunt, respectively ▪ Compared to group B, group C had a 19.4 times greater likelihood of receiving a shunt ▪ There were no significant differences in the incidences of postoperative neurologic deficits among the three groups
Conclusions	NIRS monitoring technology paired with the use of a specific treatment algorithm may aid in the selective shunting of patients undergoing CEA.

A selection of clinical studies using INVOS™ technology intraoperatively in carotid endarterectomy surgery

- Cho JW, Jang JS. Near-infrared spectroscopy versus transcranial Doppler-based monitoring in carotid endarterectomy. *Korean J Thorac Cardiovasc Surg.* 2017;50:448-452.
- Kamenskaya OV, Loginova IY, Lomivorotov VV. Brain oxygen supply parameters in the risk assessment of cerebral complications during carotid endarterectomy. *J Cardiothorac Vasc Anesth.* 2017;31(3):944-949.
- Öz K, Güler S, Kahraman Z, Özarslan M, Ertürk M. Evaluation of shunt use with near-infrared spectroscopy during carotid endarterectomy: an observational clinical study. *Turkiye Klinikleri J Cardiovasc Sci.* 2017;29(2):48-54.
- Pennekamp CWA, Immink RV, den Ruijter HM, et al. Near-infrared spectroscopy to indicate selective shunt use during carotid endarterectomy. *Eur J Vasc Endovasc Surg.* 2013;46(4):397-403.
- Uchino H, Nakamura T, Kuroda S, Houkin K, Murata J, Saito H. Intraoperative dual monitoring during carotid endarterectomy using motor evoked potentials and near-infrared spectroscopy. *World Neurosurg.* 2012;78(6):651-657.
- Picton P, Ramachandran SK, Shanks A, Dorje P, Mashour GA. Optimizing cerebral oxygenation in anaesthetized patients with carotid artery stenosis: the influence of inspired oxygen fraction. *Br J Anaesth.* 2012;108(2):326-327.
- Pennekamp CW, Immink RV, den Ruijter HM, et al. Near-infrared spectroscopy can predict the onset of cerebral hyperperfusion syndrome after carotid endarterectomy. *Cerebrovasc Dis.* 2012;34(4):314-321.
- Zogogiannis ID, Iatrou CA, Lazarides MK, et al. Evaluation of an intraoperative algorithm based on near-infrared refracted spectroscopy monitoring, in the intraoperative decision for shunt placement, in patients undergoing carotid endarterectomy. *Middle East J Anaesthesiol.* 2011;21(3):367-373.
- Tambakis CL, Papadopoulos G, Sergentanis TN, et al. Cerebral oximetry and stump pressure as indicators for shunting during carotid endarterectomy: comparative evaluation. *Vascular.* 2011;19(4):187-194.
- Ritter JC, Green D, Slim H, Tiwari A, Brown J, Rashid H. The role of cerebral oximetry in combination with awake testing in patients undergoing carotid endarterectomy under local anaesthesia. *Eur J Vasc Endovasc Surg.* 2011;41(5):599-605.
- Ali AM, Green D, Zayed H, Halawa M, El-Sakka K, Rashid HI. Cerebral monitoring in patients undergoing carotid endarterectomy using a triple assessment technique. *Interact Cardiovasc Thorac Surg.* 2011;12(3):454-457.
- Moritz S, Schmidt C, Bucher M, et al. Neuromonitoring in carotid surgery: are the results obtained in awake patients transferable to patients under sevoflurane/fentanyl anesthesia? *J Neurosurg Anesthesiol.* 2010;22(4):288-295.
- Giustiniano E, Alfano A, Battistini GM, Gavazzeni V, Spoto MR, Cancellieri F. Cerebral oximetry during carotid clamping: is blood pressure raising necessary? *J Cardiovasc Med (Hagerstown).* 2010;11(7):522-528.
- Baikoussis NG, Karanikolas M, Siminelakis S, Matsagas M, Papadopoulos G. Baseline cerebral oximetry values in cardiac and vascular surgery patients: a prospective observational study. *J Cardiothorac Surg.* 2010;5:41.
- Picton P, Chambers J, Shanks A, Dorje P. The influence of inspired oxygen fraction and end-tidal carbon dioxide on post-cross-clamp cerebral oxygenation during carotid endarterectomy under general anesthesia. *Anesth Analg.* 2010;110(2):581-587.
- Pugliese F, Ruberto F, Tosi A, et al. Regional cerebral saturation versus transcranial Doppler during carotid endarterectomy under regional anaesthesia. *Eur J Anaesthesiol.* 2009;26(8):643-647.
- Matsumoto S, Nakahara I, Higashi T, et al. Near-infrared spectroscopy in carotid artery stenting predicts cerebral hyperperfusion syndrome. *Neurology.* 2009;72(17):1512-1518.
- Stoneham MD, Lodi O, de Beer TC, Sear JW. Increased oxygen administration improves cerebral oxygenation in patients undergoing awake carotid surgery. *Anesth Analg.* 2008;107(5):1670-1675.
- Lee TS, Hines GL, Feuerman M. Significant correlation between cerebral oximetry and carotid stump pressure during carotid endarterectomy. *Ann Vasc Surg.* 2008;22(1):58-62.
- Moritz S, Kasprzak P, Arlt M, Taeger K, Metz C. Accuracy of cerebral monitoring in detecting cerebral ischemia during carotid endarterectomy: a comparison of transcranial Doppler sonography, near-infrared spectroscopy, stump pressure, and somatosensory evoked potentials. *Anesthesiology.* 2007;107(4):563-569.
- Botes K, Le Roux DA, Van Marle J. Cerebral monitoring during carotid endarterectomy – a comparison between electroencephalography, transcranial cerebral oximetry and carotid stump pressure. *S Afr J Surg.* 2007;45(2):43-46.

Fassiadis N, Zayed H, Rashid H, Green DW. Invos cerebral oximeter compared with the transcranial Doppler for monitoring adequacy of cerebral perfusion in patients undergoing carotid endarterectomy. *Int Angiol.* 2006;25(4):401-406.

Rigamonti A, Scandroglio M, Minicucci F, Magrin S, Carozzo A, Casati A. A clinical evaluation of near-infrared cerebral oximetry in the awake patient to monitor cerebral perfusion during carotid endarterectomy. *J Clin Anesth.* 2005;17(6):426-430.

El-Dawlatly AA. Carotid endarterectomy: a study of cerebral oxygen saturation vs stump pressure. *Middle East J Anaesthesiol.* 2004;17(6):1037-1044.

Mille T, Tachimiri ME, Klersy C, et al. Near infrared spectroscopy monitoring during carotid endarterectomy: which threshold value is critical? *Eur J Vasc Endovasc Surg.* 2004;27(6):646-650.

Kragsterman B, Pärsson H, Bergqvist D. Local haemodynamic changes during carotid endarterectomy: the influence on cerebral oxygenation. *Eur J Vasc Endovasc Surg.* 2004;27(4):398-402.

Cuadra SA, Zwerling JS, Feuerman M, Gasparis AP, Hines GL. Cerebral oximetry monitoring during carotid endarterectomy: effect of carotid clamping and shunting. *Vasc Endovascular Surg.* 2003;37(6):407-413.

Hirofumi O, Otone E, Hiroshi I, et al. The effectiveness of regional cerebral oxygen saturation monitoring using near-infrared spectroscopy in carotid endarterectomy. *J Clin Neurosci.* 2003;10(1):79-83.

Takeda N, Fujita K, Katayama S, Tamaki N. Cerebral oximetry for the detection of cerebral ischemia during temporary carotid artery occlusion. *Neurol Med Chir (Tokyo).* 2000;40(11):557-562.

Samra SK, Dy EA, Welch K, Dorje P, Zelenock GB, Stanley JC. Evaluation of a cerebral oximeter as a monitor of cerebral ischemia during carotid endarterectomy. *Anesthesiology.* 2000;93(4):964-970.

Fearn SJ, Picton AJ, Mortimer AJ, Parry AD, McCollum CN. The contribution of the external carotid artery to cerebral perfusion in carotid disease. *J Vasc Surg.* 2000;31(5):989-993.

de Letter JA, Sie TH, Moll FL, Algra A, Eikelboom BC, Ackerstaff GA. Transcranial cerebral oximetry during carotid endarterectomy: agreement between frontal and lateral probe measurements as compared with an electroencephalogram. *Cardiovasc Surg.* 1998;6(4):373-377.

Carlin RE, McGraw DJ, Calimlim JR, Mascia MF. The use of near-infrared cerebral oximetry in awake carotid endarterectomy. *J Clin Anesth.* 1998;10(2):109-113.

de Letter JA, Sie HT, Thomas BM, et al. Near-infrared reflected spectroscopy and electroencephalography during carotid endarterectomy: in search of a new shunt criterion. *Neurol Res.* 1998;20(suppl 1):S23-S27.

Duffy CM, Manninen PH, Chan A, Kearns CF. Comparison of cerebral oximeter and evoked potential monitoring in carotid endarterectomy. *Can J Anaesth.* 1997;44(10):1077-1081.

Samra SK, Dorje P, Zelenock GB, Stanley JC. Cerebral oximetry in patients undergoing carotid endarterectomy under regional anesthesia. *Stroke.* 1996;27(1):49-55.

Williams IM, Vohra R, Farrell A, Picton AJ, Mortimer AJ, McCollum CN. Cerebral oxygen saturation, transcranial Doppler ultrasonography and stump pressure in carotid surgery. *Br J Surg.* 1994;81(7):960-964.