
Microstream[®]

ETCO₂ SAMPLING

CapnoLine[®] H O₂ For Non-intubated Patients Infant and Neonatal Directions for Use

CapnoLine[®] H O₂

Length=3m/10ft

US Patents: www.covidien.com/patents.

© 2015 Oridion Medical 1987 Ltd.

All rights reserved.

PN 012345E




Oridion[®]

CapnoLine® H O₂

For Non-intubated Infant and Neonatal Patients

Instructions for use:

1. Connect the connector of the CapnoLine® H O₂ sampling line to the CO₂ port on the Microstream monitor.
① Screw the connector into the monitor clockwise until it can no longer be turned.
2. Attach the O₂ tubing ② to an O₂ delivery system.
3. Connect the sampling line to the patient by gently inserting the nare sections of the patient face piece into the nares of the patient. ③ The tubing may be looped around the patient's ears (as seen in figure 3), if desired, to hold the sampling line in place better.
4. Ensure that the tubing is not twisted or crimped.

CapnoLine® H O₂ for Non-intubated Patients Infant and Neonatal

For use with all Microstream-enabled monitors

Single patient use, not to be re-processed

Non-sterile / non-pyrogenic

Does not contain natural rubber latex

- Warning:** Loose or damaged connections may compromise ventilation or cause an inaccurate measurement of respiratory gases. Securely connect all components and check connections for leaks according to standard clinical procedures.
- Warning:** Do not cut or remove any part of the sample line. Cutting the sample line could lead to erroneous readings.
- Warning:** Do not cover the Nafion®** ** on the CO₂ tubing. Covering the Nafion®** will limit its effectiveness in reducing moisture and occlusions in the sampling line.
- Warning:** Check CO₂ and O₂ tubing regularly during use to ensure that no kinks are present. Kinked tubing may cause inaccurate CO₂ sampling or affect O₂ delivery to patient.
- Caution:** Ensure that tubing is not stretched during use.
- Caution:** Clear the patient's nose of mucus and other obstructions before use to optimize duration of use, CO₂ reading accuracy and O₂ delivery. Duration of use, CO₂ readings as well as O₂ inhalation may be affected if patient nares are not patent.
- Caution:** Do not attempt to clean, disinfect, sterilize or flush any part of the sampling line as this can cause damage to the monitor.
- Caution:** Dispose of sampling lines according to standard operating procedures or local regulations for the disposal of contaminated medical waste.
- Caution:** Use of a CO₂ sampling line with H in its name (indicating that it is for use in humidified environments) during MRI scanning may cause interference. These sampling lines include CapnoLine H O₂. The use of non H sampling lines is advised.
- Note:** Replace the sampling line according to hospital protocol or when a blockage is indicated by the device. Excessive patient secretions or a build-up of liquids in the airway tubing may occlude the sampling line, requiring more frequent replacement.
- Note:** When connecting a sampling line to the monitor, insert the sampling line connector clockwise into the monitor CO₂ port and turn the connector until it can no longer be turned, to ensure that it is connected securely to the monitor. This will assure that there is no leak of gases during measurement at the connection point and that measurement accuracy is not compromised.
- Note:** Following connection of the CO₂ sampling line, check that CO₂ values appear on the monitor display.
- Note:** Expect delay times for the following products to be longer than the standard by the given times: Long (4m) CapnoLines - 3 seconds; 3m CapnoLines - 1.5 seconds.
- Note:** Note: Sampling lines with H in their names include a moisture reduction component (Nafion®** or its equivalent) for use in higher humidity environments where long duration use of CO₂ sampling is required.
- Note:** Note: Sampling lines are recommended for use with oxygen provided at up to 3l/min. At higher levels of oxygen provision, dilution of CO₂ readings may occur, leading to lower CO₂ values.

Federal law (U.S.) restricts this device to sale by or on the order of a physician.

** or alternate moisture reduction technology

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 Single patient use	For use with a single patient
 This product does not contain natural rubber latex	Does not contain natural rubber latex
 -4°F -20°C Temperature limitation	Storage temperature
Rx ONLY	Rx Only
 Consult instructions for use	Consult instructions for use
 Caution, consult accompanying documents	Caution, consult accompanying documents
 MR Unsafe	Not safe for use in MR environments
 Fire hazard	Fire Hazard
 0482	CE Mark

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Oridion®, Microstream®, FilterLine®, Smart CapnoLine®, CapnoLine®, and Smart CapnoLine Guardian™ are trademarks or registered trademarks of Oridion Medical 1987 Ltd.™* and ®* are trademarks of their respective owners.

Oridion Medical 1987 Ltd.
7 Hamarpe St.
P.O. Box 45025
Jerusalem 9777407, Israel
US Toll Free: 1-888-674-3466
1-888-ORIDION

Obelis S.A.
Bd Général Wahis 53
1030 Brussels
Belgium
www.oridion.com