

## References (continued)

- <sup>17</sup> Fantoni C, Raffa S, Regoli F, et al. Cardiac resynchronization therapy improves heart rate profile and heart rate variability of patients with moderate to severe heart failure. *J Am Coll Cardiol*. November 15, 2005;46(10):1875-1882.
- <sup>18</sup> Abraham WT. Rationale and design of a randomized clinical trial to assess the safety and efficacy of cardiac resynchronization therapy in patients with advanced heart failure: the Multicenter InSync Randomized Clinical Evaluation (MIRACLE). *J Card Fail*. December 2000; 6(4):369-380.
- <sup>19</sup> Ousdigian KT, Johnson WB, Koehler JL, Wilkoff BL. Very poor rate control during atrial tachyarrhythmias predicts cardiovascular hospitalizations. *J Card Fail*. August 2006; Volume 12, Issue 6(suppl 1):S57-S58.
- <sup>20</sup> Ousdigian KT, Thompson D, Qi D, et al. A single day of AF with poor rate control predicts a subsequent week of reduced CRT. [abstract PO6-68]. *Heart Rhythm*. 2007; 4 (Suppl 5S): S379.
- <sup>21</sup> Andriulli J, Coles J, Hettrick DA, et al. Association between decreased intra-thoracic impedance and ventricular tachyarrhythmias. *Int J Cardiol*. 2007. doi:10.1016/j.ijcard.2006.11.153.
- <sup>22</sup> Gasparini M, Auricchio A, Regoli F, et al. Four-year efficacy of cardiac resynchronization therapy on exercise tolerance and disease progression: the importance of performing atrioventricular junction ablation in patients with atrial fibrillation. *J Am Coll Cardiol*. August 15, 2006; 48(4):734-743.
- <sup>23</sup> La Rovere MT, Pinna GD, Maestri R, et al. Short-term heart rate variability strongly predicts sudden cardiac death in chronic heart failure patients. *Circulation*. February 4, 2003;107(4):565-570.
- <sup>24</sup> Dialysis Patient Case. Courtesy of David Martin, MD, Cardiology Consultants, Medford, OR. 2006 Medtronic Inc., data on file, UC200601357 EN.
- <sup>25</sup> Hunt SA, Abraham WT, Chin MH, et al. ACC/AHA 2005 Guideline Update for the Diagnosis and Management of Chronic Heart Failure in the Adult: a report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines (Writing Committee to Update the 2001 Guidelines for the Evaluation and Management of Heart Failure): developed in collaboration with the American College of Chest Physicians and the International Society for Heart and Lung Transplantation: endorsed by the Heart Rhythm Society. *Circulation*. September 20, 2005;112(12):e154-e235.
- <sup>26</sup> Vinson JM, Rich MW, Sperry JC, Shah AS, McNamara T. Early readmission of elderly patients with congestive heart failure. *J Am Geriatr Soc*. December 1990;38(12):1290-1295.
- <sup>27</sup> Rathman L. Use of device diagnostics as an educational tool to improve patient adherence. *Am J Cardiol*. May 21, 2007;99(10A):29G-33G.

Considerations in this Clinicians Practical Guide are courtesy of Lisa Rathman, MSN, CRNP; Roy Small, MD, FACC; and Jill Repoley, MSN, CRNP. The Heart Group and Lancaster General Hospital, Lancaster, PA.

## References

- <sup>1</sup> Medtronic Concerto and InSync Sentry CRT-Ds and Virtuoso DR/VR ICD System Reference Guides.
- <sup>2</sup> Yu CM, Wang L, Chau E, et al. Intrathoracic impedance monitoring in patients with heart failure: correlation with fluid status and feasibility of early warning preceding hospitalization. *Circulation*. August 9, 2005;112(6):841-848.
- <sup>3</sup> Wang L. Fundamentals of intrathoracic impedance monitoring in heart failure. *Am J Cardiol*. May 21, 2007;99(10A):3G-10G.
- <sup>4</sup> Vollmann D, Nägele H, Schauerer P, et al. Clinical utility of intrathoracic impedance monitoring to alert patients with an implanted device of deteriorating chronic heart failure. *Eur Heart J*. August 2007;28(15):1835-1840.
- <sup>5</sup> Repoley J, Dukes-Graves D, Kiser C. Surgical revision of implantable device "pocket" or lead can alter thoracic impedance-based indices of worsening heart failure. *J Card Fail*. 2006;12(6) (suppl):S66.
- <sup>6</sup> Patient Case: LV Lead Dislodgement and Pocket Revision. Courtesy of R. Ward Pulliam and Jill Repoley, NP. The Heart Group and Lancaster General Hospital, Lancaster, PA. 2006 Medtronic, Inc., data on file, UC200602091 EN.
- <sup>7</sup> Small R, Tang W, Wickemeyer W, et al. Managing heart failure patients with intra-thoracic impedance monitoring: a multi-center US evaluation. *J Card Fail*. August 2007;13(6)(suppl 2):S113-S114.
- <sup>8</sup> Yu CM, Wang L, Chau E, et al. Correlation of device-based intra-thoracic impedance and patient fluid status during intravenous diuretic therapy in acute CHF. *Eur Heart J*. 2002;23(abstract suppl):158.
- <sup>9</sup> Tang W, Germany R, Andriulli J, et al. Weekly and seasonal variations in daily activity and intrathoracic impedance in heart failure patients. *J Card Fail*. August 2007;13(6)(suppl 2):S182.
- <sup>10</sup> Braunschweig F, Mortensen PT, Gras D, et al. Monitoring of physical activity and heart rate variability in patients with chronic heart failure using cardiac resynchronization devices. *Am J Cardiol*. May 1, 2005;95(9):1104-1107.
- <sup>11</sup> Fuster V, Rydén LE, Cannom DS, et al. ACC/AHA/ESC 2006 Guidelines for the Management of Patients with Atrial Fibrillation: Executive Summary: A report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines and the European Society of Cardiology Committee for Practice Guidelines (Writing Committee to Revise the 2001 Guidelines for the Management of Patients with Atrial Fibrillation): developed in collaboration with the European Heart Rhythm Association and the Heart Rhythm Society. *Circulation*. August 15, 2006;114(7):e257-e354.
- <sup>12</sup> Adamson PB, Smith AL, Abraham WT, et al. Continuous autonomic assessment in patients with symptomatic heart failure: prognostic value of heart rate variability measured by an implanted cardiac resynchronization device. *Circulation*. October 19, 2004;110(16):2389-2394.
- <sup>13</sup> Casolo GC, Stroder P, Sulla A, Chelucci A, Freni A, Zeraushek M. Heart rate variability and functional severity of congestive heart failure secondary to coronary artery disease. *Eur Heart J*. March 1995;16(3):360-367.
- <sup>14</sup> Sweeney MO, Hellkamp AS. Heart failure during cardiac pacing. *Circulation*. May 2, 2006;113(17):2082-2088.
- <sup>15</sup> Smit MD, Van Dessel PF, Nieuwland W, et al. Right ventricular pacing and the risk of heart failure in implantable cardioverter-defibrillator patients. *Heart Rhythm*. December 2006;3(12):1397-1403.
- <sup>16</sup> Adamson PB, Kleckner KJ, VanHout WL, Srinivasan S, Abraham WT. Cardiac resynchronization therapy improves heart rate variability in patients with symptomatic heart failure. *Circulation*. July 22, 2003;108(3):266-269.

### World Headquarters

Medtronic, Inc.  
710 Medtronic Parkway  
Minneapolis, MN 55432-5604  
USA  
Tel: (763) 514-4000  
Fax: (763) 514-4879  
www.medtronic.com

Medtronic USA, Inc.  
Toll-free: 1 (800) 328-2518  
(24-hour technical support  
for physicians and  
medical professionals)

### Europe

Medtronic International Trading Sàrl  
Route du Molliu 31  
CH-1131 Tolochenaz  
Switzerland  
Tel: (41 21) 802 7000  
Fax: (41 21) 802 7900  
www.medtronic.com

### Canada

Medtronic of Canada Ltd.  
6733 Kitimat Road  
Mississauga, Ontario L5N 1W3  
Canada  
Tel: (905) 826-6020  
Fax: (905) 826-6620  
Toll-free: 1 (800) 268-5346

### Asia Pacific

Medtronic International, Ltd.  
16/F Manulife Plaza  
The Lee Gardens, 33 Hysan Avenue  
Causeway Bay  
Hong Kong  
Tel: (852) 2891 4456  
Fax: (852) 2891 6830  
enquiryap@medtronic.com  
www.medtronic.com

### Latin America

Medtronic USA, Inc.  
Doral Corporate Center II  
3750 NW 87th Avenue Suite 700  
Miami, FL 33178  
USA  
Tel: (305) 500-9328  
Fax: (786) 709-4244  
www.medtronic.com

UC200802511a EN  
© Medtronic, Inc. 2008  
All Rights Reserved  
Printed in USA  
April 2008



## Clinical Application of the HF Diagnostic Trends A three-step approach to patient management



Trends:<sup>1</sup>

- OptiVol<sup>®</sup> Impedance
- AT/AF Diagnostics
- Heart Rate Variability
- Percent Pacing






# 1 REVIEW IMPEDANCE TREND AND OPTIVOL TREND DATA<sup>1-3\*</sup>

Accessible on Medtronic CareLink® Programmer and CardioSight®


**Start with the thoracic impedance trend and evaluate the patient's Daily impedance.**



The thoracic impedance trend is an accurate representation of the patient's actual Daily impedance level. The Daily impedance trend may change as the patient's status changes.

**Reference impedance initializes 34 days post-implant.<sup>3</sup>** It adapts slowly to changes in Daily impedance and acts as the patient's own control.<sup>3</sup>

**Review the OptiVol fluid index trend next . . .**



The OptiVol fluid index trend is a graphical representation of the consecutive day-to-day difference between the Daily and Reference impedance. The OptiVol fluid index will rise as intrathoracic fluid level increases. Generally, the OptiVol trends will demonstrate the following patterns as the patient's status changes.<sup>2-8</sup>

<sup>1</sup>OptiVol Fluid Management is a feature in Concerto® and InSync Sentry® CRT-Ds and Virtuoso® DR/VR ICDs.



# 1a IMPEDANCE AND OPTIVOL TREND DATA<sup>1</sup>

- Consider reasons for changes in Daily impedance – possible elements that could impact impedance including but not limited to:<sup>\*</sup>
  - Pulmonary edema<sup>2,4</sup>
  - Volume retention<sup>2,4</sup>
  - Pulmonary congestion<sup>2,4</sup>
- Consider non-heart failure explanations for changes in impedance:<sup>\*</sup>
  - Pocket or lead revision<sup>4-7</sup>
  - Pocket infection<sup>2</sup>
  - Respiratory infection<sup>2</sup>
  - Pleural/pericardial effusion<sup>3</sup>
  - Diuretic therapy changes<sup>2,8</sup>
  - Anemia<sup>7</sup>

<sup>\*</sup>Based on clinician experience using OptiVol trend data in practice and during clinical trials.<sup>2-8</sup>




Considerations in this document do not replace a physician's expert judgment. The physician's knowledge of the patient's medical condition should be considered, and clinical and device intervention considerations may be tailored to fit the patient.



## EVALUATE OTHER DEVICE TREND DATA<sup>1</sup>

Evaluate trend for:



Evaluate all trends for changes over time; typically interrelationships can be seen across multiple trends as the patient's condition changes. The lists on the next page reflect some of the more common trend changes that may be observed as a patient's status changes.



## CHECK FOR TREND CHANGES OVER TIME<sup>1</sup>

Condition improving:

Condition worsening:

Patient activity increasing<sup>7,9,10</sup>




Atrial fibrillation burden increasing<sup>11</sup>

HR variability increasing<sup>10,12,13</sup>




Average day/night HR increasing<sup>12,13</sup>

Separation between day/night HR increased<sup>13</sup>



High RV pacing in ICD patient<sup>14,15</sup>

CRT pacing maximized<sup>16-18</sup>



Ventricular rate with AF uncontrolled<sup>19,20</sup>

Ventricular rate with AF controlled<sup>21</sup>



CRT pacing below 100%<sup>18,22</sup>

Average day/night HR decreasing<sup>12,13</sup>



Patient activity declining<sup>7,9,10</sup>

Atrial fibrillation burden decreasing



HR variability decreasing<sup>10,12,13,23</sup>



Increasing



Neutral



Decreasing



## CLINICAL AND DEVICE INTERVENTION CONSIDERATIONS\*

### Daily impedance is changed, consider possible cause:

Although the most common cause, thoracic fluid accumulation is not the only reason thoracic impedance might change. Nevertheless, a thoracic impedance change may indicate a clinically relevant event that may require medical attention.<sup>2-9</sup>

#### Heart failure-related event:

- Pulmonary edema<sup>2,4</sup>
- Volume retention<sup>2,4</sup>
- Pulmonary congestion<sup>2,4</sup>
- Pocket or lead revision<sup>4-7</sup>

#### Clinically relevant event:

- Pocket infection<sup>2</sup>
- Pleural/pericardial effusion<sup>3</sup>
- Diuretic therapy changes<sup>2,8</sup>
- Anemia<sup>7</sup>
- Dialysis<sup>24</sup>
- Respiratory infection<sup>2</sup>

#### Consider the following actions:

- Corroborate with other diagnostic trend findings<sup>1</sup>
- Educate patient, i.e., dietary and/or fluid, medication compliance<sup>25-27</sup>
- Adjust HF medications to achieve optimal management<sup>25</sup>
- Titrate diuretics and appropriately monitor electrolytes<sup>25</sup>
- Refer to electrophysiology for system-related issue<sup>18,25</sup>

#### AT/AF observed, consider:<sup>11</sup>

- Corroborating with other diagnostic trend findings<sup>1</sup>
- Cardioversion
- Instituting of anticoagulation therapy
- Adjusting rhythm control therapy
- Surgical procedures

#### AF with rapid ventricular response observed, consider:<sup>11</sup>

- Corroborating with other diagnostic trend findings<sup>1</sup>
- Cardioversion
- Program Conducted AF Response "ON" in CRT-D devices (Concerto)<sup>1</sup>
- Program Mode Switch "ON" – to a nontracking pacing mode<sup>1</sup>
- Adjusting rate control therapy

#### Decrease in patient activity, consider:

- Corroborating with other diagnostic trend findings<sup>1</sup>
- Promoting regular exercise as indicated<sup>25-27</sup>

#### Increasing day/night HR with poor separation, consider:<sup>12,13,23</sup>

- Corroborating with other diagnostic trend findings<sup>1</sup>
- Atrial fibrillation
- Adjusting HF medications<sup>25</sup>
- Sleep apnea or nocturnal dyspnea

#### Persistently low HRV (< 50 ms) or declining values, consider:<sup>12,13,23</sup>

- Corroborating with other diagnostic trend findings<sup>1</sup>
- Adjusting HF medications<sup>25</sup>
- Consider if patient is in AF or atrially paced
- Not measured by device if patient is in AF<sup>1</sup>

#### Percent ventricular pacing below 100% in CRT patient, consider:<sup>18,22,25</sup>

- Corroborating with other diagnostic trend findings<sup>1</sup>
- Atrial fibrillation with rapid ventricular response<sup>18,19</sup>
- Evaluating device settings may necessitate electrophysiology consult<sup>18</sup>
- AV optimization
- Lead dislodgement<sup>5,6</sup>

#### Percent pacing higher than desired in dual chamber ICD patient, consider:<sup>14,15</sup>

- Corroborating with other diagnostic trend findings<sup>1</sup>
- AV delay
- Program MVP® Mode "ON"<sup>1</sup>
- Atrial fibrillation
- Reprogram rate response settings<sup>1</sup>

\*Courtesy of The Heart Group and Lancaster General Hospital, Lancaster, PA.

Considerations in this document do not replace a physician's expert judgment. The physician's knowledge of the patient's medical condition should be considered, and clinical and device intervention considerations may be tailored to fit the patient.



## EVALUATE THE PATIENT AND CORRELATE FINDINGS TO MAKE CLINICAL DECISION<sup>25\*</sup>

Upon completing review of the device-based diagnostics, the next step is to evaluate the patient. The patient evaluation should consist of your standard approach to patient assessment in the clinic. This step might include, but is not limited to, the following measures or order:

#### Review medication list

- Medications prescribed by HF clinicians
- Medications prescribed by other: GP, IM, Nephrologist, etc.
- Over-the-counter medications

#### Assess patient's current status

- Review chart for past medical history
- Interview patient
- Perform physical examination

#### Assess patient compliance

- Medications
- Diet
- Exercise

#### Review other diagnostic studies

- Labs
- X-rays
- Echocardiography
- Perfusion studies
- Educate patient concerning salt and/or fluid restriction and compliance with medications
- Adjust HF medications (e.g., ACE-I, BB, etc.) to achieve optimal management
- Add/increase diuretic and appropriate monitoring of electrolytes



\*Courtesy of The Heart Group and Lancaster General Hospital, Lancaster, PA.

Considerations in this document do not replace a physician's expert judgment. The physician's knowledge of the patient's medical condition should be considered, and clinical and device intervention considerations may be tailored to fit the patient.

### Brief Statement

InSync Sentry® Models 7297/7299, Concerto® Model C154DWK, and Virtuoso® Models D154AWG/D154VWC (DR/VR)

**Indications for InSync Sentry and Concerto:** The InSync Sentry and Concerto are indicated for ventricular antitachycardia pacing and ventricular defibrillation for automated treatment of life-threatening ventricular arrhythmias, and for the reduction of the symptoms of moderate to severe heart failure (NYHA Functional Class III or IV) in those patients who remain symptomatic despite stable, optimal medical therapy, and have a left ventricular ejection fraction less than or equal to 35% and a prolonged QRS duration.

**Indications for Virtuoso:** Virtuoso DR/VR devices are indicated to provide ventricular antitachycardia pacing and ventricular defibrillation for automated treatment of life-threatening ventricular arrhythmias in patients with NYHA Functional Class II/III heart failure. The Virtuoso DR device is also indicated for use in the above patients with atrial tachyarrhythmias, or those patients who are at significant risk of developing atrial tachyarrhythmias. Atrial rhythm management features, available on the Virtuoso DR, such as Atrial Rate Stabilization (ARS), Atrial Preference Pacing (APP), and Post Mode Switch Overdrive Pacing (PMOP) are indicated for the suppression of atrial tachyarrhythmias in ICD-indicated patients with atrial septal lead placement and an ICD indication. Due to the addition of the OptiVol® diagnostic feature, the Virtuoso indication is limited to NYHA Functional Class II/III heart failure patients who are indicated for an ICD. The clinical value of the OptiVol fluid monitoring diagnostic feature has not been assessed in those patients who do not have fluid retention related symptoms due to heart failure.

**Contraindications:** The InSync Sentry, Concerto, and Virtuoso are contraindicated in patients whose ventricular tachyarrhythmias may have transient or reversible causes; patients with incessant VT or VF; and patients who have a unipolar pacemaker.

**Warnings and Precautions:** Changes in a patient's disease and/or medications may alter the efficacy of the device's programmed parameters. Patients should avoid sources of magnetic and electromagnetic radiation to avoid possible underdetection, inappropriate sensing and/or therapy delivery, tissue damage, induction of an arrhythmia, device electrical reset, or device damage. Do not place transthoracic defibrillation paddles directly over the device. Certain programming and device operations may not provide cardiac resynchronization.

**Potential Complications:** Potential complications include, but are not limited to, rejection phenomena, erosion through the skin, muscle or nerve stimulation, oversensing, failure to detect and/or terminate tachyarrhythmia episodes, acceleration of ventricular tachycardia, and surgical complications such as hematoma, infection, inflammation, and thrombosis.

*See the device manual for detailed information regarding the implant procedure, indications, contraindications, warnings, precautions, and potential complications/adverse events. For further information, please call Medtronic at 1 (800) 328-2518 and/or consult Medtronic's website at [www.medtronic.com](http://www.medtronic.com).*

**Caution:** Federal law (USA) restricts these devices to sale by or on the order of a physician.