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表示径 (mm)	表示長 (mm)					
	16 mm	18 mm	20 mm	25 mm	30 mm	35 mm
3.00 mm	PED2-300-16	PED2-300-18	PED2-300-20	—	—	—
3.25 mm	PED2-325-16	PED2-325-18	PED2-325-20	—	—	—
3.50 mm	PED2-350-16	PED2-350-18	PED2-350-20	PED2-350-25	PED2-350-30	PED2-350-35
3.75 mm	PED2-375-16	PED2-375-18	PED2-375-20	PED2-375-25	PED2-375-30	PED2-375-35
4.00 mm	PED2-400-16	PED2-400-18	PED2-400-20	PED2-400-25	PED2-400-30	PED2-400-35
4.25 mm	PED2-425-16	PED2-425-18	PED2-425-20	PED2-425-25	PED2-425-30	PED2-425-35
4.50 mm	PED2-450-16	PED2-450-18	PED2-450-20	PED2-450-25	PED2-450-30	PED2-450-35
4.75 mm	PED2-475-16	PED2-475-18	PED2-475-20	PED2-475-25	PED2-475-30	PED2-475-35
5.00 mm	PED2-500-16	PED2-500-18	PED2-500-20	PED2-500-25	PED2-500-30	PED2-500-35

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
# Shield Technology™

## PIPELINE™ Flex Flow diverter system

## PIPELINE™ Flex with Shield Technology™

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

Individual Value Plot of Max Resheathing Force (Newton)  
(N=160 for each sample groupe)




Re-sheathing Force Measured in Simulated Use Testing of all time points, including as manufactured, 1 year accelerated aging and 3 years accelerated aging(N=160)

Individual Value Plot of Max Delivery Force (Newton)  
(N=183 for each sample groupe)

Delivery Force Measured in Simulated Use Testing of all time points, including as manufactured, 1 year accelerated aging and 3 years accelerated aging(N=183)

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



(Left) Flow-loop construction with leak-free seals. (Middle) Flow-loop with Pipeline Shield device. (Right) Pulsatile drive system with flow-loops mounted.

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
## Analysis of neointima development in flow diverters using optical coherence tomography imaging

Neointimal ratio on day 21 in each stent. The Pipeline Flex embolization device (PED Flex) and PED Shield stent had a significantly higher neointimal ratio than the Solitaire stent (p<0.01 and p<0.05, respectively). There is no difference between the PED Shield and PED Flex. (B) Neointimal thickness ratio on day 21 in each stent. The PED Shield stent shows a significantly higher neointimal thickness ratio than the PED Flex and Solitaire stents (p<0.05 and p<0.01, respectively). There is no difference between the PED Flex and Solitaire stents. Bars are +SEs.

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# Clinical evidence – 臨床成績<sup>5</sup>

プライマリーエンドポイントは a major stroke in the territory supplied by the treated artery or neurologic death at 1-year post-procedure

Treatment of intracranial aneurysms using the pipeline flex embolization device with shield technology: angiographic and safety outcomes at 1-year follow-up

## Primary and secondary safety outcomes at 1-year postprocedure

Endpoint	Event rate n/N (%) (95% CI)
Major stroke in the territory supplied by the treated artery or neurological death	0/50 (0%) (0%, 5.8%)
Major stroke in the territory supplied by the treated artery	0/50 (0%) (0%, 5.8%)
Neurological death	0/50 (0%) (0%, 5.8%)
Pipeline Flex embolization device with shield technology-related neurologic adverse event rate	7/50 (14.0%) (5.8%, 26.7%)

Reproduced from J Neurointerv Surg. Mario Martínez-Galdámez, et al. 11. 396-399. 2019 with permission from BMJ Publishing Group Ltd.

## Treatment of intracranial aneurysms using the pipeline flex embolization device with shield technology: angiographic and safety outcomes at 1-year follow-up (PFLEX Study)概要

目的: The Pipeline Flex Embolization Device with Shield Technology (Pipeline Shield) is an updated version of the PED which has been modified to include a surface phosphorylcholine biocompatible polymer. Its early technical success and efficacy of the Pipeline Shield for the treatment of IAs.

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Primary endpoint: The primary endpoint was a major stroke in the territory supplied by the treated artery or neurologic death at 1-year post-procedure. Angiographic outcomes were also assessed by an independent radiology laboratory at 6 months and 1 year.

Secondary endpoint: The secondary endpoint was the rate of Pipeline Flex Embolization Device with Shield Technology-related neurologic adverse events within 1 year of the procedure.

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©2019 Medtronic Japan Co., Ltd. All Rights Reserved.