

U.S. MEDTRONIC DISSECTION TRIAL 3 YEAR RESULTS


PROSPECTIVE, NON-RANDOMIZED, MULTI-CENTER TRIAL

TEVAR in Acute, Complicated Type B Aortic Dissection: Results from the Valiant Captivia IDE Trial (n=50)[†]

CLINICAL STATUS AT ONSET

Malperfusion: **86%**

Rupture: **20%**



- The Valiant™ Captivia™ stent graft system demonstrates safety and effectiveness in the treatment of acute, complicated Type B aortic dissections
- The Valiant™ Captivia™ stent graft system is stabilizing the aorta

3 YEAR AORTIC REMODELING

92.3% True Lumen diameter increase/stable

69.2% False Lumen diameter decrease/stable

75.0% Partial/complete False Lumen thrombosis

KEY OUTCOMES THROUGH 3 YEARS

0% Ruptures or Conversions

3 Patients with Dissection-related Reinterventions

RETROGRADE TYPE A AORTIC DISSECTION POST-TEVAR FOR TYPE B AORTIC DISSECTION:

- Based on the published IDE data in the chart below, Retrograde Type A Dissections (RTADs) can occur with any device.

STUDY	DEVICES STUDIED	PATIENTS	RETROGRADE TYPE A	1Y RTAD (%)
Medtronic Dissection IDE ¹	Medtronic Valiant™ Captivia™	50 acute	1 RTAD <30d 1 RTAD <1yr	4%
Gore Dissection IDE ²	GORE cTAG™*	50 acute	3 RTAD <30d	6%
Cook STABLE IDE ³	Cook Zenith™* TX2™*	24 acute 16 sub acute	2 RTAD <30d 1 RTAD <1yr	7.5%

NOTE: These trials are not powered to be compared directly.

- According to the Systematic Review by Canaud et al., **pathology, landing zone, and oversizing are significantly associated with a risk of RTAD development regardless of device**, whereas proximal device configuration is not⁴

[†]Azzizadeh A. Endovascular Repair in Acute, Complicated Type B Aortic Dissection: 3-Year Results from the Valiant US-IDE Study, VIVA 2016

[‡]Dissection-related mortality

[§]All-cause mortality

¹Valiant Captivia IFU

²Gore cTAG IFU

³Lombardi JV, Cambria RP, Nienaber CA, et al. "Prospective multicenter clinical trial (STABLE) on the endovascular treatment of complicated Type B aortic dissection using a composite device design." J Vasc Surg. 2012 Mar; 55(3):629-640.e2. doi: 10.1016/j.jvs.2011.10.022. Epub 2011 Dec 9.

⁴Canaud, L, et al. "Retrograde Aortic Dissection After Thoracic Endovascular Aortic Repair." Annals of Surgery. 2014 Aug;260(2):389-95. doi: 10.1097/SLA.0000000000000585.

OBJECTIVE[†]

Evaluate the clinical performance of the Valiant™ thoracic stent graft with the Captivia™ delivery system for the treatment of acute, complicated Type B aortic dissections.

STUDY DESIGN

- 16 U.S. centers; enrollment 2010–2012
- N = 50 patients with acute, complicated (malperfusion, rupture) Type B aortic dissection
- 1° endpoint: all-cause mortality within 30 days of index procedure

MAIN CLINICAL FINDINGS THROUGH 3 YEARS

Treatment of trial patients resulted in:

- Primary endpoint achieved with 8% all-cause mortality at 30 days
- No conversions to open repair and three patients with dissection-related reinterventions
- No incidents of post-operative rupture
- 100% delivery and deployment success
- 100% coverage of primary entry tear
- 3 CVA / 1 CVI
- Favorable remodeling over stented segment

[†]Azzadeh A. Endovascular Repair in Acute, Complicated Type B Aortic Dissection: 3-Year Results from the Valiant US-IDE Study, VIVA 2016

Medtronic

Invatec S.p.A.
Via Martiri della Libertà 7
25030 Roncadelle (BS)
Italy
Tel: +39.030.2589311

Medtronic
710 Medtronic Parkway NE
Minneapolis, MN 55432
USA
Tel: +1.763.514.4000

**Medtronic International
Trading Sàrl**
Route du Molliau 31
CH-1131 Tolochenaz
Switzerland
Tel: +41.21.802.7000

Medtronic of Canada Ltd.
99 Hereford Street
Brampton, Ontario L6Y 0R3
Canada
Tel: +1.905.460.3800

Medtronic Latin America
3750 NW 87th Avenue, Suite 700
Miami, FL 33178
USA
Tel: +1.786.709.4200

Medtronic International Ltd.
49 Changi South Avenue 2
Singapore 486056
Tel: +65.6436.5000

Medtronic Australasia Pty Ltd.
97 Waterloo Road
North Ryde, NSW 2113
Australia
Tel: +61.29857.9000

Medtronic New Zealand
Unit N16, Mezzanine Level 5
Gloucester Park Road
Onehunga
Auckland
New Zealand

Medtronic Korea Co., Ltd.
5F, Sajo Building
1001 Daechi-dong, Kangnam-ku
Seoul, 135-280
Korea
Tel: +82.2.3404.3600

Medtronic METAFZ-LLC
Office Park, Block D, 2nd Floor
PO Box 500638
Dubai Internet City | Dubai,
United Arab Emirates
Tel: +971.4.818.2666

medtronic.com/aortic

UC201608794a EE © 2017 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.™™ Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. Not for distribution in the USA or France. Printed in Europe. 07/17