

**DRUG-COATED BALLOONS SHOW  
SUPERIOR THREE-YEAR OUTCOMES  
VS. ANGIOPLASTY: RESULTS FROM  
IN.PACT SFA RANDOMIZED TRIAL**

**PRAKASH KRISHNAN, MD, FACS**

**NEW YORK, UNITED STATES**

**Medtronic  
Further. Together**

# Drug-coated Balloons Show Superior Three-Year Outcomes versus Angioplasty: Results from the IN.PACT SFA Randomized Trial

*Prakash Krishnan, MD, FACS*

*Associate Professor of Medicine*

*Icahn School of Medicine at Mount Sinai*

*New York, New York USA*

# Background

- Drug-coated balloons have significantly changed the treatment landscape for femoropoplial artery disease
- IN.PACT™ Admiral™ DCB has demonstrated best in class patency results at 1 and 2 years

		1-Year	2-Year	3-Year
Primary Patency	IN.PACT DCB <sup>1</sup>	87.5%	78.9%	?
	LUTONIX DCB <sup>2,3</sup>	73.5%	58.6%	?
CD-TLR	IN.PACT DCB <sup>1</sup>	2.4%	9.1%	?
	LUTONIX DCB <sup>2,3</sup>	12.3%	18.0%	?

- Unique DCB designs and variable clinical outcomes across technology platforms support there being no class effect
- Presentation of the 3-year results will answer the question of continued durability of the IN.PACT Admiral DCB

1. Laird JR, et al. 24-Month Results of IN.PACT SFA. J Am Coll Cardiol 2015;66:2329-38.

2. Rosenfield K, et al. 12-Month Results of LEVANT2 Trial. New Engl J Med 2015; 373:145-53.

3. LEVANT 2 2-Year Results. Presented by Laurich C, SVS Chicago 2015.

# IN.PACT SFA Trial Overview

**Objective: Assess the safety and efficacy of IN.PACT Admiral DCB vs. standard PTA for the treatment of superficial femoral and proximal popliteal artery disease due to claudication and rest pain**



## IN.PACT SFA I

150 subjects enrolled at 13 EU sites  
Sep 2010-Apr 2011



## IN.PACT SFA II

181 subjects enrolled at 44 US sites  
Apr 2012-Jan 2013

### Robust Level 1 Evidence

- Prospective, multicenter EU and US, randomized (2:1), single-blinded trial
- 331 patients enrolled:  
IN.PACT DCB (n = 220) vs. PTA (n = 111)

### Rigorous and Unbiased

- Independent and blinded Duplex Ultrasound Core Lab,[1] Angiographic Core Lab,[2] and Clinical Events Committee[3]
- Independent Safety Monitoring Board
- External monitoring with 100% source data verification

### Durability of Outcomes

- Subjects followed up to 5 years

1. VasCore DUS Core Laboratory, Boston, MA, US;
2. SynvaCor Angiographic Core Laboratory, Springfield, IL, US;
3. Clinical Events Committee and Data Safety Monitoring services provided by HCRI, Boston, MA, US

# IN.PACT SFA: Investigators and Sites



## IN.PACT SFA I

150 subjects enrolled at 13 EU sites Sep 2010-Apr 2011

M. Brodmann, Graz, Austria	G. Sorropago, Mercogliano, Italy
G. Tepe, Rosenheim, Germany	P. Peeters, Bonheiden, Belgium
T. Zeller, Bad Krozingen, Germany	F. Vermassen, Gent, Belgium
D. Scheinert, Leipzig, Germany	C. Trani, Rome, Italy
A. Micari, Palermo, Italy	M. Bosiers, Dendermonde, Belgium
I. Baumgartner, Bern, Switzerland	J. Van den Berg, Lugano, Switzerland
S. Sixt, Hamburg, Germany	



## IN.PACT SFA II

181 subjects enrolled at 44 US sites Apr 2012-Jan 2013

P. Krishnan, New York, NY, USA	C. Walker, Houma, LA, USA
C. Metzger, Kingsport, TN, USA	N. Strickman, Houston, TX, USA
A. Jain, Fremont, CA, USA	R. Fairman, Philadelphia, PA, USA
R. Sachar, Raleigh, NC, USA	S. Laster, Kansas City, MO, USA
N. Farhat, Elyria, OH, USA	W. Gray, New York, NY, USA
L. Garcia, Boston, MA, USA	V. Ramaiah, Phoenix, AZ, USA
R. Malhotra, Glendale, AZ, USA	P. Alden, Minneapolis, MN, USA
S. Germanwala, Longview, TX, USA	C. Stinis, La Jolla, CA, USA
A. Pershad, Phoenix, AZ, USA	R. Dave, Camp Hill, PA, USA
B. Bigelow, Indianapolis, IN, USA	R. Gallino, Washington, DC, USA
J. Zidar, Raleigh, NC, USA	G. Ansel, Columbus, OH, USA
S. Ahanchi, Norfolk, VA, USA	M. Schermerhorn, Boston, MA, USA
R. Feldman, Ocala, FL, USA	M. Hunter, Cincinnati, OH, USA
R. Kovach, Brown Mills, NJ, USA	M. Dake, Stanford, CA, USA
M. Goodwin, Naperville, IL, USA	J. Benenati, Miami, FL, USA
L. Marone, Pittsburgh, PA, USA	P. Schneider, Honolulu, HI, USA
M. Shishehbor, Cleveland, OH, USA	R. Serry, Poway, CA, USA
D. Chew, Des Moines, IA, USA	J. Angle, Charlottesville, VA, USA
P. Soukas, Providence, RI, USA	K. Gupta, Kansas City, KS, USA
M. Garcia, Newark, DE, USA	P. Jones, Chicago, IL, USA
M. Mewissen, Milwaukee, WI, USA	G. Petrossian, Roslyn, NY, USA
R. Brown, Waco, TX, USA	A. Patel, Morristown, NJ, USA

# IN.PACT SFA Trial

## Blinded, Independently Assessed Outcomes

### Primary Efficacy Endpoint

Primary patency within 12 months, defined as freedom from clinically-driven TLR and DUS-derived restenosis (PSVR  $\leq$  2.4)

### Primary Safety Endpoint

Freedom from device- and procedure-related death through 30 days, and freedom from target limb major amputation and clinically-driven TVR within 12 months

- *MAEs (including all individual components of the primary endpoints and key secondary endpoints) are adjudicated by the blinded CEC through 5 years*
- *Restenosis is assessed by the blinded Duplex and Angiographic Core Labs through the 3-year follow-up visits*

# IN.PACT SFA Trial

## Baseline Clinical Characteristics

	IN.PACT n = 220 subjects	PTA n = 111 subjects	P-value
<b>Age, Y ± SD</b>	67.5 ± 9.5	68.0 ± 9.2	0.612
<b>Male, % (n)</b>	65.0% (143/220)	67.6% (75/111)	0.713
<b>Diabetes, % (n)</b>	40.5% (89/220)	48.6% (54/111)	0.161
<b>Hypertension, % (n)</b>	91.4% (201/220)	88.3% (98/111)	0.431
<b>Current smoker, % (n)</b>	38.6% (85/220)	36.0% (40/111)	0.719
<b>Rutherford class, % (n)</b>			
<b>2</b>	37.7% (83/220)	37.8% (42/111)	0.898
<b>3</b>	57.3% (126/220)	55.9% (62/111)	
<b>4</b>	5.0% (11/220)	5.4% (6/111)	
<b>5</b>	0.0% (0/220)	0.9% (1/111)	
<b>ABI / TBI, ± SD <sup>[1]</sup></b>	0.769 ± 0.228	0.744 ± 0.189	0.308

1. TBI allowed in cases of incompressible vessels in IN.PACT SFA II phase

# IN.PACT SFA Trial


## Baseline Lesion Characteristics

	<b>IN.PACT</b> n = 220 Subjects, n = 221 Lesions	<b>PTA</b> n = 111 Subjects, n = 113 Lesions	<b>P-value</b>
<b>Lesion length (cm ± SD)</b>	8.94 ± 4.89	8.81 ± 5.12	0.815
<b>Total occlusions, % (n)</b>	25.8% (57/221)	19.5% (22/113)	0.222
<b>Calcification, % (n)</b>	59.3% (131/221)	58.4% (66/113)	0.907
<b>Severe calcification, % (n)</b>	8.1% (18/221)	6.2% (7/113)	0.662
<b>Provisional stenting, % (n)</b>	7.3% (16/220)	12.6% (14/111)	0.110



# IN.PACT SFA Trial

## Primary Patency<sup>1</sup> Through 3 Years




<b>Number<sup>2</sup> DCB</b>	220	213	192	149	121
<b>at risk PTA</b>	111	108	69	52	41

- Freedom from core laboratory-assessed restenosis (duplex ultrasound PSVR  $\leq 2.4$ ) and clinically-driven target lesion revascularization through 36 months (adjudicated by a Clinical Events Committee blinded to the assigned treatment)
- Number at risk represents the number of evaluable subjects at the beginning of the each 30-day window

# IN.PACT SFA Trial

## Freedom from CD-TLR through 3 Years



<b>Number<sup>1</sup></b>	<b>DCB</b>	220	215	205	175	153
<b>at risk</b>	<b>PTA</b>	111	108	93	78	70

1. Number at risk represents the number of evaluable subjects at the beginning of the each 30-day window

# IN.PACT SFA Trial

## Effectiveness Outcomes through 3 Years

	IN.PACT DCB (N=220)	PTA (N=111)	P-value†
<b>Clinically-driven TLR</b> [1]	15.2% (30/197)	31.1% (32/103)	0.002
<b>All TLR</b> [2]	16.2% (32/197)	34.0% (35/103)	< 0.001
<b>Time to First CD-TLR</b>	542.9 ± 278.2	302.9 ± 213.0	< 0.001
<b>Primary Sustained Clinical Improvement</b> [3]	68.7% (114/166)	52.6% (51/97)	0.012
<b>ABI / TBI</b> [4]	0.917 ± 0.231	0.894 ± 0.194	0.429

1. Clinically-driven TLR adjudicated by an independent Clinical Event Committee, blinded to the assigned treatment based on any re-intervention at the target lesion due to symptoms or drop of ABI of  $\geq 20\%$  or  $> 0.15$  when compared to post-procedure baseline ABI
2. Any TLR includes clinically-driven and incidental or duplex driven TLR
3. Freedom from target limb amputation, target vessel revascularization (TVR), and increase in Rutherford class
4. TBI allowed / used in case of incompressible vessels in IN.PACT SFA II phase

† Unless otherwise indicated, all tests were for superiority using the Fisher's exact test for binary variables and t-test for continuous variables.

# IN.PACT SFA Trial

## Safety Outcomes through 3 Years

3-Year Outcomes	IN.PACT DCB N=220	PTA N=111	P-value†
<b>Primary Safety Composite</b> <sup>[1]</sup>	81.2% (160/197)	64.1% (66/103)	< 0.001‡
<b>Major Adverse Events</b> <sup>[2]</sup>	27.9% (55/197)	37.9% (39/103)	0.089
<b>All-cause Death</b> *	10.7% (21/197)	1.9% (2/103)	0.006
<b>Device- or Procedure-related Death</b>	0.0% (0/197)	0.0% (0/103)	N/A
<b>Clinically-driven TVR</b>	18.8% (37/197)	35.9% (37/103)	0.002
<b>Target Limb Major Amputation</b>	0.0% (0/197)	0.0% (0/103)	N/A
<b>Thrombosis</b>	2.0% (4/197)	4.9% (5/103)	0.283

1. Freedom from 30-day device and procedure-related death and target limb major amputation and clinically-driven TVR within 36 months

2. Composite of death, clinically-driven TVR, target limb major amputation, and thrombosis

\* No deaths were adjudicated as device- or procedure-related by the CEC; Median post-index days to death: 610 days in DCB vs. 637 days in PTA

† P-values are based on Fisher's exact test for superiority with significance level of 0.05

‡ 10% Non-inferiority Test Margin with one-sided 97.5005% CL


# All-Cause Mortality through 3 Years

Causes of Death Through 3 Years	Treatment Group	Days to Death	CEC Adjudication	
			Procedure-Related (Y/N)	Device-Related (Y/N)
<b>Cardiac-related</b>				
Acute Diastolic Congestive Heart Failure	DCB	540	NO	NO
Cardiac Arrest	DCB	568	NO	NO
Cardiac Arrest	DCB	610	NO	NO
CAD	DCB	615	NO	NO
Ischemic Cardiomyopathy	DCB	699	NO	NO
Cardiac Arrest	DCB	907	NO	NO
<b>Malignancy</b>				
Metastatic Colon Cancer	PTA	397	NO	NO
GI Cancer	DCB	561	NO	NO
Worsening Pre-existing Carcinoma	PTA	877	NO	NO
<b>Respiratory-related</b>				
Acute Respiratory Failure	DCB	657	NO	NO
Hypoxic Respiratory Failure	DCB	681	NO	NO
Acute Respiratory Failure	DCB	1051	NO	NO
<b>Other</b>				
Infarction of Right Cerebral Hemisphere In Anterior & Medial Flow Region	DCB	127	NO	NO
Biliary Sepsis	DCB	168	NO	NO
Perforated Transverse Colon Secondary to Cecal Volvulus	DCB	314	NO	NO
Sepsis	DCB	374	NO	NO
Deterioration of General Condition	DCB	603	NO	NO
Dementia	DCB	679	NO	NO
Septic Shock	DCB	756	NO	NO
Hemorrhagic Stroke	DCB	788	NO	NO
Gastric Perforation	DCB	1055	NO	NO
<b>Unknown</b>				
Sudden Death	DCB	287	NO	NO
Unknown	DCB	541	NO	NO

# IN.PACT SFA Trial

## 3-Year Functional Outcomes

### 6-Minute Walk Test\*




\*Data collected in IN.PACT SFA II phase only

# IN.PACT SFA Trial

## 3-Year Functional Outcomes

### Walking Impairment



# IN.PACT SFA Trial Summary

**Superior 3-year patency and low intervention rates with the  
IN.PACT™ Admiral™ Drug-Coated Balloon**

## **Durability**

- First independently adjudicated, randomized pivotal IDE trial to show durable treatment effect with a DCB through three years

## **Safety**

- Results demonstrate long-term safety of IN.PACT Admiral DCB

## **Paradigm Shift**

- Data support IN.PACT™ Admiral™ as a first-line treatment for symptomatic femoropopliteal disease



# IN.PACT™ Admiral™ Drug-Coated PTA Balloon Catheter

## Brief Statement

### Indications for Use:

The IN.PACT™ Admiral™ Paclitaxel-Coated PTA Balloon catheter is indicated for percutaneous transluminal angioplasty, after appropriate vessel preparation, of de novo, restenotic, or in-stent restenotic lesions with lengths up to 180 mm in superficial femoral or popliteal arteries with reference vessel diameters of 4-7 mm.

### Contraindications

The IN.PACT Admiral DCB is contraindicated for use in:

- Coronary arteries, renal arteries, and supra-aortic/cerebrovascular arteries
- Patients who cannot receive recommended antiplatelet and/or anticoagulant therapy
- Patients judged to have a lesion that prevents complete inflation of an angioplasty balloon or proper placement of the delivery system
- Patients with known allergies or sensitivities to paclitaxel
- Women who are breastfeeding, pregnant or are intending to become pregnant or men intending to father children. It is unknown whether paclitaxel will be excreted in human milk and whether there is a potential for adverse reaction in nursing infants from paclitaxel exposure.

### Warnings

- Use the product prior to the Use-by Date specified on the package.
- Contents are supplied sterile. Do not use the product if the inner packaging is damaged or opened.
- Do not use air or any gaseous medium to inflate the balloon. Use only the recommended inflation medium (equal parts contrast medium and saline solution).
- Do not move the guidewire during inflation of the IN.PACT Admiral DCB.
- Do not exceed the rated burst pressure (RBP). The RBP (14 atm [1419 kPa]) is based on the results of in vitro testing. Use of pressures higher than RBP may result in a ruptured balloon with possible intimal damage and dissection.
- The safety and effectiveness of using multiple IN.PACT Admiral DCBs with a total drug dosage exceeding 20,691 µg of paclitaxel in a patient has not been clinically evaluated in the IN.PACT SFA Trial.

# IN.PACT™ Admiral™ Drug-Coated PTA Balloon Catheter

## Brief Statement

### Precautions

- This product should only be used by physicians trained in percutaneous transluminal angioplasty (PTA).
- This product is designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Assess risks and benefits before treating patients with a history of severe reaction to contrast agents.
- The safety and effectiveness of the IN.PACT Admiral DCB used in conjunction with other drug-eluting stents or drug-coated balloons in the same procedure or following treatment failure has not been evaluated.
- The extent of the patient's exposure to the drug coating is directly related to the number of balloons used. Refer to the *Instructions for Use* (IFU) for details regarding the use of multiple balloons and paclitaxel content.
- The use of this product carries the risks associated with percutaneous transluminal angioplasty, including thrombosis, vascular complications, and/or bleeding events
- Vessel preparation using only pre-dilatation was studied in the clinical study. Other methods of vessel preparation, such as atherectomy, have not been studied clinically with IN.PACT Admiral DCB.
- This product is not intended for the expansion or delivery of a stent

### Potential Adverse Effects

The potential adverse effects (e.g. complications) associated with the use of the device are: abrupt vessel closure; access site pain; allergic reaction to contrast medium, antiplatelet therapy, or catheter system components (materials, drugs, and excipients); amputation/loss of limb; arrhythmias; arterial aneurysm; arterial thrombosis; arteriovenous (AV) fistula; death; dissection; embolization; fever; hematoma; hemorrhage; hypotension/hypertension; inflammation; ischemia or infarction of tissue/organ; local infection at access site; local or distal embolic events; perforation or rupture of the artery; pseudoaneurysm; renal insufficiency or failure; restenosis of the dilated artery; sepsis or systemic infection; shock; stroke; systemic embolization; vessel spasms or recoil; vessel trauma which requires surgical repair.

Potential complications of peripheral balloon catheterization include, but are not limited to the following: balloon rupture; detachment of a component of the balloon and/or catheter system; failure of the balloon to perform as intended; failure to cross the lesion.

Although systemic effects are not anticipated, potential adverse events that may be unique to the paclitaxel drug coating include, but are not limited to: allergic/immunologic reaction; alopecia; anemia; gastrointestinal symptoms; hematologic dyscrasia (including leucopenia, neutropenia, thrombocytopenia); hepatic enzyme changes; histologic changes in vessel wall, including inflammation, cellular damage, or necrosis; myalgia/arthralgia; myelosuppression; peripheral neuropathy.

Refer to the Physician's Desk Reference for more information on the potential adverse effects observed with paclitaxel. There may be other potential adverse effects that are unforeseen at this time.

Please reference appropriate product Instructions for Use for a detailed list of indications, warnings, precautions and potential adverse effects. This content is available electronically at [www.manuals.medtronic.com](http://www.manuals.medtronic.com).

**CAUTION:** Federal law (USA) restricts this device to sale by or on the order of a physician.

[www.peripheral.medtronicendovascular.com/international](http://www.peripheral.medtronicendovascular.com/international)

➤ **[www.aortic.medtronicendovascular.com](http://www.aortic.medtronicendovascular.com)**

**Medtronic**

710 Medtronic Parkway NE  
Minneapolis, MN 55432  
USA  
Tel: +1.763.514.4000

**Medtronic International Trading Sàrl**

Route du Molliau 31  
CH-1131 Tolochenaz  
Switzerland  
Tel: +41.21.802.7000

**Medtronic International Ltd.**

49 Changi South  
Avenue 2  
Singapore 486056  
Tel: +65.6436.5000

**Medtronic of Canada Ltd.**

99 Hereford Street  
Brampton, Ontario L6Y 0R3  
Canada  
Tel: +1.905.460.3800

**Medtronic Nederland**

Earl Bakkenstraat 10  
6422 PJ Heerlen  
The Netherlands  
Tel: +31.45.566.8000

**Medtronic Australasia Pty Ltd.**

97 Waterloo Road  
North Ryde, NSW 2113  
Australia  
Tel: +61.29857.9000

**Medtronic Latin America**

3750 NW 87 Avenue, Suite 700  
Miami, FL 33178  
USA  
Tel: +1.305.500.9328

**Medtronic Ireland**

Parkmore Business Park West  
Galway  
Ireland  
Tel: +353.91.708000



**Medtronic**