

EXTEND-IA

KEY FINDINGS

The EXTEND-IA trial was conducted to test the hypothesis that anterior circulation ischemic stroke patients, selected with a "dual target" of vessel occlusion and evidence of salvageable tissue on perfusion imaging within 4.5h of onset, would have improved reperfusion and early neurological improvement when treated with endovascular thrombectomy using the Solitaire™ stent thrombectomy device after intravenous (IV) alteplase, compared to alteplase alone.



IV t-PA
RANDOMIZED



LESS THAN
5% LOSS TO
FOLLOW UP



INDEPENDENT
DATA REVIEW



An investigator-initiated, multi-center, prospective, randomized, open-label, blinded-endpoint (PROBE) study in ischemic stroke patients receiving intravenous alteplase within 4.5h of stroke onset.

100%

SOLITAIRE™ DEVICE*

0%

100%

PRIMARY ENDOVASCULAR TREATMENT

Medtronic

NO SIGNIFICANT DIFFERENCE IN STUDY POPULATIONS:

CHARACTERISTIC

	INTRAVENOUS t-PA ALONE (n=35)	COMBINED INTRAVENOUS t-PA AND ENDOVASCULAR THERAPY (n=35)*
Number of randomized subjects	35	35
Age (years) - Mean (SD)	70.2 ± 11.8	68.6 ± 12.3
Male sex - no. (%)	17 (49)	17 (49)
BL NIHSS score, Median (IQR)	13 (9-19)	17 (13-20)
History of hypertension - no. (%)	23 (66)	21 (60)
History of diabetes - no. (%)	8 (23)	2 (6)
Time from stroke onset to initiation of alteplase (min) - Median (IQR)	145 (105-180)	127 (93-162)
Site of vessel occlusion - no. (%)		
▪ Internal carotid artery (ICA)	11 (31)	11 (31)
▪ First segment of middle cerebral artery (M1)	18 (51)	20 (57)
▪ Second segment of middle cerebral artery (M2)	6 (17)	4 (11)

STUDY RESULTS:

OUTCOME


	ALTEPLASE- ONLY GROUP (n=35)	ENDOVASCULAR THERAPY GROUP (n=35)*	UNADJUSTED EFFECT SIZE (95% CI)	P VALUE
Primary outcomes:				
Median Reperfusion at 24 hr - (IQR)	37 (-0.5 to 96)	100 (100 to 100)	4.9 (2.5 to 9.5)	<0.001
Early neurologic improvement: Reduction ≥ 8 points on NIHSS or a score of 0-1 at 3 days - no. (%)	13 (37)	28 (80)	6.8 (2.3 to 20)	<0.001
Secondary outcomes:				
Independent outcome (mRS 0-2) at 90 days - no. (%)	14 (40)	25 (71)	3.8 (1.4 to 10.0)	0.009
Final mTICI 2B-3 (in patients who had an initial occlusion angiography - no./total no.%)	N/A	25/29 (86)	N/A	N/A
Death - no. (%)	7 (20)	3 (9)	0.38 (0.1 to 1.6)	0.18
Symptomatic intracerebral hemorrhage - no. (%)	2 (6)	0	-6 (-13 to 2)*	0.49
Home time (number of days spent at home during the first 90 days after the diagnosis of stroke) - median (IQR)	15 (0 to 69)	73 (47 to 86)	58 (17 to 90)	0.006

* The effect size in this category is a risk difference, as measured in percentage points for symptomatic intracerebral hemorrhage and parenchymal hematoma.

Statistically significant improvement in the rate of functional outcomes at 90 days with intervention using the Solitaire™ device as primary treatment.¹

¹ Generalized OR 2.0 (95% CI 1.2-3.8; p=0.006)

■ mRS 0-2 ■ mRS 3-6



STUDY CONCLUSION:

In ischemic stroke patients with a proximal cerebral arterial occlusion and salvageable tissue on CT perfusion imaging, early thrombectomy with the Solitaire™ stent thrombectomy device improves reperfusion, early neurologic recovery and functional outcome compared with alteplase alone.

Medtronic

9775 Toledo Way
Irvine, CA 92618
USA
Tel 877.526.7890
Fax 763.526.7888

medtronic.com/SoliFacts

* Solitaire™ FR Revascularization Device used in Endovascular therapy

SOURCE: Campbell BC, Mitchell PJ, Kleinig TJ, et al. Endovascular therapy for ischemic stroke with perfusion-imaging selection. N. Engl. J. Med. Mar 12 2015;372(11):1009-1018.

CAUTION: Federal (USA) law restricts these devices to sale distribution and use by or on order of a physician. Indications, contraindications, warnings and instructions for use for Solitaire™ X Revascularization Device can be viewed at www.medtronic.com/manuals

© 2019 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. UC201908416 EN