

# 2018 ESC SYNCOPE GUIDELINES

## SUMMARY


# NEW GUIDELINES

OVERVIEW OF UPDATED  
RECOMMENDATIONS SINCE 2009

# 2018 EUROPEAN SOCIETY OF CARDIOLOGY SYNCOPE GUIDELINES

## Goals of 2018 Task Force<sup>1</sup>

- **Reducing Cost & Admissions:** Focus of recommendations is on the need to reduce hospital admissions & unnecessary system costs.
- **Providing Practical Guidance:** Guidelines include clinical recommendations and a separate *practical instructions guide* to drive implementation.
- **Driving Multidisciplinary Approach:** Most multidisciplinary syncope guidelines task force to date – must engage across specialties to better manage syncope patients (ED, cardiology, neurology, internal medicine, geriatrics, nursing).

 ESC European Heart Journal (2018) 00, 1–69  
European Society of Cardiology doi:10.1093/eurheartj/ehy037

ESC GUIDELINES

### 2018 ESC Guidelines for the diagnosis and management of syncope

The Task Force for the diagnosis and management of syncope of the European Society of Cardiology (ESC)

Developed with the special contribution of the European Heart Rhythm Association (EHRA)

Endorsed by: European Academy of Neurology (EAN), European Federation of Autonomic Societies (EFAS), European Federation of Internal Medicine (EFIM), European Union Geriatric Medicine Society (EUGMS), European Society of Emergency Medicine (EuSEM)

Authors/Task Force Members: Michele Brignole\* (Chairperson) (Italy), Angel Moya\* (Co-chairperson) (Spain), Frederik J. de Lange (The Netherlands), Jean-Claude Deharo (France), Perry M. Elliott (UK), Alessandra Fanciulli (Austria), Artur Fedorowski (Sweden), Raffaello Furlan (Italy), Rose Anne Kenny (Ireland), Alfonso Martin (Spain), Vincent Probst (France), Matthew J. Reed (UK), Ciara P. Rice (Ireland), Richard Sutton (Monaco), Andrea Ungar (Italy), and J. Gert van Dijk (The Netherlands)

\* Corresponding authors: Michele Brignole, Department of Cardiology, Ospedale Del Tigullio, Via Don Bobbio 25, IT-16033 Lavagna, (GE) Italy. Tel: +39 0185 329 567, Fax: +39 0185 306 506, Email: mbrignole@univgda.it; Angel Moya, Arrhythmia Unit, Hospital Vall d'Hebron, P Vall d'Hebron 119-129, ES-08035 Barcelona, Spain. Tel: +34 93 2746166, Fax: +34 93 2746002, Email: amoya@vhebrons.com

ESC Committee for Practice Guidelines (CPG) and National Cardiac Societies document reviewers listed in the Appendix.

<sup>1</sup> Representing the European Academy of Neurology (EAN)

<sup>2</sup> Representing the European Federation of Internal Medicine (EFIM)

<sup>3</sup> Representing the European Society of Emergency Medicine (EuSEM)

ESC articles having participated in the development of this document:  
Associations: European Heart Rhythm Association (EHRA)  
Councils: Council on Cardiovascular Nursing and Allied Professions, Council for Cardiology Practice, Council on Cardiovascular Primary Care Working Group: Myocardial and Pericardial Diseases

The content of these European Society of Cardiology (ESC) Guidelines has been published for personal and educational use only. No commercial use is authorized. No part of the ESC Guidelines may be translated or reproduced in any form without written permission from the ESC. Permission can be obtained upon submission of a written request to: Oxford University Press, the publisher of the European Heart Journal and the party authorized to handle such permissions on behalf of the ESC (journalspermissions@oxfordjournals.org).


**Disclaimer.** The ESC Guidelines represent the views of the ESC and were produced after careful consideration of the scientific and medical knowledge and the evidence available at the time of their publication. The ESC is not responsible in the event of any contradiction, discrepancy and/or ambiguity between the ESC Guidelines and any other official recommendations or guidelines issued by the relevant public health authorities, in particular in relation to good use of healthcare or therapeutic strategies. Health professionals are encouraged to take the ESC Guidelines fully into account when exercising their clinical judgment, as well as in the determination and the implementation of preventive, diagnostic or therapeutic medical strategies; however, the ESC Guidelines do not override, in any way whatsoever, the individual responsibility of health professionals to make appropriate and accurate decisions in consideration of each patient's health condition and in consultation with that patient and, where appropriate and/or necessary, the patient's caregiver. Nor do the ESC Guidelines exempt health professionals from taking into full and careful consideration the relevant official updated recommendations or guidelines issued by the competent public health authorities, in order to manage each patient's case in light of the scientifically accepted data pursuant to their respective ethical and professional obligations. It is also the health professional's responsibility to verify the applicable rules and regulations relating to drugs and medical devices at the time of prescription.

© The European Society of Cardiology 2018. All rights reserved. For permissions please email: journalspermissions@oxfordjournals.org

1. Brignole et al Eur Heart J 2018

# 2018 EUROPEAN SOCIETY OF CARDIOLOGY SYNCOPE GUIDELINES

## New Concepts in 2018 Guidelines<sup>1</sup>



### Extensive changes based on new evidence:

- Increased role of prolonged ECG monitoring with ILR
- Revised recommendation from Class I to Class II for Holter Monitors & Tilt Tests
- Increased importance of risk stratification from ED and referral to syncope specialist

1. Brignole et al Eur Heart J 2018, Figure I in Guidelines Introduction

# GUIDELINES RECOMMENDATION CLASSIFICATIONS<sup>1</sup>

Classes of recommendations	Definition	Suggested wording to use
<b>Class I</b>	Evidence and/or general agreement that a given treatment or procedure is beneficial, useful, effective.	Is recommended/Is indicated
<b>Class II</b>	Conflicting evidence and/or a divergence of opinion about the usefulness/efficacy of the given treatment or procedure.	
<i>Class IIa</i>	<i>Weight of evidence/opinion is in favour of usefulness/efficacy.</i>	Should be considered
<i>Class IIb</i>	<i>Usefulness/efficacy is less well established by evidence/opinion.</i>	May be considered
<b>Class III</b>	Evidence or general agreement that the given treatment or procedure is not useful/effective; and in some cases may be harmful.	Is not recommended

<b>Level of evidence A</b>	Data derived from multiple randomized clinical trials or meta-analyses.
<b>Level of evidence B</b>	Data derived from a single randomized clinical trial or large non-randomized studies.
<b>Level of evidence C</b>	Consensus of opinion of the experts and/or small studies, retrospective studies, registries.

©ESC 2018

**Class I, Level of Evidence “A” is HIGHEST Recommendation**

1. Brignole et al Eur Heart J 2018, Table 1 and Table 2 in Guidelines Introduction

# COMPARISON OF CHANGES IN ILR RECOMMENDATIONS

## 2009 Guidelines<sup>1</sup>

### Class I Level B

- Indicated in early phase evaluation in patients with recurrent syncope of uncertain origin, absence of high risk criteria and a high likelihood of recurrence within battery longevity of the device
- Indicated in high risk individuals in whom comprehensive evaluation did not demonstrate a cause of syncope or lead to a specific treatment

### Class IIa Level B

- Considered in patients with suspected or certain reflex syncope presenting with frequent or traumatic syncopal episodes

**NEW!**

**NEW!**

1. Moya A. et al Eur Heart J 2009

2. Brignole et al Eur Heart J 2018

## 2018 Guidelines<sup>2</sup>

### Class I Level A

- Indicated in early phase evaluation in patients with recurrent syncope of uncertain origin, absence of high risk criteria and a high likelihood of recurrence within battery longevity of the device
- Indicated in patients with high-risk criteria in whom a comprehensive evaluation did not demonstrate a cause of syncope or lead to a specific treatment, and who do not have conventional indications for primary prevention ICD or pacemaker

### Class IIa Level B

- Considered in patients with suspected or certain reflex syncope presenting with frequent or severe syncopal episodes

### Class IIb Level B


Rule in or out arrhythmias in:

- Patients in whom epilepsy is suspected
- Patients with unexplained falls

### Class IIb Level C

ILR may be considered in patients with recurrent episodes of unexplained syncope who are at low risk for SCD and without a current indication for ICD

# T-LOC MONITORING RECOMMENDATIONS



Adapted from Ungar A: ECG monitoring Role in Syncope 2018 ESC Guidelines, EHRA 2018  
\* High Risk & Low Risk Recommendations Summarized on Slides 18 – 19

# EVIDENCE SUPPORT

RECOMMENDATIONS FOR DIAGNOSING  
UNEXPLAINED SYNCOPE



# EVIDENCE SUPPORT IN GUIDELINES

*ILR is indicated in an early phase of evaluation in patients with recurrent syncope of uncertain origin, absence of high-risk criteria, and a high likelihood of recurrence within the battery life of the device.*

**2009**

**Class I  
LOE: B**



**2018**

**Class I  
LOE: A**


- **NEW – Meta-analysis:** 5 randomized clinical trials comparing diagnostic yield of ICM monitoring to standard of care<sup>1</sup>

Study	Diagnostic yield		Relative probability	95% CI	P value
	ICM group n/N (%)	Control group n/N (%)			
RAST 2001	14/27 (52)	6/30 (20)	2.6	1.2-5.8	<b>0.01</b>
EaSyAS 2006	43/101 (43)	7/97 (7)	5.9	2.8-12	<b>0.001</b>
Da Costa 2013	15/41 (37)	4/37 (11)	3.4	1.2-9.3	<b>0.01</b>
FRESH 2014	18/39 (46)	2/39 (5)	9.0	2.2-3.6	<b>0.001</b>
EaSyAS II 2016	62/125 (50)	21/121 (17)	2.9	1.9-4.4	<b>0.001</b>
<b>Total</b>	<b>152/333 (46)</b>	<b>40/324 (12)</b>	<b>3.6</b>	<b>2.4-5.3</b>	<b>0.001</b>

**3.6X** more likely to reach a diagnosis with ICM vs standard of care

- **NEW – Large non-randomized clinical trials**

## PICTURE<sup>2</sup>



**570 patients**

## SPANISH REGISTRY<sup>3</sup>





**680 patients**

**Strong evidence supports efficacy of early use of ICM in unexplained, recurrent syncope**

1. Brignole et al Eur Heart J 2018  
 2. Edvardsson et al Clin Cardiol, 2011  
 3. Lacumza-Ruiz et al Circ J, 2013

# EVIDENCE SUPPORT IN GUIDELINES

*ILR is indicated in patients with high-risk criteria in whom a comprehensive evaluation did not demonstrate a cause of syncope or lead to a specific treatment, and who do not have conventional indications for primary prevention ICD or pacemaker indication.*



Pooled analysis of 9 studies (n=509) – diagnostic yield after “extensive” work up<sup>1</sup>

- Diagnostic yield **35%** (176/506)
- **67%** had arrhythmia associated with syncopal event

ILR monitoring in BBB patients w/ negative Electrophysiologic Study (EPS) to detect paroxysmal AVB

- NEW: DaCosta '13 – RCT
- NEW: Moya '11 – Large prospective trial (n=323)

ILR diagnostic yield in syncope patients with BBB<sup>2</sup>

Study	N	ILR documented event	ILR documented arrhythmia	ILR documented AVB
Brignole 2001	52	24	22	12
Moya 2011	108	52	45	36
Da Costa 2013	41	15	15	11
<b>Total</b>	<b>201</b>	<b>91 (45%)</b>	<b>82 (41%)</b>	<b>59 (29%)</b>

**Strong evidence to support ICM monitoring if comprehensive testing is inconclusive or does not meet current indications for a pacemaker**

# EVIDENCE SUPPORT IN GUIDELINES

*ILRs should be considered in patients with suspected or certain reflex syncope presenting with frequent or severe syncopal episodes.*

**2009**

Class IIa  
LOE: B




**2018**

Class IIa  
LOE: B

## ILR monitoring in reflex syncope patients to identify those that benefit from cardiac pacing<sup>1-3</sup>

- NEW: Brignole '12 (ISSUE – 3) – RCT - Pacing on/off in reflex syncope patients with asystole (ICM detected)

**ISSUE-3: Reduction in recurrent syncope**



- Dual-chamber pacemakers implanted in reflex syncope patients with ILR detected asystole (randomized to “on” or “off” pacing)
- 57% reduction in risk of recurrent syncope with pacing

- NEW: Brignole '15 (SUP2) – single arm – ICM detection as part of a diagnostic algorithm (asystole detected in reflex syncope patients response to pacing)
- Brignole '06 ISSUE2 – Therapy based on ICM monitoring effective in Neurally Mediated Syncope)

**Evidence is building for the need to identify asystole in reflex syncope patients as pacing therapy may be effective**

# EVIDENCE SUPPORT IN GUIDELINES

*ILR is indicated in an early phase of evaluation in patients with recurrent syncope of uncertain origin, absence of high-risk criteria, and a high likelihood of recurrence...*

*Holter monitoring should be considered in patients who have frequent syncope or presyncope (>1 episode per week).*

**2009**

**Class I  
LOE: B**

**2018**

**Class I  
LOE: A**



**Class I  
LOE: B**


**Class IIa  
LOE: B**



## 2018 ESC Syncope Guidelines Task Force:

- “Holter monitoring in syncope is inexpensive in terms of set-up costs, but **expensive in terms of cost per diagnosis.**”<sup>1</sup>
- “ILR was **more cost-effective** than a conventional [monitoring] strategy”<sup>2-5</sup>
- NEW: Edvardsson '11 and '15 PICTURE (n=570) – prospective registry study

### Cost of tests performed per patient prior to ICM implant (£)



- Prior to ILR implant, patients saw an average of 3 specialists and received a median 13 tests (mean cost of testing: £1,613)
  - 10% of patients had tests exceeding £3,540
- ILR guided diagnosis in 78% of patients with syncopal recurrence

- Krahn '03 RAST – Cost per diagnosis was lower in the ILR group (\$5,852 CAD) vs. Conventional testing (\$8,414 CAD)
- Farwell '04 EasyAS – Earlier diagnosis and treatment in the ILR group resulted in 67% lower costs for subsequent diagnostic testing and hospitalization (ILR: £406 vs. Conventional testing: £1,210)

# EVIDENCE SUPPORT

EXPANDED RECOMMENDATIONS FOR  
DIAGNOSING T-LOC

# EVIDENCE SUPPORT IN GUIDELINES


**NEW 2018**

**Class IIb  
LOE: B**

*ILR may be considered in patients with **unexplained falls**.*

- “Despite the lack of controlled trials and an overall modest quality of studies, there is strong consensus that the management of unexplained falls should be the same as that for unexplained syncope”<sup>1</sup>

**Suggested flow diagram to identify unexplained falls**



- Pooled analysis of 5 studies with an ILR in the unexplained falls population<sup>1</sup>

Study	ILR patients	ILR documented event	ILR documented arrhythmia
Armstrong 2003	6	3 (50%)	1 (15%)
Ryan 2010	71	48 (68%)	3 (4%)
Maggi 2014	29	16 (55%)	7 (24%)
Bhangu 2016	70	56 (80%)	14 (20%)
<b>Total</b>	<b>176</b>	<b>123 (70%)</b>	<b>25 (14%)</b>

- 70% of patients had a recurrent fall while monitored with an ILR**
- 14% of patients had an arrhythmia response for fall**

**Evidence is growing to support prolonged monitoring in patients with unexplained falls**

1. Brignole et al Eur Heart J 2018

# EVIDENCE SUPPORT IN GUIDELINES

**NEW 2018**

*ILR may be considered in patients in whom **epilepsy was suspected** but the treatment has proven ineffective.*

**Class IIb  
LOE: B**

- “Syncope and Epilepsy have similar clinical features, including abrupt T-LOC without warning, abnormal limb movements such as myoclonic jerks or tonic-clonic activity”<sup>1</sup>
- Pooled analysis of 6 studies with an ILR in the unconfirmed epilepsy population<sup>1</sup>

Study	ILR patients	ILR documented event	ILR documented arrhythmia
Simpson 2000	1	1 (100%)	0 (0%)
Kanjwal 2009	3	3 (100%)	3 (100%)
Zaidi 2000	10	2 (20%)	2 (20%)
Ho 2006	14	6 (43%)	0 (0%)
Petkar 2012	103	69 (67%)	28 (27%)
Maggi 2014	28	17 (61%)	8 (29%)
<b>Total</b>	<b>201</b>	<b>98 (62%)</b>	<b>41 (26%)</b>

- **62% of patients had a T-LOC event captured during ILR monitoring**
- **1 in 4 (26%) of patients had an arrhythmia response for T-LOC**

**ILR monitoring in the unconfirmed epileptic population can lead to appropriate diagnosis and therapy for T-LOC**

1. Brignole et al Eur Heart J 2018

# EVIDENCE SUPPORT IN GUIDELINES

**NEW 2018**

*Instead of an ICD, an ILR may be considered in patients with recurrent episodes of unexplained syncope who are at **low risk for SCD** and without current indication for ICD.*

**Class IIb  
LOE: C**

- Unexplained syncope in this population is defined as “suspected arrhythmic syncope”
- No studies are available to support this specific guideline (LOE: C – *expert opinion*)
- ILR are recommended to complete the diagnostic investigation in those not indicated for therapy devices and at “low risk” for SCD<sup>1</sup>
  - Structural heart disease is the most important predictive factor for a cardiac cause of syncope, with a sensibility of >90% and a specificity of ~50%
  - In many patients with CMP/IAD, individual risk stratification can be difficult due to atypical symptoms

**While there is an absence of specific trials to support long term cardiac monitoring in this population, the expert opinion is to continue cardiac monitoring in low-risk, non-ICD indicated patients**

1. Brignole et al Eur Heart J 2018



# 2018 ESC SYNCOPE GUIDELINES: CONCLUSIONS

- The increasing need to reduce inappropriate admissions & tests, and new evidence led to creation of the 2018 ESC Syncope Guidelines<sup>1</sup>
- ILR recommendation upgraded from a Class I/Level B to a Class I/Level A for unexplained syncope as meta-analysis shows 3.6X more patients were diagnosed compared to conventional tests<sup>1</sup>
- Holter monitors and tilt testing recommendations were revised from Class II to Class I due to lack of new evidence and cost implications<sup>1-5</sup>
- ILR recommendation expanded to patients with unexplained falls, unproven epilepsy, and patients at low-risk for SCD<sup>1</sup>

1. Brignole et al Eur Heart J 2018
2. Edvardsson N et al Europace 2011; 13:262-269,
3. Farwell DJ et al. Eur, Heart J 2004; 25:1287-1263
4. Krahn AD et al. J Am Coll Cardiol 2003; 42:495-501
5. Edvardsson N et al. Europace 2015; 17:1141-1148

# \*RISK ASSESSMENT RECOMMENDATIONS<sup>1</sup>

## Assessment of Syncopal Event

Low-risk	High-risk (red flag)
<ol style="list-style-type: none"> <li>1. Associated with prodrome typical of reflex syncope (e.g. light-headedness, feeling of warmth, sweating, nausea, vomiting)</li> <li>2. After unexpected unpleasant sight, sound, smell, or pain</li> <li>3. After prolonged standing or crowded, hot places</li> <li>4. During a meal or postprandial</li> <li>5. Triggered by cough, defecation, or micturition</li> <li>6. With head rotation or pressure on carotid sinus</li> <li>7. Standing from supine/sitting position</li> </ol>	<p><b>Major</b></p> <ol style="list-style-type: none"> <li>1. New onset of chest discomfort, breathlessness, abdominal pain, or headache</li> <li>2. Syncope during exertion or when supine</li> <li>3. Sudden onset palpitation immediately followed by syncope</li> </ol> <p><b>Minor</b> (high risk only if associated with SHD or abnormal ECG):</p> <ol style="list-style-type: none"> <li>1. No warning symptoms or short (&lt;10 s) prodrome</li> <li>2. Family history of SCD at young age</li> <li>3. Syncope in the sitting position</li> </ol>

## Past Medical History

Low-risk	High-risk (red flag)
<ol style="list-style-type: none"> <li>1. Long history of recurrent syncope with low-risk features with the same characteristics of the current episode</li> <li>2. Absence of structural heart disease</li> </ol>	<p><b>Major</b></p> <ol style="list-style-type: none"> <li>1. Severe structural or coronary artery disease (heart failure, low LVEF or previous myocardial infarction)</li> </ol>

1. Brignole et al Eur Heart J 2018

# \*RISK ASSESSMENT RECOMMENDATIONS<sup>1</sup>

## Physical Examination

Low-risk	High-risk (red flag)
1. Normal examination	<p><b>Major</b></p> <ol style="list-style-type: none"> <li>1. Unexplained systolic BP in the ED &lt;90 mmHg</li> <li>2. Suggestion of gastrointestinal bleed on rectal examination</li> <li>3. Persistent bradycardia (&lt;40 b.p.m.) in awake state and in absence of physical training</li> <li>4. Undiagnosed systolic murmur</li> </ol>
1. Normal ECG	<p><b>Major</b></p> <ol style="list-style-type: none"> <li>1. ECG changes consistent with acute ischaemia</li> <li>2. Mobitz II second- and third-degree AV block</li> <li>3. Slow AF (&lt;40 b.p.m.)</li> <li>4. Persistent sinus bradycardia (&lt;40 b.p.m.)</li> <li>5. Bundle branch block or IVCD</li> <li>6. Q waves consistent with CAD or cardiomyopathy</li> <li>7. Sustained and non-sustained VT</li> <li>8. Dysfunction of a pacemaker or ICD</li> <li>9. Type 1 Brugada pattern</li> <li>10. Long QT</li> </ol> <p><b>Minor</b></p> <ol style="list-style-type: none"> <li>1. ECG changes consistent with acute ischaemia</li> <li>2. Mobitz II second- and third-degree AV block</li> <li>3. Slow AF (&lt;40 b.p.m.)</li> <li>4. Persistent sinus bradycardia (&lt;40 b.p.m.)</li> <li>5. Bundle branch block or IVCD</li> <li>6. Q waves consistent with CAD or cardiomyopathy</li> <li>7. Sustained and non-sustained VT</li> <li>8. Dysfunction of a pacemaker or ICD</li> <li>9. Type 1 Brugada pattern</li> <li>10. Long QT</li> </ol>

1. Brignole et al Eur Heart J 2018

## Brief statement

See the device manual for detailed information regarding the instructions for use, the implant procedure, indications, contraindications, warnings, precautions, and potential adverse events. If using an MRI SureScan® device, see the MRI SureScan® technical manual before performing an MRI. For further information, contact your local Medtronic representative and / or consult the Medtronic website at [www.medtronic.com](http://www.medtronic.com).



[www.medtronic.com/manuals](http://www.medtronic.com/manuals)

Consult instructions for use at this website. Manuals can be viewed using a current version of any major Internet browser. For best results, use Adobe Acrobat Reader® with the browser.

# Medtronic

## Europe

Medtronic International Trading  
Sàrl.  
Route du Molliau 31  
Case postale  
CH-1131 Tolochenaz  
[www.medtronic.eu](http://www.medtronic.eu)  
Tel: +41 0 21 802 70 00  
Fax: +41 0 21 802 79 00

## United Kingdom/Ireland

Medtronic Limited  
Building 9  
Croxley Park  
Hatters Lane  
Watford  
Herts WD18 8WW  
[www.medtronic.co.uk](http://www.medtronic.co.uk)  
Tel: +44 0 1923 212213  
Fax: +44 0 1923 241004

[medtronic.eu](http://medtronic.eu)

UC201708879EE ©2017 Medtronic.  
All Rights Reserved.

Medtronic, Medtronic logo, and Further, Together are trademarks of Medtronic.  
™Third party brands are trademarks of their respective owners.  
All other brands are trademarks of a Medtronic company.