

Your clinic would like to do a check of your device the morning after your implant.



**Tomorrow morning** \_\_\_\_\_

### **FOLLOW THESE STEPS:**

- 1.** Open your MyCareLink Heart™ app and select the "My Transmissions" tile.
- 2.** Select "SEND TRANSMISSION."
- 3.** Select "YES" when asked "Does your clinic know that you are going to send a transmission?"
- 4.** Select "OK" when shown "Transmission in progress" and "Your transmission is being sent to your clinic."

Go about your day, and simply keep your phone or tablet with you until you receive a notification that your transmission was successfully received.

**YOU'RE ALL SET!**

**Medtronic**

## Your MyCareLink Heart app keeps you connected to your clinic going forward. Here are some helpful tips for remaining connected:

- ✓ Stay connected to internet via Wi-Fi or cellular data.
- ✓ Keep your smartphone or tablet near you for a minimum of three hours per day.
- ✓ Keep your phone on, the app running in the background. If you turn your phone off, simply tap on the MyCareLink Heart app icon to relaunch it.
- ✓ Be sure your mobile device and software are current with the latest updates.

**Note:** All future transmissions will be sent automatically. No need to send another manual transmission unless asked to do so by your clinic.

**Call Medtronic Patient Services for any questions regarding your MyCareLink Heart app:  
1-888-660-4616**