

PROVEN TO **PERFORM**

Five years of
clinical data for the
InterStim™ system



Medtronic

BACKED BY SCIENCE, **DRIVEN TO INNOVATION**

Over twenty years ago, Medtronic introduced Sacral Neuromodulation (SNM) therapy and transformed the treatment options for incontinence. Today, we continue our commitment to Pelvic Health and improving patients lives.

Partnering with Medtronic provides your patients with...

RELIABILITY



250,000
INTERSTIM™
PATIENTS
TREATED
WORLDWIDE

CONFIDENCE




GLOBAL
PRESENCE
TO SUPPORT
INTERSTIM™
THERAPY

FREEDOM



MORE THAN 80%
SATISFIED WITH
INTERSTIM™
THERAPY^{1,2}



BLADDER CONTROL PROVEN AT 5 YEARS³

82%

of Overactive Bladder
(OAB) patients
achieved clinical
success* at 5 years.†

45%


of Urinary Incontinence
patients were completely
continent at 5 years.†

The most common
adverse events were:
undesirable change in
stimulation, implant site
pain, and therapeutic
product ineffective.

* Clinical success was defined as $\geq 50\%$
improvement in average leaks or voids/day
or a return to normal voiding (<8 voids/day).


† Completers analysis

INTERSTIM™ DEMONSTRATED SUSTAINED LONG-TERM EFFICACY†³



† Numbers reflect completers analysis, which includes all implanted subjects with diary data at baseline and at each follow-up. Clinical success was 83% at 12 months, 76% at 36 months and 67% at 5 years using the modified completers analysis (subjects who either had baseline and follow-up evaluation or withdrew early due to device-related reasons and are considered failures).

INTERSTIM™ DEMONSTRATED SUSTAINED QUALITY OF LIFE (QOL) IMPROVEMENTS AT 5 YEARS³




*Coyne KS, et al. 2006. Determining the importance of change in the OAB-q. J Urol 176:627-32. All paired tests comparing follow-up to baseline had a p-value <0.0001.

PROVEN SAFE IN A WIDE RANGE OF PATIENTS^{4,5}

- SNM using the InterStim™ system is a safe and effective treatment for patients with OAB and non-obstructive urinary retention who have failed or are not candidates for more conservative treatments.
- InterStim™ demonstrated sustained and significant improvements in quality of life for patients with OAB and non-obstructive urinary retention at 5 years.

The most common adverse events experienced during clinical studies included pain at implant sites, new pain, lead migration, infection, technical or device problems, adverse change in bowel or voiding function, and undesirable stimulation or sensations. Any of these may require additional surgery or cause return of symptoms.

Non-obstructive Urinary Retention




Clinical success was defined as 50% or greater improvement from baseline in primary voiding diary variable(s) or subjective patient satisfaction.

**SIGNIFICANT
SYMPTOM
IMPROVEMENT
DEMONSTRATED
AT 5 YEARS**




SNM with InterStim™ therapy⁸: established by international guidelines

SPECIALIZED MANAGEMENT OF URINARY INCONTINENCE IN WOMEN: SURGICAL TREATMENT



ICI* 2017: Surgical management of fecal incontinence⁹



ACE: Antegrade Continence Enema, **BI:** Biomaterial injection
SNS: Sacral Nerve Stimulation, *Not widely available.

Novel therapies

- Magnetic anal sphincter
- Puborectal sling
- Radiofrequency energy treatment
- Stem cell therapy
- Vaginal pessary - Eclipse™

BOWEL CONTROL PROVEN AT 5 YEARS¹⁰

The most common adverse events through 5 years ($\geq 5\%$ of patients, n=120) are implant site pain, parasthesia, change in sensation of stimulation, implant site infection, urinary incontinence, neurostimulator battery depletion, diarrhea, pain in extremity, undesirable change in stimulation, and buttock pain.

89%

of bowel patients achieved clinical success* at 5 years.^{10†}


36%

of bowel patients achieved complete continence.^{10†}

* Clinical success was defined as $\geq 50\%$ reduction of FI episodes/week.


† Completers analysis

INTERSTIM™ DEMONSTRATED SUSTAINED LONG-TERM EFFICACY FOR FECAL INCONTINENCE¹⁰



Completers analysis with the proportion of patients with therapeutic success ($\geq 50\%$ improvement)

INTERSTIM™ DEMONSTRATED SUSTAINED QOL IMPROVEMENTS AT 5 YEARS¹⁰



Study data demonstrate significant improvement in patient quality of life for patients with chronic FI, including physical and psychological well-being, as determined by a variety of accepted measures.³ Mean FIQOL component scales at each visit for patients with at least 5 years of follow-up (higher score indicates better quality of life).

* Fecal Incontinence Quality of Life (FIQOL)

More than
250,000
patients treated
globally.

See the device manual for detailed information regarding the instructions for use, implant procedure, indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative or consult the Medtronic website at www.medtronic.ca

References

1. Foster et al. In patients undergoing neuromodulation for intractable urge incontinence a reduction in 24-hr pad weight after the initial test stimulation best predicts long-term patient satisfaction. *Neurourol Urodyn.* 2007;26(2):213-7.
2. van Wunnik BP, Govaert B, Leong R, Nieman FH, Baeten CG. Patient experience and satisfaction with sacral neuromodulation: results of a single-center sample survey. *Dis Colon Rectum.* 2011 Jan;54(1):95-100.
3. Siegel S, Noblett K, Mangel J, et al. Five Year Follow-up Results of a Prospective, Multicenter Study in Overactive Bladder Subjects Treated with Sacral Neuromodulation. *J Urol.* 2017. doi: 10.1016/j.juro.2017.07.010
4. van Voskuilen AC, Oerlemans DJ, Weil EH, de Bie RA, van Kerrebroeck PE. Long term results of neuromodulation by sacral nerve stimulation for lower urinary tract symptoms: a retrospective single center study. *Eur Urol.* 2006 Feb;49(2):366-72.
5. van Kerrebroeck PE, van Voskuilen AC, Heesakkers JP et al. Results of sacral neuromodulation therapy for urinary voiding dysfunction: outcomes of a prospective, worldwide clinical study. *J Urol.* 2007 Nov;178(5):2029-34.
6. White WM, Dobmeyer-Dittrich C, Klein FA, Wallace LS. Sacral nerve stimulation for treatment of refractory urinary retention: long-term efficacy and durability. *Urology.* 2008 Jan;71(1):71-4.
7. Peeters K, Sahai A, De Ridder D, Van Der Aa F. Long-term follow-up of sacral neuromodulation for lower urinary tract dysfunction. *BJU Int.* 2014 May;113(5):789-94.
8. EAU guidelines 2018 :uroweb.org/guideline/urinary-incontinence/#4
9. Herold A., Lehur P.A., Matzel K.E., O'Connell P.R. (eds.). *Coloproctology*, Springer 2017, DOI 10.1007/978-3-662-53210-2_9.
10. Hull T, Giese C, Wexner SD, et al. Long-term durability of sacral nerve stimulation therapy for chronic fecal incontinence. *Dis Colon Rectum.* 2013;56(2):234-245.

Medtronic

Medtronic Canada

16720 Trans-Canada, Suite 200
Kirkland, QC, H9H 4M7
877.664.8926 [t]
800.567.1939 [f]

medtronic.ca

UC201804916 EC © 2021 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.™* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. CA-SI-0455-E Rev. 2021/08