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# Exacta™ Medical Device Correction Letter Training

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### Exacta™ Medical Device Correction Letter Exacta™ 医疗器械纠正措施函

#### Issue Description 问题描述

Medtronic has received customer complaints reporting cracks and/or leaks on the stopcock of Becker™ and Exacta™ EDMS devices.

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On these devices, the stopcocks may be in **three** different locations in the EDMS, depending on the specific system configuration.

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

Cracks or leaks of the stopcock can compromise line integrity and, as a result, present a potential for infection.

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


## Exacta™ Medical Device Correction Letter    Exacta™ 医疗器械纠正措施函

### Recommended Mitigations    建议的缓解措施

1. Prior to use, inspect all stopcocks and connections to ensure that connections are secure and that there are no visible cracks in the stopcocks.  
使用前，检查所有三通阀及其连接，确保连接牢固且三通阀无可见裂痕。

2. If cracks or leaks are identified, do not use the device and return it to Medtronic.

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\*All images are examples of potential cracks.  
所有图像均为潜在裂痕的示例

Images provided by Medtronic Quality Engineering. 美敦力质量工程团队提供的图像

3




Medtronic

## Exacta™ Medical Device Correction Letter    Exacta™ 医疗器械纠正措施函

### Recommended Mitigations    建议的缓解措施

3. The system must be pre-filled with sterile isotonic saline solution prior to connecting to the patient.  
在将系统连接至患者之前，系统必须先用无菌等渗生理盐水预充。

4. Once filled, **check all connections** to ensure that fittings are connected and that there is no leak-free.  
预充后，请检查所有连接，确保连接紧固且无泄漏。

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## Exacta™ Medical Device Correction Letter Exacta™ 医疗器械纠正措施函

### Recommended Mitigations 建议的缓解措施

5. All connections should be **finger tightened**. Do not use any instruments to tighten.

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
Over tightening can cause cracks and leaks to occur.

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6. Whenever cleaning with alcohol, or a disinfectant containing alcohol, **allow cleaned components to air dry completely** prior to connecting the system.

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## Exacta™ Medical Device Correction Letter Exacta™ 医疗器械纠正措施函

### Patient Management Recommendations 患者管理建议

- If the system develops cracks or leaks during use, replace it using sterile technique and return the damaged system to Medtronic. If the leak occurs with the Patient Line Stopcock (see image on right), consider application of a hemostat or other clamp to the **proximal** patient line while coordinating a replacement if occlusion of the patient line does not create a risk to the patient.

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- As stated in the product labeling, all patients with EVDs should be monitored for evidence of infection. If a system is found to have cracks or leaks, the system should be replaced, and the patient should continue to be monitored for evidence of infection.

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## Exacta™ Medical Device Correction Letter Exacta™ 医疗器械纠正措施函

### Customer Required Actions 客户行动要求

- Ensure customers complete and return the customer confirmation form to acknowledge that they have read and understand the letter.  
确保客户完成并返回客户确认函，以确认他们已经阅读并理解客户通知函的内容。  
o Form within FCA letter 主动召回通知函随附的确认函
- Ensure customers post a copy of this notification near impacted product as a reminder of the issue and recommended mitigating actions.  
确保客户在受影响产品附近张贴本通知函的复印件以提醒该问题的存在及建议的缓解措施。
- Ensure customers retain a copy of this letter and completed response form for their records.  
确保客户保留一份此通知函及填写完整的确认函以供记录。

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