

Medtronic

Cryptogenic Stroke and Underlying Atrial Fibrillation (CRYSTAL-AF)¹



Objectives


- To assess whether a long-term cardiac monitoring strategy with an insertable cardiac monitor (ICM) is superior to standard monitoring for the detection of atrial fibrillation (AF) in patients with cryptogenic stroke at six months (primary endpoint) and 12 months follow-up (secondary endpoint).
- Determine actions taken after patient is diagnosed with AF.

Study overview

- 441 patients
- 1:1 randomized trial comparing the yield of AF detection through continuous monitoring with Reveal LINQ[™] ICM versus standard medical care in cryptogenic or transient ischemic attack (TIA) patients.

Results

- The Reveal LINQ ICM is superior to standard medical care for the detection of AF in patients with a cryptogenic stroke
- At three years, AF was detected at a rate of 30% in the ICM arm versus 3% in the standard follow-up arm
- 30-day monitoring would not be sufficient in this patient population
 - Median time to AF detection was 84 days over 12 months of follow-up
- Patients had sustained periods of AF and physicians took action
 - 92% of patients in the ICM arm had a longest daily burden of AF of > six minutes
 - Vast majority of patients (97%) who had AF detected were prescribed OAC



8.8x more AF detected at 36 months: 30% in ICM group versus 3% in control

See back side for definitions, and inclusion and exclusion criteria.

Definitions

AF: An episode of irregular heart rhythm, without detectable P-waves, of at least 30 seconds duration. AF episodes that qualified for analysis were adjudicated by an independent committee.

Patient inclusion criteria

- ≥ 40 years of age
- Cryptogenic stroke (or clinical TIA), with infarct seen on MRI or CT, within the previous 90 days; and no mechanism (including AF) determined after:
 - 12-lead ECG
 - 24-hour ECG monitoring (e.g., Holter)
 - Transesophageal echocardiography (TEE)
 - CTA or MRA of head and neck to rule out arterial source
 - Screening for hypercoagulable states in patients < 55 years old

Patient exclusion criteria

- History of AF or atrial flutter
- Permanent indication or contraindication for anticoagulation
- Indication for pacemaker or ICD

Continuous monitoring with ICM is **superior to standard medical care** for the detection of AF after cryptogenic stroke.

Conclusion:

“Atrial fibrillation was more frequently detected with an ICM than with conventional follow-up in patients with a recent cryptogenic stroke. Atrial fibrillation after cryptogenic stroke was most often asymptomatic and paroxysmal and thus unlikely to be detected by strategies based on symptom-driven monitoring or intermittent short-term recordings.”

Reference

¹ Sanna T, Diener HC, Passman RS, et al. Cryptogenic Stroke and Underlying Atrial Fibrillation (CRYSTAL AF). *N Engl J Med.* June 26, 2014;370(26):2478-2486.

Reveal LINQ™ Insertable Cardiac Monitor

Indications

The Reveal LINQ Insertable Cardiac Monitor (ICM) is an implantable patient-activated and automatically-activated monitoring system that records subcutaneous ECG and is indicated in the following cases:

- patients with clinical syndromes or situations at increased risk of cardiac arrhythmias
- patients who experience transient symptoms such as dizziness, palpitation, syncope, and chest pain, that may suggest a cardiac arrhythmia

The device has not been tested specifically for pediatric use.

Contraindications

There are no known contraindications for the implant of the Reveal LINQ ICM. However, the patient's particular medical condition may dictate whether or not a subcutaneous, chronically implanted device can be tolerated.

Warnings/Precautions

Patients with the Reveal LINQ ICM should avoid sources of diathermy, high sources of

radiation, electrosurgical cautery, external defibrillation, lithotripsy, therapeutic ultrasound, and radiofrequency ablation to avoid electrical reset of the device, and/or inappropriate sensing as described in the Medical procedure and EMI precautions manual. MRI scans should be performed only in a specified MR environment under specified conditions as described in the Reveal LINQ MRI Technical Manual.

Potential Complications

Potential complications of the Reveal LINQ device include, but are not limited to, device rejection phenomena (including local tissue reaction), device migration, infection, and erosion through the skin.

See the device manuals for detailed information regarding the implant procedure, indications, contraindications, warnings, precautions, and potential complications/adverse events. For further information, please call Medtronic at 1-800-328-2518 and/or consult the Medtronic website at medtronic.com.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Medtronic
710 Medtronic Parkway
Minneapolis, MN 55432-5604
USA

Toll-free in USA: 800.633.8766
Worldwide: +1.763.514.4000

UC201405904c EN ©2023 Medtronic.
Minneapolis, MN. All Rights Reserved.
Printed in USA. 01/2023

medtronic.com

Medtronic