

# Medtronic

# Accuracy matters.

Experience  
the difference



AccuRhythm™ AI algorithms for  
the Reveal LINQ™ insertable  
cardiac monitor (ICM)



# What can AccuRhythm AI do for your Reveal LINQ ICM clinic?

Experience the difference of Reveal LINQ ICM now with AccuRhythm AI algorithms for AF and pause. Reveal LINQ, the most studied insertable cardiac monitor (ICM) in the market – trusted in over one million implants<sup>1</sup> – is now enhanced with proven artificial intelligence (AI), providing data that's more actionable – without compromise.<sup>2</sup> Preserving over 99.9% true pause and 98.2% true AF alerts relative to Reveal LINQ ICM, AccuRhythm AI algorithms significantly improve the clinic experience.<sup>2</sup>

Validation with data from real-world ICM patients demonstrated:

## AF algorithm

**89.5%** **98.2%**

Reduced false alerts<sup>2</sup>

Preserved true alerts<sup>2</sup>

## Pause algorithm

**80.2%** **99.9%**

Reduced false alerts<sup>2</sup>

Preserved true alerts<sup>2</sup>



To alleviate clinician review burden, these algorithms address the two most common sources of ICM false alerts – atrial fibrillation (AF) and Pause.<sup>2-4</sup>

The AccuRhythm AI platform applies deep learning AI algorithms to Reveal LINQ ICM data flowing into the CareLink™ network.


## How does AccuRhythm AI work?

AccuRhythm AI algorithms were rigorously trained and developed based on over one million professionally adjudicated ECGs to preserve sensitivity and provide data-driven insights free from bias. Once performance was validated, the AI algorithms were locked to ensure an unbiased lens.<sup>2,4-7</sup>

Bypass logic ensures that clinically relevant events are still sent to you for review.

## Greater confidence


How do AccuRhythm AI algorithms work?



Data transmission through the patient's MyCareLink™ monitor.

# Reduce false alerts

AF and Pause generate **~90% of false alerts** in the ICM space<sup>3</sup>



The AccuRhythm AI algorithms can save clinicians approximately **186 hours** of false alert review yearly for every 200 Reveal LINQ ICM patients.<sup>†‡</sup>

**85% cumulative reduction in Reveal LINQ ICM false alerts<sup>2</sup>**

<sup>†</sup>The validation study performance and time study results were projected onto 104,315 Reveal LINQ patients to calculate false alert reduction per year in 200 Reveal LINQ ICM patients.

<sup>‡</sup>Estimated time savings was calculated based on 11.3 minutes per non-actionable ICM transmission.<sup>8</sup>

# Seamless experience

Moving product updates into the cloud enables application of algorithms to Reveal LINQ devices already implanted. This allows clinics to experience immediate benefits from AccuRhythm AI algorithms.

Allow us to seamlessly transform your experience by leveraging Medtronic history, expertise, and data-driven insights in the cardiac monitoring space.

**Talk to your Medtronic sales representative to learn more about AccuRhythm AI algorithms.**



## References

- <sup>1</sup> Medtronic Reveal™ ICM family data. Data on file. 2024.
- <sup>2</sup> Radtke A, Hall M. AccuRhythm AI AF & Pause Algorithms White Paper. April 2023. Medtronic data on file.
- <sup>3</sup> O'Shea CJ, Middeldorp ME, Hendriks JM, et al. Remote Monitoring of Implantable Loop Recorders: False-Positive Alert Episode Burden. *Circ Arrhythm Electrophysiol*. November 2021;14(11):e009635.
- <sup>4</sup> AccuRhythm™ AI Clinician Manual Supplements.
- <sup>5</sup> Sachin Khane R, Surdi AD. Gender differences in the prevalence of electrocardiogram abnormalities in the elderly: a population survey in India. *Iran J Med Sci*. June 2012;37(2):92-99.
- <sup>6</sup> Mansi IA, Nash IS. Ethnic differences in electrocardiographic amplitude measurements. *Ann Saudi Med*. November–December 2004;24(6):459-464.
- <sup>7</sup> Santhanakrishnan R, Wang N, Larson MG, et al. Racial Differences in Electrocardiographic Characteristics and Prognostic Significance in Whites Versus Asians. *J Am Heart Assoc*. March 25, 2016;5(3):e002956. Erratum in: *J Am Heart Assoc*. July 2016;5(7).
- <sup>8</sup> Seiler A, Biundo E, Di Bacco M, et al. Clinic Time Required for Remote and In-Person Management of Patients With Cardiac Devices: Time and Motion Workflow Evaluation. *JMIR Cardio*. October 15, 2021;5(2):e27720.

## Brief Statement

### Reveal LINQ™ Insertable Cardiac Monitor (ICM)

#### System and Accessories

#### Indications

The Reveal LINQ Insertable Cardiac Monitor (ICM) is an implantable patient-activated and automatically-activated monitoring system that records subcutaneous ECG and is indicated in the following cases:

- patients with clinical syndromes or situations at increased risk of cardiac arrhythmias
- patients who experience transient symptoms such as dizziness, palpitation, syncope and chest pain, that may suggest a cardiac arrhythmia.

The device has not been tested specifically for pediatric use.

#### Contraindications

There are no known contraindications for the implant of the Reveal LINQ ICM or its accessories. However, the patient's particular medical condition may dictate whether or not a subcutaneous, chronically implanted device can be tolerated.

#### Warnings/Precautions

Patients with the Reveal LINQ ICM should avoid sources of diathermy, high sources of radiation, electrosurgical cautery, external defibrillation, lithotripsy, therapeutic ultrasound and radiofrequency ablation to avoid electrical reset of the device, and/or inappropriate sensing as described in the Medical procedure and EMI precautions manual. MRI scans should be performed only in a specified MR environment under specified conditions as described in the Reveal LINQ MRI Technical Manual.

Accessories available for use with Reveal LINQ ICM may experience connectivity or performance issues. See product manuals for details and troubleshooting instructions.

#### Potential Adverse Events or Potential Complications

Potential complications of the Reveal LINQ device include, but are not limited to, device rejection phenomena (including local tissue reaction), device migration, infection, and erosion through the skin.

There are no known adverse events associated with the use of any Reveal LINQ ICM accessory.

See the device manuals for detailed information regarding the implant procedure, indications, contraindications, warnings, precautions, and potential complications/adverse events. For further information, please call Medtronic at (800) 328-2518 (Technical Services), (800) 551-5544 (Patient Services) and/or consult Medtronic's website at [www.medtronic.com](http://www.medtronic.com).

**Caution:** Federal law (USA) restricts these devices to sale by or on the order of a physician.

#### AccuRhythm AI ECG Classification System

**Intended Use:** The intended use of the system is to reduce false positive cardiac arrhythmia episodes.

**Contraindications:** There are no known contraindications for AccuRhythm AI Models ZA400, ZA410, or ZA420.

**Precaution:** The AccuRhythm AI ECG classification system may incorrectly adjudicate a true positive episode as an AI false episode, causing that episode to be suppressed in the remote monitoring system.

See the device manual for detailed information regarding the intended use, contraindications, warnings, precautions, and potential complications / adverse events. For further information, call Medtronic Technical Services at (800) 328-2518 and/or consult Medtronic's website at [www.medtronic.com](http://www.medtronic.com).

**Caution:** Federal law (USA) restricts these devices to sale by or on the order of a physician.

# Medtronic

Medtronic  
710 Medtronic Parkway  
Minneapolis, MN 55432-5604 USA

Toll-free in USA: 800.633.8766  
Worldwide: +1.763.514.4000

[medtronic.com](http://medtronic.com)

UC202310336 EN ©2024 Medtronic.  
Minneapolis, MN.  
All Rights Reserved.  
Printed in USA. 02/2024

Medtronic and the Medtronic logo are trademarks of Medtronic. Adobe and Acrobat Reader are registered trademarks of Adobe Systems incorporated in the United States and/or other countries. ™™Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.