

In-Vitro Thrombogenicity Assessment of Polymer Filament Modified and Native Platinum Embolic Coils

Clinical Paper Review

Prepared by
PVH Scientific
Communications
for Health Care
Professionals

Authors Gaurav Girdhar; Megan Read; JiHae Sohn, Chirag Shah; Sanjay Shrivastava

Journal *Journal of the Neurological Sciences* 339 (2014) 97-101,
[https://www.jns-journal.com/article/S0022-510X\(14\)00054-9/fulltext](https://www.jns-journal.com/article/S0022-510X(14)00054-9/fulltext)

Purpose

To evaluate the prothrombogenic effects of platinum embolic coils with filaments and platinum embolic coils without filaments.

Methods

- In-vitro study utilizing a quantitative method that tracks the formation of thrombin upon exposure of the test samples to human platelet rich plasma.
- Thrombin generation was measured via a slow binding fluorogenic substrate. Fluorescence was measured at one-minute intervals and representative thrombograms were generated for each device tested
- Devices: Axiom™ MicroFX™ nylon detachable coil system; Axiom™ MicroFX™ PGLA detachable coil system; Axiom™ detachable coil system (sample size of 2 coils in each of 9 unique batches of human platelet rich plasma = 18 total tests per device type)
- Study endpoints evaluated:
 - Peak Thrombin[†]
 - Rate of Thrombin Generation[‡]
 - Endogenous Thrombin Potential (ETP)[§]

Results

Axiom MicroFX coils exhibited significantly higher peak thrombin, rate of thrombin, and Endogenous Thrombin Potential (ETP) as compared to bare platinum Axiom MicroFX coils.

Figure 1: Representative thrombograms for Axiom, Axiom MicroFX PGLA, and Axiom MicroFX Nylon samples.





Figure 2: Graphical summary of test result values normalized against bare Axiom coil results (mean +/- S.E.) for the three parameters over 18 tests per device type.



Discussion

- Coils with nylon filaments exhibited significantly higher thrombogenicity than both coils with PGLA filaments and bare platinum coils. Coils with PGLA filaments significantly enhanced thrombogenicity relative to bare platinum coils.
- Because these tests were performed in microplates, the findings provide a quantitative measure of thrombus growth around the coil in the absence of blood flow. This provides a direct thrombogenicity comparison between coils composed of different materials independent of dynamic physiological conditions.
- The only difference between the test sample groups was the presence or absence of microfilaments, and therefore, the improved thrombin generation seen with Axiom MicroFX Nylon or Axiom MicroFX PGLA is attributed to the presence of the nylon or PGLA fibers.

Author's Conclusions

Coils with added filaments produce a stronger thrombogenic response as compared to bare platinum coils. The impact on thrombogenicity is dependent on the filament material (nylon or PGLA).

Funding Source

Medtronic, formerly Covidien

[†]Peak Thrombin – maximum thrombin generated in the experiment

[‡]Rate of Thrombin – gradient between the max amount of thrombin generation and the first inflection point

[§]Total Amount of Thrombin Generated – defined as ETP or endogenous thrombin potential (area under the curve)

If you are located in the United States, please refer to the brief statement(s) below to review applicable indications, safety and warning information. If you are located outside the United States, see the device manual for detailed information regarding instructions for use, the implant procedure, indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative and/or consult the Medtronic website at www.medtronic.eu.

Axiom™ and Axiom™ Prime detachable coils

The Axiom™ and Axiom™ Prime detachable coils are not intended for all patients and may not be the appropriate treatment for all clinical scenarios. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Axiom™ and Axiom™ Prime detachable coils are intended for the endovascular embolization of intracranial aneurysms. Axiom™ and Axiom™ Prime detachable coils are also intended for the embolization of other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistulae.