

# Medtronic

Engineering the extraordinary

SINUVA™ (mometasone furoate) sinus implant

## RESOLVE II phase 3 clinical trial

Mometasone furoate sinus implants for patients  
with chronic rhinosinusitis with nasal polyps



### INDICATION

SINUVA™ sinus implant is a corticosteroid-eluting implant indicated for the treatment of chronic rhinosinusitis with nasal polyps in adult patients  $\geq 18$  years of age who have had ethmoid sinus surgery.

### IMPORTANT SAFETY INFORMATION

#### CONTRAINDICATIONS

Patients with known hypersensitivity to mometasone furoate and any of the ingredients of the SINUVA™ sinus implant.

**Please see additional important safety information on the following page.**

A phase 3 clinical trial -  
endpoints

A phase 3 clinical trial -  
methods

Path to improvement -  
co-primary efficacy

Path to improvement -  
secondary efficacy

Safety information

# A phase 3 clinical trial of mometasone furoate sinus implants for patients with CRSwNP

## Objective

To evaluate the efficacy and safety of in-office bilateral placement of SINUVA™ (mometasone furoate) sinus implant, a bioabsorbable corticosteroid-eluting sinus implant containing mometasone furoate, vs mometasone furoate nasal spray alone in patients  $\geq 18$  years of age with nasal polyps and a history of ethmoid sinus surgery<sup>1</sup>

## Endpoints

### Co-primary efficacy:

- Change from baseline to Day 90 in bilateral polyp grade, as determined by an independent blinded panel on a scale of 0-8 with each side grading from 0 (no visible nasal polyps) to 4 (nasal polyps completely obstructing nasal cavity)<sup>1</sup>
- Change from baseline to Day 30 in nasal obstruction/congestion score, as determined by eDiary, on a scale of 0 (no symptoms) to 3 (severe symptoms)<sup>1</sup>

### Secondary efficacy endpoints:

- Proportion of patients still indicated for repeat endoscopic sinus surgery (ESS) at Day 90 based on clinical investigator assessment using study-specific criteria<sup>1</sup>

- Change from baseline to Day 90 in nasal obstruction/congestion score, as determined by patients, on a scale of 0 (no symptoms) to 3 (severe symptoms)<sup>1</sup>
- Change in percent ethmoid sinus obstruction at Day 90, as determined by the independent, blinded panel on a 100-mm visual analogue scale<sup>1</sup>
- Decreased sense of smell score change from baseline to Day 90, as determined by patients on a six point Likert scale of 0 (absent) to 5 (very severe)<sup>1</sup>
- Decreased facial pain/pressure score change from baseline to Day 90, as determined by patients on a 6-point Likert scale of 0 (absent) to 5 (very severe)<sup>1</sup>

## IMPORTANT SAFETY INFORMATION (continued)

### WARNINGS AND PRECAUTIONS

**Local Nasal Adverse Reactions:** Monitor nasal mucosa adjacent to the SINUVA™ sinus implant for any signs of bleeding (epistaxis), irritation, infection, or perforation. Avoid use in patients with nasal ulcers or trauma.

**Please see additional important safety information on the following page.**

A phase 3 clinical trial -  
endpoints

A phase 3 clinical trial -  
methods

Path to improvement -  
co-primary efficacy

Path to improvement -  
secondary efficacy

Safety information

# A phase 3 clinical trial of mometasone furoate sinus implants for patients with nasal polyps

## Objective

To evaluate the efficacy and safety of in-office bilateral placement of SINUVA™ (mometasone furoate) sinus implant, a bioabsorbable corticosteroid-eluting sinus implant containing mometasone furoate, vs mometasone furoate nasal spray alone in patients  $\geq 18$  years of age with nasal polyps and a history of ethmoid sinus surgery<sup>1</sup>

## Methods

### Study design:

Randomized, controlled, double-blind, multicenter study with 300 patients followed for 90 days<sup>1</sup>

### Study population:

- $\geq 18$  years of age
- Diagnosed with CRS
- Had undergone prior bilateral total ethmoidectomy
- Indicated for revision ESS
- Presented with moderate to severe nasal obstruction/congestion symptoms and recurrent bilateral sinus obstruction due to sinonasal polyposis despite the use of intranasal cortico-steroid sprays and recent high dose steroids

### Treatment:

- 201 patients were randomized to the SINUVA™ treatment arm where they underwent bilateral placement of the SINUVA™ in the ethmoid sinuses
- 99 patients were randomized to the control arm where they received a sham procedure
- Patients in both study arms received once daily mometasone furoate nasal spray (200  $\mu$ g) through Day 90

## IMPORTANT SAFETY INFORMATION (continued)

### WARNINGS AND PRECAUTIONS (continued)

**Glaucoma and Cataracts:** Nasal steroids may result in development of glaucoma and/or cataracts. Glaucoma, cataracts, and clinically significant elevation of intraocular pressure were not observed in patients from the treatment group of one randomized controlled clinical study (N = 53) who underwent bilateral placement of SINUVA™ sinus implants. Close monitoring is warranted in patients with a change in vision or with a history of increased intraocular pressure, glaucoma, and/or cataracts.

**Please see additional important safety information on the following page.**

A phase 3 clinical trial -  
endpoints

A phase 3 clinical trial -  
methods

Path to improvement -  
co-primary efficacy

Path to improvement -  
secondary efficacy

Safety information

# Discover the path to improvement

## Co-primary efficacy outcomes<sup>†,‡,§,1,3</sup>

In RESOLVE II (n=300), SINUVA™ (mometasone furoate) sinus implant patients experienced statistically significant reductions in bilateral polyp grade and nasal obstruction/congestion score

### 2.3x improvement<sup>§</sup> in bilateral polyp grade

from baseline to day 90 compared to once-daily mometasone furoate nasal spray alone

Mean change (SD) for SINUVA™ implant -0.56 (1.06) vs. -0.15 (0.91) with control, mean group difference (95% CI) - 0.35 (-0.60, -0.09); (P=0.0073)<sup>1,3</sup>

### Significant improvement<sup>§</sup> in nasal obstruction/congestion score

from baseline to day 30 compared to once-daily mometasone furoate nasal spray alone

Mean change (SD) for SINUVA™ implant -0.80 (0.73) vs. -0.56 (0.62) with control, mean group difference (95% CI) -0.23 (-0.39, -0.06); (P=0.0074)<sup>1,3</sup>

SD, standard deviation.  
† Patients in the SINUVA™ group also received once daily mometasone furoate nasal spray daily.  
‡ Co-primary endpoints: Change from baseline to day 90 in bilateral polyp grade, as determined by an independent panel on a scale of 0 (no visible nasal polyps) to 4 (nasal polyps completely obstructing nasal cavity). Change from baseline to day 30 in nasal obstruction/congestion score, as determined by patients on a scale of 0 (no symptoms) to 3 (severe symptoms).  
§ Improvement calculated in post-hoc analysis as [(T-C)/C] \*100 using values reported in the SINUVA™ prescribing information and Kern (2018).  
Study design: The RESOLVE II study was a randomized, controlled, double-blind, multicenter study with 300 patients. Patients were ≥18 years of age with chronic sinusitis who had prior bilateral total ethmoidectomy, but were indicated for revision endoscopic surgery because they presented with moderate to severe nasal obstruction/congestion symptoms and recurrent bilateral sinus obstruction due to sinonasal polyposis despite the use of intranasal corticosteroid sprays and recent high dose steroids. 201 patients were randomized to the SINUVA™ implant treatment arm where they underwent bilateral placement of SINUVA™ implant in the ethmoid sinuses. 99 patients were randomized to the control arm where they received a sham procedure. Patients in both study arms received once-daily mometasone furoate nasal spray (200 µg) through day 90.

### IMPORTANT SAFETY INFORMATION (continued)

### WARNINGS AND PRECAUTIONS (continued)

**Hypersensitivity Reactions:** Hypersensitivity reactions, including rash, pruritus, and angioedema have been reported with the use of corticosteroids. Please see additional important safety information on the following page.

A phase 3 clinical trial - endpoints

A phase 3 clinical trial - methods


Path to improvement - co-primary efficacy

Path to improvement - secondary efficacy

Safety information

# Discover the path to improvement

Secondary efficacy endpoints (baseline to day 90)<sup>†,1,2</sup>



Reduced eligibility for repeat surgery



Delivered up to 90 days of sustained symptom relief



Significantly reduced sinus obstruction



Improved patients' sense of smell

SINUVA™ (mometasone furoate) sinus implant placement resulted in a **61% reduction in the proportion of patients still indicated for repeat ESS** vs. 37% in the mometasone furoate nasal-spray-only control arm at day 90. Odds ratio (95% CI) 2.69 (1.63, 4.44); (P=0.0004)

**Patients treated with SINUVA™ did not experience a significant improvement in their self-reported facial pain/pressure score.** The mean change (SD) for SINUVA™ implant was -0.77 (1.21) vs. -0.90 (1.27) with control (P=0.9130)

SD, standard deviation.

† Secondary endpoints: Proportion of patients still indicated for repeat ESS at day 90 despite ongoing intranasal steroid use based on clinical investigator assessment using study-specific criteria. Change from baseline to day 90 in instantaneous nasal obstruction/congestion score (daily diary), as determined by patients on a scale of 0 (no symptoms) to 3 (severe symptoms). Change in percent ethmoid sinus obstruction at day 90, as determined by the independent, blinded panel on a 100 mm visual analogue scale (VAS). Decreased sense of smell and facial pain/pressure score change from baseline to day 90, as determined by patients on a six point Likert scale of 0 (absent) to 5 (very severe). P-values for secondary endpoints were prespecified and adjusted for multiplicity.<sup>12</sup>

## IMPORTANT SAFETY INFORMATION (continued)

### WARNINGS AND PRECAUTIONS (continued)

**Immunosuppression and Risk of Infections:** Persons who are using drugs that suppress the immune system, such as corticosteroids, including SINUVA™ sinus implants are more susceptible to infections than healthy individuals. The safety and effectiveness of SINUVA™ sinus implant have not been established in pediatric patients less than 18 years of age and SINUVA™ is not indicated for use in this population. Corticosteroids should be used with caution, if at all, in patients with active or quiescent tuberculosis infection of the respiratory tract; untreated systemic fungal, bacterial, viral, or parasitic infections; or ocular herpes simplex.

**Please see additional important safety information on the following page.**

A phase 3 clinical trial - endpoints

A phase 3 clinical trial - methods


Path to improvement - co-primary efficacy

Path to improvement - secondary efficacy

Safety information

# Discover the path to improvement

Secondary efficacy endpoints (baseline to day 90)<sup>†,1,2</sup>



Reduced eligibility for repeat surgery



Delivered up to 90 days of sustained symptom relief



Significantly reduced sinus obstruction



Improved patients' sense of smell

- Patients treated with SINUVA™ (mometasone furoate) sinus implant experienced sustained symptomatic **improvements in nasal obstruction/congestion** at day 90
- Mean change (SD) for SINUVA™ implant was -0.93 (0.80) vs. -0.69 (0.79) with control, mean group difference (95% CI) -0.27 (-0.48, -0.07); (P=0.0248)

**Patients treated with SINUVA™ did not experience a significant improvement in their self-reported facial pain/pressure score.** The mean change (SD) for SINUVA™ implant was -0.77 (1.21) vs. -0.90 (1.27) with control (P=0.9130)

SD, standard deviation.

† Secondary endpoints: Proportion of patients still indicated for repeat ESS at day 90 despite ongoing intranasal steroid use based on clinical investigator assessment using study-specific criteria. Change from baseline to day 90 in instantaneous nasal obstruction/congestion score (daily diary), as determined by patients on a scale of 0 (no symptoms) to 3 (severe symptoms). Change in percent ethmoid sinus obstruction at day 90, as determined by the independent, blinded panel on a 100 mm visual analogue scale (VAS). Decreased sense of smell and facial pain/pressure score change from baseline to day 90, as determined by patients on a six point Likert scale of 0 (absent) to 5 (very severe). P-values for secondary endpoints were prespecified and adjusted for multiplicity.

## IMPORTANT SAFETY INFORMATION (continued)

### WARNINGS AND PRECAUTIONS (continued)

**Hypercorticism and Adrenal Suppression:** If corticosteroid effects such as hypercorticism and adrenal suppression appear in patients, consider sinus implant removal.

**Please see additional important safety information on the following page.**

A phase 3 clinical trial - endpoints

A phase 3 clinical trial - methods


Path to improvement - co-primary efficacy

Path to improvement - secondary efficacy

Safety information

# Discover the path to improvement

Secondary efficacy endpoints (baseline to day 90)<sup>†,2</sup>



Reduced eligibility for repeat surgery



Delivered up to 90 days of sustained symptom relief



Significantly reduced sinus obstruction



Improved patients' sense of smell

- Patients treated with SINUVA™ (mometasone furoate) sinus implant had a significantly greater **decrease in percent ethmoid sinus obstruction** at day 90
- Mean change (SD) for SINUVA™ implant was -11.3 (18.1) vs. -1.9 (14.4) with control, mean group difference (95% CI) -7.96 (-12.10, -3.83); (P=0.0007)

**Patients treated with SINUVA™ did not experience a significant improvement in their self-reported facial pain/pressure score.** The mean change (SD) for SINUVA™ implant was -0.77 (1.21) vs. -0.90 (1.27) with control (P=0.9130)

SD, standard deviation.

† Secondary endpoints: Proportion of patients still indicated for repeat ESS at day 90 despite ongoing intranasal steroid use based on clinical investigator assessment using study-specific criteria. Change from baseline to day 90 in instantaneous nasal obstruction/congestion score (daily diary), as determined by patients on a scale of 0 (no symptoms) to 3 (severe symptoms). Change in percent ethmoid sinus obstruction at day 90, as determined by the independent, blinded panel on a 100 mm visual analogue scale (VAS). Decreased sense of smell and facial pain/pressure score change from baseline to day 90, as determined by patients on a six point Likert scale of 0 (absent) to 5 (very severe). P-values for secondary endpoints were prespecified and adjusted for multiplicity.

## IMPORTANT SAFETY INFORMATION (continued)

### ADVERSE REACTIONS

The most common adverse reactions observed (> 1% of subjects) in clinical studies were asthma, headache, epistaxis, presyncope, bronchitis, otitis media, and nasopharyngitis.

**Please see additional important safety information on the following page.**

A phase 3 clinical trial - endpoints

A phase 3 clinical trial - methods

Path to improvement - co-primary efficacy


Path to improvement - secondary efficacy

Safety information



# Discover the path to improvement

Secondary efficacy endpoints (baseline to day 90)<sup>†,1,3</sup>



Reduced eligibility for repeat surgery



Delivered up to 90 days of sustained symptom relief



Significantly reduced sinus obstruction



Improved patients' sense of smell

- Patients treated with SINUVA™ (mometasone furoate) sinus implant experienced a **significant improvement in their self-reported decreased sense of smell score** at day 90 on a six-point Likert scale
- Mean change (SD) for SINUVA™ implant was -1.20 (1.66) vs. -0.76 (1.60) with control, mean group difference (95% CI) -0.46 (-0.85 -0.06); (P=0.0470)

**Patients treated with SINUVA™ did not experience a significant improvement in their self-reported facial pain/pressure score.** The mean change (SD) for SINUVA™ implant was -0.77 (1.21) vs. -0.90 (1.27) with control (P=0.9130)

SD, standard deviation.

† Secondary endpoints: Proportion of patients still indicated for repeat ESS at day 90 despite ongoing intranasal steroid use based on clinical investigator assessment using study-specific criteria. Change from baseline to day 90 in instantaneous nasal obstruction/congestion score (daily diary), as determined by patients on a scale of 0 (no symptoms) to 3 (severe symptoms). Change in percent ethmoid sinus obstruction at day 90, as determined by the independent, blinded panel on a 100 mm visual analogue scale (VAS). Decreased sense of smell and facial pain/pressure score change from baseline to day 90, as determined by patients on a six point Likert scale of 0 (absent) to 5 (very severe). P-values for secondary endpoints were prespecified and adjusted for multiplicity.

## IMPORTANT SAFETY INFORMATION (continued)

### POSTMARKETING EXPERIENCE

The following adverse reactions have been identified during post-approval use of the SINUVA™ sinus implant. These events include implant migration, lack of efficacy, nasal pain, headache, epistaxis.

Please see additional important safety information on the following page.

A phase 3 clinical trial - endpoints

A phase 3 clinical trial - methods

Path to improvement - co-primary efficacy

Path to improvement - secondary efficacy

Safety information



# SINUVA™ (mometasone furoate) sinus implant safety data

Adverse reactions with > 1% incidence rate and more common in treatment than the control group in combined data from RESOLVE I and RESOLVE II

| Adverse reactions | Treatment<br>(n=254) n (%) | Control<br>(MF nasal spray)<br>(n=146) n (%) |
|-------------------|----------------------------|--|
| Asthma            | 12 (4.7)                   | 6 (4.1)                                      |
| Headache          | 9 (3.5)                    | 5 (3.4)                                      |
| Epistaxis         | 6 (2.4)                    | 2 (1.4)                                      |
| Presyncope        | 6 (2.4)                    | 3 (2.1)                                      |
| Bronchitis        | 5 (2.0)                    | 2 (1.4)                                      |
| Otitis media      | 5 (2.0)                    | 2 (1.4)                                      |
| Nasopharyngitis   | 3 (1.2)                    | 1 (0.7)                                      |

SINUVA™ Indication & Important Safety Information

INDICATION

SINUVA™ sinus implant is a corticosteroid-eluting implant indicated for the treatment of chronic rhinosinusitis with nasal polyps in adult patients ≥ 18 years of age who have had ethmoid sinus surgery.

IMPORTANT SAFETY INFORMATION

CONTRAINDICATIONS

Patients with known hypersensitivity to mometasone furoate and any of the ingredients of the SINUVA™ sinus implant.

WARNINGS AND PRECAUTIONS

**Local Nasal Adverse Reactions:** Monitor nasal mucosa adjacent to the SINUVA™ sinus implant for any signs of bleeding (epistaxis), irritation, infection, or perforation. Avoid use in patients with nasal ulcers or trauma.

**Glaucoma and Cataracts:** Nasal steroids may result in development of glaucoma and/or cataracts. Glaucoma, cataracts, and clinically significant elevation of intraocular pressure were not observed in patients from the treatment group of one randomized controlled clinical study (N = 53) who underwent bilateral placement of SINUVA™ sinus implants. Close monitoring is warranted in patients with a change in vision or with a history of increased intraocular pressure, glaucoma, and/or cataracts.

**Hypersensitivity Reactions:** Hypersensitivity reactions, including rash, pruritus, and angioedema have been reported with the use of corticosteroids.

**Immunosuppression and Risk of Infections:** Persons who are using drugs that suppress the immune system, such as corticosteroids, including SINUVA™ sinus implants are more susceptible to infections than healthy individuals. The safety and effectiveness of SINUVA™ sinus implant have not been established in pediatric patients less than 18 years of age and SINUVA™ is not indicated for use in this population. Corticosteroids should be used with caution, if at all, in patients with active or quiescent tuberculosis infection of the respiratory tract; untreated systemic fungal, bacterial, viral, or parasitic infections; or ocular herpes simplex.

**Hypercorticism and Adrenal Suppression:** If corticosteroid effects such as hypercorticism and adrenal suppression appear in patients, consider sinus implant removal.

ADVERSE REACTIONS

The most common adverse reactions observed (> 1% of subjects) in clinical studies were asthma, headache, epistaxis, presyncope, bronchitis, otitis media, and nasopharyngitis.

POSTMARKETING EXPERIENCE

The following adverse reactions have been identified during post-approval use of the SINUVA™ sinus implant. These events include implant migration, lack of efficacy, nasal pain, headache, epistaxis.

Rx only. Please see Full Prescribing Information for SINUVA™ attached and at [SINUVA.com/PI](https://www.sinuva.com/PI).

A phase 3 clinical trial - endpoints

A phase 3 clinical trial - methods

Path to improvement - co-primary efficacy

Path to improvement - secondary efficacy

Safety information

References

- 1. Kern RC, Stolovitzky JP, Silvers SL, et al. A phase 3 trial of mometasone furoate sinus implants for chronic sinusitis with recurrent nasal polyps. Int Forum Allergy Rhinol. 2018;8:471-481.
- 2. Data on file, Intersect ENT. RESOLVE II CSR R 28017 Rev. 3.0; February 18, 2017.
- 3. Data on file, Intersect ENT, Inc. CR-00014; Supplemental Statistical Analyses Rev 1.0; January 17, 2018.
- 4. SINUVA™ [Prescribing Information]. Menlo Park, CA: Intersect ENT; 2023.

For more information, visit [SINUVA.com/hcp](https://sinuva.com/hcp).

Medtronic

ENT

6743 Southpoint Drive N  
Jacksonville, FL 32216  
USA  
Toll free: (800) 874-5797  
Telephone: (904) 296-9600  
Fax: (800) 678-3995

©2024 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.  
10/2024 - US-ENT-2400414 v 1.0 - [WF#14792182]

A phase 3 clinical trial - endpoints

A phase 3 clinical trial - methods

Path to improvement - co-primary efficacy

Path to improvement - secondary efficacy

Safety information