

Quick reference guide: Integrated nebulizer




The Puritan Bennett™ 980 ventilator includes an integrated nebulizer option – designed to power and control the Aerogen™ nebulizer to deliver prescribed medication to patients.

Connect the nebulizer

1. Open the connector panel door
2. Firmly insert the sensor cable into the connector port labeled **Aerogen Nebulizer**
3. Connect the other end to the Aerogen™ Solo or Pro nebulizer electrical connection port

Assemble the nebulizer

1. Fit the Aerogen™ Solo or Pro nebulizer firmly into the appropriate T-piece adapter
 - a. Adult T-adapter
 - b. Pediatric T-adapter
 - c. Neonatal T-adapter with adapters



Start and stop a timed nebulization cycle

1. Add medication to the Aerogen™ Solo or Pro nebulizer
2. Touch **Nebulizer** icon on the GUI
3. Choose **Solo or Pro** button
4. Touch **Time** button and turn the dial to the desired nebulization time
5. Touch **START**
6. Allow nebulizer to deliver the medication for the allotted time
7. Touch **Extend** button to extend nebulization time by time chosen on step 4
8. Touch **Stop** button to stop the nebulizer

Start and stop a continuous nebulization cycle

1. Prime the syringe and connect it to the nebulizer
2. Touch **Nebulizer** icon on the GUI
3. Touch **Solo** button
4. Touch **Continuous Nebulization** (∞) button
5. Touch **START** button
6. Allow the nebulizer to deliver the medication
7. Touch **Stop** button to stop the nebulizer

Note: The nebulizer does not auto-shutoff. When the medication reservoir is empty, manually stop the nebulizer by touching the **Stop** button.



Consult the operator manual for a full set of instructions for use and indications for use, including warnings and contraindications.