


## Puritan Bennett™ 980 ventilator

# IE Sync™ software quick reference guide

The IE Sync™ software option for the Puritan Bennett™ 980 (PB980) ventilator provides a noninvasive method of breath triggering and cycling that may enhance patient - ventilator synchrony compared to conventional flow triggering and cycling for adult patients with airflow obstruction and weak inspiratory efforts.<sup>1</sup>

When the IE Sync trigger type is selected while in SPONT mode with the pressure support (PS) or volume support (VS) spontaneous breath type selected,  $I_{SYNC}$  is the setting for trigger sensitivity and  $E_{SYNC}$  is the setting for cycle sensitivity on the ventilator.



### IE Sync™ setup and use

1. Select **Invasive Ventilation Type**.
2. Select **SPONT Mode**.
3. Select **PS** or **VS Spontaneous Type**.
4. Select **IE Sync Trigger Type**.
5. At start up, both  $I_{SYNC}$  and  $E_{SYNC}$  settings default to the mid-point "0".
6. Monitor data values for  $f_{TOT}$ ,  $V_T$ , and  $T_{I,SPONT}$  and pressure and flow waveforms to ensure patient comfort and synchrony.


If needed, adjust  $I_{SYNC}$  and  $E_{SYNC}$  settings according to the instructions described later in this guide.

**Note:** If settings adjustment is needed, adjust  $E_{SYNC}$  first because the timing for breath cycling may impact triggering for the breath that follows.



## E\_SYNC cycling adjustment

1. E\_SYNC settings are mapped to patient lung mechanics (obstructive, restrictive, and normal).



- Patients with obstructive lung conditions
- Patients with normal lung conditions
- Patients with restrictive lung conditions

2. Examine  $V_T$  and  $T_{I\text{SPONT}}$  data values and the flow and pressure waveforms to ensure that PS/VS breaths cycle to exhalation when the patient stops inhaling.
3. If breath cycling-off is delayed past the end of patient inhalation ( $V_T$  and  $T_{I\text{SPONT}}$  data values are elevated/waveforms show late cycling), adjust E\_SYNC to a lower (earlier) setting.

**Note:** This may be more common on patients with airflow limitation/obstructive lung conditions.

4. If breath cycling-off occurs before the patient stops inhaling ( $V_T$  and  $T_{I\text{SPONT}}$  data values are reduced /waveforms show early cycling), adjust E\_SYNC to a higher (later) setting.

**Note:** This may be more common on patients with restrictive lung conditions.

1. Internal engineering performance studies.

## I\_SYNC triggering adjustment

1. Compare the  $f_{TOT}$  data value with:

- Number of observed patient efforts
- Respiratory rate value shown on the patient bedside monitor and
- Flow waveform to ensure every patient effort is matched with a triggered breath from the PB980 ventilator.

2. Optimize the I\_SYNC setting to the lowest setting that allows for comfortable, reliable breath triggering without auto-triggering.

**Note:** Patients with cardiac oscillations or shivering may require a higher (less sensitive) I\_SYNC setting.

©2022 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 09/2022-US-RE-2100666-[WF#3554454]

[medtronic.com/covidien](https://www.medtronic.com/covidien)

**Medtronic**