

Improve endotracheal tube management in neonatal patients

The presence of an endotracheal tube (ETT) places patients at risk for several complications. Among neonatal patients, invasive mechanical ventilation tends to be required in the most premature, vulnerable patients who can suffer significant harm from endotracheal tube-related issues, such as unplanned extubations, tube misplacement, or ETT suctioning-related stress.¹⁻⁶

Utilize this guide for background on the types of endotracheal tube-related issues and stressful interventions that may be mitigated with additional monitoring via the SonarMed™ airway monitoring system.



TABLE OF CONTENTS



Managing endotracheal tubes in neonatal patients	3
SonarMed™ airway monitoring system overview	4
Endotracheal tube related adverse events	5
Unplanned Extubations	6
Endotracheal tube malposition	7
Tube movement in neonates	8
Suctioning-related harm	9
Minimizing stress	10
Chest radiographs	11
Safe kangaroo care	12
References	13

MANAGING ENDOTRACHEAL TUBES IN NEONATAL PATIENTS

OVERVIEW

Keeping an endotracheal tube (ETT) in place can be challenging in neonates as they are not typically paralyzed.⁷ Movement can lead to migration of the ETT, which can lead to dislodgement or an unplanned extubation (UE).^{8,9}

UEs are a significant safety concern for neonates. UEs are the most common adverse event during mechanical ventilation in the neonatal intensive care unit (NICU). The risk of UEs in neonates is higher when compared to other populations. UEs often result in an emergent reintubation and may cause cardiovascular collapse, leading to an increase in hospital length of stay and costs.^{4,10}

Learn more:

[Unplanned extubations](#)[Endotracheal tube malposition](#)[Tube movement in neonates](#)

The frequency and timing of ETT suctioning, one of the most common procedures in the NICU, is controversial.^{11,12} Though it is essential to prevent life threatening obstruction of the endotracheal tube, the procedure is not without complications and adverse effects.¹¹ Therefore, it should be used as infrequently as possible while still preventing ETT obstruction. While clinical assessment is recommended to inform whether suctioning is indicated, there is limited evidence evaluating which clinical signs are most predictive of suctioning need.^{11,13}

Learn more:

[Suctioning-related harm](#)

Evidence has demonstrated that over-stimulation and caregiving interventions associated with NICU care may be detrimental to neurobehavioral outcomes.¹⁴ A review of the frequency of stressful or painful procedures average that neonates undergo 16 such procedures on average per day.¹⁵

Learn more:

[Chest radiographs](#)[Safe kangaroo care](#)

SONARMED™ AIRWAY MONITORING SYSTEM OVERVIEW

SONARMED™ AIRWAY MONITORING SYSTEM OVERVIEW

A minimally disruptive solution that provides actionable airway management information

The SonarMed™ airway monitoring system may help improve a clinician's ability to manage a patient's airway — offering peace of mind to caregivers — by providing precise, continuous, real-time monitoring of the ETT position and patency.

- Helps reduce unplanned extubations
- Optimizes suctioning
- Monitors ETT movement



Throughout the duration of an intubation, the SonarMed™ airway monitoring system:

- Measures the location of the ETT tip within the trachea, and may assist in detection of any movement, to help reduce unplanned extubations
- Assists in identifying the location of any obstructions within the ETT, and the percentage obstructed, for optimized suctioning and removal
- Measures the circumference of the patient's trachea at the tip of the ETT to observe any movement toward a smaller or larger passageway

ENDOTRACHEAL TUBE RELATED **ADVERSE EVENTS**

UNPLANNED EXTUBATIONS

UEs are a frequent and harmful adverse event in mechanically ventilated neonates. UE and subsequent re-intubation is associated with several adverse effects, including cardiopulmonary deterioration, ventilator-associated pneumonia, hypoxemia, hypercarbia, laryngeal trauma, and subglottic stenosis.¹⁸ There are multiple patient-specific and treatment-specific factors that contribute to the risk of UEs.

Frequent:

- 14-41% of infants experience a UE during their NICU hospitalization⁴
- UE rates vary by institution and from 0.54 to 16.1 UE per 100 ventilator days⁴

Harmful and costly:

- Compared to matched patients, patients who suffer UEs have:⁴
 - Excess MV LOS of 6.5 days
 - Excess inflation-adjusted total hospital costs of \$49 587

Multifactorial underlying causes:

- Risk factors for UEs include:¹⁹
 - Age (younger patients)
 - Inadequate tube fixation
 - Agitation
 - Copious secretions
 - Performance of patient procedures
 - Nursing workload
 - Manipulation at the time of UE

ENDOTRACHEAL TUBE RELATED ADVERSE EVENTS

UNPLANNED EXTUBATIONS (CONT'D)

Klugman, D. et al.²⁰

Assessment of an Unplanned Extubation Bundle to Reduce Unplanned Extubations in Critically Ill Neonates, Infants, and Children


JAMA Pediatr 2020 Vol. 174 Issue 6 Pages e200268

Key takeaways:

- UEs rate are higher in the NICU
- Interventions were less effective in the NICU (17.6% absolute reduction) compared to the PICU (20.6% absolute reduction)
- Even after implementation of a quality improvement bundle, the rate of UEs in the NICU was still higher than the goal of less than 0.95 UEs rate for 100 vent days set by Children's Hospitals' Solutions for Patient Safety initiative²¹

STUDY INFORMATION	
STUDY DESIGN	Quality improvement historical controlled study
METHODS	Forty-three children's hospital implemented an unplanned intubation bundle consisting of standardized procedures for: <ul style="list-style-type: none"> ▪ Securing, assessing, repositioning, and manipulating ETTs ▪ Managing tubes during high-risk situations
RESULTS	<ul style="list-style-type: none"> ▪ Overall UE rate declined from 1.135 to 0.862 per 100 ventilator days ▪ UE rates varied significantly by unit type (figure 1) ▪ Nearly 60 of UEs required reintubation ▪ UE-related cardiovascular collapse declined from 0.041 to 0.025 per 100 ventilators days

UE rate by unit before and after implementation of quality improvement bundle



*No change during study period

ENDOTRACHEAL TUBE RELATED **ADVERSE EVENTS**


35%
OF NICU
PATIENTS
experience ETT
misplacement²⁶

ENDOTRACHEAL TUBE MALPOSITION

In preterm neonates, the difference between bronchial intubation and risk for unplanned extubation is only a few centimeters.²²

Inappropriate positioning of the endotracheal tube is associated with a number of adverse events, including atelectasis, hyperinflation, pneumothorax, cardiac arrhythmias, or UE.^{23,24}


- In a prospective review of neonatal patients whose endotracheal tube depth was set according to the AHA/AAP Pediatric Advanced Life Support guidelines, Volsko et al. found that 69% of endotracheal tubes were malpositioned, with 11.8% positioned above the thoracic inlet and 88.2% positioned < 1 cm above the carina.²²



**ETT
RIDING LOW**
Increases the risk of the tube
entering a bronchial stem.⁴

Consequences of placing the tube too low:

- In patients where bronchial intubation occurs, ventilation is impaired, and the patient is at-risk for ventilator induced lung injury by excessively ventilating one lung and inadequately ventilating the other.^{24,27,28}
- Premature patients with immature lungs are particularly vulnerable in this circumstance due to low compliance associated with surfactant deficiency and fluid-filled alveoli, which increases risk of pneumothorax in one lung and atelectasis in the other.^{23,29}



**ETT
RIDING HIGH**
Increases the risk of tube
dislodgement from the trachea.⁴

Consequences of placing the tube too high:

- In a pre/post intervention study investigating the influence of an unplanned intubation prevention bundle that emphasized placing the tube tip below T1, Morris et al. found that high endotracheal tubes are associated with a greater risk of UE.⁸

ENDOTRACHEAL TUBE RELATED **ADVERSE EVENTS**



TUBE MOVEMENT IN NEONATES

Patient movement is associated with migration of the endotracheal tube, which even after initial proper placement may lead to bronchial intubation or unplanned extubation.⁹ Extension of the cervical spine will cause the tube to move towards the head. When the head is in a flexed position, the tube will move away from the head.⁹

**Movement from patient care and repositioning may contribute to UEs.¹⁸
In neonates, movement is typically not restricted in order to limit sedative exposure.¹⁷**

In a trial of an unplanned extubation quality improvement bundle, Klugman et al.²⁰ identified the following high-risk procedures as requiring two clinicians for monitoring tube position and preventing movement:

- Bedside imaging procedures
- Bedside invasive procedures
- Kangaroo care/parent holding
- Routing repositioning
- Switching beds
- Early mobility

Consequently, many unplanned extubation prevention programs include assessment of tube positioning at regular periods as well as before and after high-risk procedures.^{18, 20, 30, 31}

ENDOTRACHEAL TUBE RELATED **ADVERSE EVENTS**

ETT SUCTIONING-RELATED HARM

The frequency and timing of suctioning ETTs in neonatal patients is controversial. Though it is essential to prevent life threatening complete obstruction of the endotracheal tube, the procedure is not without complications and adverse effects.¹¹ Therefore, it should be used as infrequently as possible while still preventing obstruction. While clinical assessment is recommended to inform whether suctioning is indicated,¹¹ there is limited evidence evaluating which clinical signs are most predictive of suctioning need.¹³

Consequences of under-suctioning:

- **Partial occlusion:** Increased airway resistance, which may lead to increased work of breathing and potentially delayed weaning.³³
- **Complete occlusion:** Medical emergency requiring immediate intervention³³

Adverse consequences of suctioning:

- Greater risk of ventilator associated pneumonia³⁴
- Pain and distress³⁵
- Negative end-expiratory pressure,³⁶ leading to:
 - Reduction in end expiratory lung volume regardless of suctioning method³⁷⁻³⁹
 - Alveolar derecruitment^{38, 39}
 - Hypoxemia³⁹
- Excessively high expiratory tidal volume and airway pressures following suctioning in patients receiving volume guaranteed ventilation⁴⁰
- Hemodynamic instability
 - Decreases in respiratory rate⁴¹
 - Lower heart rate⁴¹
 - Prolonged increases of mean cerebral blood flow velocity⁴²

CHEST RADIOGRAPHS

The American Academy of Pediatrics lists daily chest x-rays as one of top five tests or treatments that lack evidence of efficacy and therefore contribute to “unnecessary utilization of staffing or material resources.”⁴³ Though chest x-rays are considered the “gold standard” for confirming tube placement or misplacement, they contribute to radiation exposure, can be delayed, and require manipulation of the patients, which may move the tube.⁴⁴

- A quality initiative published by Ridore et al. to restrict the frequency of chest x-rays found that they could reduce the number of x-rays by nearly half, saving almost \$1.6m annually, without an increase in the number unplanned extubations.⁴⁵
- As mentioned above, Klugmen et al. considers chest x-rays to one of the procedures that is associated with a high risk of unplanned extubations.²⁰
- Finally, a prospective study documenting the frequency that neonatal patients undergo painful or stressful procedures found that x-rays account for 6% of all stressful procedures.¹⁵

MINIMIZING STRESS

SAFE KANGAROO CARE

Kangaroo mother care, which is defined as skin-to-skin contact and frequent breastfeeding, has been demonstrated to be associated with multiple positive outcomes in low birthweight infants, including reduced mortality and nosocomial infections as well as increased weight gain and growth.⁴⁶ Unfortunately, kangaroo care, due to significant patient repositioning, is also considered a risk factor for UEs.^{19,20}

In fact, a survey of 1,122 nurse managers in hospitals in the United States found that 77% of respondents consider fear of UE to be a significant barrier to kangaroo care implementation.⁴⁷



REFERENCES

1. Carvalho FL, Mezzacappa MA, Calil R, Machado Hda C. Incidence and risk factors of accidental extubation in a neonatal intensive care unit. *J Pediatr (Rio J)*. May-Jun 2010;86(3):189-95. doi:10.2223/jped.1999
2. Dargaville PA, Gerber A, Johansson S, et al. Incidence and Outcome of CPAP Failure in Preterm Infants. *Pediatrics*. Jul 2016;138(1)doi:10.1542/peds.2015-3985
3. Fanaroff AA, Stoll BJ, Wright LL, et al. Trends in neonatal morbidity and mortality for very low birthweight infants. *Am J Obstet Gynecol*. Feb 2007;196(2):147.e1-8. doi:10.1016/j.ajog.2006.09.014
4. Hatch LD, 3rd, Scott TA, Slaughter JC, et al. Outcomes, Resource Use, and Financial Costs of Unplanned Extubations in Preterm Infants. *Pediatrics*. Jun 2020;145(6)doi:10.1542/peds.2019-2819
5. Moya FR, Mazela J, Shore PM, et al. Prospective observational study of early respiratory management in preterm neonates less than 35 weeks of gestation. *BMC Pediatr*. May 11 2019;19(1):147. doi:10.1186/s12887-019-1518-3
6. Stoll BJ, Hansen NI, Bell EF, et al. Neonatal outcomes of extremely preterm infants from the NICHD Neonatal Research Network. *Pediatrics*. Sep 2010;126(3):443-56. doi:10.1542/peds.2009-2959
7. Wei JL, Bond J. Management and prevention of endotracheal intubation injury in neonates. *Curr Opin Otolaryngol Head Neck Surg*. Dec 2011;19(6):474-7. doi:10.1097/MOO.0b013e32834c7b5c
8. Morris HF, Schuller L, Archer J, et al. Decreasing Unplanned Extubation in the Neonatal ICU With a Focus on Endotracheal Tube Tip Position. *Respir Care*. Nov 2020;65(11):1648-1654. doi:10.4187/respcare.07446
9. Lange M, Jonat S, Nikischin W. Detection and correction of endotracheal-tube position in premature neonates. *Pediatr Pulmonol*. Dec 2002;34(6):455-61. doi:10.1002/ppul.10200
10. Pavlek LR, Dillard J, Ryszen G, Hone E, Shepherd EG, Moallem M. Short-term complications and long-term morbidities associated with repeated unplanned extubations. *J Perinatol*. Mar 2021;41(3):562-570. doi:10.1038/s41372-021-00927-9
11. Gonçalves RL, Tsuzuki LM, Carvalho MG. Endotracheal suctioning in intubated newborns: an integrative literature review. *Rev Bras Ter Intensiva*. Jul-Sep 2015;27(3):284-92. doi:10.5935/0103-507x.20150048
12. Hough JL, Barton J, Jardine LA. A quality appraisal using the AGREE II instrument of endotracheal tube suction guidelines in neonatal intensive care units. *Aust Crit Care*. Mar 19 2021;doi:10.1016/j.aucc.2021.02.001
13. Tume LN, Copnell B. Endotracheal Suctioning of the Critically Ill Child. *J Pediatr Intensive Care*. Jun 2015;4(2):56-63. doi:10.1055/s-0035-1556747
14. Cong X, Wu J, Vittner D, et al. The impact of cumulative pain/stress on neurobehavioral development of preterm infants in the NICU. *Early Hum Dev*. May 2017;108:9-16. doi:10.1016/j.earlhumdev.2017.03.003
15. Carbajal R, Rousset A, Danan C, et al. Epidemiology and treatment of painful procedures in neonates in intensive care units. *Jama*. Jul 2 2008;300(1):60-70. doi:10.1001/jama.300.1.60
16. Hatfield LA, Murphy N, Karp K, Polomano RC. A Systematic Review of Behavioral and Environmental Interventions for Procedural Pain Management in Preterm Infants. *J Pediatr Nurs*. Jan-Feb 2019;44:22-30. doi:10.1016/j.pedn.2018.10.004
17. McPherson C, Ortinau CM, Vesoulis Z. Practical approaches to sedation and analgesia in the newborn. *J Perinatol*. Mar 2021;41(3):383-395. doi:10.1038/s41372-020-00878-7
18. Powell BM, Gilbert E, Volsko TA. Reducing Unplanned Extubations in the NICU Using Lean Methodology. *Respir Care*. Dec 2016;61(12):1567-1572. doi:10.4187/respcare.04540
19. Lucas da Silva PS, de Carvalho WB. Unplanned extubation in pediatric critically ill patients: a systematic review and best practice recommendations. *Pediatr Crit Care Med*. Mar 2010;11(2):287-94. doi:10.1097/PCC.0b013e3181b80951
20. Klugman D, Melton K, Maynard PO, et al. Assessment of an Unplanned Extubation Bundle to Reduce Unplanned Extubations in Critically Ill Neonates, Infants, and Children. *JAMA Pediatr*. Jun 1 2020;174(6):e200268. doi:10.1001/jamapediatrics.2020.0268
21. <https://www.solutionsforpatientsafety.org/about-us/our-goals/>.
22. Volsko TA, McNinch NL, Prough DS, Bigham MT. Adherence to Endotracheal Tube Depth Guidelines and Incidence of Malposition in Infants and Children. *Respir Care*. Sep 2018;63(9):1111-1117. doi:10.4187/respcare.06024
23. Chung HW, Lee WT, Chen HL. Reexamining the ideal depth of endotracheal tube in neonates. *Pediatr Neonatol*. Jun 2018;59(3):258-262. doi:10.1016/j.pedneo.2017.10.001
24. Mahajan A, Hoftman N, Hsu A, Schroeder R, Wald S. Continuous monitoring of dynamic pulmonary compliance enables detection of endobronchial intubation in infants and children. *Anesth Analg*. Jul 2007;105(1):51-6. doi:10.1213/01.ane.0000268119.55909.b4
25. Razak A, Faden M. Methods for Estimating Endotracheal Tube Insertion Depth in Neonates: A Systematic Review and Meta-Analysis. *Am J Perinatol*. Jan 30 2020;doi:10.1055/s-0039-3402747
26. Chowdhry R, Dangman B, Pinheiro JM. The concordance of ultrasound technique versus X-ray to confirm endotracheal tube position in neonates. *J Perinatol*. Jul 2015;35(7):481-4. doi:10.1038/jp.2014.240
27. Pinheiro JM, Munshi UK. Factors contributing to endobronchial intubation in neonates. *Pediatr Crit Care Med*. Jan 2015;16(1):54-8. doi:10.1097/pcc.0000000000000270
28. Thayyil S, Nagakumar P, Gowers H, Sinha A. Optimal endotracheal tube tip position in extremely premature infants. *Am J Perinatol*. Jan 2008;25(1):13-6. doi:10.1055/s-2007-995221
29. Bartle RM, Miller AG, Diez AJ, Smith PB, Gentile MA, Puia-Dumitrescu M. Evaluating Endotracheal Tube Depth in Infants Weighing Less Than 1 Kilogram. *Respir Care*. Mar 2019;64(3):243-247. doi:10.4187/respcare.06202
30. Hu X, Zhang Y, Cao Y, Huang G, Hu Y, McArthur A. Prevention of neonatal unplanned extubations in the neonatal intensive care unit: a best practice implementation project. *JBI Database System Rev Implement Rep*. Nov 2017;15(11):2789-2798. doi:10.11124/jbisrir-2016-003249
31. Kandil SB, Emerson BL, Hooper M, et al. Reducing Unplanned Extubations Across a Children's Hospital Using Quality Improvement Methods. *Pediatr Qual Saf*. Nov-Dec 2018;3(6):e114. doi:10.1097/pq9.0000000000000114
32. Boqué MC, Gualis B, Sandiumenge A, Rello J. Endotracheal tube intraluminal diameter narrowing after mechanical ventilation: use of acoustic reflectometry. *Intensive Care Med*. Dec 2004;30(12):2204-9. doi:10.1007/s00134-004-2465-4
33. Shah C, Kolfel MH. Endotracheal tube intraluminal volume loss among mechanically ventilated patients. *Crit Care Med*. Jan 2004;32(1):120-5. doi:10.1097/01.Ccm.0000104205.96219.D6
34. Yuan TM, Chen LH, Yu HM. Risk factors and outcomes for ventilator-associated pneumonia in neonatal intensive care unit patients. *J Perinat Med*. 2007;35(4):334-8. doi:10.1515/jpm.2007.065
35. Väilitalo PAJ, van Dijk M, Krekels EJJ, et al. Pain and distress caused by endotracheal suctioning in neonates is better quantified by behavioural than physiological items: a comparison based on item response theory modelling. *Pain*. Aug 2016;157(8):1611-1617. doi:10.1097/j.pain.0000000000000485
36. Nakstad ER, Opdahl H, Heyerdahl F, Borchsenius F, Skjønberg OH. Manual ventilation and open suction procedures contribute to negative pressures in a mechanical lung model. *BMJ Open Respir Res*. 2017;4(1):e000176. doi:10.1136/bmjresp-2016-000176
37. Copnell B, Dargaville PA, Ryan EM, et al. The effect of suction method, catheter size, and suction pressure on lung volume changes during endotracheal suction in piglets. *Pediatr Res*. Oct 2009;66(4):405-10. doi:10.1203/PDR.0b013e3181b337b9
38. Lindgren S, Odenstedt H, Erlandsson K, Grivans C, Lundin S, Stenqvist O. Bronchoscopic suctioning may cause lung collapse: a lung model and clinical evaluation. *Acta Anaesthesiol Scand*. Feb 2008;52(2):209-18. doi:10.1111/j.1399-6576.2007.01499.x
39. Maggiore SM, Lellouche F, Pigeot J, et al. Prevention of endotracheal suctioning-induced alveolar derecruitment in acute lung injury. *Am J Respir Crit Care Med*. May 1 2003;167(9):1215-24. doi:10.1164/rccm.200203-1950C
40. Kiraly NJ, Tingay DG, Mills JF, Dargaville PA, Copnell B. Volume not guaranteed: closed endotracheal suction compromises ventilation in volume-targeted mode. *Neonatology*. 2011;99(1):78-82. doi:10.1159/000316854
41. Barbosa AL, Cardoso MV, Brasil TB, Scocchi CG. Endotracheal and upper airways suctioning: changes in newborns' physiological parameters. *Rev Lat Am Enfermagem*. Nov-Dec 2011;19(6):1369-76. doi:10.1590/s0104-11692011000600013
42. Kaiser JR, Gauss CH, Williams DK. Tracheal suctioning is associated with prolonged disturbances of cerebral hemodynamics in very low birth weight infants. *J Perinatol*. Jan 2008;28(1):34-41. doi:10.1038/sj.jp.7211848
43. Ho T, Dukhovny D, Zupancic JA, Goldmann DA, Horbar JD, Pursley DM. Choosing Wisely in Newborn Medicine: Five Opportunities to Increase Value. *Pediatrics*. Aug 2015;136(2):e482-9. doi:10.1542/peds.2015-0737
44. Zaytseva A, Kurepa D, Ahn S, Weinberger B. Determination of optimal endotracheal tube tip depth from the gum in neonates by X-ray and ultrasound. *J Matern Fetal Neonatal Med*. Jun 2020;33(12):2075-2080. doi:10.1080/14767058.2018.1538350
45. Ridore M, Pastor W, Bulas D, McKenney S, Soghier L, Short B. Initiation and Compliance with a Chest X-Ray (CXR) Reduction Protocol in the Children's National Medical Center Neonatal Intensive Care Unit (NICU). *Pediatrics*. 2018;142(1 MeetingAbstract):221.
46. Conde-Agudelo A, Diaz-Rossello JL. Kangaroo mother care to reduce morbidity and mortality in low birthweight infants. *Cochrane Database Syst Rev*. Aug 23 2016;2016(8):Cd002771. doi:10.1002/14651858.CD002771.pub4
47. Engler AJ, Ludington-Hoe SM, Cusson RM, et al. Kangaroo care: national survey of practice, knowledge, barriers, and perceptions. *MCN Am J Matern Child Nurs*. May-Jun 2002;27(3):146-53. doi:10.1097/00005721-200205000-00004

© 2021 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 07/2021-US-RE-2100209-[WF#5304052]

6135 Gunbarrel | Avenue Boulder, CO 80301 | 800.635.5267 [medtronic.com/covidien](https://www.medtronic.com/covidien)

Medtronic

Table of contents

Managing endotracheal tubes in neonatal patients

SonarMed™ airway monitoring system overview

Endotracheal tube related adverse events

Minimizing stress



12