

Clinical data “simplified”

The data proving significant, safe, and sustained blood pressure reductions.<sup>1-4</sup>

The SPYRAL clinical program

**SPYRAL HTN-OFF MED**

Pivotal, randomized, sham-controlled trial to show safety and efficacy of the Symplcity Spyrall RDN system in uncontrolled hypertension patients in the absence of antihypertensive medications (n = 331).<sup>3</sup>

**SPYRAL HTN-ON MED**

Randomized, sham-controlled trial to show safety and efficacy of the Symplcity Spyrall RDN system in uncontrolled hypertension patients in the presence of antihypertensive medications (n = 337).<sup>4</sup>


**Global SYMPLICTY Registry**

Prospective, single-arm, all-comer, real-world registry to demonstrate long-term durability and safety in patients with uncontrolled hypertension treated with the Symplcity RDN system (n = 3,300).<sup>1,2</sup>



9 mmHg mean reduction


in office systolic blood pressure (SBP) in patients on and off medications at primary endpoint follow-up.



Results may vary.



Excellent safety profile with low adverse event rates and stable kidney function




Results may vary.




17 mmHg mean reduction

in office SBP real-world patients at 3 years<sup>1,2</sup>



Blood pressure reductions across high-risk subgroups at 3 years<sup>†</sup>



Results may vary.

† Includes Symplcity Spyrall and Flex catheters.  
‡ Resistant hypertension defined as OSBP >150 mmHg, ≥3 antihypertensive medications.  
§ CKD defined as eGFR < 60 mL/min/1.73m<sup>2</sup>.  
¶ Elderly defined as 65 years or older.  
1. Mahfoud F, Kandzari DE, Kario K, et al. Long-term efficacy and safety of renal denervation in the presence of antihypertensive drugs (SPYRAL HTN-ON MED): a randomised, sham-controlled trial. *Lancet*. April 9, 2022;399(10234):1401-1410.  
2. Mahfoud F, et al. Outcomes following radiofrequency renal denervation according to antihypertensive medications: subgroup analysis of the Global SYMPLICTY Registry DEFINE. *EuroPCR* 2023.  
3. Böhm M, Kario K, Kandzari DE, et al. Efficacy of catheter-based renal denervation in the absence of antihypertensive medications (SPYRAL HTN-OFF MED Pivotal): a multicentre, randomised, sham-controlled trial. *Lancet*. May 2, 2020;395(10234):1444-1451.  
4. Kandzari D, Townsend R, Kario K, et al. Safety and Efficacy of Renal Denervation in Patients Taking Antihypertensive Medications. *J Am Coll Cardiol*. November 7, 2023;82(19):1809-1823.  
5. Schlaich M, et al. Long-term safety and efficacy of renal denervation with the Symplcity Spyrall catheter in the Global SYMPLICTY Registry. Presented at American Society of Nephrology Kidney Week, San Diego, CA, November 4-7, 2021.  
6. Mahfoud F, Mancia G, Schmieder R, et al. Renal Denervation in high-risk patients with hypertension. *J Am Coll Cardiol*. June 16, 2020;75(23):2879-2888.  
7. Mahfoud, F, et al. Blood pressure reduction after catheter-based renal denervation in patients with cardiovascular disease in the Global SYMPLICTY Registry. *ESH* 2022, Athens.

**Indications**  
The Symplcity Spyrall™ renal denervation system is indicated to reduce blood pressure as an adjunctive treatment in patients with hypertension in whom lifestyle modifications and antihypertensive medications do not adequately control blood pressure.

**Contraindications**  
The Symplcity Spyrall system is contraindicated in patients with any of the following conditions: • Renal artery diameter < 3 mm or > 8 mm • Renal artery fibromuscular dysplasia (FMD) • Stented renal artery (< 3 months prior to RDN procedure) • Renal artery aneurysm • Renal artery diameter stenosis > 50% • Pregnancy • Presence of abnormal kidney (or secreting adrenal) tumor • Iliac/femoral artery stenosis precluding insertion of the catheter.

**Warnings and Precautions**  
A thorough understanding of the technical principles, clinical applications, and risks associated with vascular access techniques and percutaneous transluminal catheterization in renal arteries is necessary before using this device. The safety and efficacy of the Symplcity Spyrall system has not been established in patients with isolated systolic hypertension or in patients with prior renal artery interventions including renal stents, renal angioplasty, or prior renal denervation. The Symplcity Spyrall system has not yet been studied in patients who are breastfeeding, under the age of 18, or with secondary hypertension. • Avoid treatment with the Symplcity Spyrall™ catheter within 5 mm of any diseased area or stent. • Implantable pacemakers (IPGs) and implantable cardioverter defibrillators (ICDs) or other active implants may be adversely affected by RF ablation. Refer to the implantable device's Instructions for Use. • The patient's heart rate may drop during the ablation procedure. • Proper pain medication should be administered at least 10 min before ablating renal nerves.

**Potential Adverse Events**  
Potential adverse events associated with use of the renal denervation device or the interventional procedures include, but are not limited to, the following conditions: • Allergic reaction to contrast • Arterial damage, including injury from energy application, dissection, or perforation • Arterial spasm or stenosis • Arterio-enteric fistula • AV fistula • Bleeding • Blood clots or embolism • Bruising • Cardiopulmonary arrest • Complications associated with medications commonly utilized during the procedure, such as narcotics, anxiolytics, or other pain or anti-spasmodic medications • Death • Deep vein thrombosis • Edema • Electrolyte imbalance • Heart rhythm disturbances, including bradycardia • Hematoma • Hematoma – retroperitoneal • Hematuria • Hypertension • Hypotension (may cause end organ hypoperfusion) • Infection • Kidney damage including renal failure or perforation • Myocardial infarction • Nausea or vomiting • Pain or discomfort • Peripheral ischemia • Pulmonary embolism • Proteinuria • Pseudoaneurysm • Radiocontrast nephropathy • Renal artery aneurysm • Skin burns from failure of the dispersive electrode pad • Stroke • Other potential adverse events that are unforeseen at this time.

Please reference appropriate product Instructions for Use and User Manual for more information regarding indications, contraindications, warnings, precautions, and potential adverse events.

**CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician.

For further information, please call and/or consult Medtronic at 800-633-8766 or the Medtronic website at medtronic.com.