

Medtronic

Global value dossier for
LigaSure™ technology



Preclinical data and summary

Prepared by: **Coreva Scientific**
Im Mühlenbruch 1
53227 Königswinter, Germany
Phone: +49 (0) 2223 781 8010, ext. 53
E-mail: jh@coreva-scientific.com

Version No. 7.0

Date: November 2022

©2022 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.
11/2022 - US-SE-2200137

Contents

1	Disclaimer and scope.....	4
1.1	Products and aliases	4
1.2	Data sources	4
1.3	Analysis.....	5
2	Introduction to LigaSure™ technology.....	6
2.1	Overview	6
2.2	Prelude	7
2.3	The surgical application: Hemostasis	7
2.3.1	Various methods of controlling blood loss.....	7
2.3.2	Electrosurgery.....	8
2.4	LigaSure™ technology	9
2.4.1	Energy platform	10
2.4.2	Preclinical and laboratory data	10
2.4.3	Nano-coating.....	13
2.4.4	A history of innovation	14
2.4.5	Surgical indications	15
2.5	Extensive worldwide usage.....	15
2.6	Dossiers regarding LigaSure™ devices in procedures	16
3	Key messages for LigaSure™ devices in clinical applications	17
4	Source data tables	19
4.1	Data table summary: LigaSure™ technology preclinical data.....	19
5	Structured literature search details.....	23
5.1	Searches performed	23
5.2	Screening literature results	27
5.2.1	Screening general search (EMBASE) results	27
5.2.2	Screening returned results for evidence (PubMed)	27
6	Overview: Indications of LigaSure™ devices	29

7	References.....	33
8	Geographical references.....	42

List of Tables

Table 3-1	General LigaSure™ technology outcomes across five indications (colectomy, hysterectomy, thyroidectomy, breast surgery and gastric surgery).....	17
Table 4-1	Publication data on LigaSure™ technology in preclinical settings.....	19
Table 5-1	Structured searched in PubMed to identify relevant LigaSure™ technology data (original searches in 2017).....	23
Table 5-2	Structured searched in PubMed to identify relevant LigaSure™ technology data (searches for new indications in 2022).....	25
Table 5-3	Structured search in EMBASE to identify use of LigaSure™ devices	27
Table 5-4	Description of exclusion criteria.....	28
Table 6-1	LigaSure™ devices and their indications.....	29

List of Figures

Figure 2-1	Methods of surgical hemostasis.....	8
Figure 2-2	Monopolar and bipolar electrosurgery.....	9
Figure 2-3	LigaSure™ device supraphysiological burst pressure measurements	11
Figure 2-4	Vessel sealing times in experimental models	12
Figure 2-5	Histological measures of extent of effect for LigaSure™ devices.....	13
Figure 2-6	Timeline of introduction of LigaSure™ devices and supporting technology	14
Figure 2-7	LigaSure™ devices have been applied in clinical settings globally	16

1 Disclaimer and scope

1.1 Products and aliases

This document refers to devices making use of LigaSure™ technology. This technology includes, for laparoscopic procedures, the LigaSure™ 5mm blunt tip device (37 cm and 44 cm), the LigaSure™ Maryland jaw device (37 cm and 44 cm), the LigaSure Advance™ Pistol Grip device, the LigaSure Atlas™ device and the LigaSure™ dolphin tip device. For open surgeries, applicable devices include the LigaSure Precise™ device, the LigaSure Impact™ device, the LigaSure™ Small Jaw device, the LigaSure™ 5mm blunt tip device (23 cm), the LigaSure™ Maryland jaw device (23 cm), the LigaSure Atlas™ device (20 cm), the LigaSure™ dolphin tip device (20 cm), the LigaSure™ curved jaw device, the LigaSure™ Exact device, and other reusable LigaSure™-branded instruments. LigaSure™ devices are available with uncoated jaws or with nonstick jaws coated with hexamethyldisiloxane (HMDSO).

LigaSure™ Exact is the only LigaSure™ device cleared for breast procedures, and only in the USA. Indications, contraindications, warnings, precautions and procedure steps may vary between products and models, and availability may vary by jurisdiction. Please always refer to indication labeling for your jurisdiction and read all applicable instructions for use provided with the products. An overview of all LigaSure™ devices, their indications and the countries/regions where these indications apply is provided in chapter 6 Overview: Indications of LigaSure™ devices.

1.2 Data sources

Data regarding the clinical application of LigaSure™ devices were derived from searches of literature published within the last 14 years in PubMed (last searches performed in September 2021 and January 2022) and EMBASE (last search performed in March 2022). General text searches for LigaSure™ device references were performed in EMBASE to include all published instances of the use of the technology without restriction of time of publication. These results provide an overview of the extent of usage of the technology.

For clinical evidence data of LigaSure™ technology usage, structured searches were performed using PubMed for peer-reviewed literature applying consistent exclusion criteria across searches for specific surgical indications (see Section 5, Structured literature search details). Results were restricted to publications based on data obtained within the most recent 14 years of the search (2007 and onwards), and excluded editorials/commentaries, articles with no abstract, those that did not report

relevant clinical data (such as animal or ex vivo studies), those that reported data on fewer than 20 patients, those which were not focused on outcomes related to the technology (that is, no mention of vessel sealing or hemostasis), and articles which did not reference LigaSure™ devices or generic LigaSure™ technology (electrosurgical or radiofrequency bipolar vessel sealing).

1.3 Analysis

Clinical results from individual studies are presented as reported (including indication of statistical significance where determined). Where data are amalgamated from multiple reports, please note that the individual studies will vary in terms of design, protocol, surgical technique and patient population, which may limit conclusions drawn from direct comparison and relevant analysis of statistical significance. The resulting figures, however, provide insight into clinical outcomes that have been achieved using LigaSure™ devices in vessel sealing during surgical procedures.

2 Introduction to LigaSure™ technology

2.1 Overview

Summary and key messages

1. **Answers unmet surgical need:** Blood loss is still a major concern during surgical procedures^{1,2} and LigaSure™ technology enables efficient, electrosurgical division and sealing of tissues where the control of blood loss is required.³
2. **History of development:** LigaSure™ devices span more than 20 years of continuous development and improvement to meet current healthcare needs.
3. **Nano-coating:** Nonstick coating on the jaws of the LigaSure™ device reduces eschar build-up and sticking, and increases cleanability compared to uncoated or reprocessed LigaSure™ devices.^a Nano-coated LigaSure™ device jaws also show a significantly lower sticking rate than similar competitive devices.^b
4. **Extensive worldwide usage:** LigaSure™ technology has over 1,700 published reports of its use in a broad range of surgical procedures spanning 66 countries around the world.^{G1-G401}
5. **Excellent preclinical performance data:** LigaSure™ devices have been tested in preclinical animal⁴⁻⁹ and ex vivo human tissue^{10,11} models, demonstrating sealing times of 10 seconds or faster with the LVS generator,^{3,12-14} or 2-4 seconds with the ForceTriad™ energy platform.^{13,14} The applicator provides seal quality resistant to supraphysiological burst pressures for vessels ranging from small (2-3 mm)^{7,10,12,14} to large (up to and including 7 mm).^{4,5,7,8,14}

^a Based on internal test report #RE00065124, LF4418 Report Marketing Claims for Reprocessed Jaws. Oct. 4 to Oct. 7, 2016 and internal test report #RE00065374, LF18XX Product Claims Report for Reprocessed Jaws. Sep. 22 to Sep. 29, 2016. Bench test data may not be indicative of clinical results.

^b Based on internal test report #RE00073194, Market Research - Report - Tissue Sticking Comparison of the Ethicon® G2, Voyant® 5 mm Fusion, LigaSure™ LF1737, and LigaSure™ LF1937 Devices. Jan. 18, 2017.

2.2 Prelude

LigaSure™ technology is used in surgical procedures to divide and seal vessels up to and including 7 mm in diameter. The devices are electrosurgical in nature, using current delivered to patient tissues to effect tissue sealing. The complete LigaSure™ vessel sealing system comprises the vessel sealing device (the LigaSure™ sealer/divider), and the energy platform (such as the ForceTriad™ energy platform, the Valleylab™ LS10 generator or the Valleylab™ FT10 energy platform). The vessel sealing device delivers a combination of pressure and electrical current to tissues, the current is provided in a smart algorithm by the energy platform, using tissue-sensing technology (TissueFect™ technology^c) to control energy delivery. The technology is suitable for use wherever the division and ligation of vessels is desired during general open or minimally invasive surgical procedures. This document presents background information, and associated dossiers provide clinical evidence demonstrating the application of LigaSure™ technology in improving surgical procedures and patient outcomes.

2.3 The surgical application: Hemostasis


In surgery, for the patient's health, as well as to ensure maintenance of adequate visualization, the operative field must stay free of excess blood and other fluids. This is accomplished by hemostatic techniques, whereby dissection is only carried out in blood-vessel-free planes, and when necessary, hemostatic devices are used.

2.3.1 Various methods of controlling blood loss

Methods for controlling bleeding during surgery vary in terms of the immediacy and desired length of action of their effect. A general, non-comprehensive overview of some hemostatic methods is shown in Figure 2-1. The more durable means of preventing excessive blood loss include making use of intrinsic components from the patient's own tissues, from providing a matrix to support the natural clotting process, to electrosurgical techniques which alter the vessel collagen to collapse or fuse the vessel lumen to create a permanent seal. In the case of some electrosurgical instruments, tissue/vessel division is performed simultaneously to ligation, thereby preventing blood loss from the separated vessel.

^c White paper, Valleylab™ FT10 Energy Platform TissueFect™ Technology. McHenry J, Dunning J and Waggle K (2015) 10/2015 US150755 [REF#479324]. Bench test data may not be indicative of clinical results.

Figure 2-1 Methods of surgical hemostasis



Shown are several options (a non-comprehensive list) for hemostasis that can be intraoperatively deployed. Placement along the timeline indicates a typical expected duration of effect, assuming no failure. The intrinsic methods make use of body components, such as the natural clotting process captured in the matrix of the biological agents indicated, or the fusion of tissues/structures for the mechanical. Extrinsic methods rely on external structures to provide the seal, which may be subsequently removed (such as pads) or may be absorbed (as in the case of some biologics). Adapted from references.^{1,2}

2.3.2 Electrosurgery

As mentioned (section 2.3.1), electrosurgery is one method of hemostasis. These technologies create an electrical current through the patient to generate localized bursts of heat to coagulate tissue. The electricity oscillates at radiofrequencies (hence sometimes referenced as radiofrequency or RF electrocoagulation), which are well above frequencies for normal nerve function, preventing devices from triggering muscle stimulation.

The waveform of energy delivery determines the function of the device:


- Cut mode: dissects tissue and uses a low voltage with a long duty cycle
- Coagulation mode: cauterizes tissue using a higher voltage with a shorter duty cycle
- Blend mode: uses intermediate voltage and duty cycle length to achieve both dissection and coagulation of the tissues simultaneously with dissection.

The energy from these waveforms is delivered to patient tissues via:

- Monopolar delivery: the electric circuit passes from the generator through the patient, and back to the source via a separately positioned return electrode (Figure 2-2A)

- Bipolar delivery: the electric circuit passes only through tissue positioned between graspers of the device, obviating potential issues of return electrode misplacement or detachment (Figure 2-2B)

Figure 2-2 Monopolar and bipolar electrosurgery



Monopolar and bipolar forms of electrosurgery differ in the paths taken by the current. In the schematic on the left (A) monopolar electrosurgery with monitoring of the return electrode is shown. Current flows along the magenta arrows from the source (lower left) and through the patient, returning via the patient return electrode (PRE, blue rectangle). In contrast, for bipolar surgery (B), the current passes only through tissue grasped between the forceps, reducing risks to the patient, including from misplaced or detached PREs. Illustration adapted from Golpaygani et al. 2016.¹⁵

2.4 LigaSure™ technology

LigaSure™ technology is based on bipolar electrosurgery, commonly referenced as electrothermal bipolar vessel sealing (EBVS). The complete system for the sealing of vessels during surgery (previously known as the LigaSure™ Vessel Sealing System, or LVSS) comprises the bipolar device, the energy platform (including the ForceTriad™ energy platform, the Valleylab™ LS10 generator, or the Valleylab™ FT10 energy platform), and the TissueFect™ sensing system, which, using data from the applicator device analyzed at 3,333 to 434,000 data points per second (dependent on the generator), senses the properties of the tissue being grasped to optimize vessel sealing performance^d. The LS10 generator and the ForceTriad™ and FT10 energy platforms include the TissueFect™ technology. The sealing occurs through the application of a potential difference between the electrodes of the grasping device, which denatures collagen and elastin within the vessel wall and surrounding connective tissue. The addition of extreme pressure applied by the instrument causes

^d White paper, Valleylab™ FT10 Energy Platform TissueFect™ Technology. McHenry J, Dunning J and Wagle K (2015) 10/2015 US150755 [REF#479324]. Bench test data may not be indicative of clinical results.

the denatured and reformed collagen and elastin in the tissue bundle or blood vessel to form a hemostatic plug or seal. As a method engaged for hemostasis (section 2.3.1), use of the system is intended to provide an alternative to clips or suture for vessel closure during and after surgery.

2.4.1 Energy platform

The current required for the tissue dividing and vessel sealing functions of LigaSure™ devices is provided by generators designed to provide multiple waveforms for any type of electrosurgery. LigaSure™ device-compatible systems include the ForceTriad™ energy platform and the newer generation Valleylab™ LS10 generator and Valleylab™ FT10 energy platform. These devices (ForceTriad™ and FT10 energy platforms) provide modes to deliver monopolar, bipolar, and LigaSure™ technology-specific waveforms (the LS10 generator provides only the LigaSure™ technology waveform). In its contact with tissue, the LigaSure™ device provides energy delivery and acts as a sensor. The device provides data to the energy delivery platform and is subjected to feedback regulation via the TissueFect™ system. The system provides data to the instrument regarding the resulting effects on the tissue. In LigaSure™ technology mode, the system can be set to automatically complete the sealing, registering an auditory signal when the cycle is complete. A preclinical comparison of a LigaSure™ device as part of the older generation LVSS and the same device with the ForceTriad™ energy platform in blood vessels of various type and size from a porcine model found the use of the ForceTriad™ energy platform to be associated with statistically significantly shorter sealing times ($p < 0.05$) and higher burst pressures ($p < 0.01$), although for both power sources, the LigaSure™ device achieved supraphysiological burst strengths.¹³ The latest version of the power generator, the FT10 energy platform, demonstrates still further improvements on the ForceTriad™.^d energy platform.

2.4.2 Preclinical and laboratory data


To minimize patient risk, and to enable analysis of parameters that would not be possible in the live patient, experimental studies to assess instrument performance have been performed. These consist of animal studies (mostly porcine) of various types of vessels, and ex vivo human tissue, immediately after tissue removal or post-fixation.^e Details of studies can be found in Table 4-1.

^e These studies may not be indicative of clinical performance.

2.4.2.1 Burst pressure

A key parameter not readily measurable in vivo in patients is the pressure to which newly divided and sealed vessels can be subjected before seal failure, or the burst pressure Figure 2-3. Across numerous analyses, in vessels extracted from human patients perioperatively (mesenteric vessels)¹⁰ and in animal studies of vessels (porcine,^{4,5,7-9,14,16} bovine³ and ovine¹⁷), LigaSure™ devices have achieved burst pressures well above normal physiological pressures.

Figure 2-3 LigaSure™ device supraphysiological burst pressure measurements




Data are compiled from several comparable studies in animal and tissue and one study in human tissue removed during surgery (mesenteric vessels¹⁰) and examined ex vivo. Animal data for small vessels are all from porcine models,^{5,7-9,14} while animal large-vessel data consist of porcine,^{4,5,7,8,14} bovine³ and ovine.¹⁷ Bars represent weighted means across studies of the ratio of the mean reported burst pressure as a multiple of physiological pressure (taken as 250 mm for blood vessels), N is shown for the total number of individual measurements, and error is indicated as standard error of the mean (SEM).

2.4.2.2 Sealing time

In addition to the security of the created seal, the average time for the device to complete the sealing cycle is another operative parameter for a user to consider. The latter is in part dependent on the power generator used, where across different sizes of vessels tested in animal models from 2-7 mm, shorter sealing times have been observed using the ForceTriad™ energy platform^{13,14} with LigaSure™ devices in comparison to the original LigaSure™ Vessel Sealing System Generator^{3,9,13,14,17} (Figure 2-4). Direct comparisons in the same study of the ForceTriad™ energy platform and

the earlier version generator, statistically significantly shorter ($p < 0.05$) seal times were observed. Even given the performance improvement observed with the ForceTriad™ energy platform over the original generator, the newer version FT10 demonstrates still further reductions in seal time over the ForceTriad™ energy platform, illustrating the continued improvement in performance in the LigaSure™ technology energy platform^f.

Figure 2-4 Vessel sealing times in experimental models



Although multiple factors will influence the vessel seal time, the power generator used can also have a significant impact on the time to complete the vessel seal using LigaSure™ devices. In animal models, the ForceTriad™ energy platform^{13,14} has demonstrated shorter sealing times than with the original version LigaSure™ Vessel Sealing System generator (LVSGen)^{3,9,13,14,17} ($p < 0.001$, two-tailed t-test), although both have mean times under 10 seconds. Data aggregated from animal studies (LVSGen: porcine,^{9,12-14} ovine,¹⁷ and bovine³; ForceTriad™ technology, porcine^{13,14}) only where the generator used was described. In direct head-to-head comparisons of the LVSGen and ForceTriad™ technology in porcine vessels of various sizes, a significant reduction in seal time was also observed^{13,14} ($p < 0.05$). Total counts of individual measurements are indicated on the bars, and error bars represent the standard error of the mean (SEM).


2.4.2.3 Thermal damage

A further consideration for use of LigaSure™ devices in clinical applications is the extent of thermal changes caused by the device to assess the potential for unintended damage that may occur in patient tissues. In the preclinical setting, the temperatures generated by the grasping jaws of the LigaSure™ device and the extent of tissue necrosis have been measured. At the point of sealing, mean maximum temperatures (T_{max}) of the instrument (measured by thermal camera imaging of the jaw surface at $85 \pm 3^\circ\text{C}$ ⁸) and of the cutting blade (103.1°C ¹²) have been measured in porcine models, while 2 mm from the device, a mean T_{max} of $55.5 \pm 2.4^\circ\text{C}$ has been determined in bovine tissue.³ Separately from temperatures, histological examination of tissue samples have been used to measure the range of effect (Figure 2-5) in human

^f White paper, Valleylab™ FT10 Energy Platform TissueFect™ Technology. McHenry J, Dunning J and Waggle K (2015) 10/2015 US150755 [REF#479324]. Bench test data may not be indicative of clinical results.

(mesenteric vessels from intestine¹⁰) and animal models (porcine carotid artery segments⁹ and multiple vessel types¹²). The weighted mean distance of tissue changes from the sealing point across these studies is 2.9 ± 0.3 mm.

Figure 2-5 Histological measures of extent of effect for LigaSure™ devices



Reported tissue effects (lateral surface spread and full thickness through all vessel layers) from histological examinations from different studies. Study identifier is indicated in parentheses beneath the model system: Hruby 2007,¹² Noble 2011,¹⁰ and Mantke 2011.⁹ Error bars correspond to standard deviation (no error value reported for Hruby 2007)

2.4.3 Nano-coating

LigaSure™ devices are available with uncoated jaws or with hexamethyldisiloxane (HMDSO) coating applied to the jaw and seal plate surfaces. This nonstick nano-coating enhances the performance of the LigaSure™ device.

Eschar build-up causes the jaws of the device to stick to the tissue. It can affect the performance of the device and require multiple cleanings of the device jaws, which affects procedural flow. LigaSure™ devices with nonstick coating help to prevent eschar build-up: the average eschar build-up is reduced by 15% in nano-coated devices compared to uncoated devices, and nonstick-coated device jaws can also be cleaned more effectively.⁹ The improved cleanability and ease of cleaning enable greater procedural efficacy.⁹

The nano-coated LigaSure™ devices are intended for single-use only, and the results in the previous paragraph are based on tests according to the expected conditions of a single surgical procedure. However, it is conceivable that third-party reprocessed devices have possibly been subjected to process steps, tools, and/or chemicals that may negatively affect the device's HMDSO nano-coating, which may result in an increase in tissue sticking force in reprocessed jaws compared to nano-coated jaws. In


⁹ Based on internal test report #RE00071599, LF19XX MJC marketing claims report conducted on porcine tissue. Feb. 7 to Feb. 22, 2017. Bench test data may not be indicative of clinical results.

fact, the HMDSO nano-coating of nonstick-coated LigaSure™ device jaws is degraded during standard reprocessing procedures, resulting in a significant increase in eschar build-up and sticking as well as reduced cleanability compared to LigaSure™ device jaws with intact nonstick coating.^h Finally, it has been shown that the rate of sticking is significantly lower in LigaSure™ devices with HMDSO coating, not just compared to uncoated LigaSure™ devices, but also compared to similar competing devices such as the Curved EnSeal® G2 and the Voyant™ 5mm Fusion device.ⁱ

2.4.4 A history of innovation

Since the introduction of the first LigaSure™ device in 1998, continual development has seen the release of multiple devices employing LigaSure™ technology for open and laparoscopic electro-surgical applications, spanning a wide array of indicated surgeries (Figure 2-6).

Figure 2-6 Timeline of introduction of LigaSure™ devices and supporting technology



^h Based on internal test report #RE00065124, LF4418 Report Marketing Claims for Reprocessed Jaws. Oct. 4 to Oct. 7, 2016 and internal test report #RE00065374, LF18XX Product Claims Report for Reprocessed Jaws. Sep. 22 to Sep. 29, 2016. Bench test data may not be indicative of clinical results.

ⁱ Based on internal test report #RE00073194, Market Research - Report - Tissue Sticking Comparison of the Ethicon® G2, Voyant® 5 mm Fusion, LigaSure™ LF1737, and LigaSure™ LF1937 Devices. Jan. 18, 2017.

2.4.5 Surgical indications

The LigaSure™ vessel sealing system (including the devices and energy delivery) is designed for use in minimally invasive (laparoscopic) or open procedures where ligation and division of vessels, tissue bundles, and lymphatics is desired. Note that different devices have different specific surgical indications, and these should always be confirmed in the instructions for use accompanying the device according to the local jurisdiction. Please refer to chapter 6 for an overview of all LigaSure™ devices and their indications.


2.5 Extensive worldwide usage

Publications referencing LigaSure™ devices demonstrating good performance in surgical settings originate from a wide range of countries (Figure 2-7). LigaSure™ technology was also included among medical supplies (“ligasure vessel sealing equipment, a steam sterilizer, obstetric labor table and intensive care ward beds”) provided by the World Health Organization to a camp in Damascus, Syrian Arab Republic in 2013 as part of its health response to the situation at the time^j. Overall, there is published evidence of LigaSure™ device use in:

Argentina,^{G1,G2} Australia,^{G3-G18} Austria,^{G19-G27} Bahamas,^{G28} Belgium,^{G29-G34} Brazil,^{G35-G37} Bulgaria,^{G38} Canada,^{G39-G47} Chile,^{G48,G49} China,^{G50-G64} Costa Rica,^{G65} Croatia,^{G66,G67} Czech Republic,^{G68-G71} Denmark,^{G72-G76} Egypt,^{G77-G92} Estonia,^{G93} Finland,^{G94-G96} France,^{G97-G115} Germany,^{G116-G127} Greece,^{G128-G133} Hong Kong,^{G134,G135} Hungary,^{G136} India,^{G137-G143} Iran,^{G144-G153} Iraq,^{G154,G155} Ireland,^{G156-G160} Israel,^{G161-G167} Italy,^{G113,G168-G187} Jamaica,^{G188} Japan,^{G189-G209} Jordan,^{G210-G212} Kuwait,^{G213} Lebanon,^{G34,G214-G216} Libya,^{G217} Malaysia,^{G218} Mexico,^{G219-G224} Mongolia,^{G225} Morocco,^{G226} Netherlands,^{G227-G232} New Zealand,^{G233} Norway,^{G234-G237} Oman,^{G238} Pakistan,^{G239-G241} Peru,^{G242} Poland,^{G243-G250} Portugal,^{G251,G252} Romania,^{G253-G260} Russia,^{G261-G267} Saudi Arabia,^{G268-G271} Serbia,^{G272-G274} Singapore,^{G275-G280} Slovakia,^{G281} Slovenia,^{G67} South Africa,^{G282,G283} South Korea,^{G284-G298} Spain,^{G299-G313} Sweden,^{G314,G315} Switzerland,^{G316,G317} Taiwan,^{G318-G330} Thailand,^{G331,G332} Trinidad and Tobago,^{G333} Turkey,^{G334-G362} Ukraine,^{G363,G364} United Kingdom,^{G18,G365-G377} United States of America^{G16,G378-G398} and Venezuela.^{G399-G401}

^j World Health Organization, Health Situation Report, Syrian Arab Republic, Jordan, Lebanon, Iraq, Issue 16, reporting period 21 June – 9 July, 2013; available at http://www.who.int/hac/crises/syr/sitreps/syria_regional_sitrep_21june_9july2013.pdf

Figure 2-7 LigaSure™ devices have been applied in clinical settings globally



Highlighted in blue are countries of origin for publications demonstrating the successful use of LigaSure™ technology in clinical settings across all surgeries for which the devices have been indicated.

2.6 Dossiers regarding LigaSure™ devices in procedures

In this series of dossiers, use of LigaSure™ devices in hysterectomy, colectomy, thyroidectomy, bariatric surgery, and breast surgery is discussed. An overview of the indications of the different LigaSure™ devices is detailed in chapter 6. It should be noted that breast surgery is only indicated for LigaSure™ Exact and only in the USA.

3 Key messages for LigaSure™ devices in clinical applications^k

LigaSure™ technology has also demonstrated excellent outcomes for patients in clinical settings. Key, common findings from analyses across hysterectomy, colorectal procedures (colectomy), thyroidectomy, breast surgery and gastric surgery are indicated below (Table 3-1). Individual dossiers for each surgery are available for further details.

Table 3-1 General LigaSure™ technology outcomes across five indications (colectomy, hysterectomy, thyroidectomy, breast surgery and gastric surgery)

General Extensive worldwide usage: LigaSure™ technology has over 1,300 published reports of its use in a broad range of surgical procedures spanning 66 countries around the world.^{G1-G401} Compared with conventional ligation: Decreased operative time compared with conventional ligation: Across studies reporting operating time (total 35 in the search window for the five procedures), LigaSure™ devices have on average reduced procedure time by 16.6% ($\pm 3.0\%$ SEM^l) compared to conventional ligation.¹⁸⁻⁵² Time differences (difference in mean times within these studies) range from longer procedures by 17.9 minutes (of 85.6 min axillary lymph node dissection procedures)⁵⁰ to time reductions of 106 min (of 349 min radical hysterectomy procedures).³⁶ Reduced intraoperative blood loss compared with conventional ligation: Across studies reporting volumetric blood loss (total 22 in the search window for the five procedures), the use of LigaSure™ devices has reduced these losses by an average of 28.9% ($\pm 3.9\%$ SEM^l) compared to conventional ligation.^{18,19,25,27-39,41,44-46,48,52} Difference in mean blood loss within these studies have ranged from a non-significant increase of 28 ml (of 260 ml in open radical gastrectomy)²⁸ to a reduction of 1,025 ml (of 4,223 ml in radical hysterectomy procedures).³⁹

^k Always refer to product labeling for indications for use of the associated LigaSure™ device (see section 1.1 for more details and chapter 6 for an overview of indications of the different LigaSure™ devices).

^l SEM, standard error of the mean. Source data used in calculation available on request  MD21154
Ligasure, background

- **Rates of operative complications are similar or reduced with LigaSure™ technology:** LigaSure™ devices have reduced or had equivalent risks of complications compared to conventional ligation in several procedures. These include statistically significant reductions in infections after hysterectomy,⁴³ and transient hypocalcemia⁵³ and hypoparathyroidism³⁵ after thyroidectomy. In breast surgery, LigaSure™ devices have achieved a significant reduction in both drainage duration^{23,26,34,40,42,44,50} and volume.^{23,34,44,50,54}
- **Potential to reduce costs due to complications:** Interventions and additional operating room usage required to treat complications can considerably increase costs.^{38,55-57} Improvement of patient safety endpoints could potentially decrease complication-related cost burdens.
- **Substantial cost savings reported with LigaSure™ devices:** In comparison to suture ligation, LigaSure™ devices have been reported to reduce hospital,³⁸ staffing,³⁸ and per-procedure costs.^{18,48,58,59}

Compared with other energy devices:

- **Faster or equivalent operative time compared to other energy devices:** In comparison to other energy devices, LigaSure™ technology has shown similar^{29,60-65} and significantly reduced operative time for monopolar scissors^{48,66} and ultrasonic technology.⁶⁷⁻⁶⁹
- **Reduced or equivalent intraoperative blood loss:** Compared to other energy devices, LigaSure™ technology has achieved similar^{29,60,61,64} and significantly lower blood loss.⁶⁷⁻⁶⁹
- **Equivalent complication rates:** There have been no increased risks of complications with LigaSure™ technology compared to other energy-based methods of hemostasis.^{29,47,48,53,62,63,66,69-74}
- **Lowers the cost of care:** Substantial savings with LigaSure™ technology have been achieved when compared to reusable monopolar scissors,^{66,75} even when specifically including costs of disposal of the single-use device,⁷⁵ and costs were comparable to ultrasonic technology device costs.^{58,73,76}

4 Source data tables

Refer to section 1.2 for scope of literature presented in the following data tables.

4.1 Data table summary: LigaSure™ technology preclinical data™

Table 4-1 Publication data on LigaSure™ technology in preclinical settings

Source	Model	Device	Generator	Measure	Detail	N	Value	Error
Chekan EG, et al., 2015 ⁴	porcine	LS Blunt tip	NR	BP	5-7mm	66	823.05	NR
				Sticking	5-7mm	66	0.133	NR
Timm RW, et al., 2014 ⁵	porcine	LS Blunt tip	NR	BP	<5mm	40	836	NR
					5-7mm	40	591	NR
Lesser TG, et al., 2013 ⁷⁷	human lung	LS Blunt tip	ForceTriad™ generator	BP	small arteries	15	315	213.1
					thick arteries	15	156	42.5
Takebayashi K, et al., 2013 ⁶	porcine lymph	NR	NR	BP	NR	NR	290	NR
				Leakage	NR	NR	0.133	NR
Milsom J, et al., 2012 ⁷	porcine	LS V	ForceTriad™ generator	BP	2-3mm	20	822.6	131
					4-5mm	20	740.9	124
					6-7mm	20	788.9	492
				TS (mm)	2-7mm	41	1.3	0.4
					4-5mm	20	2.1	1
					6-7mm	20	2.1	0.7

m Data from preclinical studies may not be indicative of clinical performance.

Source	Model	Device	Generator	Measure	Detail	N	Value	Error
Seehofer D, et al., 2012⁸	porcine	LS V	ForceTriad™ generator	BP	2-4mm	50	800	NR
					5-7mm	50	615	NR
				TS (degC), thermosensor	max T	8	86	NR
				TS (degC), thermocamera	max T	8	85	NR
Klar M, et al., 2011¹⁷	ovine	LS Atlas	LVSGen	BP	5.35mm	40	484	229
				SF	5.35mm	44	0.091	NA
				Time (s)	5.35mm	40	9.1	2.5
Mantke R, et al., 2011⁹	porcine	LS Dolphin	LVSGen	BP	<5mm	37	965	364
				Time (s)	<5mm	37	11.7	2
				TS (mm)	depth	37	2.5	0.59
Noble EJ, et al., 2011¹⁰	human intestine	LS Atlas	NR	BP	arteries	9	2005	NR
					veins	21	1295	NR
				TS (mm)	depth	30	3.37	1.44
Pick DL, et al., 2011¹⁶	porcine	NR	NR	BP	NR	24	518.2	244.7
Katsuno G, et al., 2010¹³	porcine	LS Atlas	ForceTriad™ generator	Time (s)	3.8mm	20	3.5	0.7
			LVSGen	Time (s)	3.7mm	20	7.6	1.9
Tsunezuka Y, et al., 2010¹¹	human lung	LS V	NR	BP	thick PA (>5mm)	12	399	147
					thin PA (<5mm)	13	604	82

Source	Model	Device	Generator	Measure	Detail	N	Value	Error
Newcomb WL, et al., 2009 ¹⁴	porcine	LS V	ForceTriad™ generator	BP	2-3mm	13	999.43	NR
					4-5mm	13	884.67	NR
					6-7mm	13	386.53	NR
				Device failure (%)	2-3mm	13	0	NR
					4-5mm	13	0	NR
					6-7mm	13	0	NR
				Time (s)	2-3mm	13	2.85	1.54
					4-5mm	13	3	2.36
					6-7mm	13	3.54	2.62
			LVSGen	BP	2-3mm	13	744.43	NR
					4-5mm	13	1261.45	NR
					6-7mm	13	645.16	NR
				Device failure (%)	2-3mm	13	13.3	NR
					4-5mm	13	0	NR
					6-7mm	13	0	NR
				Time (s)	2-3mm	13	4.07	1.56
					4-5mm	13	7.2	2.38
					6-7mm	13	6.77	2.62
Lamberton GR, et al., 2008 ³	bovine	LS V	LVSGen	BP	5mm	10	385	76
					Time (s)	20	10	0.9
				TS (degC)	max T, 5mm vessels	10	55.5	2.4

Source	Model	Device	Generator	Measure	Detail	N	Value	Error
Hruby GW, et al., 2007 ¹²	porcine	LS V	LVSGen	BP	< 3mm	73	569.4	NR
					3-5mm	38	557.8	NR
					5-7mm	3	533.7	NR
				Time (s)	< 3mm	73	4.1	NR
					3-5mm	38	5.7	NR
					5-7mm	3	8	NR
				TS (degC)	max T	10	103.14	NR
				TS (mm)	depth	2	4.5 (arterial)	NR

degC, degrees Celsius; NR, not reported; LVSGen, LigaSure™ vessel sealing system generator; PA, pulmonary artery; s, seconds. Measures: BP, burst pressure; SF, seal failure rate; Sticking, tissue sticking rate; Time, cycle completion time seconds; TS, thermal spread.

5 Structured literature search details

5.1 Searches performed

Structured searches were performed to identify literature reporting on clinical applications of LigaSure™ technology. The searches were divided into two streams: one to identify the most recent clinical evidence of the use of LigaSure™ devices (within the last 10 years), and a second parallel search to identify all clinical applications of LigaSure™ technology, regardless of time, to identify settings in which the devices have been used in patient care. The searches for hysterectomy and colorectal procedures were originally performed on February 27, 2017 (Table 5-1). Specifically, for thyroid surgery, the search was composed to capture any reference to LigaSure™ technology. On full text screening, only those publications were retained that applied to the LigaSure™ Small Jaw device. The PubMed search for peer-reviewed clinical evidence in thyroidectomy was originally performed on March 17, 2017 (Table 5-1). New searches for hysterectomy, colorectal and thyroidectomy procedures were undertaken in September 2021, using the same inclusion and exclusion criteria and search strings to identify any publications since the original search up to August 31, 2021. Further PubMed searches for clinical evidence on LigaSure™ technology in two additional procedures (breast surgery and gastric surgery) were performed on January 7, 2022 (Table 5-2).

Table 5-1 Structured searched in PubMed to identify relevant LigaSure™ technology data (original searches in 2017)

Index	Aim	Search string	Hits
#1	LigaSure™ devices by product name	ligasure*[tiab] OR LigaSure™[tiab] OR ligasurev*[tiab]	614
#2	Generic names for LigaSure™ devices	EBVS[tiab] OR BVSS[tiab] OR ((bipolar[tiab] OR ((high[tiab] OR advanced[tiab]) AND energy[tiab]) OR radiofrequency[tiab]) AND (diathermy[tiab] OR cautery[tiab] OR electrocautery[tiab] OR electrocauterization[tiab] OR electrocauterisation[tiab] OR cauterization[tiab] OR cauterisation[tiab] OR coagulation[tiab] OR ((vessel[tiab] OR tissue[tiab]) AND (sealing[tiab] OR sealer[tiab]))) OR electrosurgical[tiab] OR electrocoagulation[tiab] OR electrocoagulation[MeSH]	16,742
#3	Publication years of interest - original search	"2007/01/01"[PDat]:"2018/12/31"[PDat]	9,418,446

Index	Aim	Search string	Hits
#4	Invalid publication types	"Case Reports" [ptyp] OR "Clinical Conference" [ptyp] OR "Comment" [ptyp] OR "Editorial" [ptyp] OR "Letter"[ptyp] OR "Retracted Publication"[ptyp] OR "Congresses"[ptyp] OR "Duplicate Publication"[ptyp]	3,240,780
#5	Non-clinical data	"ex vivo"[tiab] OR cadaver[tiab] OR "deceased donor"[tiab]	78,271
#6	Animal data	Search terms for animal studies ⁷⁸	6,291,401
#7	LigaSure™ devices (by name or generic), restricted to publication years, excluding animal, non-clinical, and non-primary data	(#1 OR #2) AND #3 NOT (#4 OR #5 OR #6)	3,656
#8	All hysterectomy procedures	(hysterectomy[tiab] OR hysterectomy[MeSH] OR "hysterectomy,vaginal"[MeSH]) AND (surgery[tiab] OR "surgical procedure"[tiab] OR "General Surgery"[Mesh] OR "Surgical Procedures, Operative"[Mesh])	34,056
#9	All colorectal procedures	colectomy[tiab] OR sigmoidectomy[tiab] OR hemicolectomy[tiab] OR "colon resection"[tiab] OR "bowel resection"[tiab] OR colectomy[MeSH] OR ((colon[tiab] OR colorectal[tiab]) AND (polyp[tiab] OR polyps[tiab] OR polypectomy[tiab])) OR hemorrhoidectomy[tiab] OR haemorrhoidectomy[tiab] OR Hemorrhoids[MeSH] OR hemorrhoidectomy[MeSH] OR colostomy[MeSH] OR "Colorectal Surgery"[MeSH] OR "colorectal cancer"[tiab]	120,805
#10	All thyroidectomy procedures	Thyroidectomy[MeSH] OR thyroidectomy[tiab]	26,518
#11	LigaSure™ device + hysterectomy	#7 AND #8	136
#12	LigaSure™ device + colorectal	#7 AND #9	164
#13	LigaSure™ device + thyroidectomy	#7 AND #10	102

The original searches designed above revealed few relevant comparative studies reporting on LigaSure™ device usage in colectomy, although broader searches indicated wider use of the technology. For colectomy surgeries only then, the above searches were supplemented with hand searches, manual screening of article contents (to confirm usage of a LigaSure™ device when not otherwise mentioned in the abstract) and the inclusion of publications after 2007, which may also have referenced pre-2007 data.

Table 5-2 Structured searched in PubMed to identify relevant LigaSure™ technology data (searches for new indications in 2022)

Index	Aim	Search string	Hits
#1	LigaSure™ devices by product name	ligasure*[tiab] OR LigaSure™[tiab] OR ligasurev*[tiab]	802
#2	Generic names for LigaSure™ devices	EBVS[tiab] OR BVSS[tiab] OR ((bipolar[tiab] OR ((high[tiab] OR advanced[tiab]) AND energy[tiab]) OR radiofrequency[tiab]) AND (diathermy[tiab] OR cautery[tiab] OR electrocautery[tiab] OR electrocauterization[tiab] OR electrocauterisation[tiab] OR cauterization[tiab] OR cauterisation[tiab] OR coagulation[tiab] OR ((vessel[tiab] OR tissue[tiab]) AND (sealing[tiab] OR sealer[tiab]))) OR electrosurgical[tiab] OR electrocoagulation[tiab] OR electrocoagulation[MeSH]	19,818
#3	Publication years of interest - original search	"2010/12/31"[PDat]:"2022/01/04"[PDat]	12,752,634
#4	Invalid publication types	"Case Reports" [ptyp] OR "Clinical Conference" [ptyp] OR "Comment" [ptyp] OR "Editorial" [ptyp] OR "Letter"[ptyp] OR "Retracted Publication"[ptyp] OR "Congresses"[ptyp] OR "Duplicate Publication"[ptyp]	4,049,793
#5	Non-clinical data	"ex vivo"[tiab] OR cadaver[tiab] OR "deceased donor"[tiab]	112,219
#6	Animal data	Search terms for animal studies ⁷⁸	7,468,331
#7	LigaSure™ devices (by name or generic), restricted to publication years, excluding animal, non-clinical, and non-primary data	(#1 OR #2) AND #3 NOT (#4 OR #5 OR #6)	4,369

Index	Aim	Search string	Hits
#8	All breast surgery procedures	"Breast surgery"[tiab] OR (breast[tiab] AND (surgery[tiab] OR reduction[tiab] OR augmentation[tiab] OR implant[tiab] OR reconstruction[tiab] OR DIEP[tiab] OR flap[tiab] OR lymphadenectomy[tiab] OR "Lymph Node Excision"[MeSH] OR (lymph[tiab] AND node[tiab] AND (dissection[tiab] OR excision[tiab] OR biopsy[tiab]))) OR mastectomy[MeSH] OR mastectomy[tiab] OR lumpectomy[tiab]	103,885
#9	All gastric surgery procedures	Gastrectomy[tiab] OR gastrectomies[tiab] OR Gastrectomy[MeSH] OR Pyloromyotomy[MeSH] OR Pyloromyotomy[tiab] OR pyloroplasty[tiab] OR ((stomach[tiab] OR gastric [tiab] OR bariatric[tiab]) AND surgery[tiab]) OR gastroplasty[tiab] OR "Bariatric surgery"[MeSH] OR "Gastric Bypass"[MeSH] OR "gastric bypass"[tiab] OR Roux-en-Y[tiab] OR RYGB[tiab] OR "sleeve gastrectomy"[tiab] OR Billroth[tiab] OR BPD/DS[tiab] OR (BPD[tiab] AND DS[tiab]) OR ((biliopancreatic[tiab] OR bilio-pancreatic[tiab]) AND diversion[tiab] AND duodenal[tiab] AND switch[tiab])	109,521
#10	LigaSure™ device + breast surgery	#7 AND #8	63
#11	LigaSure™ device + gastric surgery	#7 AND #9	76

The second search, to reveal all clinical instances of LigaSure™ device use, including incidental mentions of use of the technology even if not the focus of the study, was performed using EMBASE. This database includes coverage of non-PubMed-indexed journals and congress reports to thus provide a comprehensive survey of LigaSure™ device use (Table 5-3). A differential search was performed May 2017 to retrieve new records from 2017 and updated or added from 2016, yielding 166 publications to screen for additional settings where LigaSure™ technology has been used. An additional differential search was performed in March 2022 to retrieve records that were added since the last search, yielding 810 publications to screen for additional settings where LigaSure™ technology has been used.

Table 5-3 Structured search in EMBASE to identify use of LigaSure™ devices

Index	Aim	Search string	Hits
1	LigaSure™ technology (referenced by name) excluding ex vivo data and select animals	((ligasure or ligasurev or LigaSure™) not (preclinical or "ex vivo" or cadaver or animal or dog or cat or pig or veterinary or veterinarian)).mp.	1,390

Results from the various searches were exported from their respective databases as search strings (PubMed) and .RIS files (EMBASE) for integration into the Sourcerer software utility for literature screening and review. Duplicate articles returned from the different sources were automatically removed.

5.2 Screening literature results

5.2.1 Screening general search (EMBASE) results

The purpose of the EMBASE search as described (Section 5.1) was to capture as many references (by name) of LigaSure™ device use in both peer-reviewed and non-peer-reviewed literature. As such, the results were not screened for exclusion criteria, but for surgery type and geographical location to verify clinical application of the technology.

5.2.2 Screening returned results for evidence (PubMed)

The articles returned from the PubMed searches were subjected to screening to identify articles relevant to the targeted indication (hysterectomy) for further deeper analysis. The criteria, and rationale are presented in Table 5-4.

Table 5-4 Description of exclusion criteria

Criterion	Explanation
Data pre-2007	Although the article was published after 2007 according to the search terms, the data referenced within cover a range prior to 2007
Articles with no abstract	At the level of top level screening, no informed decision regarding evidence or quality can be made without an abstract
Editorial/commentary	Articles that are commentaries or letters responding to other articles are not included for evidence recovery
Non-clinical	Articles which do not present any relevant clinical data, including patient surveys and experimental/ex vivo studies which were captured
Fewer than 20 patients	For higher quality evidence, studies of fewer than 20 patients are excluded
Non-targeted surgery	For a given surgical search area, if the focus of surgery of the article is for another, the study is excluded. Examples include the capture of mention of “bariatric procedures” for a study of appendectomy
Not vessel sealing focused	The mention of LigaSure™ technology or related technology is incidental and not the focus of the study with no data relevant to the performance of the technology
Not LigaSure™ technology	The reference by generic name to the technology cannot be conclusively identified as LigaSure™ technology at abstract level, or the reference to generic terms in the search such as electrocautery catches technology not relevant to bipolar vessel sealing

6 Overview: Indications of LigaSure™ devices

Table 6-1 gives an overview of the indications for use of all LigaSure™ devices as well as related BiZact™ bipolar devices and in which regions/countries these indications are applicable. This table was accurate at the time of writing. Indications, contraindications, warnings, precautions, and procedure steps may vary between products and models, and availability may vary by jurisdiction. Please always refer to indication labeling for your jurisdiction and read all applicable instructions for use provided with the products.

Table 6-1 LigaSure™ devices and their indications

Device Name	Reorder Code	Indications for Use	Region
BiZact™ Tonsillectomy Device	BZ4212	The BiZact™ device is a bipolar instrument intended for use in open surgical procedures where ligation and division of vessels, tissue bundles, and lymphatics is desired. The tissue fusion function of the device can be used on vessels (arteries and veins) and lymphatics up to and including 3 mm diameter. The BiZact™ device is indicated for use in general open surgical procedures. It is also indicated for ENT procedures, including adult, children and adolescent tonsillectomies for the ligation and division of vessels, tissue bundles and lymphatics 2-3 mm away from unintended thermally-sensitive structures. The BiZact™ device has not been shown to be effective for tubal sterilization or tubal coagulation for sterilization procedures. Do not use for these procedures.	OUS
BiZact™ Tonsillectomy Device	BZ4212A	The BiZact™ device is a bipolar instrument intended for use in open surgical procedures where ligation and division of vessels, tissue bundles, and lymphatics is desired. The tissue fusion function of the device can be used on vessels (arteries and veins) and lymphatics up to and including 3 mm diameter. The BiZact™ device is indicated for use in open general surgical procedures. It is also indicated for adult, children, and adolescent ENT procedures (3 years of age and above), including tonsillectomy, for the ligation and division of vessels, tissue bundles, and lymphatics 2-3 mm away from unintended thermally-sensitive structures. The BiZact™ device has not been shown to be effective for tubal sterilization or tubal coagulation for sterilization procedures. Do not use for these procedures.	US Only

Device Name	Reorder Code	Indications for Use	Region
LigaSure Impact™ Instrument	LF4418	The LigaSure™ Sealer/Divider is a bipolar electrosurgical instrument intended for use in open surgical procedures where ligation and division of vessels, tissue bundles, and lymphatics is desired. The LigaSure™ Sealer/Divider can be used on vessels (arteries, veins, pulmonary arteries, pulmonary veins) up to and including 7 mm. It is indicated for use in general surgery and such surgical specialties as urologic, vascular, thoracic, and gynecologic. Procedures may include, but are not limited to, Nissen fundoplication, colectomy, cholecystectomy, adhesiolysis, hysterectomy, oophorectomy, etc. The LigaSure™ system has not been shown to be effective for tubal sterilization or tubal coagulation for sterilization procedures. Do not use the LigaSure™ system for these procedures.	Global
LigaSure™ Curved, Small Jaw	LF1212	The LigaSure™ Curved, Small Jaw, Open Sealer/Divider is a bipolar electrosurgical instrument intended for use in open general surgical procedures where ligation and division of vessels (up to 7 mm in diameter), tissue bundles, and lymphatics is performed, such as urologic, thoracic, plastic, and reconstructive. Procedures may include, but are not limited to, bowel resections, gall bladder procedures, Nissen fundoplication, and adhesiolysis. The LigaSure™ system has not been shown to be effective for tubal sterilization or tubal coagulation for sterilization procedures. Do not use the LigaSure™ system for these procedures.	OUS
LigaSure™ Curved, Small Jaw	LF1212A	The LigaSure™ Curved, Small Jaw, Open Sealer/Divider is a bipolar electrosurgical instrument intended for use in open general surgical procedures where ligation and division of vessels (up to 7 mm in diameter), tissue bundles, and lymphatics is performed, such as urologic, thoracic, plastic, and reconstructive. Procedures may include, but are not limited to, bowel resections, gall bladder procedures, Nissen fundoplication, and adhesiolysis. The instrument is also indicated for open ENT procedures in adults (thyroidectomy, radical neck dissection, parotidectomy, and tonsillectomy) for ligation and division of vessels, lymphatics and tissue bundles 2-3 mm away from unintended thermally-sensitive structures such as nerves and parathyroid glands. The LigaSure™ system has not been shown to be effective for tubal sterilization or tubal coagulation for sterilization procedures. Do not use the LigaSure™ system for these procedures.	US Only

Device Name	Reorder Code	Indications for Use	Region
LigaSure™ Exact Dissector	LF2019	The LigaSure™ Sealer/Divider is a bipolar electrosurgical instrument intended for use in open surgical procedures where ligation and division of vessels, tissue bundles, and lymphatics is desired. The LigaSure™ Sealer/Divider can be used on vessels (arteries and veins) up to and including 7 mm. It is indicated for use in general surgery and in such surgical specialties as urologic, thoracic, plastic, and reconstructive. Procedures may include, but are not limited to, bowel resections, gall bladder procedures, Nissen fundoplication, and adhesiolysis. The instrument is also indicated for open ENT procedures in adults (thyroidectomy, radical neck dissection, parotidectomy, and tonsillectomy) for ligation and division of vessels, lymphatics and tissue bundles 2-3 mm away from unintended thermally-sensitive structures such as nerves and parathyroid glands. The LigaSure™ system has not been shown to be effective for tubal sterilization or tubal coagulation for sterilization procedures. Do not use the LigaSure™ system for these procedures.	OUS
LigaSure™ Exact Dissector	LF2019	The LigaSure™ Sealer/Divider is a bipolar electrosurgical instrument intended for use in open surgical procedures where ligation and division of vessels, tissue bundles, and lymphatics is desired. The LigaSure™ Sealer/ Divider can be used on vessels (arteries and veins) up to and including 7 mm. It is indicated for use in general surgery and in such surgical specialties as urologic, thoracic, plastic, and reconstructive. Procedures may include, but are not limited to, bowel resections, gall bladder procedures, Nissen fundoplication, breast procedures, and adhesiolysis. The instrument is also indicated for open ENT procedures in adults (thyroidectomy, radical neck dissection, parotidectomy, and tonsillectomy) for ligation and division of vessels, lymphatics and tissue bundles 2-3 mm away from unintended thermally-sensitive structures such as nerves and parathyroid glands. The LigaSure™ system has not been shown to be effective for tubal sterilization or tubal coagulation for sterilization procedures. Do not use the LigaSure™ system for these procedures.	US only
LigaSure™ Maryland Jaw	LF1923, LF1937, LF1944	The LigaSure™ Sealer/Divider is a bipolar electrosurgical instrument intended for use in minimally invasive or open surgical procedures where ligation and division of vessels, tissue bundles, and lymphatics is desired. The LigaSure™ Sealer/Divider can be used on vessels (arteries and veins) up to and including 7 mm. It is indicated for use in general surgery and such surgical specialties as urologic, vascular, thoracic, and gynecologic. These may include, but are not limited to, such procedures a Nissen fundoplication, colectomy, cholecystectomy, adhesiolysis, hysterectomy, oophorectomy, etc. The LigaSure™ system has not been shown to be effective for tubal sterilization or tubal coagulation for sterilization procedures. Do not use the LigaSure™ system for these procedures.	Global

Device Name	Reorder Code	Indications for Use	Region
LigaSure™ Maryland Jaw Thoracic	LF1930T	The LigaSure™ sealer/divider is a bipolar electrosurgical instrument intended for use in minimally invasive or open surgical procedures where ligation and division of vessels, tissue bundles, and lymphatics is desired. The LigaSure™ sealer/divider can be used on vessels (arteries, veins, pulmonary arteries, pulmonary veins) up to and including 7 mm. It is indicated for use in general surgery and such surgical specialties as urologic, vascular, thoracic, and gynecologic. Procedures may include, but are not limited to, Nissen fundoplication, colectomy, cholecystectomy, adhesiolysis, hysterectomy, oophorectomy, etc. The LigaSure™ system has not been shown to be effective for tubal sterilization or tubal coagulation for sterilization procedures. Do not use the LigaSure™ system for these procedures.	Global
LigaSure™ Blunt Tip	LF1823, LF1837, LF1844	The LigaSure™ Sealer/Divider is a bipolar electrosurgical instrument intended for use in minimally invasive or open surgical procedures where ligation and division of vessels, tissue bundles, and lymphatics is desired. The LigaSure™ Sealer/Divider can be used on vessels (arteries and veins) up to and including 7 mm. It is indicated for use in general surgery and such surgical specialties as urologic, vascular, thoracic, and gynecologic. Procedures may include, but are not limited to, Nissen fundoplication, colectomy, cholecystectomy, adhesiolysis, hysterectomy, oophorectomy, etc. The LigaSure™ system has not been shown to be effective for tubal sterilization or tubal coagulation for sterilization procedures. Do not use the LigaSure™ system for these procedures.	Global

OUS, Outside of the US

This table was accurate at the time of writing. Indications, contraindications, warnings, precautions, and procedure steps may vary between products and models, and availability may vary by jurisdiction. Please always refer to indication labeling for your jurisdiction and read all applicable instructions for use provided with the products.

7 References

1. Samudrala S. Topical Hemostatic Agents in Surgery: A Surgeon's Perspective. *AORN Journal*. 2008;88(3):S2-S11. doi:10.1016/S0001-2092(08)00586-3.
2. Sileshi B, Achneck HE, Lawson JH. Management of surgical hemostasis: topical agents. *Vascular*. 2008;16 Suppl 1:S22-8.
3. Lamberton GR, Hsi RS, Jin DH, Lindler TU, Jellison FC, Baldwin DD. Prospective comparison of four laparoscopic vessel ligation devices. *J Endourol*. 2008;22(10):2307-2312. doi:10.1089/end.2008.9715.
4. Chekan EG, Davison MA, Singleton DW, Mennone JZ, Hinoul P. Consistency and sealing of advanced bipolar tissue sealers. *Med Devices (Auckl)*. 2015;8:193-199. doi:10.2147/MDER.S79642.
5. Timm RW, Asher RM, Tellio KR, Welling AL, Clymer JW, Amaral JF. Sealing vessels up to 7 mm in diameter solely with ultrasonic technology. *Med Devices (Auckl)*. 2014;7:263-271. doi:10.2147/MDER.S66848.
6. Takebayashi K, Shiomi H, Naka S, et al. Utility of a microwave surgical instrument in sealing lymphatic vessels. *American journal of surgery*. 2013;206(2):229-233. doi:10.1016/j.amjsurg.2012.07.048.
7. Milsom J, Trencheva K, Monette S, et al. Evaluation of the safety, efficacy, and versatility of a new surgical energy device (THUNDERBEAT) in comparison with Harmonic ACE, LigaSure V, and EnSeal devices in a porcine model. *Journal of Laparoendoscopic and Advanced Surgical Techniques*. 2012;22(4):378-386. doi:10.1089/lap.2011.0420.
8. Seehofer D, Mogl M, Boas-Knoop S, et al. Safety and efficacy of new integrated bipolar and ultrasonic scissors compared to conventional laparoscopic 5-mm sealing and cutting instruments. *Surg Endosc*. 2012;26(9):2541-2549. doi:10.1007/s00464-012-2229-0.
9. Mantke R, Halangk W, Habermann A, et al. Efficacy and safety of 5-mm-diameter bipolar and ultrasonic shears for cutting carotid arteries of the hybrid pig. *Surg Endosc*. 2011;25(2):577-585. doi:10.1007/s00464-010-1224-6.
10. Noble EJ, Smart NJ, Challand C, Sleigh K, Oriolowo A, Hosie KB. Experimental comparison of mesenteric vessel sealing and thermal damage between one

bipolar and two ultrasonic shears devices. *Br J Surg*. 2011;98(6):797-800. doi:10.1002/bjs.7433.

11. Tsunetzuka Y, Waseda R, Yachi T. Electrothermal bipolar vessel sealing device LigaSureV for pulmonary artery ligation--burst pressure and clinical experiences in complete video-assisted thoracoscopic major lung resection for lung cancer. *Interact Cardiovasc Thorac Surg*. 2010;11(3):229-233. doi:10.1510/icvts.2010.239087.
12. Hrubby GW, Marruffo FC, Durak E, et al. Evaluation of surgical energy devices for vessel sealing and peripheral energy spread in a porcine model. *J Urol*. 2007;178(6):2689-2693. doi:10.1016/j.juro.2007.07.121.
13. Katsuno G, Nagakari K, Fukunaga M. Comparison of two different energy-based vascular sealing systems for the hemostasis of various types of arteries: a porcine model-evaluation of LigaSure ForceTriad™. *Journal of Laparoendoscopic and Advanced Surgical Techniques*. 2010;20(9):747-751. doi:10.1089/lap.2009.0438.
14. Newcomb WL, Hope WW, Schmelzer TM, et al. Comparison of blood vessel sealing among new electrosurgical and ultrasonic devices. *Surg Endosc*. 2009;23(1):90-96. doi:10.1007/s00464-008-9932-x.
15. Tavakoli Golpaygani A, Movahedi MM, Reza M. A Study on Performance and Safety Tests of Electrosurgical Equipment. *J Biomed Phys Eng*. 2016;6(3):175-182. Published September 1, 2016.
16. Pick DL, Kolla SB, Mucksavage P, et al. Sprayed fibrin sealant as the sole hemostatic agent for porcine laparoscopic partial nephrectomy. *J Urol*. 2011;185(1):291-297. doi:10.1016/j.juro.2010.09.002.
17. Klar M, Haberstroh J, Timme S, Fritzsche G, Gitsch G, Denschlag D. Comparison of a reusable with a disposable vessel-sealing device in a sheep model: efficacy and costs. *Fertil Steril*. 2011;95(2):795-798. doi:10.1016/j.fertnstert.2010.09.014.
18. Adamina M, Champagne BJ, Hoffman L, Ermlich MB, Delaney CP. Randomized clinical trial comparing the cost and effectiveness of bipolar vessel sealers versus clips and vascular staplers for laparoscopic colorectal resection. *Br J Surg*. 2011;98(12):1703-1712. doi:10.1002/bjs.7679.
19. Arowolo OA, Olasehinde O, Adisa AO, Adeyemo A, Alatise OI, Wuraola F. Early Experience with LigaSure Thyroidectomy in a Nigeria Teaching Hospital. *Niger J Surg*. 2019;25(1):64-69. doi:10.4103/njs.NJS_40_18.

20. Aydin C, Yildiz A, Kasap B, Yetimalar H, Kucuk I, Soylu F. Efficacy of electrosurgical bipolar vessel sealing for abdominal hysterectomy with uterine myomas more than 14 weeks in size: a randomized controlled trial. *Gynecol Obstet Invest.* 2012;73(4):326-329. doi:10.1159/000336400.
21. Bakacak Z, Bakacak M, Uzkar A, et al. The efficacy of LigaSure™ open instruments in cases of cesarean hysterectomy due to placenta percreta: a retrospective, record-based, comparative study. *J Matern Fetal Neonatal Med.* 2021;34(6):960-965. doi:10.1080/14767058.2020.1846177.
22. Bircan HY, Inal A, Ozcelik U, et al. LigaSure® versus clamp tie technique for thyroid surgery; decreased operative time versus increased inflammatory effect: a prospective randomized study. *European Review for Medical and Pharmacological Sciences.* 2014;18(14):1997-2005.
23. Chang YW, Kim HS, Jung SP, et al. Comparison of skin-sparing mastectomy using LigaSure™ Small Jaw and electrocautery. *World J Surg Oncol.* 2017;15(1):129. doi:10.1186/s12957-017-1199-z.
24. Coerper S, Stengl W. Einführung der Gefäßversiegelung in einem Schilddrüsenzentrüm: Risiken und wirtschaftlicher Nutzen. *Zentralblatt für Chirurgie.* 2018;143(5):543-549. doi:10.1055/s-0043-102570.
25. Coiro S, Frattaroli FM, Lucia F de, et al. A comparison of the outcome using Ligasure™ small jaw and clamp-and-tie technique in thyroidectomy: a randomized single center study. *Langenbecks Arch Surg.* 2015;400(2):247-252. doi:10.1007/s00423-014-1270-y.
26. Cortadellas T, Córdoba O, Espinosa-Bravo M, et al. Electrothermal bipolar vessel sealing system in axillary dissection: a prospective randomized clinical study. *Int J Surg.* 2011;9(8):636-640. doi:10.1016/j.ijso.2011.08.002.
27. Elhao M, Abdallah K, Serag I, El-Laithy M, Agur W. Efficacy of using electrosurgical bipolar vessel sealing during vaginal hysterectomy in patients with different degrees of operative difficulty: a randomised controlled trial. *Eur J Obstet Gynecol Reprod Biol.* 2009;147(1):86-90. doi:10.1016/j.ejogrb.2009.07.011.
28. Fujita J, Takiguchi S, Nishikawa K, et al. Randomized controlled trial of the LigaSure vessel sealing system versus conventional open gastrectomy for gastric cancer. *Surg Today.* 2014;44(9):1723-1729. doi:10.1007/s00595-014-0930-4.

29. Gambardella C, Clarizia G, Patrone R, et al. Advanced hemostasis in axillary lymph node dissection for locally advanced breast cancer: new technology devices compared in the prevention of seroma formation. *BMC Surg.* 2019;18(Suppl 1):125. doi:10.1186/s12893-018-0454-8.
30. Gizzo S, Burul G, Di Gangi S, et al. LigaSure vessel sealing system in vaginal hysterectomy: safety, efficacy and limitations. *Arch Gynecol Obstet.* 2013;288(5):1067-1074. doi:10.1007/s00404-013-2857-1.
31. Grøndal AYK, Høgsbro M, Pryds K, Pedersen HB, Jacobsen H. Intra- and postoperative complications using LigaSure™ Small Jaw in patients undergoing thyroidectomy: a register-based study. *Eur Arch Otorhinolaryngol.* 2021;278(11):4491-4500. doi:10.1007/s00405-021-06685-w.
32. Hirano J, Kinoshita S, Kazama T, et al. Value of Vessel Sealing System in Skin-Sparing Mastectomy in Patients with Early Breast Cancer. *Clinics in Surgery.* 2018;3(1).
33. Hirunwiwatkul P, Tungkavivachagul S. A multicenter, randomized, controlled clinical trial of LigaSure small jaw vessel sealing system versus conventional technique in thyroidectomy. *Eur Arch Otorhinolaryngol.* 2013;270(7):2109-2114. doi:10.1007/s00405-012-2289-8.
34. Inoue Y, Yamashita N, Ueo H, et al. The Clinical Usefulness of the LigaSure™ Small Jaw in Axillary Lymph Node Dissection in Patients with Breast Cancer. *Anticancer Res.* 2018;38(4):2359-2362. doi:10.21873/anticancer.12483.
35. Kuboki A, Nakayama T, Konno W, et al. New technique using an energy-based device versus conventional technique in open thyroidectomy. *Auris Nasus Larynx.* 2013;40(6):558-562. doi:10.1016/j.anl.2013.04.005.
36. Kyo S, Mizumoto Y, Takakura M, et al. Experience and efficacy of a bipolar vessel sealing system for radical abdominal hysterectomy. *Int J Gynecol Cancer.* 2009;19(9):1658-1661. doi:10.1111/IGC.0b013e3181a840d1.
37. Lakeman M, Kruitwagen RF, Vos MC, Roovers J-PWR. Electrosurgical bipolar vessel sealing versus conventional clamping and suturing for total abdominal hysterectomy: a randomized trial. *Journal of Minimally Invasive Gynecology.* 2008;15(5):547-553. doi:10.1016/j.jmig.2008.05.011.
38. Lakeman MME, The S, Schellart RP, et al. Electrosurgical bipolar vessel sealing versus conventional clamping and suturing for vaginal hysterectomy: a randomised

controlled trial. *BJOG*. 2012;119(12):1473-1482. doi:10.1111/j.1471-0528.2012.03484.x.

39. Lauroy A, Verhaeghe C, Vidal F, Parant O, Legendre G, Guerby P. Perioperative outcomes using LigaSure compared with conventional technique in peripartum hysterectomy. *Arch Gynecol Obstet*. 2020;301(1):229-234. doi:10.1007/s00404-019-05398-0.
40. Lee JS, Kim DG, Lee JW, et al. Usefulness of the LigaSure™ small jaw sealing device for breast reconstruction with a latissimus dorsi flap. *J Plast Surg Hand Surg*. 2019;53(5):295-300. doi:10.1080/2000656X.2019.1612753.
41. Nakajima K, Nezu R, Ito T, Nishida T. Hand-assisted laparoscopic restorative proctocolectomy for ulcerative colitis: the optimization of instrumentation toward standardization. *Surg Today*. 2010;40(9):840-844. doi:10.1007/s00595-009-4157-8.
42. Nespoli L, Antolini L, Stucchi C, Nespoli A, Valsecchi MG, Gianotti L. Axillary lymphadenectomy for breast cancer. A randomized controlled trial comparing a bipolar vessel sealing system to the conventional technique. *Breast*. 2012;21(6):739-745. doi:10.1016/j.breast.2012.08.003.
43. Nouri K, Ott J, Demmel M, Promberger R, Huber JC, Mayerhofer K. Bipolar vessel sealing increases operative safety in laparoscopic-assisted vaginal hysterectomy. *Arch Gynecol Obstet*. 2011;283(1):91-95. doi:10.1007/s00404-009-1333-4.
44. Seki T, Hayashida T, Takahashi M, Jinno H, Kitagawa Y. A randomized controlled study comparing a vessel sealing system with the conventional technique in axillary lymph node dissection for primary breast cancer. *Springerplus*. 2016;5(1):1004. doi:10.1186/s40064-016-2710-7.
45. Silva-Filho AL, Rodrigues AM, Vale de Castro Monteiro M, et al. Randomized study of bipolar vessel sealing system versus conventional suture ligature for vaginal hysterectomy. *Eur J Obstet Gynecol Reprod Biol*. 2009;146(2):200-203. doi:10.1016/j.ejogrb.2009.03.014.
46. Takiguchi N, Nagata M, Soda H, et al. Multicenter randomized comparison of LigaSure versus conventional surgery for gastrointestinal carcinoma. *Surg Today*. 2010;40(11):1050-1054. doi:10.1007/s00595-009-4234-z.
47. Testini M, Pasculli A, Di Meo G, et al. Advanced vessel sealing devices in total thyroidectomy for substernal goitre: A retrospective cohort study. *Int J Surg*. 2016;35:160-164. doi:10.1016/j.ijsu.2016.09.092.

48. Tou S, Malik AI, Wexner SD, Nelson RL. Energy source instruments for laparoscopic colectomy. *Cochrane Database Syst Rev*. 2011;(5):CD007886. doi:10.1002/14651858.CD007886.pub2.
49. Trilling B, Riboud R, Abba J, Girard E, Faucheron JL. Energy vessel sealing systems versus mechanical ligature of the inferior mesenteric artery in laparoscopic sigmoidectomy. *Int J Colorectal Dis*. 2016;31(4):903-908. doi:10.1007/s00384-016-2508-7.
50. Tukenmez M, Agcaoglu O, Aksakal N, et al. The use of Ligasure vessel sealing system in axillary dissection; effect on seroma formation. *Chirurgia (Bucur)*. 2014;109(5):620-625.
51. Uludag SS, Teksoz S, Arian AE, et al. Effect of energy-based devices on voice quality after total thyroidectomy. *Eur Arch Otorhinolaryngol*. 2017;274(5):2295-2302. doi:10.1007/s00405-016-4444-0.
52. Zhou H, Yi W, Zhang J, et al. Short- and long-term outcomes of LigaSure versus conventional surgery for curative gastric cancer resection: a matched pair analysis. *Gastric Cancer*. 2015;18(4):843-849. doi:10.1007/s10120-014-0418-3.
53. Zhang L, Li N, Yang X, Chen J. A meta-analysis comparing the outcomes of LigaSure Small Jaw versus clamp-and-tie technique or Harmonic Focus Scalpel in thyroidectomy. *Medicine*. 2017;96(11):e6141. doi:10.1097/MD.00000000000006141.
54. Kim MK, Lee S. LigaSure for the Creation of Bloodless Breast Pockets in Patients Undergoing Transaxillary Breast Augmentation. *Plast Reconstr Surg Glob Open*. 2020;8(12):e3295. doi:10.1097/GOX.00000000000003295.
55. Al-Qurayshi Z, Robins R, Hauch A, Randolph GW, Kandil E. Association of Surgeon Volume With Outcomes and Cost Savings Following Thyroidectomy: A National Forecast. *JAMA Otolaryngol Head Neck Surg*. 2016;142(1):32-39. doi:10.1001/jamaoto.2015.2503.
56. Selby LV, Gennarelli RL, Schnorr GC, et al. Association of Hospital Costs With Complications Following Total Gastrectomy for Gastric Adenocarcinoma. *JAMA Surg*. 2017;152(10):953-958. doi:10.1001/jamasurg.2017.1718.
57. Sivakumar J, Alnimri F, Johnson MA, Ward S, Chong L, Hii MW. Health economic analysis of curative-intent gastrectomy for gastric carcinoma and the costs related

to post-operative complications. *ANZ J Surg.* 2021;91(1-2):E1-E6.
doi:10.1111/ans.16014.

58. Pons Y, Gauthier J, Ukkola-Pons E, et al. Comparison of LigaSure vessel sealing system, harmonic scalpel, and conventional hemostasis in total thyroidectomy. *Otolaryngol Head Neck Surg.* 2009;141(4):496-501.
doi:10.1016/j.otohns.2009.06.745.
59. Chang L-Y, O'Neill C, Suliburk J, Sidhu S, Delbridge L, Sywak M. Sutureless total thyroidectomy: a safe and cost-effective alternative. *ANZ J Surg.* 2011;81(7-8):510-514. doi:10.1111/j.1445-2197.2010.05492.x.
60. Janssen PF, Brölmann HAM, van Kesteren PJM, et al. Perioperative outcomes using LigaSure compared with conventional bipolar instruments in laparoscopic hysterectomy: a randomised controlled trial. *BJOG.* 2011;118(13):1568-1575.
doi:10.1111/j.1471-0528.2011.03089.x.
61. Aytan H, Nazik H, Narin R, Api M, Tok EC. Comparison of the use of LigaSure, HALO PKS cutting forceps, and ENSEAL tissue sealer in total laparoscopic hysterectomy: a randomized trial. *Journal of Minimally Invasive Gynecology.* 2014;21(4):650-655. doi:10.1016/j.jmig.2014.01.010.
62. Velotti N, Manigrasso M, Di Lauro K, et al. Comparison between LigaSure™ and Harmonic® in Laparoscopic Sleeve Gastrectomy: A Single-Center Experience on 422 Patients. *J Obes.* 2019;2019:3402137. doi:10.1155/2019/3402137.
63. Tsamis D, Natoudi M, Arapaki A, et al. Using Ligasure™ or Harmonic Ace® in Laparoscopic Sleeve Gastrectomies? A Prospective Randomized Study. *Obes Surg.* 2015;25(8):1454-1457. doi:10.1007/s11695-014-1551-0.
64. Kirmizi S, Kayaalp C, Karagul S, et al. Comparison of Harmonic scalpel and Ligasure devices in laparoscopic Roux-en-Y gastric bypass. *Wideochir Inne Tech Maloinwazyjne.* 2017;12(1):28-31. doi:10.5114/wiitm.2017.66641.
65. Lee CM, Da Park W, Park S, Kim J-H, Park S-H, Kim C-S. Lymph Node Dissection Using Bipolar Vessel-Sealing Device During Reduced Port Laparoscopic Distal Gastrectomy for Gastric Cancer: Result of a Pilot Study from a Single Institute. *J Laparoendosc Adv Surg Tech A.* 2017;27(11):1101-1108.
doi:10.1089/lap.2016.0685.
66. Hubner M, Demartines N, Muller S, Dindo D, Clavien P-A, Hahnloser D. Prospective randomized study of monopolar scissors, bipolar vessel sealer and ultrasonic

shears in laparoscopic colorectal surgery. *Br J Surg*. 2008;95(9):1098-1104. doi:10.1002/bjs.6321.

67. Demirturk F, Aytan H, Caliskan AC. Comparison of the use of electrothermal bipolar vessel sealer with harmonic scalpel in total laparoscopic hysterectomy. *J Obstet Gynaecol Res*. 2007;33(3):341-345. doi:10.1111/j.1447-0756.2007.00533.x.
68. Campagnacci R, Sanctis A de, Baldarelli M, Rimini M, Lezoche G, Guerrieri M. Electrothermal bipolar vessel sealing device vs. ultrasonic coagulating shears in laparoscopic colectomies: a comparative study. *Surg Endosc*. 2007;21(9):1526-1531. doi:10.1007/s00464-006-9143-2.
69. Guerrieri M, Sanctis A, Baldarelli M, et al. Electrothermal bipolar energy-based device in laparoscopic right colectomy: our experience. *Minerva Chir*. 2008;63(6):455-460.
70. Rimonda R, Arezzo A, Garrone C, Allaix ME, Giraudo G, Morino M. Electrothermal bipolar vessel sealing system vs. harmonic scalpel in colorectal laparoscopic surgery: a prospective, randomized study. *Dis Colon Rectum*. 2009;52(4):657-661. doi:10.1007/DCR.0b013e3181a0a70a.
71. Dionigi G, Boni L, Rausei S, et al. The safety of energy-based devices in open thyroidectomy: a prospective, randomised study comparing the LigaSure™ (LF1212) and the Harmonic® FOCUS. *Langenbecks Arch Surg*. 2012;397(5):817-823. doi:10.1007/s00423-011-0898-0.
72. Hammad AY, Deniwar A, Al-Qurayshi Z, Mohamed HE, Rizwan A, Kandil E. A Prospective Study Comparing the Efficacy and Surgical Outcomes of Harmonic Focus Scalpel Versus LigaSure Small Jaw in Thyroid and Parathyroid Surgery. *Surg Innov*. 2016;23(5):486-489. doi:10.1177/1553350616639143.
73. Hwang SO, Jung JH, Park HY, Kim WW. A Prospective, Randomized Study between the Small Jaw® and the Harmonic Focus® in Open Thyroidectomy. *Otolaryngol Head Neck Surg*. 2014;150(6):943-948. doi:10.1177/0194599814527730.
74. Simillis C, Hompes R, Penna M, Rasheed S, Tekkis PP. A systematic review of transanal total mesorectal excision: is this the future of rectal cancer surgery? *Colorectal Disease*. 2016;18(1):19-36. doi:10.1111/codi.13151.
75. Holloran-Schwartz MB, Gavard JA, Martin JC, Blaskiewicz RJ, Yeung PP. Single-Use Energy Sources and Operating Room Time for Laparoscopic Hysterectomy: A

Randomized Controlled Trial. *Journal of Minimally Invasive Gynecology*. 2016;23(1):72-77. doi:10.1016/j.jmig.2015.08.881.

76. Rahbari R, Mathur A, Kitano M, et al. Prospective randomized trial of ligasure versus harmonic hemostasis technique in thyroidectomy. *Ann Surg Oncol*. 2011;18(4):1023-1027. doi:10.1245/s10434-010-1251-5.
77. Lesser TG, Wolfram F, Boltze C. Sealing of pulmonary arteries with LigaSure: in vivo and ex vivo examinations. *J Thorac Cardiovasc Surg*. 2013;145(6):1525-1528. doi:10.1016/j.jtcvs.2012.11.009.
78. Hooijmans CR, Tillema A, Leenaars M, Ritskes-Hoitinga M. Enhancing search efficiency by means of a search filter for finding all studies on animal experimentation in PubMed. *Lab Anim*. 2010;44(3):170-175. doi:10.1258/la.2010.009117.

8 Geographical references

- G1. Bignon H, Buela E, Martinez-Ferro M. Which Is the best vessel-sealing method for pediatric thoracoscopic lobectomy? *J Laparoendosc Adv Surg Tech A*. 2010;20(4):395-398.
- G2. Ferro MM, Vidaurreta S, Buela E, Bignon H. Thoracoscopic lung lobectomy: Preliminary results in pediatric patients. *Archivos Argentinos de Pediatría*. 2010;108(2):153-157.
- G3. Chang LY, O'Neill C, Suliburk J, Sidhu S, Delbridge L, Sywak M. Sutureless total thyroidectomy: a safe and cost-effective alternative. *ANZ J Surg*. 2011;81(G7-8):510-514.
- G4. Chiow AK, Lewin J, Manoharan B, Cavallucci D, Bryant R, O'Rourke N. Intercostal and transthoracic trocars enable easier laparoscopic resection of dome liver lesions. *HPB (Oxford)*. 2015;17(4):299-303. doi:10.1111/hpb.12336.
- G5. Ding Z, Wable M, Rane A. Use of Ligasure bipolar diathermy system in vaginal hysterectomy. *Journal of Obstetrics and Gynaecology*. 2005;25(1):49-51. doi:10.1080/01443610400024609.
- G6. Dong L, Suh H, Karantanis W, Jia S, Yang Y, Loi K. Evaluation of Micronutrient Status Post Laparoscopic Sleeve Gastrectomy: an Australian Perspective. *Obesity Surgery*. 2021;31(3):1099-1104. doi:10.1007/s11695-020-05089-9.
- G7. Dudi-Venkata NN, Houli N, Weinberg L, Nikfarjam M. Laparoscopic partial splenectomy performed by monopolar saline-cooled radiofrequency coagulation. *J Laparoendosc Adv Surg Tech A*. 2014;24(7):502-505. doi:10.1089/lap.2014.0143.
- G8. O'Neill CJ, Chang LY, Suliburk JW, Sidhu SB, Delbridge LW, Sywak MS. Sutureless thyroidectomy: surgical technique. *ANZ J Surg*. 2011;81(G7-8):515-518.
- G9. Perera M, Anabell L, Page D, Harding T, Gnaneswaran N, Chan S. Risk factors for post-thyroidectomy haematoma. *J Laryngol Otol*. 2016;130 Suppl 1(GS1):S20-5. doi:10.1017/S0022215115003199.
- G10. Rahman AA, Riaz AA, Kuzinkovas V. Laparoscopic sleeve gastrectomy and the "sydney Patch": Triple reinforcement of staple line: A study with 546 patients. *Bariatric Surgical Practice and Patient Care*. 2018;13(4):145-150. doi:10.1089/bari.2018.0038.
- G11. Rhou, Y. J. J., Pather S, Loadsman JA, Campbell N, Philp S, Carter J. Direct hospital costs of total laparoscopic hysterectomy compared with fast-track open hysterectomy at a tertiary hospital: A retrospective case-controlled study. *Australian and New Zealand Journal of Obstetrics and Gynaecology*. 2015;55(6):584-587.

- G12. Rieger NA, Lam FF. Single-incision laparoscopically assisted colectomy using standard laparoscopic instrumentation. *Surgical Endoscopy and Other Interventional Techniques*. 2010;24(4):888-890.
- G13. Soon PS, Yeh MW, Sywak MS, Sidhu SB. Use of the ligaSure vessel sealing system in laparoscopic adrenalectomy. *ANZ J Surg*. 2006;76(9):850-852. doi:10.1111/j.1445-2197.2006.03881.x.
- G14. Suhardja TS, Norhadi S, Ee E, Hodgkins B. Comparison of the Thunderbeat and Other Energy Devices in Laparoscopic Colorectal Resection: A Single-Center Experience. *Journal of Laparoendoscopic and Advanced Surgical Techniques*. 2018;28(12):1417-1421. doi:10.1089/lap.2018.0208.
- G15. Surve A, Rao R, Cottam D, et al. Early Outcomes of Primary SADI-S: an Australian Experience. *Obesity Surgery*. 2020;30(4):1429-1436. doi:10.1007/s11695-019-04312-6.
- G16. Surve A, Rao R, Cottam D. Laparoscopic Single Anastomosis Duodeno-Ileal Bypass with Sleeve Gastrectomy: Surgical Technique. *Obesity Surgery*. 2020;30(11):4684-4686. doi:10.1007/s11695-020-04847-z.
- G17. Wong C, Goh A, Merkur H. Comparison of surgical outcomes using Gyrus PKS™ vs LigaSure™ in total laparoscopic hysterectomy: A randomised controlled trial. *Aust N Z J Obstet Gynaecol*. 2020;60(5):790-796. doi:10.1111/ajo.13217.
- G18. Page PL, Martin D, Taylor C, et al. Does hiatal repair affect gastroesophageal reflux symptoms in patients undergoing laparoscopic sleeve gastrectomy? *Surgical Endoscopy and Other Interventional Techniques*. 2018;32(5):2373-2380. doi:10.1007/s00464-017-5935-9.
- G19. Anegg U, Lindenmann J, Maier A, Renner H, Matzi V, Smolle-Juttner FM. Does ligasure in mediastinal lymph node dissection impact drainage duration? Preliminary results. *Interactive Cardiovascular and Thoracic Surgery*. 2012;15:S33-S33.
- G20. Brunner W, Schirnhöfer J, Frass R, et al. Single-access laparoscopic sigmoid resections. *Surgical Endoscopy and Other Interventional Techniques*. 2011;25:S28-S28.
- G21. Gnant M. Challenges and Controversies in Breast Surgery. *Breast Care (Basel)*. 2019;14(4):185-186. doi:10.1159/000502134.
- G22. Hochstätter R, Bertholin Galvez MY, Hartleb R, et al. Voiding dysfunction after colectomy with neo-urethra preformation in transgender men. *International Urogynecology Journal*. 2020;31(GSUPPL 1):S219-S220. doi:10.1007/s00192-020-04555-3.

- G23. Panhofer P, Rothe S, Schutz M, et al. Morbidity reduction using the vessel sealing device LigaSure™ in breast cancer surgery. *European Surgery - Acta Chirurgica Austriaca*. 2015;47(4):150-156.
- G24. Perathoner A, Weissenbacher A, Sucher R, Laimer E, Pratschke J, Mittermair R. Significant weight loss and rapid resolution of diabetes and dyslipidemia during short-term follow-up after laparoscopic sleeve gastrectomy. *Obes Surg*. 2013;23(12):1966-1972. doi:10.1007/s11695-013-1038-4.
- G25. Pimpl K, Brunner WB, Schirnhöfer JS, et al. First case of transumbilical single port laparoscopic pancreas resection. *Surgical Endoscopy and Other Interventional Techniques*. 2011;25:S134-S134.
- G26. Sucher R, Resch T, Cakar-Beck F, Pratschke J, Mittermair R. Single port access sleeve gastrectomy is reasonable! *European Surgery - Acta Chirurgica Austriaca*. 2012;44:5.
- G27. Tamussino K, Afschar P, Reuss J, Perschler M, Ralph G, Winter R. Electrosurgical bipolar vessel sealing for radical abdominal hysterectomy. *Gynecol Oncol*. 2005;96(2):320-322. doi:10.1016/j.ygyno.2004.09.021.
- G28. Downes RO, McFarlane M, Diggiss C, Iferenta J. Single incision cholecystectomy using a clipless technique with LigaSure in a resource limited environment: The Bahamas experience. *Int J Surg Case Rep*. 2015;11:104-109. doi:10.1016/j.ijscr.2015.04.034.
- G29. Dapri G, Donckier V, Himpens J, Cadiere GB. Single-incision transumbilical laparoscopic left hepatic lobectomy and Nissen fundoplication. *Surgical Endoscopy and Other Interventional Techniques*. 2012;26:S448-S448.
- G30. Dapri G, Donckier V. Transumbilical single-incision laparoscopic liver surgery. *Surgical Endoscopy and Other Interventional Techniques*. 2014;28:S193-S193.
- G31. Hassayoune N, Mhallem Gziri M, Lentini A, et al. Severe Gestational Gigantomastia: from Mastectomy to Staged Autologous Breast Reconstruction. A Case Report. *JPRAS Open*. 2021;29(18):65-70. doi:10.1016/j.jpra.2021.04.011.
- G32. Khelif K, Maassarani F, Dassonville M, Laet MH. Thoracoscopic thoracic duct sealing with LigaSure in two children with refractory postoperative chylothorax. *J Laparoendosc Adv Surg Tech A*. 2007;17(1):137-139. doi:10.1089/lap.2006.0537.
- G33. Lafosse A, Vandeputte C, Sabor I, Mahaudens P, Denoel C. About an extreme case of giant lymphoedema of the upper limb: Ligasure may be useful and a multidisciplinary approach is mandatory. *Annales de Chirurgie Plastique et Esthétique*. 2011;56(4):325-328.

- G34. Chelala E, El Hajj Moussa W, Rizk S, Assaker N. Consecutive Versus Selective Primary and Revisional Single Incision Laparoscopic Bariatric Surgery: Personal Experience in 330 Cases. *Obesity Surgery*. 2020;30(4):1515-1526. doi:10.1007/s11695-019-04356-8.
- G35. Ettinger JE, Ramos AC, Azaro E, et al. Staplerless laparoscopic gastric bypass: a new option in bariatric surgery. *Obes Surg*. 2006;16(5):638-645. doi:10.1381/096089206776944896.
- G36. Sakamoto E, Dias AR, Ramos M, et al. Laparoscopic Completion Total Gastrectomy for Remnant Gastric Cancer. *Journal of Laparoendoscopic and Advanced Surgical Techniques*. 2021;31(7):803-807. doi:10.1089/lap.2020.0569.
- G37. Silva-Filho AL, Rodrigues AM, Vale de Castro Monteiro, M., et al. Randomized study of bipolar vessel sealing system versus conventional suture ligation for vaginal hysterectomy. *Eur J Obstet Gynecol Reprod Biol*. 2009;146(2):200-203. doi:10.1016/j.ejogrb.2009.03.014.
- G38. Dimov P, Kjoshev K. Vessel sealing device in thyroid surgery-sharing experience. *European Surgical Research*. 2015;55:97-98.
- G39. Cheang T, Hanna SS, Wright FC, Law, C. H. L. Use of a collagen-sealing device in hepatic resection: A comparative analysis to standard resection technique. *HPB (Oxford)*. 2006;8(3):194-199.
- G40. Dubuc-Lissoir J. Use of a new energy-based vessel ligation device during laparoscopic gynecologic oncologic surgery. *Surg Endosc*. 2003;17(3):466-468. doi:10.1007/s00464-002-8536-0.
- G41. Garneau PY, Abouzahr O, Garofalo F, et al. Decreasing complication rates for one-stage conversion band to laparoscopic sleeve gastrectomy: A retrospective cohort study. *Journal of Minimal Access Surgery*. 2020;16(3):264-268. doi:10.4103/jmas.JMAS_86_18.
- G42. Helpman L, Covens A. Laparoscopic appendectomy using the LigaSure vessel sealing device at the time of gynecologic surgery. *Journal of Gynecologic Surgery*. 2012;28(3):188-191.
- G43. Iwanicki S, Cenaiko D, Robert M. Retrospective review of vaginal salpingectomies: Success rates and complications. *International Urogynecology Journal and Pelvic Floor Dysfunction*. 2014;25(G1 SUPPL. 1):S165-S166.

- G44. Kim DY, Boushey RP. The Use of Alternative Energy Sources and New Techniques for the Treatment of Hemorrhoidal Disease. *Seminars in Colon and Rectal Surgery*. 2007;18(3):181-186.
- G45. Lennox G, Pantazi S, Keunen J, van Mieghem T, Allen L. Minimally invasive surgical management of a second trimester pregnancy in a rudimentary uterine horn. *J Obstet Gynaecol Can*. 2013;35(5):468-472.
- G46. Palter VN, MacLellan S, Ashamalla S. Laparoscopic translevator approach to abdominoperineal resection for rectal adenocarcinoma: feasibility and short-term oncologic outcomes. *Surg Endosc*. 2016;30(7):3001-3006. doi:10.1007/s00464-015-4589-8.
- G47. Singh P, O'Connell D, Langille M, Dziegielewski P, Allegretto M, Harris J. LigaSure versus conventional hemostasis in thyroid surgery: prospective randomized controlled trial. *J Otolaryngol Head Neck Surg*. 2010;39(4):378-384.
- G48. Burgos AM, Braghetto I, Csendes A, et al. Gastric leak after laparoscopic-sleeve gastrectomy for obesity. *Obes Surg*. 2009;19(12):1672-1677. doi:10.1007/s11695-009-9884-9.
- G49. Farias C, Fernandez JI, Ovalle C, Cabrera C, La Maza J de. Simplified transumbilical sleeve gastrectomy technique using rigid instruments. *Obesity Surgery*. 2011;21(8):1009-1010.
- G50. Chen HL, Woo XB, Cui J, Chen CQ, Peng JS. Ligasure versus stapled hemorrhoidectomy in the treatment of hemorrhoids: a meta-analysis of randomized control trials. *Surg Laparosc Endosc Percutan Tech*. 2014;24(4):285-289. doi:10.1097/SLE.0000000000000009.
- G51. Gao Y, Duan Y, Zheng B, Pan W, Wang H, Li P. Laparoscopic liver resection in children: A report of 4 cases. *J Laparoendosc Adv Surg Tech A*. 2013;23(12):A78-A79.
- G52. Hong D, Cheng J, Wang Z, et al. Comparison of two laparoscopic splenectomy plus pericardial devascularization techniques for management of portal hypertension and hypersplenism. *Surg Endosc*. 2015;29(12):3819-3826. doi:10.1007/s00464-015-4147-4.
- G53. Hong J, Zhu X, Yao L, et al. Ligating of periphery vessels could prevent staple line bleeding effectively in laparoscopic sleeve gastrectomy. *Obesity Surgery*. 2019;29(3):149. doi:10.1007/s11695-019-03899-0.
- G54. Huanwei C, Feiwen D. Pure laparoscopic right hemihepatectomy via anterior approach. *Surg Endosc*. 2016;30(12):5621. doi:10.1007/s00464-016-4903-0.

- G55. Li H, Wei Y, Li B, Peng B. Modified Thoracoscopic Hepatectomy For Segment VIII: A Case Report. *Medicine (Baltimore)*. 2016;95(22):e3801. doi:10.1097/MD.0000000000003801.
- G56. Li Z, Chen L, Wang J, et al. A single institution experience using the LigaSure vessel sealing system in video-assisted thoracoscopic surgery for primary spontaneous pneumothorax. *J Biomed Res*. 2014;28(6):494-497. doi:10.7555/JBR.28.20130098.
- G57. Mao X-C, Chen C, Wang K-J. Efficacy and safety of LigaSure™ small jaw instrument in thyroidectomy: a 1-year prospective observational study. *Eur Arch Otorhinolaryngol*. 2018;275(5):1257-1263. doi:10.1007/s00405-018-4912-9.
- G58. Wang M, Wang X, Zhang H, Cai Y, Peng B. Laparoscopic splenectomy for portal hypertension secondary to liver cirrhosis: Ligasure combined with ultrasound scalpel versus ultrasound scalpel. *Surgical Endoscopy and Other Interventional Techniques*. 2015;29:S520-S520.
- G59. Wu S, Fan Y, Kong J, Su Y, Tian Y. Transumbilical single-incision laparoscopic splenectomy using conventional instruments: Initial single center experience. *Surgical Endoscopy and Other Interventional Techniques*. 2013;27:S499-S499.
- G60. Wu J, Tang Z, Zhao G, et al. Incidence and risk factors for postoperative pancreatic fistula in 2089 patients treated by radical gastrectomy: A prospective multicenter cohort study in China. *Int J Surg*. 2022;98:106219. doi:10.1016/j.ijsu.2021.106219.
- G61. Wu S, Chen Y, Tian Y, Jing K. Transumbilical single-incision laparoscopic multiple organ procedures: initial experience of 20 cases. *J Laparoendosc Adv Surg Tech A*. 2013;23(1):56-59. doi:10.1089/lap.2012.0317.
- G62. Yuan CH, Xiu DR, Jia YM, Xiong JW, Zhang TL. [Laparoscopic liver tumor resection of clinical experience in 126 patients]. *Zhonghua Wai Ke Za Zhi*. 2013;51(9):776-779.
- G63. Zhan XL, Ji Y, Wang YD. Laparoscopic splenectomy for hypersplenism secondary to liver cirrhosis and portal hypertension. *World J Gastroenterol*. 2014;20(19):5794-5800. doi:10.3748/wjg.v20.i19.5794.
- G64. Zhou H, Yi W, Zhang J, et al. Short- and long-term outcomes of LigaSure versus conventional surgery for curative gastric cancer resection: a matched pair analysis. *Gastric Cancer*. 2015;18(4):843-849. doi:10.1007/s10120-014-0418-3.
- G65. Esmeral J, Carrasquilla A, Zafrani J. Laparoscopic sleeve gastrectomy, v-loc staple line reinforcement, 480 cases without leak or bleeding. *Obesity Surgery*. 2014;24(8):1338.

- G66. Druzijanic N, Pogorelic Z, Perko Z, Mrklic I, Tomic S. Comparison of lateral thermal damage of the human peritoneum using monopolar diathermy, Harmonic scalpel and LigaSure. *Can J Surg*. 2012;55(5):317-321. doi:10.1503/cjs.000711.
- G67. Arslani N, Patrlj L, Rajkovic Z, Papes D, Altarac S. A randomized clinical trial comparing Ligasure versus stapled hemorrhoidectomy. *Surg Laparosc Endosc Percutan Tech*. 2012;22(1):58-61. doi:10.1097/SLE.0b013e318247d966.
- G68. Brodak M, Kosina J, Holub L, Husek P, Pacovsky J. Comparison between harmonic scalpel and ligasure vessel sealing system in laparoscopic radical prostatectomy. *European Urology, Supplements*. 2010;9(6):635.
- G69. Cecka F, Ferko A, Jon B, Subrt Z, Kasparova P, Repak R. Pancreatic Castleman disease treated with laparoscopic distal pancreatectomy. *Hepatobiliary Pancreat Dis Int*. 2013;12(3):332-334.
- G70. Hora M, Eret V, Urge T, et al. Antegrade mini-invasive nephroureterectomy: laparoscopic nephrectomy, transurethral excision of ureterovesical junction and lower abdominal incision. *Urol Int*. 2009;83(3):264-270. doi:10.1159/000241664.
- G71. Hora M, Stransky P, Urge T, et al. Transperitoneal less (laparoendoscopic single-site surgery) adrenalectomy (AE)-matched case-control study with standard laparoscopic adrenalectomy (SLAE). *Journal of Endourology*. 2014;28:A247-A248.
- G72. Grøndal AYK, Høgsbro M, Pryds K, Pedersen HB, Jacobsen H. Intra- and postoperative complications using LigaSure™ Small Jaw in patients undergoing thyroidectomy: a register-based study. *Eur Arch Otorhinolaryngol*. 2021;278(11):4491-4500. doi:10.1007/s00405-021-06685-w.
- G73. Hahn CH, Sorensen CH. LigaSure small jaws versus cold knife dissection in superficial parotidectomy. *Eur Arch Otorhinolaryngol*. 2013;270(4):1489-1492. doi:10.1007/s00405-012-2204-3.
- G74. Jensen JS, Antonsen HK, Durup J. Two years of experience with robot-assisted anti-reflux surgery: A retrospective cohort study. *International Journal of Surgery*. 2017;39:260-266. doi:10.1016/j.ijsu.2017.02.014.
- G75. Nielsen, H. U. K., Trolle W, Rubek N, Homoe P. New technique using LigaSure for endoscopic mucomyotomy of Zenker's diverticulum: Diverticulotomy made easier. *Laryngoscope*. 2014;124(9):2039-2042.

- G76. Schulze S, Damgaard B, Jorgensen LN, Larsen SS, Kristiansen VB. Cystic duct closure by sealing with bipolar electrocoagulation. *JSLS*. 2010;14(1):20-22. doi:10.4293/108680810X12674612014347.
- G77. Abd Ellatif, M. E., Abdallah E, Askar W, et al. Long term predictors of success after laparoscopic sleeve gastrectomy. *Int J Surg*. 2014;12(5):504-508. doi:10.1016/j.ijssu.2014.02.008.
- G78. Abdallah A, Abdel-Hakim M, El-Feel A. Laparo-endoscopic single-site radical cystectomy with orthotopic urinary diversion: Technique, feasibility, and the 3-year follow-up. *Arab J Urol*. 2014;12(3):229-233. doi:10.1016/j.aju.2014.05.003.
- G79. Al Anwar A, Youssef T, El Shinawi M, Abdel Moniem M. Stapled versus stapleless laparoscopic splenectomy in the management of hematological disorders of the spleen 'building our experience in egypt'. *Surgical Endoscopy and Other Interventional Techniques*. 2013;27:S44-S44.
- G80. Allam IS, Makled AK, Gomaa IA, El Bishry, G. M., Bayoumy HA, Ali DF. Total laparoscopic hysterectomy, vaginal hysterectomy and total abdominal hysterectomy using electrosurgical bipolar vessel sealing technique: a randomized controlled trial. *Archives of Gynecology and Obstetrics*. 2014:no pagination-no pagination.
- G81. Bessa SS. LigaSure™ vs. conventional diathermy in excisional hemorrhoidectomy: A prospective, randomized study. *Diseases of the Colon and Rectum*. 2008;51(6):940-944.
- G82. El Ghazeery M, Elsayaf M, Ashour M, Metwaly M, Hashish M. Effect of laparoscopic sleeve gastrectomy on metabolic syndrome in adolescent obesity. *Journal of Pediatric Endoscopic Surgery*. 2020;2(4):175-181. doi:10.1007/s42804-020-00075-5.
- G83. Elbanna ME, Mashaal A, El-Anwar A, Fouad O. Stapleless vs. stapled laparoscopic splenectomy. *Surgical Endoscopy and Other Interventional Techniques*. 2011;25:S141-S141.
- G84. Elhao M, Abdallah K, Serag I, El-Laithy M, Agur W. Efficacy of using electrosurgical bipolar vessel sealing during vaginal hysterectomy in patients with different degrees of operative difficulty: a randomised controlled trial. *Eur J Obstet Gynecol Reprod Biol*. 2009;147(1):86-90. doi:10.1016/j.ejogrb.2009.07.011.
- G85. Fareed M, El-Awady S, Abd-El monaem H, Aly A. Randomized trial comparing LigaSure to closed Ferguson hemorrhoidectomy. *Tech Coloproctol*. 2009;13(3):243-246. doi:10.1007/s10151-009-0520-4.

- G86. Labib K, El-Mansy H, Ibrahim M. Laparoscopic Hysterectomy of a 20 weeks uterus weighting 1 Kg with extensive adhesions. *Gynecological Surgery*. 2016;13(G1 Supplement 1):S173-S173.
- G87. Nabil TM, Moustafa KG, Balamoun HA. Comparison of Laparoscopic Sleeve Gastrectomy and Laparoscopic Gastric Plication: Two-Year Follow-Up Results. *Bariatric Surgical Practice and Patient Care*. 2019;14(1):3-6. doi:10.1089/bari.2018.0028.
- G88. Osman A, Helmy AS, Mikhail S, AlAyat AA, Serour DK, Ibrahim MY. Early Effects of Laparoscopic Sleeve Gastrectomy and Laparoscopic One-Anastomosis Gastric Bypass on Portal Venous Flow: a Prospective Cohort Study. *Obesity Surgery*. 2021;31(6):2410-2418. doi:10.1007/s11695-021-05267-3.
- G89. Sakr MF, Moussa MM, Elserafy M. Ligasure hemorrhoidectomy versus Stapled hemorrhoidopexy: a prospective randomized clinical trial. *Minerva Chir*. 2010;65(3):251-258.
- G90. Shady NW, Farouk HA, Sallam HF. Perioperative Outcomes of LigaSure Versus Standard Ligature Technique among Overweight and Obese Women Undergoing Abdominal Hysterectomy: A Randomized Clinical Trial. *Journal of Gynecologic Surgery*. 2021;37(4):345-351. doi:10.1089/gyn.2020.0158.
- G91. Youssef T, Mahdy T, Farid M, Latif AA. Thyroid surgery: Use of the LigaSure Vessel Sealing System versus conventional knot tying. *International Journal of Surgery*. 2008;6(4):323-327.
- G92. Youssef T, Abdalla E, El-Alfy K, Dawoud I, Morshed M, Farid M. Impact of Botulinum Neurotoxin Pyloric Injection During Laparoscopic Sleeve Gastrectomy on Postoperative Gastric Leak: a Clinical Randomized Study. *Obes Surg*. 2016;26(3):494-504. doi:10.1007/s11695-015-1794-4.
- G93. Lepner U, Vaasna T. Ligasure vessel sealing system versus conventional vessel ligation in thyroidectomy. *Scand J Surg*. 2007;96(1):31-34.
- G94. Kiiski J, Kaartinen I, Kotaluoto S, Kuokkanen H. Modified approach for endoscopic harvest of the latissimus dorsi free flap with CO2 insufflation and standard laparoscopic equipment. *Microsurgery*. 2017;37(5):383-387. doi:10.1002/micr.30127.
- G95. Kossi J, Luostarinen M, Kontula I, Laato M. Laparoscopic sigmoid and rectal resection using an electrothermal bipolar vessel sealing device. *J Laparoendosc Adv Surg Tech A*. 2007;17(6):719-722. doi:10.1089/lap.2006.0238.

- G96. Luoto TT, Pakarinen MP, Koivusalo A. Long-term outcomes after pediatric splenectomy. *Surgery*. 2016;159(6):1583-1590. doi:10.1016/j.surg.2015.12.014.
- G97. Blanc P, Pradat C, Breton C, Kassir R. Comparison Between Caïman® and Ligasure® in Laparoscopic Sleeve Gastrectomy: a Retrospective Study of 200 Patients. *Obesity Surgery*. 2020;30(7):2804-2806. doi:10.1007/s11695-020-04588-z.
- G98. Catheline JM, Fysekidis M, Dbouk R, et al. Weight loss after sleeve gastrectomy in super superobesity. *J Obes*. 2012;2012:959260. doi:10.1155/2012/959260.
- G99. Evrard S, Becouarn Y, Brunet R, Fonck M, Larrue C, Mathoulin-Pelissier S. Could bipolar vessel sealers prevent bile leaks after hepatectomy? *Langenbeck's Archives of Surgery*. 2007;392(1):41-44.
- G100. Fakhry N, Michel J, Lagier A, Turner F, Dessi P, Giovanni A. Usefulness of the LigaSure™ vessel sealing system during major head and neck cancer surgery: A prospective study. *Oral Oncology*. 2011;47:S144-S144.
- G101. Fakhry N, Michel J, Santini L, et al. Utility of the LigaSure vessel sealing system during major head and neck cancer surgery. *J Laryngol Otol*. 2012;126(8):815-817. doi:10.1017/S0022215112001168.
- G102. Gaillard M, Tranchart H, Maitre S, Perlemuter G, Lainas P, Dagher I. Preoperative Detection of Sarcopenic Obesity Helps to Predict the Occurrence of Gastric Leak After Sleeve Gastrectomy. *Obesity Surgery*. 2018;28(8):2379-2385. doi:10.1007/s11695-018-3169-0.
- G103. Gaillard M, Tranchart H, Lainas P, Ferretti S, Perlemuter G, Dagher I. Single-port laparoscopic sleeve gastrectomy as a routine procedure in 1000 patients. *Surgery for obesity and related diseases : official journal of the American Society for Bariatric Surgery*. 2016;12(7):1270-1277. doi:10.1016/j.soard.2016.01.011.
- G104. Guédon ACP, Meij SEP, Osman KNMMH, et al. Deep learning for surgical phase recognition using endoscopic videos. *Surg Endosc*. 2021;35(11):6150-6157. doi:10.1007/s00464-020-08110-5.
- G105. Kassir R, Blanc P, Lointier P, Breton C, Tiffet O. Using Ligasure™ or Harmonic Ace® in Laparoscopic Sleeve Gastrectomies? A Prospective Randomised Study. *Obesity Surgery*. 2015;25(10):1944-1946. doi:10.1007/s11695-015-1809-1.
- G106. Kassir R, Blanc P, Lointier P, et al. Four tips and tricks for performing laparoscopic gastric bypass: A tip for patient positioning for performing laparoscopic gastric bypass: French position. *American Surgeon*. 2015;81(11):E388-E389.

<https://www.embase.com/search/results?subaction=viewrecord&id=L607351507&from=export>.

- G107. Lainas P, Filippo G de, Di Giuro G, Mikhael R, Bougneres P, Dagher I. Laparoscopic Sleeve Gastrectomy for Adolescents Under 18 Years Old with Severe Obesity. *Obesity Surgery*. 2020;30(1):267-273. doi:10.1007/s11695-019-04150-6.
- G108. Lainas P, Dupond-Athenor A, Tranchart H, Dagher I. Safety and Feasibility of Single-Port Sleeve Gastrectomy Following Liver Transplantation. *Obesity Surgery*. 2018;28(3):874-876. doi:10.1007/s11695-017-3080-0.
- G109. Lainas P, Tranchart H, Tzanis D, Dagher I. Single-incision laparoscopic rectal resection for cancer: A video. *Surgical Endoscopy and Other Interventional Techniques*. 2016;30:S170-S170.
- G110. Lainas P, Derienne J, Zervaki S, et al. Left Hypochondrium or Transumbilical Single-Incision Laparoscopic Sleeve Gastrectomy for the Treatment of Severe Obesity: Surgical Technique and Results of a Tertiary Referral Bariatric Center. *Obesity Surgery*. 2021;31(11):5063-5070. doi:10.1007/s11695-021-05563-y.
- G111. Lauroy A, Verhaeghe C, Vidal F, Parant O, Legendre G, Guerby P. Perioperative outcomes using LigaSure compared with conventional technique in peripartum hysterectomy. *Arch Gynecol Obstet*. 2020;301(1):229-234. doi:10.1007/s00404-019-05398-0.
- G112. Leroy J, Cahill RA, Asakuma M, Dallemagne B, Marescaux J. Single-access laparoscopic sigmoidectomy as definitive surgical management of prior diverticulitis in a human patient. *Archives of Surgery*. 2009;144(2):173-179.
- G113. Mauriello C, Chouillard E, D'alessandro A, et al. Retrospective comparison of single-port sleeve gastrectomy versus three-port laparoscopic sleeve gastrectomy: A propensity score adjustment analysis. *Obesity Surgery*. 2018;28(7):2105-2112. doi:10.1007/s11695-018-3244-6.
- G114. Pourcher G, Filippo G de, Ferretti S, Piquard C, Dagher I, Bougneres P. Short-term results of single-port sleeve gastrectomy in adolescents with severe obesity. *Surg Obes Relat Dis*. 2015;11(1):65-69. doi:10.1016/j.soard.2014.05.029.
- G115. Tranchart H, Rebibo L, Gaillard M, et al. Short-term outcomes of single-port versus conventional laparoscopic sleeve gastrectomy: a propensity score matched analysis. *Surgical Endoscopy and Other Interventional Techniques*. 2020;34(9):3978-3985. doi:10.1007/s00464-019-07175-1.

- G116. Dingemann C, Petersen C, Kuebler JF, Ure BM, Lacher M. Laparoscopic transperitoneal heminephrectomy for duplex kidney in infants and children: a comparative study. *J Laparoendosc Adv Surg Tech A*. 2013;23(10):889-893. doi:10.1089/lap.2013.0194.
- G117. Gehrig T, Muller-Stich BP, Kenngott H, et al. LigaSure versus conventional dissection technique in pancreatoduodenectomy: a pilot study. *Am J Surg*. 2011;201(2):166-170. doi:10.1016/j.amjsurg.2010.02.023.
- G118. Goretzki PE, Schwarz K, Lammers BJ. Implementing the general use of dissection devices in thyroid surgery from prospective randomized trial to daily use. *Surg Technol Int*. 2009;18:86-92.
- G119. Hasanov M, Denschlag D, Seemann E, Gitsch G, Woll J, Klar M. Bipolar vessel-sealing devices in laparoscopic hysterectomies: a multicenter randomized controlled clinical trial. *Arch Gynecol Obstet*. 2018;297(2):409-414. doi:10.1007/s00404-017-4599-y.
- G120. Karcz WK, Kuesters S, Marjanovic G, Grueneberger JM. Duodeno-enteral omega switches - more physiological techniques in metabolic surgery. *Wideochir Inne Tech Maloinwazyjne*. 2013;8(4):273-279. doi:10.5114/witm.2013.39647.
- G121. Küsters S, Marjanovic G, Baumann T, Kulemann B, Holzner P, Karcz WK. The laparoscopic banded gastric bypass - Operation technique. *Wideochir Inne Tech Maloinwazyjne*. 2010;5(1):38-41. doi:10.5114/witm.2010.13606.
- G122. Lacher M, Dingemann C, Petersen C, Kuebler JF, Ure BM. Laparoscopic transperitoneal heminephrectomy for duplex kidney above and below one year of age: Feasibility and long-term results. *J Laparoendosc Adv Surg Tech A*. 2013;23(12):A-141.
- G123. Lieber J, Urla CI, Baden W, Schafer J, Kirschner HJ, Fuchs J. Experiences and challenges of thoroscopic lung surgery in the pediatric age group. *Int J Surg*. 2015;23(GPt A):169-175. doi:10.1016/j.ijsu.2015.10.005.
- G124. Patrzyk M, Glitsch A, Hoene A, Bernstorff W von, Heidecke CD. Laparoscopic partial splenectomy using a detachable clamp with and without partial splenic embolisation. *Langenbeck's Archives of Surgery*. 2011;396(3):397-402.
- G125. Turial S, Engel V, Sultan T, Schier F. Closure of the cystic duct during laparoscopic cholecystectomy in children using the LigaSure Vessel Sealing System. *World journal of surgery*. 2011;35(1):212-216.
- G126. Turial S, Trendafilow M, Schier F. LigaSure™ vessel sealing system (LVSS) used for appendix transection during laparoscopic appendectomy for acute appendicitis in children. *European Surgery - Acta Chirurgica Austriaca*. 2011;43:5.

- G127. Uzunoglu FG, Bockhorn M, Fink JA, et al. LigaSure™ Vs. Conventional Dissection Techniques in Pancreatic Surgery-a Prospective Randomised Single-Centre Trial. *J Gastrointest Surg*. 2013;17(3):494-500.
- G128. Daskalopoulos G, Karyotis I, Heretis I, Delakas D. Electrothermal bipolar coagulation for radical prostatectomies and cystectomies: a preliminary case-controlled study. *Int Urol Nephrol*. 2004;36(2):181-185.
- G129. Manouras A, Markogiannakis H, Koutras AS, et al. Thyroid surgery: comparison between the electrothermal bipolar vessel sealing system, harmonic scalpel, and classic suture ligation. *American Journal of Surgery*. 2008;195(1):48-52.
- G130. Prokopakis EP, Lachanas VA, Helidonis ES, Velegrakis GA. The use of LigaSure vessel sealing system in patients undergoing parotid gland surgery. *Otolaryngology - Head and Neck Surgery*. 2006;135(1):172-173.
- G131. Sdralis E, Argentou M, Mead N, Kehagias I, Alexandridis T, Kalfarentzos F. A prospective randomized study comparing patients with morbid obesity submitted to sleeve gastrectomy with or without omentectomy. *Obesity Surgery*. 2013;23(7):965-971.
- G132. Tepetes K, Christodoulidis G, Spyridakis EM, Chatzitheofilou C. Tissue preserving hepatectomy by a vessel sealing device. *J Surg Oncol*. 2008;97(2):165-168. doi:10.1002/jso.20889.
- G133. Tsamis D, Natoudi M, Arapaki A, et al. Using Ligasure™ or Harmonic Ace® in Laparoscopic Sleeve Gastrectomies? A Prospective Randomized Study. *Obesity Surgery*. 2015;25(8):1454-1457. doi:10.1007/s11695-014-1551-0.
- G134. Ho NS, Keung HY, Leung YW, Wong CY. Comparison of surgical outcome on thyroidectomy using the LigaSure Small Jaw and the Harmonic Focus+. *Surgical Practice*. 2019;23:15.
<https://www.embase.com/search/results?subaction=viewrecord&id=L632062431&from=export>.
- G135. Lai C-S, Liu S-W, Wong J-T, Wong S-H, Ng E-W. Portomesenteric vein thrombosis following laparoscopic greater curve plication for morbid obesity. *Obesity Research and Clinical Practice*. 2018;12(6):578-582. doi:10.1016/j.orcp.2018.10.001.
- G136. Sztipits T, Meszaros P, Duboczki Z, Tarpay A, Mersich T. Laparoscopic subtotal gastric resection with Roux-en-Y reconstruction for antral adenocarcinoma – Case report. *Magy Seb*. 2015;68(6):235-238. doi:10.1556/1046.68.2015.6.4.

- G137. Gadkari Y, Singbal N, Kumbhare S, Joshi S, Puntambekar SP. Duplication of Steps of Laparoscopic Radical Hysterectomy By Gynaecology Fellows. *J Minim Invasive Gynecol*. 2015;22(G6S):S120. doi:10.1016/j.jmig.2015.08.344.
- G138. Gupta G, Madnani M, Patel T, et al. Pancreatic transaction using ligasure system in laparoscopic distal pancreatectomy: Case report. *HPB (Oxford)*. 2014;16:728.
- G139. Kriplani A, Garg P, Sharma M, Lal S, Agarwal N. A review of total laparoscopic hysterectomy using ligasure uterine artery-sealing device: AIIMS experience. *J Laparoendosc Adv Surg Tech A*. 2008;18(6):825-829.
- G140. Kumar S, Khandelwal N, Kumar A, et al. Simultaneous living donor liver transplant with sleeve gastrectomy for metabolic syndrome and NASH-related ESLD—First report from India. *Indian Journal of Gastroenterology*. 2017;36(3):243-247. doi:10.1007/s12664-017-0753-5.
- G141. Misra MC, Aggarwal S, Guleria S, Seenu V, Bhalla AP. Clipless and sutureless laparoscopic surgery for adrenal and extra-adrenal tumors. *JSLs*. 2008;12(3):252-255.
- G142. Roy KK, GC N, Singhal S, et al. Impact of energy devices on the post-operative systemic immune response in women undergoing total laparoscopic hysterectomy for benign disease of the uterus. *J Turk Ger Gynecol Assoc*. 2018;19(1):1-6. doi:10.4274/jtgga.2017.0076.
- G143. Sugoor D, Joshi G, Puntambekar SP, Kumbhare S, Sharma V, Parikh H. Single institutional experience of 410 cases of type B & type C (querleu morrow classification) laparoscopic radical hysterectomy. *Journal of Minimally Invasive Gynecology*. 2015;22(G6 SUPPL. 1):S91-S92.
- G144. Alibakhshi A, Khashayar P, Mahmudzade H, et al. Low lying rectal anastomosis with stapler in proctectomy patients: The experience of a new technique. *Annals of Oncology*. 2010;21:vi116-vi117.
- G145. Bouini NR, Pour MA, Hadian HS, Pour, R. A. M. Randomised clinical trial comparing ligasure hemorrhoidectomy with conventional hemorrhoidectomy. *Journal of Mazandaran University of Medical Sciences*. 2012;22(93):66-73.
- G146. Fazeli MS, Safari S, Kazemeini A, et al. A prospective study comparing ligasure and open hemorrhoidectomy. *Tehran University Medical Journal*. 2011;69(8):495-501.
- G147. Foroutan H, Banani SA, Hoseini, S. M. V., Ali B, Soltan N. Laparoscopic splenectomy without stapler. *J Laparoendosc Adv Surg Tech A*. 2011;21(4):A42-A42.

- G148. Hosseini S, Hosseini SA, Al-Hurry A, Khazraei H, Ganji F, Sadeghi F. Comparison of early results and complications between multi-and single-port sleeve gastrectomy: A randomized clinical study. *Iranian Journal of Medical Sciences*. 2017;42(3):251-257. <https://www.embase.com/search/results?subaction=viewrecord&id=L616099269&from=export>.
- G149. Niazi M, Talebpour M, Maleki A, Nejhad K. Short-term outcome of laparoscopic total vertical gastric plication (TVGP); Interim report of a long-term prospective study. *Obesity Surgery*. 2011;21(8):1119.
- G150. Ramouz A, Rasihashemi SZ, Safaeiyan A, Hosseini M. Comparing postoperative complication of LigaSure Small Jaw instrument with clamp and tie method in thyroidectomy patients: a randomized controlled trial IRCT2014010516077N1. *World J Surg Oncol*. 2018;16(1):154. doi:10.1186/s12957-018-1448-9.
- G151. Shabahang H, Maddah G, Tavassoli A, et al. Laparoscopic splenectomy: ligasure or clip ligation? *Surg Laparosc Endosc Percutan Tech*. 2012;22(2):136-138. doi:10.1097/SLE.0b013e318246d9a4.
- G152. Tavassoli A, Abdollahi A, Golmohammadzadeh H. Resection and reconstruction of small intestine with ligasure and linear stapler: A comparative study. *Tehran University Medical Journal*. 2012;70(5):289-294.
- G153. Zeinoddini A, Heidari R, Talebpour M. Laparoscopic gastric plication in morbidly obese adolescents: A prospective study. *Surgery for Obesity and Related Diseases*. 2014;10(6):1135-1140.
- G154. Hussein KK. Use of bipolar cautery in thyroidectomy. *Indian Journal of Public Health Research and Development*. 2019;10(3):524-527. doi:10.5958/0976-5506.2019.00551.5.
- G155. Jawad AM. Comparison of Surgical Outcomes between Patients Treated with The Harmonic Scalpel and Ligasure Device During Total/ Near Total Open Thyroidectomy. *IJFMT*. 2020;14(2). doi:10.37506/ijfmt.v14i2.3068.
- G156. Browne C, Lonergan PE, Healy N, et al. Zero-ischaemia laparoscopic partial nephrectomy-a single surgeon experience in a tertiary referral hospital. *Irish Journal of Medical Science*. 2015;184(G9 SUPPL. 1):S384-S384.
- G157. Martin ST, Heeney A, Pierce C, O'Connell PR, Hyland JM, Winter DC. Use of an electrothermal bipolar sealing device in ligation of major mesenteric vessels during

laparoscopic colorectal resection. *Tech Coloproctol*. 2011;15(3):285-289.
doi:10.1007/s10151-011-0707-3.

G158. Moftah M, Nazour F, Cunningham M, Cahill RA. Single port laparoscopic surgery for patients with complex and recurrent Crohn's disease. *J Crohns Colitis*. 2014;8(9):1055-1061. doi:10.1016/j.crohns.2014.02.003.

G159. Ronan A, Narouz F, Cunningham M. Single port laparoscopic ileocolic redo resection for recurrent crohn's stricture. *Surgical Endoscopy and Other Interventional Techniques*. 2015;29:S291-S291.

G160. Smolarek S, Salih A, Kazanowski M, Walsh TN. Laparoscopic assisted total gastrectomy for gastric cancer - operative technique. *Wideochir Inne Tech Maloinwazyjne*. 2015;10(1):133-137. doi:10.5114/wiitm.2015.49092.

G161. Abu sneineh M. OAGB After LAGB Video Presentation. *Surgery for obesity and related diseases : official journal of the American Society for Bariatric Surgery*. 2019;15(10):S255-S256. doi:10.1016/j.soard.2019.09.005.

G162. Becker A, Khromov Y, Sayfan J. Delayed bleeding following ligasure hemorrhoidectomy - Is there a remedy? *Techniques in Coloproctology*. 2010;14(1):89.

G163. Pencovich N, Lahat G, Goldray O, Abu-Abeid S, Klausner JM, Meron Eldar S. Safety and Outcome of Laparoscopic Sleeve Gastrectomy Following Removal of Adjustable Gastric Banding: Lessons from 109 Patients in a Single Center and Review of the Literature. *Obesity Surgery*. 2017;27(5):1266-1270. doi:10.1007/s11695-016-2463-y.

G164. Stanevsky Y, Benjamin S, Tsivian M, Sidi A, Tsivian A. Clipless laparoscopic partial adrenalectomy. *Surgical Endoscopy and Other Interventional Techniques*. 2013;27:S209-S209.

G165. Tsivian A, Benjamin S, Sidi AA. A sealed laparoscopic nephroureterectomy: a new technique. *Eur Urol*. 2007;52(4):1015-1019. doi:10.1016/j.eururo.2006.10.036.

G166. Tsivian A, Stanevsky Y, Tsivian M, Benjamin S, Sidi A. Intramural tunnel ureterectomy for treatment of upper tract urothelial tumor. *Journal of Urology*. 2015;193(G4 SUPPL. 1):e499-e499.

G167. Yakirevitch A, Nakache G, Lipschitz N, Alon EE, Wolf M, Talmi YP. Use of the vessel sealing system in tracheostomy. *Isr Med Assoc J*. 2013;15(6):275-276.

G168. Tolone S, Bondanese M, Ruggiero R, et al. Outcomes of sutureless total thyroidectomy in elderly. *Int J Surg*. 2016;33 Suppl 1:S16-9. doi:10.1016/j.ijsu.2016.05.056.

- G169. Angrisani L, Santonicola A, Hasani A, Nosso G, Capaldo B, Iovino P. Five-year results of laparoscopic sleeve gastrectomy: effects on gastroesophageal reflux disease symptoms and co-morbidities. *Surg Obes Relat Dis*. 2016;12(5):960-968. doi:10.1016/j.soard.2015.09.014.
- G170. Cassinotti E, David G, Colombo EM, et al. Transanal minimally invasive surgery for total mesorectal excision: A new approach for low rectal cancer. *Surgical Endoscopy and Other Interventional Techniques*. 2015;29:S257-S257.
- G171. Cesana G, Cioffi S, Giorgi R, et al. Proximal Leakage After Laparoscopic Sleeve Gastrectomy: an Analysis of Preoperative and Operative Predictors on 1738 Consecutive Procedures. *Obesity Surgery*. 2018;28(3):627-635. doi:10.1007/s11695-017-2907-z.
- G172. Coiro S, Frattaroli FM, Lucia F de, et al. A comparison of the outcome using Ligasure™ small jaw and clamp-and-tie technique in thyroidectomy: a randomized single center study. *Langenbecks Arch Surg*. 2015;400(2):247-252. doi:10.1007/s00423-014-1270-y.
- G173. Divisi D, Di Leonardo G, Vico A de, Crisci R. Electrocautery versus Ultracision versus LigaSure in Surgical Management of Hyperhidrosis. *Thorac Cardiovasc Surg*. 2015;63(8):729-734. doi:10.1055/s-0035-1550339.
- G174. Gambardella C, Clarizia G, Patrone R, et al. Advanced hemostasis in axillary lymph node dissection for locally advanced breast cancer: new technology devices compared in the prevention of seroma formation. *BMC Surg*. 2019;18(GSuppl 1):125. doi:10.1186/s12893-018-0454-8.
- G175. Grieco M, Apa D, Spoletini D, Grattarola E, Carlini M. Major vessel sealing in laparoscopic surgery for colorectal cancer: a single-center experience with 759 patients. *World J Surg Oncol*. 2018;16(1):101. doi:10.1186/s12957-018-1402-x.
- G176. Martucci N, Tracey M, La Rocca A, La Manna C, Luca G de, Rocco G. A pilot prospective randomized, controlled trial comparing LigaSure™ tissue fusion technology with the ForceTriad™ energy platform to the electrosurgical pencil on rates of atrial fibrillation after pulmonary lobectomy and mediastinal lymphadenectomy. *European journal of cardio-thoracic surgery : official journal of the European Association for Cardio-thoracic Surgery*. 2015;47(1):e13-e18.
- G177. Milito G, Cadeddu F, Lisi G. Haemorrhoidectomy with LigaSure™. *Colorectal Disease*. 2014;16:157.

- G178. Nespoli L, Antolini L, Stucchi C, Nespoli A, Valsecchi MG, Gianotti L. Axillary lymphadenectomy for breast cancer. A randomized controlled trial comparing a bipolar vessel sealing system to the conventional technique. *Breast*. 2012;21(6):739-745. doi:10.1016/j.breast.2012.08.003.
- G179. Olmi S, Uccelli M, Oldani A, et al. Laparoscopic Surgery of Gastric Cancer with D2 Lymphadenectomy and Omentum Preservation: Our 10 Years Experience. *Journal of Laparoendoscopic and Advanced Surgical Techniques*. 2020;30(7):749-758. doi:10.1089/lap.2019.0781.
- G180. Olmi S, Caruso F, Uccelli M, Cioffi S, Ciccarese F, Cesana G. Laparoscopic sleeve gastrectomy combined with Rossetti fundoplication (R-Sleeve) for treatment of morbid obesity and gastroesophageal reflux. *Surgery for obesity and related diseases : official journal of the American Society for Bariatric Surgery*. 2017;13(12):1945-1950. doi:10.1016/j.soard.2017.08.017.
- G181. Petrucci F, Antoniazzi G, Valerio S. Use of the ligasure Maryland jaw in the radical prostatectomy: Our experience. *Anticancer Research*. 2015;35(6):3637-3638. <http://ar.iiarjournals.org/content/35/6/3615.full.pdf>.
- G182. Romanini MV, Vidal C, Godoy J, Morovic CG. Laparoscopically harvested omental flap for breast reconstruction in Poland syndrome. *J Plast Reconstr Aesthet Surg*. 2013;66(11):e303-9. doi:10.1016/j.bjps.2013.08.003.
- G183. Ruggiero R, Docimo G, Bosco A, et al. Update on sutureless thyroidectomy. *G Chir*. 2018;39(1):45-50. doi:10.11138/gchir/2018.39.1.045.
- G184. Testini M, Gurrado A, Lissidini G, et al. A single-surgeon randomized trial comparing conventional technique, ligasure vessel sealing system and harmonic scalpel during total thyroidectomy for substernal goiter. *European Surgical Research*. 2015;55:67.
- G185. Testini M, Pasculli A, Di Meo G, et al. Advanced vessel sealing devices in total thyroidectomy for substernal goitre: A retrospective cohort study. *Int J Surg*. 2016;35:160-164. doi:10.1016/j.ijsu.2016.09.092.
- G186. Velotti N, Manigrasso M, Di Lauro K, et al. Comparison between LigaSure™ and Harmonic® in Laparoscopic Sleeve Gastrectomy: A Single-Center Experience on 422 Patients. *J Obes*. 2019;2019:3402137. doi:10.1155/2019/3402137.
- G187. Zappa MA, Musolino C, Giorgini E, Porta A. Ligasure Maryland: A new approach in colorectal laparoscopic surgery. *Surgical Endoscopy and Other Interventional Techniques*. 2015;29:S96-S96.

- G188. Roberts PO, Mitchell, D. I. G., Toppin PJ, et al. Laparoscopic liver resection in Jamaica. *West Indian Medical Journal*. 2012;61:42-43.
- G189. Chikazawa K, Imai K, Irie Y, et al. Proper device usage in total laparoscopic hysterectomy :A retrospective analysis in single institute. *Journal of Obstetrics and Gynaecology Research*. 2018;44(8):1590-1591. doi:10.1111/jog.13762.
- G190. Fujita J, Takiguchi S, Nishikawa K, et al. Randomized controlled trial of the LigaSure vessel sealing system versus conventional open gastrectomy for gastric cancer. *Surg Today*. 2014;44(9):1723-1729. doi:10.1007/s00595-014-0930-4.
- G191. Fukunaga M, Lee Y, Nagakari K, et al. Laparoscopic total pelvic exenteration for rectal cancer. *Colorectal Disease*. 2015;17:65.
- G192. Kawahara H, Watanabe K, Hojo S, et al. Clipless laparoscopic total proctocolectomy for ulcerative colitis. *Journal of Gastroenterology and Hepatology (Australia)*. 2016;31:439.
- G193. Hirano J, Kinoshita S, Kazama T, et al. Value of Vessel Sealing System in Skin-Sparing Mastectomy in Patients with Early Breast Cancer. *Clinics in Surgery*. 2018;3(1).
- G194. Inoue Y, Yamashita N, Ueo H, et al. The Clinical Usefulness of the LigaSure™ Small Jaw in Axillary Lymph Node Dissection in Patients with Breast Cancer. *Anticancer Research*. 2018;38(4):2359-2362. doi:10.21873/anticancer.12483.
- G195. Ishikawa H, Morita Y, Shozu M, Sasaki H. Robot-assisted radical hysterectomy for early cervical cancer. *International Journal of Gynecological Cancer*. 2018;28:106. doi:10.1097/01.IGC.0000546279.09648.02.
- G196. Kanamoto M, Imura S, Morine Y, et al. Effective use of a vessel-sealing system for laparoscopic unroofing of liver cysts. *Asian journal of endoscopic surgery*. 2015;8(1):91-94.
- G197. Kobayashi M, Okamoto K, Maeda H, et al. Efficacy of hand-assisted laparoscopic splenectomy for portal hypertension patients with thrombocytopenia. *Surgical Endoscopy and Other Interventional Techniques*. 2015;29:S546-S546.
- G198. Kunisaki C, Makino H, Kimura J, et al. Application of reduced-port laparoscopic total gastrectomy in gastric cancer preserving the pancreas and spleen. *Gastric Cancer*. 2015;18(4):868-875. doi:10.1007/s10120-014-0441-4.
- G199. Kunisaki C, Makino H, Oshima T, et al. Application of the transorally inserted anvil (OrVil™) after laparoscopy-assisted total gastrectomy. *Surgical Endoscopy and Other Interventional Techniques*. 2011;25(4):1300-1305.

- G200. Kuroda H, Dejima H, Mizumo T, Sakakura N, Sakao Y. A new LigaSure technique for the formation of segmental plane by intravenous indocyanine green fluorescence during thoracoscopic anatomical segmentectomy. *J Thorac Dis*. 2016;8(6):1210-1216. doi:10.21037/jtd.2016.04.61.
- G201. Kushibiki T, Ebihara Y, Hontani K, et al. The surgical outcomes of totally laparoscopic stomach-partitioning gastrojejunostomy for gastric outlet obstruction: A retrospective, cohort study. *Surgical Laparoscopy, Endoscopy and Percutaneous Techniques*. 2018;28(2):e49-e53. doi:10.1097/SLE.0000000000000501.
- G202. Nishioka K, Iwazawa T, Watanabe R, et al. Thoracoscopic mediastinal lymph node dissection for lung cancer by the energy device (ligasure). *Surgical Endoscopy and Other Interventional Techniques*. 2015;29:S292-S292.
- G203. Ono S, Tsuji Y, Baba K, Usui Y, Yanagisawa S, Maeda K. New operative strategy for refractory microcystic lymphangioma. *Surg Today*. 2014;44(6):1184-1187. doi:10.1007/s00595-013-0806-z.
- G204. Seki T, Hayashida T, Takahashi M, Jinno H, Kitagawa Y. A randomized controlled study comparing a vessel sealing system with the conventional technique in axillary lymph node dissection for primary breast cancer. *Springerplus*. 2016;5(1):1004. doi:10.1186/s40064-016-2710-7.
- G205. Takiguchi N, Nagata M, Soda H, et al. Multicenter randomized comparison of LigaSure versus conventional surgery for gastrointestinal carcinoma. *Surg Today*. 2010;40(11):1050-1054. doi:10.1007/s00595-009-4234-z.
- G206. Toishi M, Yoshida K, Agatsuma H, et al. Usefulness of vessel-sealing devices for <7 mm diameter vessels: A randomized controlled trial for human thoracoscopic lobectomy in primary lung cancer. *Interactive Cardiovascular and Thoracic Surgery*. 2014;19(3):448-455.
- G207. Tsunetzuka Y, Tanaka N, Fujimori H. Thoracoscopic lobectomy using vessel sealing system: Comparison liga sure blunt tip versus liga sure Maryland for systematic lymph node dissection and pulmonary artery ligation. *Innovations: Technology and Techniques in Cardiothoracic and Vascular Surgery*. 2015;10:S112-S112.
- G208. Yoshimoto M, Endo K, Hanaki T, et al. Effectiveness of the LigaSure Small Jaw Vessel-Sealing System in Hepatic Resection. *Yonago Acta Med*. 2014;57(2):93-98.

- G209. Zuiki T, Hosoya Y, Sakuma Y, et al. Laparoscopic gastric devascularization without splenectomy is effective for the treatment of gastric varices. *International Journal of Surgery Case Reports*. 2016;19:119-123. doi:10.1016/j.ijscr.2015.12.036.
- G210. Bani-Irshaid I, Nussair B. Ligasure Vessel Sealing System versus Conventional Suture Ligation Method in Vaginal Hysterectomy. *JRMS*. 2014;21(4):25-31. doi:10.12816/0008062.
- G211. Obeidat F, Shanti H. Early Weight Loss as a Predictor of 2-Year Weight Loss and Resolution of Comorbidities After Sleeve Gastrectomy. *Obesity Surgery*. 2016;26(6):1173-1177. doi:10.1007/s11695-015-1903-4.
- G212. Shanti H, Almajali N, Al-Shamaileh T, Samarah W, Mismar A, Obeidat F. Helicobacter pylori Does not Affect Postoperative Outcomes After Sleeve Gastrectomy. *Obesity Surgery*. 2017;27(5):1298-1301. doi:10.1007/s11695-016-2470-z.
- G213. AlSabah S, Alsharqawi N, Almulla A, et al. Approach to Poor Weight Loss After Laparoscopic Sleeve Gastrectomy: Re-sleeve Vs. Gastric Bypass. *Obes Surg*. 2016;26(10):2302-2307. doi:10.1007/s11695-016-2119-y.
- G214. Abboud B, Sleilaty G, Braidy C, Melkane A, Nasr F. Thyroidectomy in patients at high-risk of bleeding: can it be safely performed? *Minerva Chir*. 2009;64(6):565-571.
- G215. Hussein M. Laparoscopic sleeve gastrectomy state of the Art technique. *Obesity Surgery*. 2013;23(8):1102.
- G216. Hussein MK. Laparoscopic sleeve gastrectomy-retro gastric medial to lateral approach is the state of the art technique. *Surgical Endoscopy and Other Interventional Techniques*. 2015;29:S204-S204.
- G217. Essadi F, Elmehashi M, Sharkasi A, Ataweel SM. Ligasure versus traditional suture ligation in abdominal hysterectomy. *Fertility and Sterility*. 2009;92(G3 SUPPL. 1):S124-S124.
- G218. Yeoh WS, Lee YS, Shanggar K, Selvalingam S, Razack AH, Ong TA. Audit of living related donor nephrectomy in University Malaya Medical Centre. *BJU International*. 2016;118:19-20.
<https://www.embase.com/search/results?subaction=viewrecord&id=L614675164&from=export>.
- G219. Bonfante E, Marquez G, Castillo M, Castro N, Moctezuma J. Minilaparotomy - A minimally invasive (3-4cm incision) approach for total hysterectomy, even in large

uterus. Some suggestions to the original technique. *Int J Gynaecol Obstet*. 2009;107:S131-S131.

- G220. Hurtado Estrada G. Modified radical hysterectomy for cervical cancer Ib1. *International Journal of Gynecological Cancer*. 2019;29:A80. doi:10.1136/ijgc-2019-IGCS.180.
- G221. Leal C, Ceron R, Rubio V, Unda ME. Ultrasonic Energy (Harmonic Ace) Versus Advanced Bipolar Energy (Ligasure) in a Laparoscopic Hysterectomies. *J Minim Invasive Gynecol*. 2015;22(G6S):S166. doi:10.1016/j.jmig.2015.08.620.
- G222. Lopez J, Vilallonga R, Targarona E, Balague C, Balibrea JM. Can ligasure be used to perform sleeve gastrectomy? *Surgical Endoscopy and Other Interventional Techniques*. 2015;29:S205-S205.
- G223. Luna-Ortiz K, Villavicencio-Valencia V, Carmona-Luna T, Cano-Valdez AM, Herrera Gomez A. Carotid body tumour resection with LigaSure device. *Acta Otorrinolaringol Esp*. 2010;61(1):6-11. doi:10.1016/j.otorri.2009.07.003.
- G224. Navarro Martin JM, Rodriguez B, Escalera S, Lopez JC, Gonzalez M. 2989 Uterus with Complete Double Cervix and Complete Longitudinal Vaginal Septum, A Case Report. *Journal of Minimally Invasive Gynecology*. 2019;26(7):S206. doi:10.1016/j.jmig.2019.09.422.
- G225. Gurbadam U. "Adult pancreatic hemangioma" case report and literature review. *Surgical Endoscopy and Other Interventional Techniques*. 2014;28:368.
- G226. Rabii R, El Mrini M, Elkattani Y, Aboutaieb R, Meziane F. Retroperitoneal adrenalectomy using ligasure. *Journal of Endourology*. 2013;27:A373-A373.
- G227. Janssen PF, Brolmann, H. A. M., van Kesteren PJM, et al. Perioperative outcomes using LigaSure compared with conventional bipolar instruments in laparoscopic hysterectomy: A randomised controlled trial. *BJOG: An International Journal of Obstetrics and Gynaecology*. 2011;118(13):1568-1575.
- G228. Lakeman M, Kruitwagen RF, Vos MC, Roovers JP. Electrosurgical bipolar vessel sealing versus conventional clamping and suturing for total abdominal hysterectomy: a randomized trial. *J Minim Invasive Gynecol*. 2008;15(5):547-553. doi:10.1016/j.jmig.2008.05.011.
- G229. Lakeman MM, The S, Schellart RP, et al. Electrosurgical bipolar vessel sealing versus conventional clamping and suturing for vaginal hysterectomy: a randomised controlled trial. *BJOG*. 2012;119(12):1473-1482. doi:10.1111/j.1471-0528.2012.03484.x.

- G230. Nieuwenhuis DH, Draaisma WA, Verberne, G. H. M., van Overbeeke AJ, Consten, E. C. J. Transanal endoscopic operation for rectal lesions using two-dimensional visualization and standard endoscopic instruments: A prospective cohort study and comparison with the literature. *Surgical Endoscopy and Other Interventional Techniques*. 2009;23(1):80-86.
- G231. Ramanah R, Mottet N, Chehab M, Riethmuller D. Rouhier's procedure : Modified lefort colpocleisis with hysterectomy for complete pelvic organ eversion. *Neurourology and Urodynamics*. 2015;34:S261-S261.
- G232. Schiphorst AH, Twigt BA, Elias SG, van Dalen T. Randomized clinical trial of LigaSure versus conventional suture ligation in thyroid surgery. *Head Neck Oncol*. 2012;4:2. doi:10.1186/1758-3284-4-2.
- G233. Suisted P, Chittenden B. Perioperative outcomes of total laparoscopic hysterectomy at a regional hospital in New Zealand. *Australian and New Zealand Journal of Obstetrics and Gynaecology*. 2017;57(1):81-86. doi:10.1111/ajo.12570.
- G234. Fretland AA, Rosok BI, Edwin B. Laparoscopic right hepatectomy extended to segment 4 a: Anterior approach. *HPB (Oxford)*. 2014;16:703.
- G235. Fretland AA, Sokolov A, Postriganova N, et al. Open liver resection for colorectal metastases induced a significant increase in interleukin-6 as compared to laparoscopic technique. *HPB (Oxford)*. 2014;16:250-251.
- G236. Langebrekke A, Qvigstad E. Total Laparoscopic Hysterectomy with Single-Port Access without Vaginal Surgery. *Journal of Minimally Invasive Gynecology*. 2009;16(5):609-611.
- G237. Rosok BI, Fretland A, Edwin B. Laparoscopic liver resection of lesions in the posterior segments using an anterolateral approach. *HPB (Oxford)*. 2014;16:703-704.
- G238. Machado NO, Al Kindy N, Chopra PJ. Laparoscopic Splenectomy Using LigaSure. *JSLS*. 2010;14(4):547-552. doi:10.4293/108680810X12924466008286.
- G239. Bakhtiar N, Moosa FA, Jaleel F, Qureshi NA, Jawaaid M. Comparison of hemorrhoidectomy by LigaSure with conventional Milligan Morgan's hemorrhoidectomy. *Pak J Med Sci*. 2016;32(3):657-661. doi:10.12669/pjms.323.9976.
- G240. Shaikh S, Laghari AA, Abro S, Sushel C, Ali N, Memon A. Efficacy of electrosurgical vessel sealer (LigaSure™) devices in modified radical mastectomy in primary breast carcinoma. *Pakistan Journal of Medical and Health Sciences*. 2021;15(3):1105-1107. <https://www.embase.com/search/results?subaction=viewrecord&id=L2012008967&from=export>.

- G241. Shoukat H, Iqbal M, Ullah S, Mirza A, Dar UF. Comparison of hemorrhoidectomy using bipolar diathermy vs harmonic scalpel. *Pakistan Journal of Medical and Health Sciences*. 2016;10(2):489-491.
- G242. Vinatea J de. Laparoscopic partial splenectomy. *Surgical Endoscopy and Other Interventional Techniques*. 2013;27:S43-S43.
- G243. Adamczewski Z, Krol A, Kaluzna-Markowska K, Brzezinski J, Lewinski A, Dedecjus M. Lateral spread of heat during thyroidectomy using different haemostatic devices. *Ann Agric Environ Med*. 2015;22(3):491-494. doi:10.5604/12321966.1167721.
- G244. Dziki L, Mik M, Trzcinski R, et al. Surgical treatment of haemorrhoidal disease -the current situation in Poland. *Przegląd Gastroenterologiczny*. 2016;11(2):111-114.
- G245. Jakimowicz T, Nazarewski S, Otto M, Ciacka T, Lewandowska D, Szmidt J. Hand assisted laparoscopic donor nephrectomy-single center experience. *Transplantation*. 2012;94:1100.
- G246. Lewandowski T, Giaro M, Fiedorowicz W. Application of the LigaSure™ vessel sealing system for closure of the cystic duct during laparoscopic cholecystectomy. *Polski Przegląd Chirurgiczny*. 2006;78(7):797-800.
- G247. Losin M, Golebiewski A, Czauderna P. Laparoscopic single incision approach to neuroblastoma: Early experience of a single center. *J Laparoendosc Adv Surg Tech A*. 2011;21(4):A68-A68.
- G248. Panek W, Szydelko T, Lewandowski J, et al. Simultaneous laparoscopic adrenalectomy and laparoscopic nephron-sparing surgery - new experience with port placement. *Wideochir Inne Tech Maloinwazyjne*. 2013;8(4):357-360. doi:10.5114/wiitm.2011.35351.
- G249. Rogulski R, Adamowicz-Salach A, Matysiak M, et al. Laparoscopic splenectomy for hereditary spherocytosis-preliminary report. *Eur J Haematol*. 2016;96(6):637-642. doi:10.1111/ejh.12649.
- G250. Szyca R, Leksowski K. Assessment of patients' quality of life after hemorrhoidectomy using the LigaSure device. *Wideochir Inne Tech Maloinwazyjne*. 2015;10(1):68-72.
- G251. Ferreira H, Smith AV, Vilaça J. Ghost Ileostomy in Anterior Resection for Bowel Endometriosis: Technical Description. *J Minim Invasive Gynecol*. 2020;27(5):1014-1016. doi:10.1016/j.jmig.2019.09.769.

- G252. Ferreira HC, Soares S, Caldas R, Morgado A, Tome A. Total laparoscopic hysterectomy with bilateral adnexectomy using a ligasure maryland jaw instrument. *Gynecological Surgery*. 2016;13(G1 Supplement 1):S166-S166.
- G253. Gingu C, Patrascu S, Surcel C, et al. Nonclamping open partial nephrectomy complications are decreased by "LigaSure" haemostasis and tumor bed compression with perirenal fat flap. *European Urology, Supplements*. 2011;10(9):653.
- G254. Graur F, Cote A, Szasz A, et al. [Left laparoscopic lobectomy]. *Chirurgia (Bucur)*. 2009;104(5):611-616.
- G255. Maghiar AM, Ciurtin DH, Sfirlea M, et al. Laparoscopic appendectomy: Our experience. *Surgical Endoscopy and Other Interventional Techniques*. 2013;27:S461-S461.
- G256. Maghiar AM, Sookha PR, Ciurtin DH, et al. Ovarian cysts requiring anexectomy-sils is another option-case presentation. *Surgical Endoscopy and Other Interventional Techniques*. 2012;26:S86-S86.
- G257. Martian BV, Diaconescu BI, Andrei G, Bratu MR, Craciun R, Beuran M. Wise steps in laparoscopic right hemicolectomy: Standardization of the technique in current use. *Surgical Endoscopy and Other Interventional Techniques*. 2016;30:S33-S33.
- G258. Molnar C, Botoncea M, Butiurca VO, Nicolescu CL, Molnar CV, Grigorescu BL. Monobloc Total Thyroidectomy Using Sealing Devices: Preliminary Results and Comparative Analysis. *Chirurgia (Bucur)*. 2019;114(5):668-673.
doi:10.21614/chirurgia.114.5.668.
- G259. Molnar C, Voidazan S, Rad CC, et al. Total thyroidectomy with LigaSure Small Jaw versus conventional thyroidectomy - a clinical study. *Chirurgia (Bucur)*. 2014;109(5):608-612.
- G260. Patrascu S, Georgescu EF, Bica M, et al. Retrospective analysis on the feasibility of laparoscopic approach for perforated acute appendicitis. *Surgical Endoscopy and Other Interventional Techniques*. 2016;30:S81-S81.
- G261. Gestkov K, Esakov Y, Barsky B. Staplerless lung volume reduction surgery performing by ligasure system. *Interactive Cardiovascular and Thoracic Surgery*. 2011;13:S38-S38.
- G262. Kashintcev A, Kochanenko NY, Paveletc KV, Imyaninov EN. Comparison of dissection techniques in pancreatoduodenectomy. *Pancreatology*. 2013;13(2):e42-e43.
- G263. Poddoubnyi IV, Glybina TM, Tolstov KN. Single incision laparoscopic surgery (SILS) for adnexal pathology in adolescent girls. *J Laparoendosc Adv Surg Tech A*. 2011;21(4):A64-A64.

- G264. Rasulov AO, Mamedli ZZ, Gordeyev SS, Kozlov NA, Dzhumabaev HE. Short-term outcomes after transanal and laparoscopic total mesorectal excision for rectal cancer. *Tech Coloproctol*. 2016;20(4):227-234. doi:10.1007/s10151-015-1421-3.
- G265. Sazhin AV, Mosin SV, Tjagunov AE, Lebedeva JN. Results of laparoscopic surgery in patients with colorectal cancer. *Surgical Endoscopy and Other Interventional Techniques*. 2013;27:S108-S108.
- G266. Stognii N. The combination of mobilization aorta by 'suitcase handle' technique and the operation of proximal 'local prosthesis' in surgery of type b aortic dissection. *Interactive Cardiovascular and Thoracic Surgery*. 2011;12:S44-S44.
- G267. Tulina IA, Tsarkov PV, Kravchenko A. Laparoscopic 'no touch' intersphincteric rectal resection: Three step technique. *Surgical Endoscopy and Other Interventional Techniques*. 2014;28:S175-S175.
- G268. Bamehriz F, Althuwaini S, Alobaid O, et al. Is body mass index ≥ 50 kg/m² a predictor of higher morbidity for patients who have undergone laparoscopic sleeve gastrectomy? *Saudi Journal of Anaesthesia*. 2018;12(4):556-558. doi:10.4103/sja.SJA_140_18.
- G269. Ghasoup A, Sadieh O, Al Qurashi T. Laparoscopic splenectomy, is sutures and staplers mandatory? *Surgical Endoscopy and Other Interventional Techniques*. 2015;29:S234-S234.
- G270. Troisi RI, Elsheikh YM, Shagrani MA, Broering D. First fully laparoscopic donor hepatectomy for pediatric liver transplantation using the indocyanine green near-infrared fluorescence imaging in the Middle East: A case report. *Annals of Saudi Medicine*. 2014;34(4):354-357.
- G271. Zubaidi A, AlDeghaither S, Alkhayal K, Al-Obeed O. Single-incision laparoscopic colorectal surgery: A report of 33 cases in Saudi Arabia. *Annals of Saudi Medicine*. 2016;36(4):282-287.
- G272. Doklestic K, Karamarkovic A, Stefanovic B, et al. The efficacy of three transection techniques of the liver resection: a randomized clinical trial. *Hepatogastroenterology*. 2012;59(117):1501-1506. doi:10.5754/hge11552.
- G273. Kecmanovic DM, Pavlov MJ, Ceranic MS, Kerkez MD, Rankovic VI, Masirevic VP. Bulk agent *Plantago ovata* after Milligan-Morgan hemorrhoidectomy with LigaSure™. *Phytotherapy Research*. 2006;20(8):655-658.

- G274. Sretenovic A, Perisic V, Krstic Z, et al. Warren shunt combined with partial splenectomy for children with extrahepatic portal hypertension, massive splenomegaly, and severe hypersplenism. *Surg Today*. 2013;43(5):521-525. doi:10.1007/s00595-012-0405-4.
- G275. Lee S-H, Nguyen TK, Ong W-S, et al. Comparing the Utility and Surgical Outcomes of Harmonic Focus Ultrasonic Scalpel with Ligasure Small Jaw Bipolar Device in Thyroidectomies: A Prospective Randomized Controlled Trial. *Ann Surg Oncol*. 2019;26(13):4414-4422. doi:10.1245/s10434-019-07806-w.
- G276. Loh WL, Liauw J, Ngoi SS. Video presentation of laparoscopic splenectomy for massive splenomegaly. *Annals of the Academy of Medicine Singapore*. 2012;41(G9 SUPPL. 2):S28-S29.
- G277. Oh H, Kow, A. W. C. Single incision laparoscopic lateral sectionectomy of the liver for a large, simple hepatic cyst. *HPB (Oxford)*. 2015;17:283.
- G278. Tan KY, Zin T, Sim HL, Poon PL, Cheng A, Mak K. Randomized clinical trial comparing LigaSure haemorrhoidectomy with open diathermy haemorrhoidectomy. *Tech Coloproctol*. 2008;12(2):93-97. doi:10.1007/s10151-008-0405-y.
- G279. Tan YH, Lim YM, Ng YW, Tiong HY. Taking a Step Forward in Laparoscopic Donor Nephrectomy: Transvaginal Retrieval of Donor's Kidney. *J Laparoendosc Adv Surg Tech A*. 2016;26(9):721-724. doi:10.1089/lap.2016.0284.
- G280. Wijerathne S, Kow A, Shridhar I, Madhavan K, Chang S. Experience of single-incision laparoscopic de-roofing of liver cysts combined with cholecystectomy. *HPB (Oxford)*. 2014;16:697-698.
- G281. Cingel V, Zabojsnikova L, Kurucova P, Varga I. First experience with single incision laparoscopic surgery in Slovakia: concomitant cholecystectomy and splenectomy in an 11-year-old girl with hereditary spherocytosis. *Biomed Pap Med Fac Univ Palacky Olomouc Czech Repub*. 2014;158(3):479-485. doi:10.5507/bp.2013.058.
- G282. Cronje HS, Coning EC de. Electrosurgical bipolar vessel sealing during vaginal hysterectomy. *Int J Gynaecol Obstet*. 2005;91(3):243-245. doi:10.1016/j.ijgo.2005.08.019.
- G283. van Niekerk ML. Laparoscopic repair of Morgagni diaphragmatic hernia in children. *South African Journal of Surgery*. 2009;47(1):14-16.
<https://www.embase.com/search/results?subaction=viewrecord&id=L354407989&from=export>.

- G284. Chang YW, Kim HS, Jung SP, et al. Comparison of skin-sparing mastectomy using LigaSure™ Small Jaw and electrocautery. *World J Surg Oncol*. 2017;15(1):129. doi:10.1186/s12957-017-1199-z.
- G285. Hwang SO, Jung JH, Park HY, Kim WW. A Prospective, Randomized Study between the Small Jaw® and the Harmonic Focus® in Open Thyroidectomy. *Otolaryngol Head Neck Surg*. 2014;150(6):943-948. doi:10.1177/0194599814527730.
- G286. Jang HN, Kim SH, Cho YM, Park DJ. Bariatric Surgery for Cowden Syndrome with PTEN Mutation: a Case Report. *Obesity Surgery*. 2021;31(5):2316-2318. doi:10.1007/s11695-021-05231-1.
- G287. Jeong O, Jung MR, Ryu SY. Impact of various types of comorbidities on the outcomes of laparoscopic total gastrectomy in patients with gastric carcinoma. *Journal of Gastric Cancer*. 2018;18(3):253-263. doi:10.5230/jgc.2018.18.e27.
- G288. Kim SJ, Choi BJ, Lee SC. Novel approach of single-port laparoscopic appendectomy as a solo surgery: A prospective cohort study. *International Journal of Surgery*. 2015;21:1-7.
- G289. Kim YN, Kim SY, Kim HI. A comparison of the ligasure and harmonic scalpel in laparoscopic gastrectomy. *Surgical Endoscopy and Other Interventional Techniques*. 2013;27:S418-S418.
- G290. Kim MK, Lee S. LigaSure for the Creation of Bloodless Breast Pockets in Patients Undergoing Transaxillary Breast Augmentation. *Plast Reconstr Surg Glob Open*. 2020;8(12):e3295. doi:10.1097/GOX.00000000000003295.
- G291. Lee JS, Kim DG, Lee JW, et al. Usefulness of the LigaSure™ small jaw sealing device for breast reconstruction with a latissimus dorsi flap. *J Plast Surg Hand Surg*. 2019;53(5):295-300. doi:10.1080/2000656X.2019.1612753.
- G292. Lee CM, Da Park W, Park S, Kim J-H, Park S-H, Kim C-S. Lymph Node Dissection Using Bipolar Vessel-Sealing Device During Reduced Port Laparoscopic Distal Gastrectomy for Gastric Cancer: Result of a Pilot Study from a Single Institute. *J Laparoendosc Adv Surg Tech A*. 2017;27(11):1101-1108. doi:10.1089/lap.2016.0685.
- G293. Park MW, Baek S-K, Park E-H, Jung K-Y. Long-term voice outcome after thyroidectomy using energy based devices. *Auris Nasus Larynx*. 2018;45(3):527-532. doi:10.1016/j.anl.2017.08.007.
- G294. Ryu S-W. How can we overcome the blind loop syndrome after laparoscopic sleeve gastrectomy with jejunal bypass. *Surgery for obesity and related diseases : official*

Journal of the American Society for Bariatric Surgery. 2019;15(10):S157.
doi:10.1016/j.soard.2019.08.326.

- G295. Shin H. Current Trends in and Indications for Endoscopy-Assisted Breast Surgery for Breast Cancer. *Adv Exp Med Biol*. 2021;1187:567-590. doi:10.1007/978-981-32-9620-6_30.
- G296. Kim T-J, Song T, Choi CH, Lee J-W, Bae D-S, Kim B-G. Comparison of laparoendoscopic single-site hysterectomies: laparoscopic hysterectomy with some vaginal component versus laparoscopically assisted vaginal hysterectomy. *Journal of Laparoendoscopic and Advanced Surgical Techniques*. 2014;24(4):254-259. doi:10.1089/lap.2013.0507.
- G297. Yang J, Na YJ, Song YJ, Choi OH, Lee SK, Kim HG. The effectiveness of laparoendoscopic single-site surgery (LESS) compared with conventional laparoscopic surgery for ectopic pregnancy with hemoperitoneum. *Taiwan J Obstet Gynecol*. 2016;55(1):35-39. doi:10.1016/j.tjog.2015.02.005.
- G298. Yang JD, Lee J, Lee JS, Kim E-K, Park CS, Park HY. Aesthetic scar-less mastectomy and breast reconstruction. *J Breast Cancer*. 2021;24(1):22-33. doi:10.4048/jbc.2021.24.e11.
- G299. Alcaraz A, Musquera M, Izquierdo L, Garcia-Cruz E, Ciudin A, Peri L. Simultaneous native kidney and allograft laparoscopic nephrectomy for synchronous tumor. *Journal of Urology*. 2012;187(G4 SUPPL. 1):e869-e869.
- G300. Castellvi J, Sueiras A, Espinosa J, Vallet J, Gil V, Pi F. LigaSure™ versus diathermy hemorrhoidectomy under spinal anesthesia or pudendal block with ropivacaine: A randomized prospective clinical study with 1-year follow-up. *International Journal of Colorectal Disease*. 2009;24(9):1011-1018.
- G301. Cortadellas T, Cordoba O, Espinosa-Bravo M, et al. Electrothermal bipolar vessel sealing system in axillary dissection: a prospective randomized clinical study. *Int J Surg*. 2011;9(8):636-640. doi:10.1016/j.ijsu.2011.08.002.
- G302. Diez Redondo P, Nunez Rodriguez H, Penas I, et al. Flexible endoscopic treatment of zenker's diverticulum: A safe and efficient method. *United European Gastroenterology Journal*. 2013;1(G1 SUPPL. 1):A438-A438.
- G303. Diez-Redondo P, Brullet E, Juzgado-Lucas D, et al. Endoscopic zenker's diverticulotomy using a laparoscopic tissue sealing device and diverticuloscope: A safe, effective and durable technique. Results of a Spanish multicenter study. *Gastrointest Endosc*. 2016;83(G5 SUPPL. 1):AB581-AB582.

- G304. Balen EM, Suarez J, Oronoz B, Marin G, Lera JM. Hand-assisted laparoscopic restorative total proctocolectomy (ileoanal pouch) dividing the ano-rectum with the novel radial cartridge for endo-GIA. *Surgical Endoscopy and Other Interventional Techniques*. 2014;28:S130-S130.
- G305. Ferrer-Márquez M, Rubio-Gil F, Solvas-Salmerón MJ, et al. One Anastomosis Gastric Bypass in Patient with Situs Ambiguous. *Obesity Surgery*. 2021;31(5):2346-2347. doi:10.1007/s11695-021-05265-5.
- G306. Komorowski AL, Dominguez Usero D, Ramos Curado C, Martin Hidalgo Rodil, J. R. Simple technique of liver metastasectomy without Pringle manoeuvre. *European Journal of Surgical Oncology*. 2012;38(9):747-748.
- G307. Morales-Conde S, Dominguez G, Gomez JC, et al. Magnetic-assisted single-port sleeve gastrectomy. *Surg Innov*. 2013;20(4):NP9-11. doi:10.1177/1553350611427548.
- G308. Morales-Conde S, Rubio-Manzanares Dorado M, Barranco A, et al. Single-port Roux-en-Y gastric bypass: toward a less invasive procedure with the same surgical results and a better cosmesis. *Surg Endosc*. 2013;27(11):4347-4353. doi:10.1007/s00464-013-3070-9.
- G309. Ordonez J, Vidal O, Del Diaz Gobbo G, et al. Laparoscopic approach of a infrarenal paraganglioma: Case report. *Surgical Endoscopy and Other Interventional Techniques*. 2015;29:S289-S289.
- G310. Ortiz-Oshiro E, Exposito PB, Sierra JM, Gonzalez JD, Barbosa DS, Fernandez-Represa JA. Laparoscopic and robotic distal gastrectomy for gastrointestinal stromal tumour: case report. *Int J Med Robot*. 2012;8(4):491-495. doi:10.1002/rcs.1456.
- G311. Papis D, Parodi M, Herrerias F, et al. [Ligasure vs conventional diathermy in excisional haemorrhoidectomy: a retrospective study of a single center]. *Acta Gastroenterol Latinoam*. 2013;43(4):284-287.
- G312. Ruiz-Tovar J, Zubiaga L, Llaverio C, Diez M, Arroyo A, Calpena R. Serum cholesterol by morbidly obese patients after laparoscopic sleeve gastrectomy and additional physical activity. *Obes Surg*. 2014;24(3):385-389. doi:10.1007/s11695-013-1082-0.
- G313. Sabench Perefferrer F, Molina López A, Vives Espelta M, et al. Weight Loss Analysis According to Different Formulas after Sleeve Gastrectomy With or Without Antral Preservation: a Randomised Study. *Obesity Surgery*. 2017;27(5):1254-1260. doi:10.1007/s11695-016-2454-z.

- G314. Attner P, Hemlin C, Soderman AC. Ligasure versus diathermy scissors tonsillectomy: A controlled randomized study. *Acta Otolaryngol.* 2010;130(10):1180-1184. doi:10.3109/00016481003702544.
- G315. Chiaro MD, Blomberg J, Segersvard R, Rangelova E, Ansorge C, Verbeke C. The use of Ligasure has no negative impact on histology of resection margins in pancreatoduodenectomy (PD) specimens. *Pancreatology.* 2014;14(G3 SUPPL. 1):S39-S39.
- G316. Hubner M, Demartines N, Muller S, Dindo D, Clavien PA, Hahnloser D. Prospective randomized study of monopolar scissors, bipolar vessel sealer and ultrasonic shears in laparoscopic colorectal surgery. *Br J Surg.* 2008;95(9):1098-1104. doi:10.1002/bjs.6321.
- G317. Hubner M, Hahnloser D, Hetzer F, Muller MK, Clavien PA, Demartines N. A prospective randomized comparison of two instruments for dissection and vessel sealing in laparoscopic colorectal surgery. *Surg Endosc.* 2007;21(4):592-594. doi:10.1007/s00464-006-9034-6.
- G318. Chen CW, Lai CW, Chang YJ, Chen CM, Hsiao KH. Results of 666 consecutive patients treated with LigaSure hemorrhoidectomy for symptomatic prolapsed hemorrhoids with a minimum follow-up of 2 years. *Surgery.* 2013;153(2):211-218. doi:10.1016/j.surg.2012.06.035.
- G319. Chen CW, Lai CW, Chang YJ, Hsiao KH. Modified LigaSure hemorrhoidectomy for the treatment of hemorrhoidal crisis. *Surg Today.* 2014;44(6):1056-1062. doi:10.1007/s00595-013-0650-1.
- G320. Chiang F-Y, Lee KD, Tae K, et al. Comparison of hypocalcemia rates between LigaSure and clamp-and-tie hemostatic technique in total thyroidectomies. *Head Neck.* 2019;41(10):3677-3683. doi:10.1002/hed.25884.
- G321. Hong MK, Wang JH, Chu TY, Ding DC. Laparoendoscopic single-site hysterectomy with Ligasure is better than conventional laparoscopic assisted vaginal hysterectomy. *Gynecology and Minimally Invasive Therapy.* 2014;3(3):78-81.
- G322. Lee C-L, Wu K-Y, Huang C-Y, Yen C-F. Comparison of LigaSure™ tissue fusion system and a conventional bipolar device in hysterectomy via natural orifice transluminal endoscopic surgery (NOTES): A randomized controlled trial. *Taiwan J Obstet Gynecol.* 2019;58(1):128-132. doi:10.1016/j.tjog.2018.11.024.

- G323. Lee W-J, Chen J-C, Yao W-C, Taou J-J, Lee Y-C, Ser K-H. Transumbilical 2-site laparoscopic Roux-en-Y gastric bypass: Initial results of 100 cases and comparison with traditional laparoscopic technique. *Surgery for obesity and related diseases : official journal of the American Society for Bariatric Surgery*. 2012;8(2):208-213. doi:10.1016/j.soard.2010.12.004.
- G324. Lin WJ, Wang CC, Jiang RS, Huang YC, Ho HC, Liu SA. A prospective randomised trial of LigaSure Small Jaw® versus conventional neck dissection in head and neck cancer patients. *Clin Otolaryngol*. 2016. doi:10.1111/coa.12702.
- G325. Liu C-H, Wang C-C, Wu C-W, et al. Comparison of Surgical Complications Rates Between LigaSure Small Jaw and Clamp-and-Tie Hemostatic Technique in 1,000 Neuro-Monitored Thyroidectomies. *Front Endocrinol (Lausanne)*. 2021;12:638608. doi:10.3389/fendo.2021.638608.
- G326. Liu CH, Wu C-W. Safety of sutureless thyroidectomy with the application of ligasure small jaw in 1000 thyroidectomies. *Otolaryngology - Head and Neck Surgery*. 2020;163(G1 SUPPL):P57-P58. doi:10.1177/0194599820934780.
- G327. Pan H-M, Lin C-H, Hsu F-C, et al. Revisional Single Anastomosis Sleeve Ileal (SASI) Bypass for Failed Gastric Clipping and Proximal Jejunal Bypass (GC-PJB) (Video Report). *Obesity Surgery*. 2021;31(12):5510-5513. doi:10.1007/s11695-021-05733-y.
- G328. Wang JY, Lu CY, Tsai HL, et al. Randomized controlled trial of LigaSure with submucosal dissection versus Ferguson hemorrhoidectomy for prolapsed hemorrhoids. *World J Surg*. 2006;30(3):462-466. doi:10.1007/s00268-005-0297-1.
- G329. Wang S-F, Cheng H-T, Hsu J-T, et al. Simple and reliable method for gastric subepithelial tumor localization using endoscopic tattooing before totally laparoscopic resection. *Journal of Personalized Medicine*. 2021;11(9). doi:10.3390/jpm11090855.
- G330. Zachariah SK, Chang PC, Ooi AS, Hsin MC, Kin Wat, J. Y., Huang CK. Laparoscopic sleeve gastrectomy for morbid obesity: 5 years experience from an Asian center of excellence. *Obes Surg*. 2013;23(7):939-946. doi:10.1007/s11695-013-0887-1.
- G331. Hirunwiwatkul P, Tungkavivachagul S. A multicenter, randomized, controlled clinical trial of LigaSure small jaw vessel sealing system versus conventional technique in thyroidectomy. *Eur Arch Otorhinolaryngol*. 2013;270(7):2109-2114. doi:10.1007/s00405-012-2289-8.

- G332. Pattana-Arun J, Sooriprasoet N, Sahakijrungruang C, Tantiphlachiva K, Rojanasakul A. Closed vs ligasure hemorrhoidectomy: a prospective, randomized clinical trial. *J Med Assoc Thai*. 2006;89(4):453-458.
- G333. Naraynsingh V, Cawich S, Hassranah D, Maharaj R, Islam S, Singh Y. Retrograde Thyroidectomy for preservation of the External Branch of the Superior Laryngeal Nerve: A case series. *International Journal of Surgery Case Reports*. 2018;53:517-521. doi:10.1016/j.ijscr.2017.04.030.
- G334. Acar F, Çolak B, Şahin M. Comparison of primer and revisional laparoscopic mini-gastric bypass (Mgb) for failed restrictive procedures: 2-year results at a tertiary center revisional mini-gastric bypass (rmgbp) for failed restrictive procedures: 2-year results. *Annals of Clinical and Analytical Medicine*. 2021;12(2):124-128. doi:10.4328/ACAM.20210.
- G335. Aslaner A, Ongen A, Kosar M, et al. Relation between weight loss and age after laparoscopic sleeve gastrectomy. *European Review for Medical and Pharmacological Sciences*. 2015;19(8):1398-1402.
- G336. Atabey M, Aykota MR, Özel M, Arslan G. Short-term changes and correlations of plasma spexin, kisspeptin, and galanin levels after laparoscopic sleeve gastrectomy. *Surgery Today*. 2021;51(4):651-658. doi:10.1007/s00595-021-02240-6.
- G337. Aykan Yuksel B, Karadag B, Mulayim B. Comparison of the efficacy and safety of two advanced vessel sealing technologies in total laparoscopic hysterectomy. *J Obstet Gynaecol Res*. 2019;45(11):2220-2227. doi:10.1111/jog.14096.
- G338. Aytan H, Nazik H, Narin R, Api M, Tok EC. Comparison of the Use of LigaSure, HALO PKS Cutting Forceps, and ENSEAL Tissue Sealer in Total Laparoscopic Hysterectomy: A Randomized Trial. *Journal of Minimally Invasive Gynecology*. 2014;21(4):650-655.
- G339. Aziret M, Karaman K, Ercan M, Bostancı EB, Akoğlu M. Laparoscopic sleeve gastrectomy on a morbidly obese patient with situs inversus totalis: A case study and systematic review of the literature. *Obesity Research and Clinical Practice*. 2017;11(5):144-151. doi:10.1016/j.orcp.2016.12.003.
- G340. Bakacak Z, Bakacak M, Uzkar A, et al. The efficacy of LigaSure™ open instruments in cases of cesarean hysterectomy due to placenta percreta: a retrospective, record-based, comparative study. *J Matern Fetal Neonatal Med*. 2021;34(6):960-965. doi:10.1080/14767058.2020.1846177.

- G341. Bayrak M, Altintas Y. Comparing reinforcement options for the staple line of laparoscopic sleeve gastrectomy. *Journal of the College of Physicians and Surgeons Pakistan*. 2019;29(7):661-664. doi:10.29271/jcpsp.2019.07.661.
- G342. Binboga S, Isiksacan N, Kasapoglu P, et al. Plasma presepsin in determining gastric leaks following bariatric surgery. *Turkish Journal of Biochemistry*. 2019;44(4):565-573. doi:10.1515/tjb-2018-0220.
- G343. Büber İ, Aykota MR, Sevgican C, Adalı MK. The Effects of Laparoscopic Sleeve Gastrectomy on Cardiac Diastolic Function, Aortic Elasticity, and Atrial Electromechanics Delay. *Obesity Surgery*. 2021;31(8):3571-3578. doi:10.1007/s11695-021-05431-9.
- G344. Derici S, Atila K, Bora S, Yener S. Simple, Safe, and Cost-Effective Technique for Resected Stomach Extraction in Laparoscopic Sleeve Gastrectomy. *Gastroenterol Res Pract*. 2016;2016:7090128. doi:10.1155/2016/7090128.
- G345. Deveci HS, Kule ZA, Onur E, et al. Use of a monopolar laparoscopic electrosurgical instrument for transoral laryngoscopic removal of a giant fibrovascular polyp of the esophagus inlet. *Esophagus*. 2012:1-4.
- G346. Erdem MR, Tepeler A, Gunes M, et al. Laparoscopic decortication of hilar renal cysts using LigaSure. *JSLS : Journal of the Society of Laparoendoscopic Surgeons / Society of Laparoendoscopic Surgeons*. 2014;18(2):301-307.
- G347. Erdem H, Sisik A. Could the measurement of pressure in the section of stomach excised during laparoscopic sleeve gastrectomy give an indication of potential leakage in the remnant stomach? *International Journal of Clinical and Experimental Medicine*. 2017;10(5):8492-8499.
<https://www.embase.com/search/results?subaction=viewrecord&id=L616389345&from=export>.
- G348. Eyuboglu E, Ipek T. Laparoscopic floppy Nissen fundoplication: 16 years of experience from the historical clinic of Rudolph Nissen. *Hepatogastroenterology*. 2011;58(G110-111):1607-1610. doi:10.5754/hge10654.
- G349. Kaplan K, Turgut E, Okut G, Bag YM, Sumer F, Kayaalp C. Helicobacter pylori Increases Gastric Compliance on Resected Stomach After Laparoscopic Sleeve Gastrectomy. *Obesity Surgery*. 2021;31(11):4776-4780. doi:10.1007/s11695-021-05616-2.

- G350. Kirmizi S, Kayaalp C, Karagul S, et al. Comparison of Harmonic scalpel and Ligasure devices in laparoscopic Roux-en-Y gastric bypass. *Wideochir Inne Tech Maloinwazyjne*. 2017;12(1):28-31. doi:10.5114/wiitm.2017.66641.
- G351. Lale A, Aygen E, Kirkil C, Artas H, Yur M. Efficacy of Staple Line Reinforcement with Omentopexy during Laparoscopic Sleeve Gastrectomy on Postoperative Complications: Experience of a Single Center. *Surgical Laparoscopy, Endoscopy and Percutaneous Techniques*. 2021;31(2):181-187. doi:10.1097/SLE.0000000000000863.
- G352. Ozturk K, Kaya I, Turhal G, Ozturk A, Gursan G, Akyildiz S. A comparison of electrothermal bipolar vessel sealing system and electrocautery in selective neck dissection. *Eur Arch Otorhinolaryngol*. 2016;273(11):3835-3838. doi:10.1007/s00405-016-3999-0.
- G353. Senol K, Ferhatoglu MF, Kocaeli AA, Dundar HZ, Kaya E. Clinical Features and Short-Term Outcomes of Bariatric Surgery in Morbidly Obese Patients: Institutional Experience at a Rural Hospital. *Bariatric Surgical Practice and Patient Care*. 2021;16(1):61-67. doi:10.1089/bari.2020.0110.
- G354. Sood A, Pao K, Singh K, et al. Sutureless thyroidectomy with electrosurgical devices: A tertiary care center experience. *World Journal of Endocrine Surgery*. 2018;10(1):55. <https://www.embase.com/search/results?subaction=viewrecord&id=L622019469&from=export>.
- G355. Taşkın S, Şükür YE, Altın D, et al. Bipolar Energy Instruments in Laparoscopic Uterine Cancer Surgery: A Randomized Study. *Journal of Laparoendoscopic and Advanced Surgical Techniques*. 2018;28(6):645-649. doi:10.1089/lap.2017.0639.
- G356. Tukenmez M, Agcaoglu O, Aksakal N, et al. The use of Ligasure vessel sealing system in axillary dissection; effect on seroma formation. *Chirurgia (Bucur)*. 2014;109(5):620-625.
- G357. Tukenmez M, Ozden BC, Agcaoglu O, et al. Videoendoscopic single-port nipple-sparing mastectomy and immediate reconstruction. *J Laparoendosc Adv Surg Tech A*. 2014;24(2):77-82. doi:10.1089/lap.2013.0172.
- G358. Tuluçe K, Kara C, Tuluçe SY, et al. Early Reverse Cardiac Remodeling Effect of Laparoscopic Sleeve Gastrectomy. *Obesity Surgery*. 2017;27(2):364-375. doi:10.1007/s11695-016-2301-2.
- G359. Türkan A, Akkurt G, Yalaza M, et al. Effect of LigaSure™, Monopolar Cautery, and Bipolar Cautery on Surgical Margins in Breast-Conserving Surgery. *Breast Care (Basel)*. 2019;14(4):194-199. doi:10.1159/000493985.

- G360. Ulubay M, Ozturk M, Firatligil FB, et al. Comparison of the Perioperative Outcomes of a Vessel Sealing Instrument-Assisted Technique with a Conventional Technique in Abdominal Myomectomy. *J Clin Diagn Res*. 2016;10(1):QC01-3. doi:10.7860/JCDR/2016/14825.7053.
- G361. Yalcin Y, Yalcin SE. Perioperative outcomes of bipolar energy instruments in total laparoscopic hysterectomy. *Ginekol Pol*. 2019;90(11):640-644. doi:10.5603/GP.2019.0112.
- G362. Yavuz A, Bulus H, Taş A, Aydın A. Evaluation of Stump Pressure in Three Types of Appendectomy: Harmonic Scalpel, LigaSure, and Conventional Technique. *Journal of Laparoendoscopic and Advanced Surgical Techniques*. 2016;26(12):950-953. doi:10.1089/lap.2015.0551.
- G363. Grubnik V, Iuzvak OM, Grubnik Yu V, Fomenko VA. Laparoscopic operations for cirrhosis complicated by variceal bleeding. *Surgical Endoscopy and Other Interventional Techniques*. 2016;30:S19-S19.
- G364. Havrysh YI, Pavlovskyy MP, Krsychuk IV, Shvydky YB, Havrysh RY. Laparoscopic gastric plication is the bariatric efficient operation. Analysis of our experience. *Surgical Endoscopy and Other Interventional Techniques*. 2013;27:S136-S136.
- G365. Almerie MQ, Kerrigan DD. Laparoscopic Sleeve Gastrectomy with Simultaneous Laparoscopic Cystogastrostomy in a Patient with Super Obesity and a Pancreatic Pseudocyst. *Obesity Surgery*. 2021;31(4):1859-1861. doi:10.1007/s11695-020-05135-6.
- G366. Glover AR, Gundara JS, Lee JC, Sywak MS, Delbridge LW, Sidhu SB. Thermal sealing systems with and without tissue divider for total thyroidectomy. *ANZ J Surg*. 2014;84(5):383-385. doi:10.1111/ans.12406.
- G367. Hefni MA, Bhaumik J, El-Toukhy T, et al. Safety and efficacy of using the LigaSure vessel sealing system for securing the pedicles in vaginal hysterectomy: Randomised controlled trial. *BJOG: An International Journal of Obstetrics and Gynaecology*. 2005;112(3):329-333.
- G368. Heneghan HM, Kerrigan DD. Laparoscopic SADI-S as a salvage procedure for failed gastric bypass and sleeve gastrectomy: How I do it. *Surgery for obesity and related diseases : official journal of the American Society for Bariatric Surgery*. 2018;14(5):715-718. doi:10.1016/j.soard.2018.01.043.
- G369. Kumar S, Raobiakady R, Watkins D, Terlizzo M, Bhogal RH. Radical resection of large metastatic non-functioning pancreatic neuroendocrine carcinoma complicated by

splenic vein thrombosis and sinistral portal hypertension. *International Journal of Surgery Case Reports*. 2021;81:105724. doi:10.1016/j.ijscr.2021.105724.

- G370. Lawes DA, Palazzo FF, Clifton MA. The use of LigaSure™ haemorrhoidectomy in patients taking oral anticoagulation therapy. *Colorectal Disease*. 2004;6(2):111-112.
- G371. Miyagi K, Rossi SH, Malata CM, Forouhi P. Novel use of LigaSure Impact™ electro-surgical bipolar vessel sealing system in skin-sparing mastectomy. *Journal of Plastic, Reconstructive and Aesthetic Surgery*. 2015;68(6):e126-e128.
- G372. Peters CJ, Botterill I, Ambrose NS, Hick D, Casey J, Jayne DG. LigaSure™ vs conventional diathermy haemorrhoidectomy: Long-term follow-up of a randomised clinical trial. *Colorectal Disease*. 2005;7(4):350-353.
- G373. Rossi SH, Miyagi K, Malata CM, Forouhi P. Novel use of ligasure impact™ electrothermal bipolar diathermy vessel sealing system in skin-sparing mastectomy. *European Surgical Research*. 2015;55:65.
- G374. Simillis C, Charalambides M, Ariyanayagam D, et al. A systematic review and network meta-analysis comparing energy devices used for vessel ligation in colorectal surgery. *Colorectal Disease*. 2019;21:19. doi:10.1111/codi.14774.
- G375. Sivalingam S, Misra G. A Large cervical fibroid: Overcoming a technical challenge encountered at total laparoscopic hysterectomy. *Bipolar Disorders*. 2016;18:228.
- G376. Vadamalayan B, Lakshminarayanan BB, Fitzpatrick E, et al. Laparoscopic colectomy for ulcerative colitis with associated auto-immune liver disease and portal hypertension. *Journal of Pediatric Gastroenterology and Nutrition*. 2016;62:329.
- G377. Vintila AI, Marven S. Laparoscopic splenectomy: Which technique? *J Laparoendosc Adv Surg Tech A*. 2011;21(4):A57-A57.
- G378. Benlice C, Antoine HJ, Schauer PR. Laparoscopic Conversion of a Vertical Banded Gastroplasty to a Sleeve Gastrectomy in a Morbidly Obese Patient with a Complicated Medical History. *Obesity Surgery*. 2018;28(12):4095. doi:10.1007/s11695-018-3513-4.
- G379. Deniwar A, Mohamed SE, Hammad, A. R. Y., Al-Qurayshi Z, Aslam R, Kandil E. Efficacy and surgical outcomes of ligasure small jaw versus harmonic focus scalpel in thyroid and parathyroid surgery. *Journal of the American College of Surgeons*. 2014;219(G4 SUPPL. 1):e86-e86.
- G380. Elizabeth Best CA, Johnson-Obaseki SE. Do alternative sources of cautery decrease operative time for thyroidectomy? *Otolaryngology - Head and Neck Surgery*. 2019;161(2):P161. doi:10.1177/0194599819858142.

- G381. Ferko N, Wright GWJ, Syed I, Naoumtchik E, Tommaselli GA, Gangoli G. A device category economic model of electrosurgery technologies across procedure types: a U.S. hospital budget impact analysis. *J Med Econ.* 2021;24(1):524-535. doi:10.1080/13696998.2021.1915626.
- G382. Ginn DN, Shiber LD, Jan AG, et al. Comparison of bipolar electrosurgical devices for use in gynecologic laparoscopy. *J Minim Invasive Gynecol.* 2016;23(G7 Supplement 1):S39-S40.
- G383. Guan Z, Liu J, Liew-Spilger A, Rezai S, Guan X. Transvaginal Natural Orifice Transluminal Endoscopic Surgery Hysterectomy (VNOTES): A Walkthrough. *Journal of Minimally Invasive Gynecology.* 2019;26(7):S8. doi:10.1016/j.jmig.2019.09.031.
- G384. Hammad AY, Deniwar A, Al-Qurayshi Z, Mohamed HE, Rizwan A, Kandil E. A Prospective Study Comparing the Efficacy and Surgical Outcomes of Harmonic Focus Scalpel Versus LigaSure Small Jaw in Thyroid and Parathyroid Surgery. *Surg Innov.* 2016;23(5):486-489. doi:10.1177/1553350616639143.
- G385. Holloran-Schwartz MB, Gavard JA, Martin JC, Blaskiewicz RJ, Yeung, P. P., Jr. Single-Use Energy Sources and Operating Room Time for Laparoscopic Hysterectomy: A Randomized Controlled Trial. *J Minim Invasive Gynecol.* 2016;23(1):72-77. doi:10.1016/j.jmig.2015.08.881.
- G386. Kong K, Hwang RY, Rolandelli RH, Nemeth ZH. Differences in Post-Thyroidectomy Complications in Use of Vessel Sealant Devices Compared to Traditional Methods. *Journal of the American College of Surgeons.* 2019;229(4):e15. doi:10.1016/j.jamcollsurg.2019.08.858.
- G387. Lee E, Tong JY, Pasick LJ, et al. Complications associated with energy-based devices during thyroidectomy from 2010–2020. *World Journal of Otorhinolaryngology - Head and Neck Surgery.* 2021. doi:10.1016/j.wjorl.2021.04.008.
- G388. Love BR, Campbell R, Joyner J, Floyd D. Laparoscopic 10mm ligasure resection of left cornual pregnancy. *Journal of Minimally Invasive Gynecology.* 2011;18(G6 SUPPL. 1):S113-S113.
- G389. Marks J, Salem J, Josse J, Anderson B, Schoonyoung H. Initial single port robot experience in colorectal surgery. *Surg Endosc.* 2019;33:S82. doi:10.1007/s00464-019-06706-0.
- G390. Milsom JW, Trencheva K, Momose K, Peev M, Shukla P, Garrett K. A randomized controlled trial comparing THUNDERBEAT (TB, Olympus, Japan) to the Maryland

LigaSure™ (Medtronic, USA) energy device during laparoscopic colon surgery. *Diseases of the Colon and Rectum*. 2019;62(6):e353. doi:10.1097/DCR.0000000000001415.

G391. Okoli FC, Golfeyz S, Mongiu AK, Albert JI. Colonoscopic-assisted resection of a prolapsed sigmoid colon lipoma. *American Journal of Gastroenterology*. 2020;115(GSUPPL):S912-S913. doi:10.14309/01.ajg.0000709108.36978.47.

G392. Parikh M, Somoza E, Chopra A, et al. Thrombophilia prevalence in patients seeking laparoscopic sleeve gastrectomy: extended chemoprophylaxis may decrease portal vein thrombosis rate. *Surgery for obesity and related diseases : official journal of the American Society for Bariatric Surgery*. 2020;16(7):839-843. doi:10.1016/j.soard.2020.03.008.

G393. Perez-Milicua G, Hakim J, Adeyemi O, Dietrich JE. Longitudinal Vaginal Septum Resection Using the Ligasure Device. *J Pediatr Adolesc Gynecol*. 2016;29(6):e95-e96. doi:10.1016/j.jpbg.2016.06.002.

G394. Rothenberg SS, Shipman K, Kay S, et al. Thoracoscopic segmentectomy for congenital and acquired pulmonary disease: a case for lung-sparing surgery. *J Laparoendosc Adv Surg Tech A*. 2014;24(1):50-54. doi:10.1089/lap.2013.0337.

G395. Shakir S, Spencer AB, Kozak GM, Nathan SL, Soriano IS, Kanchwala SK. Laparoscopically Assisted DIEP Flap Harvest Minimizes Fascial Incision in Autologous Breast Reconstruction. *Plastic and reconstructive surgery*. 2020;146(3):265e-275e. doi:10.1097/PRS.00000000000007048.

G396. Shakir S, Spencer AB, Piper M, Kozak GM, Soriano IS, Kanchwala SK. Laparoscopy allows the harvest of the DIEP flap with shorter fascial incisions as compared to endoscopic harvest: A single surgeon retrospective cohort study. *Journal of Plastic, Reconstructive and Aesthetic Surgery*. 2021;74(6):1203-1212. doi:10.1016/j.bjps.2020.10.098.

G397. Thompson, I. M., 3rd, Kappa SF, Morgan TM, et al. Blood loss associated with radical cystectomy: a prospective, randomized study comparing Impact LigaSure vs. stapling device. *Urol Oncol*. 2014;32(1):45 e11-5. doi:10.1016/j.urolonc.2013.06.006.

G398. Tozzi F, Lewis A, Guye M, Lee B, Singh G. LigaSure™-assisted laparoscopic left lateral sectionectomy for HCC. *Annals of Surgical Oncology*. 2016;23(G1 SUPPL. 1):S45-S45.

- G399. Castellanos Coronel, L. M., Rodriguez JI, Castellanos Pinto, L. M., Parra JC, Yaniz AR. Comparison of vessel sealing system and conventional ligation in thyroidectomy. *Salus*. 2014;18(1):13-17.
- G400. Leyba JL, Llopis SN, Aulestia SN, Ochoa R, Azuaje E. Stapleless laparoscopic sleeve gastrectomy. Preliminary report. *Surgery for obesity and related diseases : official journal of the American Society for Bariatric Surgery*. 2017;13(4):701-704. doi:10.1016/j.soard.2016.11.026.
- G401. Parilli AM, Garcia W, Galdon I, Contreras G. Laparoscopic resection of giant hepatic cystadenoma. *J Laparoendosc Adv Surg Tech A*. 2013;23(12):A-58.

©2022 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 11/2022 - US-SE-2200137