


**Medtronic**

Engineering the extraordinary

# Breathe easier



Reducing  
surgical smoke in  
the operating room

PlasmaBlade™ X 3.0S  
LIGHT dissection device



Aquamantys™ 6.0  
bipolar sealer

# Minimize smoke. Maximize performance.

## PlasmaBlade™ X Dissection Device

Precise dissection and coagulation of soft tissue

### How it works

A 99.5%-insulated electrode precisely concentrates brief pulses of radiofrequency energy, resulting in lower operating temperatures than traditional electrosurgery.<sup>24</sup>

### Key effects

- Equivalence to scalpel in healed incision strength, inflammatory response following skin incision, and healed scar width<sup>15,16,17</sup>
- 74% less thermal injury depth observed when compared to traditional electrosurgery ( $p < 0.001$ )<sup>17</sup>
- Reduced surgical smoke<sup>14\*</sup>
- Maintains cutting effectiveness and hemostatic ability even when submerged in liquified tissue or blood<sup>15</sup>

### Benefits

- Single “skin-to-skin” instrument eliminates device exchanges and, as a result, the risk of scalpel injuries incurred during exchanges<sup>18,19</sup>
- May be used in closer proximity to adjacent structures with less risk of thermal exposure<sup>15</sup>
- Improved OR efficiency due to elimination of device exchanges<sup>18</sup>

## Aquamantys™ Bipolar Sealer

Hemostatic sealing of soft tissue and bone

### How it works

Using a proprietary combination of radiofrequency energy and saline, Aquamantys™ bipolar sealers operate at nearly 200°C less than traditional electrosurgical devices<sup>23</sup>

### Key effects

- Broad-plane hemostasis
- Reduced intraoperative and total blood loss<sup>20</sup>
- Reduced surgical smoke<sup>14</sup>

### Benefits

- Reduced transfusion rates<sup>20</sup>
- Maintained hemoglobin levels<sup>20</sup>
- Reduced post-operative drainage<sup>21</sup>
- Fewer hematomas<sup>22</sup>
- Reduced post-operative incidence of hemarthrosis<sup>21</sup>
- Reduced length of stay and OR time
- Fewer discharges to home health or skilled nursing facilities<sup>22</sup>

Lose the  
smoke.

Keep the  
performance.



# The risk is clear.

The surgical community has been examining the hazards of exposure to electrosurgical smoke for more than 40 years.<sup>1,2</sup> The effluent generated when tissue is vaporized and charred during electrosurgical dissection and coagulation has been shown to contain:

- Volatile organic compounds<sup>3,4</sup>
- Known carcinogens and mutagens<sup>6,7</sup>
- Infectious viral, bacterial, and aerosolized malignant cells<sup>8,9</sup>
- A host of poorly characterized, microscopic particulates<sup>6,10,11,12</sup>

Chronic inhalation of these components by operating room personnel poses a potential long-term health hazard.<sup>5,6,9</sup>

Evidence has shown that 1 gram of tissue produces a smoke plume with an equivalent mutagenicity to six unfiltered cigarettes.<sup>13</sup> In addition, surgical masks are good at capturing larger sized particles, generally  $\geq 5$  micrometers, but do not provide adequate protection in filtering smoke.<sup>13</sup>

# So are the results.

Medtronic's advanced energy devices use proprietary technology to operate at lower temperatures than traditional electrosurgery.<sup>23,24</sup> The result is low-temperature hemostasis and dissection with technology that has been shown to produce less surgical smoke than traditional instruments.<sup>14</sup>

# A difference that's easy to see

In a benchtop study comparing PlasmaBlade™ X technology and Aquamantys™ bipolar sealers to a traditional electro-surgical instrument, researchers collected smoke produced during activation on steak with a suction device. The images shown on the right show smoke particulate captured from filters located within each device.


In addition, particulate composition collected from the study was weighed and compared in the charts below.

## PlasmaBlade™ X technology produced:

- 95% less surgical smoke in cut mode<sup>14</sup> (n=6; p=.01)
- 84% less surgical smoke in coagulation mode<sup>14</sup> (n=6; p=.01)

The Aquamantys™ bipolar sealer produced a nearly negligible level of surgical smoke compared to traditional electro-surgery.<sup>14</sup>

## CUT



PlasmaBlade™ X Technology CUT 6



Traditional Electro-surgical Instrument CUT 40 watts


## Area COAG



PlasmaBlade™ X Technology Area COAG 6




Traditional Electro-surgical Instrument Area COAG 40 watts




Aquamantys™ 6.0 Bipolar Sealer Area COAG 170 watts

## Average weight of surgical smoke particulate: CUT mode<sup>14</sup>



- Electro-surgical pencil (ES)
- PlasmaBlade™ X technology (PB)

## Average weight of surgical smoke particulate: COAG mode<sup>14</sup>



- Electro-surgical pencil (ES)
- PlasmaBlade™ X technology (PB)
- Aquamantys™ bipolar sealer (AQM)

## Ordering information

Description	Catalog number
Aquamantys™ 6.0 bipolar sealer	23-112-1
Aquamantys™ 2.3 bipolar sealer	23-113-1
Aquamantys™ MBS	23-301-1
AEX™ generator	40-405-1
PlasmaBlade™ X 4.0	PS200-040
PlasmaBlade™ X 3.0S	PS210-030S
PlasmaBlade™ X 4.0S	PS200-041S
PlasmaBlade™ X 3.0S LIGHT	PS210-030S-LIGHT

## References

1. Sebben JE. The hazards of electrosurgery. *J Am Acad Dermatol* 1987; 16:869-872.
2. Ulmer BC. The hazards of surgical smoke. *AORN J* 2008; 87:721-34; quiz 735.
3. Al Sahaf OS, Vega-Carrascal I, Cunningham FO, McGrath JP, Bloomfield FJ. Chemical composition of smoke produced by high-frequency electrosurgery. *Ir J Med Sci* 2007;176:229-232.
4. Gianella M, Sigrist MW. Chemical analysis of surgical smoke by infrared laser spectroscopy. *Appl Physics B* 2012; 109:485-496.
5. Hill DS, O'Neill JK, Powell RJ, Oliver DW: Surgical smoke—A health hazard in the operating theatre: A study to quantify exposure and a survey of the use of smoke extractor systems in UK plastic surgery units. *J Plast, Reconstr & Aesthetic Surg* 2012; 65:911-916.
6. Tseng HS, Liu SP, Uang SN et al. Cancer risk of incremental exposure to polycyclic aromatic hydrocarbons in electrocautery smoke for mastectomy personnel. *World J Surg Oncol* 2014; 12:31.
7. Krones CJ, Conze J, Hoelzl F et al. Chemical composition of surgical smoke produced by electrocautery, harmonic scalpel and argon beaming—a short study. *Eur Surg* 2007; 39:118-121.
8. Sawchuk WS, Weber PJ, Lowy DR, Dzubow LM. Infectious papillomavirus in the vapor of warts treated with carbon dioxide laser or electrocoagulation: detection and protection. *J Am Acad Dermatol* 1989; 21:41-49.
9. Mowbray N, Ansell J, Warren N, Wall P, Torkington J. Is surgical smoke harmful to theater staff? A systematic review. *Surg Endosc* 2013; 27:3100-3107.
10. Brüske-Hohlfeld I, Preissler G, Jauch K-W et al. Surgical smoke and ultrafine particles. *J Occup Med Toxicol*; 2008, 3:31.
11. DesCôteaux J-G, Picard P, Poulin EC, Baril M: Preliminary study of electrocautery smoke particles produced in vitro and during laparoscopic procedures. *Surg Endosc* 1996; 10:152-158.
12. Barrett WL, Surgical smoke: a review of the literature. *Surg Endosc* 2003; 17:979-987.
13. Alp E. et al; Surgical smoke and infection control. *The Hospital Infection Society* 2006; 62, 1-5.
14. Data on file.
15. Loh SA, Carlson GA, Chang EI, Huang E, Palanker D, Gurtner GC. Comparative healing of surgical incisions created by the PEAK PlasmaBlade, conventional electrosurgery, and a scalpel. *Plast Reconstr Surg*. 2009;124(6):1849-1859.
16. Chang EI, Carlson GA, Vose JG, Huang EJ, Yang GP. Comparative healing of rat fascia following incision with three surgical instruments. *J Surg Res*. 2011;167(1):47-54.
17. Ruidiaz ME, Messmer D, Atmodjo DY, et al. Comparative healing of human cutaneous surgical incisions created by the PEAK PlasmaBlade, conventional electrosurgery, and a standard scalpel. *Plast Reconstr Surg* 2011; Jul 128(1):104-111.
18. Data on file.
19. Vose JG, McAdara-Berkowitz J. Reducing scalpel injuries in the operating room. *AORN J*. 2009; 90(6):867-872.
20. Marulanda GA, Ulrich SD, Seyler TM et al. Reductions in blood loss with a bipolar sealer in total hip arthroplasty. *Expert Rev Med Devices* 2008; 5(2):125-131.
21. Rosenberg AG. Reducing blood loss in total joint surgery with a saline-coupled bipolar sealing technology. *J Arthrop* 2007.
22. Ackerman SJ, Tapia CI, Baik R, Pivec R, Mont MA. Use of a bipolar sealer in total hip arthroplasty: medical resource use and costs using a hospital administrative database. *Orthopedics*. 2014;37(5):e472-481.  
\*Please note: evidence derived from in vitro, in vivo, or clinical studies, and performance may not be specifically established in all procedures.
23. Geller DA, Tsung A, Maheshwari V, Rutstein LA, Fung JJ, Marsh JW. Hepatic resection in 170 patients using saline-cooled radiofrequency coagulation. *HPB (Oxford)*. 2005;7(3):208-213.
24. Data on file.

# Medtronic

## Medtronic

826 Coal Creek Circle  
Louisville, CO 80027  
USA

For further information, please call 866-777-9400.

[medtronicadvancedenergy.com](http://medtronicadvancedenergy.com)

Rx only. Refer to product instruction manual/package insert for instructions, warnings, precautions and contraindications.

© 2022 Medtronic. All rights reserved. Medtronic, Medtronic logo and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. UC202010094a EN