

INTRATHECAL BACLOFEN (ITB) THERAPY VALUE SUMMARY

For treatment of severe,
generalised spasticity



SPASTICITY IS A DISORDERED SENSORI-MOTOR CONTROL, RESULTING FROM AN UPPER MOTOR NEURON LESION, PRESENTING AS INTERMITTENT OR SUSTAINED INVOLUNTARY ACTIVATION OF MUSCLES¹

THE ORIGIN OF SPASTICITY

RECOGNIZING THE SIGNS OF SPASTICITY



Spasticity is often described as tight, stiff muscles or spasms that may make movement, posture and balance difficult.

When spasticity is perceived by the individual or caregiver as hindering body function, activities and/or participation, it is referred to as **disabling spasticity**².

Disabling spasticity is also referred to as severe spasticity by some clinicians/authors.



THE ORIGIN OF SPASTICITY

CEREBRAL ORIGIN

- STROKE
- CEREBRAL PALSY (CP)
- BRAIN INJURY

OR

SPINAL ORIGIN

- SPINAL CORD INJURY (SCI)
- MULTIPLE SCLEROSIS (MS)*

* Multiple sclerosis is routinely considered of spinal origin due to the location of the lesions in the spinal cord, although patients may also present with brain lesions³.

RECOGNIZING THE SIGNS OF SPASTICITY



SEVERAL SCALES
HAVE BEEN
DEVELOPED
AND VALIDATED
**TO ASSESS
SPASTICITY**

While some scales are indication-specific, others are commonly used across indications. The most commonly used scales include:

- (Modified) Ashworth Scale (AS/MAS)
- Functional Independence Measure (FIM)
- Electrophysiological Measures
- Gait and Movement Analysis
- Walking Tests
(e.g., 10 metres, 6 minutes)

REFERENCES

1. Pandyan AD et al. Spasticity: clinical perceptions, neurological realities and meaningful measurement. Disabil Rehabil 2005; 27: 2–6. (SPASM Consortium).
2. Burns AS, Lanig I, Grabljevec K, New PW, Bensmail D, Ertzgaard P, Nene AV. Optimizing the management of disabling spasticity following spinal cord damage The Ability Network an international initiative. Arch Phys Med Rehabil. 2016 Dec;97(12):2222-2228.
3. Saval et al. Intrathecal Baclofen for Spasticity Management: A Comparative Analysis of Spasticity of Spinal vs Cortical Origin. Journal of Spinal Cord Medicine. Volume 33. Number 1. 2010.

EPIDEMIOLOGY

THE PREVALENCE OF SEVERE SPASTICITY VARIES WITH THE UNDERLYING CONDITION

SPASTICITY AFFECTS AN ESTIMATED 12 MILLION PEOPLE WORLDWIDE¹



SPINAL ORIGIN

MULTIPLE SCLEROSIS

550,000 persons in Europe have multiple sclerosis².

Up to **30%** are living with severe spasticity³⁻⁷.



SPINAL CORD INJURY

200,000 persons in Europe are living with spinal cord injury⁸.

Over **30%** are living with severe spasticity⁹⁻¹⁴.



CEREBRAL ORIGIN

STROKE

1.1 million persons suffer a stroke in Europe each year¹⁵.

Up to **13%** are living with a more severe degree of spasticity¹⁶⁻²¹.



CEREBRAL PALSY

Close to 1 million persons in Europe have cerebral palsy²².

Over **40%** are living with severe spasticity²³⁻²⁷.



BRAIN INJURY


Prevalence is not well characterized.



MULTIPLE SCLEROSIS (MS)



Several patient survey-based studies have found that **65-84%** of patients with MS experience spasticity³⁻⁶.



17-29% of patients in these studies report severe spasticity³⁻⁷.

SPASTICITY IS A COMMON COMPLICATION IN MULTIPLE SCLEROSIS AND OFTEN AFFECTS THE LEGS



MICHELLE
ITB Therapy Patient

SPINAL CORD INJURY (SCI)



Across studies, spasticity is reported to develop in **40-68%** of patients with SCI¹⁰⁻¹³.

Patients with cervical injuries are more likely to experience spasticity than patients with thoracic or lumbosacral level injuries¹⁰⁻¹¹.



Nearly **30%** of patients with SCI may experience disabling spasticity^{10-12,14}.

INCIDENCE AND PREVELANCE OF SCI DIFFER BETWEEN DEVELOPING AND DEVELOPED COUNTRIES⁹



IAN
ITB Therapy Patient

STROKE



Several studies have observed that **17-43%** of post-stroke patients develop spasticity¹⁶⁻²¹.



One study reports that **13%** of post-stroke patients develop severe disabling spasticity²⁰⁻²¹.

SPASTICITY IS A WELL-KNOWN DISABLING CONSEQUENCE OF STROKE



PATSY
ITB Therapy Patient

CEREBRAL PALSY (CP)



A comparison of 28 CP registry studies from 8 geographic regions, reports that **90%** of patients with CP have a spastic motor type²⁵.



One study on functional status (using the definition for motor function suggested by the Surveillance of cerebral palsy in Europe, which may be a proxy for severe spasticity) reports that **42.3%** of patients have severe disability caused by the underlying disease²⁶⁻²⁷.

CEREBRAL PALSY IS THE MOST COMMON CAUSE OF SPASTICITY IN CHILDREN AND YOUNG ADULTS²³⁻²⁴



MATT
ITB Therapy Patient

BRAIN INJURY

An estimated 5% of all traumatic brain injuries are classified as severe, with contractures due to spasticity occurring in up to **85%** of patients with severe traumatic brain injuries²⁸.



EPIDEMIOLOGICAL EVIDENCE FOR SPASTICITY DEVELOPMENT AFTER BRAIN INJURY IS NOT WELL-CHARACTERIZED. INCIDENCE ESTIMATES VARY FROM 13-63%²⁹⁻³⁰

REFERENCES

1. American Association of Neurological Surgeons Website. Spasticity. <http://www.aans.org/en/Patient%20Information/Conditions%20and%20Treatments/Spasticity.aspx>. (Accessed March 2019)
2. Atlas of Multiple Sclerosis (online), MS International Federation -2013 update. Available at: <https://www.msif.org/about-us/who-we-are-and-what-we-do/advocacy/atlas/> (accessed on March 2019).
3. Oreja-Guevara C, Gonzalez-Segura D, Vila C. Spasticity in multiple sclerosis: results of a patient survey. *Int. J. Neurosci.* 2013;123(6):400-408.
4. Goodin DS. Survey of multiple sclerosis in Northern California. Northern California MS Study Group. *Mult. Scler.* 1999;5(2):78-88.
5. Barnes MP, Kent RM, Semlyen JK, McMullen KM. Spasticity in multiple sclerosis. *Neurorehabil. Neural Repair.* 2003;17(1):66-70.
6. Rizzo MA, Hadjimichael OC, Preiningerova J, Vollmer TL. Prevalence and treatment of spasticity reported by multiple sclerosis patients. *Mult. Scler.* 2004;10(5):589-595.
7. Flachenecker P, Henze T, Zettl UK. Spasticity in patients with multiple sclerosis--clinical characteristics, treatment and quality of life. *Acta Neurol. Scand.* 2014;129(3):154-162.
8. Antibodies against Nogo-A to enhance plasticity, regeneration and functional recovery after acute spinal cord injury, a multicenter European clinical proof of concept trial, European Commission <https://cordis.europa.eu/project/rcn/198795/factsheet/en> (accessed on March 2019)
9. Singh A, Tetreault L, Kalsi-Ryan S, Nouri A, Fehlings MG. Global prevalence and incidence of traumatic spinal cord injury. *Clin. Epidemiol.* 2014;6:309-331.
10. Skold C, Levi R, Seiger A. Spasticity after traumatic spinal cord injury: nature, severity, and location. *Arch. Phys. Med. Rehabil.* 1999;80(12):1548-1557.
11. Levi R, Hultling C, Seiger A. The Stockholm Spinal Cord Injury Study: 2. Associations between clinical patient characteristics and post-acute medical problems. *Paraplegia.* 1995;33(10):585-594.
12. Maynard FM, Karunas RS, Waring WP, 3rd. Epidemiology of spasticity following traumatic spinal cord injury. *Arch. Phys. Med. Rehabil.* 1990;71(8):566-569.
13. Noreau L, Proulx P, Gagnon L, Drolet M, Laramée MT. Secondary impairments after spinal cord injury: a population-based study. *Am. J. Phys. Med. Rehabil.* 2000;79(6):526-535.
14. Johnson RL, Gerhart KA, McCray J, Menconi JC, Whiteneck GG. Secondary conditions following spinal cord injury in a population-based sample. *Spinal Cord.* 1998;36(1):45-50.
15. Truelsen T, Piechowski-Jo'z'wiak C, Bonitaa R et al. Stroke incidence and prevalence in Europe: a review of available data. *European Journal of Neurology* 2006, 13: 581-598.
16. Sommerfeld DK, Gripenstedt U, Welmer AK. Spasticity after stroke: an overview of prevalence, test instruments, and treatments. *Am. J. Phys. Med. Rehabil.* 2012;91(9):814-820.
17. Lundstrom E, Terent A, Borg J. Prevalence of disabling spasticity 1 year after first-ever stroke. *Eur. J. Neurol.* 2008;15(6):533-539.
18. Welmer AK, von Arbin M, Widen Holmqvist L, Sommerfeld DK. Spasticity and its association with functioning and health-related quality of life 18 months after stroke. *Cerebrovasc. Dis.* 2006;21(4):247-253.
19. Urban PP, Wolf T, Uebele M, et al. Occurrence and clinical predictors of spasticity after ischemic stroke. *Stroke.* 2010;41(9):2016-2020.
20. Watkins CL, Leathley MJ, Gregson JM, Moore AP, Smith TL, Sharma AK. Prevalence of spasticity post stroke. *Clin. Rehabil.* 2002;16(5):515-522.
21. Dvorak EM, Ketchum NC, McGuire JR. The underutilization of intrathecal baclofen in poststroke spasticity. *Top Stroke Rehabil.* 2011;18(3):195-202.
22. Surveillance of Cerebral Palsy in Europe, Prevalence and characteristics of children with cerebral palsy in Europe, *Developmental Medicine & Child Neurology* 2002, 44: 633-640.
23. Albright AL. Spasticity and movement disorders in cerebral palsy. *J. Child Neurol.* 1996;11 Suppl 1(Suppl 1):S1-4.
24. Hutchinson R, Graham HK. Management of spasticity in children. In: Barnes MP, Johnson GR, eds. *Upper Motor Neurone Syndrome and Spasticity: Clinical Management and Neurophysiology*. 2 ed. Cambridge, UK: Cambridge University Press; 2008:214-239.
25. Reid SM, Carlin JB, Reddihough DS. Distribution of motor types in cerebral palsy: how do registry data compare? *Dev. Med. Child Neurol.* 2011;53(3):233-238.
26. Prevalence and characteristics of children with cerebral palsy in Europe. *Dev Med Child Neural.* 2002; Sept; 44:633-40.
27. Kirby RS, Wingate MS, Van Naarden Braun Ket al. Prevalence and functioning of children with CP in 4 areas of the US in 2006: a report from the Autism and Developmental Disabilities Monitoring Network. *Res in Dev Disabilities.* 2011 ;32(2):462-469.
28. Bose, P, Hou J, and Thompson FJ. Brain Neurotrauma: Molecular, Neuropsychological, and Rehabilitation Aspects. Chapter 14 Traumatic Brain Injury (TBI)-Induced Spasticity Neurobiology, Treatment, and Rehabilitation (2015) Taylor & Francis Group, LLC. PMID: 26269896.
29. Singer BJ, Jegasothy GM, Singer KP, Allison GT, Dunne JW. Incidence of ankle contracture after moderate to severe acquired brain injury. *Arch Phys Med Rehab* 2004 Sept;85(9):1465-9.
30. Williams G, Banky M, Olver J. Distribution of Lower Limb Spasticity Does Not Influence Mobility Outcome Following Traumatic Brain Injury: An Observational Study *J Head Trauma Rehabil* Vol. 30, No. 5, pp. E49-E57.

Home

What is
Spasticity?

The Burden
of Spasticity

Epidemiology

Humanistic Burden

Economic Impact

Treatment
Options

Intrathecal
Baclofen
Therapy

References

Glossary

SIGNIFICANT HUMAN IMPACT



SPASTICITY IS
A BURDENSOME
CONDITION AND
HAS BEEN SHOWN
TO IMPACT ON ALL
ASPECTS OF PATIENTS'
**HEALTH-RELATED
QUALITY OF LIFE,
ACROSS DISEASE
INDICATIONS.**

**MULTIPLE
SCLEROSIS**

**SPINAL CORD
INJURY**

STROKE

**IMPACT OF
SPASTICITY
SEVERITY**


MULTIPLE SCLEROSIS (MS)

OTHER IMPACTS OF SPASTICITY



In a 2013 Spanish study, scores on the Short Form-12 (SF-12®) quality of life survey indicated that **MS-related spasticity** had a **more pronounced effect on patients' general health and on the physical**, rather than mental, aspects of well-being¹.

MS patients with spasticity also demonstrate a quality of life decrement compared to a general population^{1,2}.




HRQoL: HEALTH-RELATED QUALITY OF LIFE

* In the absence of normative SF-12 data, normative SF-36 data from a Spanish population is used here as an approximate comparison³

ITB VALUE SUMMARY **THE BURDEN OF SPASTICITY**

OTHER IMPACTS OF SPASTICITY

A 2014 cross-sectional study from Germany, which focused on patients who had experienced multiple sclerosis-related spasticity for at least 12 months, reported general agreement in the **most disturbing symptoms associated with spasticity** according to both physicians (N=411) and patients (N=398), being **stiffness** and **mobility restriction**⁴.




SPINAL CORD INJURY (SCI)

OTHER IMPACTS OF SPASTICITY




From a 2014 UK survey of 162 persons with SCI who experienced spasticity, significant **decreases in health-related quality of life**, as measured by the SF-36 survey (version 2.0), were observed **in all health domains, when compared to a general population^{5,6}**.



OTHER IMPACTS OF SPASTICITY

A 2017 multi-centre study from the Netherlands examining community-dwelling SCI patients who used a wheelchair for everyday mobility (N=110), found that **more than half** of participants **experienced restrictions on functioning due to spasticity, persisting at 5 years post-discharge** from their first inpatient rehabilitation.⁷

Functional restrictions due to spasticity⁷



STROKE



IMPACT OF POST-STROKE SPASTICITY ON HRQOL



IMPACT OF POST-STROKE SPASTICITY ON CAREGIVERS



Post-stroke patients with spasticity suffer from a **reduction in function of 35–38%** compared to those without spasticity ($p \leq 0.001$)⁸.





BI: BARTHEL INDEX

RMI: RIVERMEAD MOBILITY INDEX

ITB VALUE SUMMARY **THE BURDEN OF SPASTICITY**

STROKE

Spasticity **negatively impacts** the health-related quality of life (HRQoL) of stroke survivors compared to those without spasticity ($p \leq 0.001$)⁹.



EQ-5D: EUROQOL-5D


SF-12 PCS: SHORT FORM-12 PHYSICAL COMPONENT SUMMARY

SS-QOL: STROKE-SPECIFIC QUALITY OF LIFE

STROKE


Caregivers of stroke survivors with spasticity suffer from **increased anxiety or depression**¹⁰.

Physician diagnosis



- None (86 caregivers)
- Depression (34 caregivers)
- Anxiety (33 caregivers)


Patient Health Questionnaire-9 (PHQ-9)
categories of depression



- None (68 caregivers)
- Mild depression (30 caregivers)
- Moderate depression (28 caregivers)
- Moderate severe depression (16 caregivers)
- Severe depression (11 caregivers)

QUALITY OF LIFE WORSENS WITH INCREASING SEVERITY OF SPASTICITY

Data from a German cross-sectional, burden-of-disease study on patients who had experienced multiple sclerosis-related spasticity for at least 12 months reported **increasing impact on health-related quality of life with increasing spasticity severity** (N=389)¹¹.



EQ-5D: EUROQOL-5D

HEALTH-RELATED
QUALITY OF LIFE
DECREASE BECOMES
MORE PROMINENT
**WITH AN INCREASING
SEVERITY OF SPASTICITY¹¹**

REFERENCES

1. Arroyo R, Massana M, Vila C. Correlation between spasticity and quality of life in patients with multiple sclerosis: the CANDLE study. *Int J Neurosci*. 2013 Dec;123(12):850-8.
2. Alonso J, Regidor E, Barrio G, Prieto L, Rodríguez C, de la Fuente L. Population reference values of the Spanish version of the Health Questionnaire SF-36. *Med Clin (Barc)*. 1998 Oct 10;111(11):410-6.
3. The SF-12 was derived from, and is a shorter version of, the SF-36 survey. Both measure eight health domains: physical functioning, role physical, bodily pain, general health, vitality, social functioning, role emotional, and mental health, score on a 0–100 scale (with higher scores indicating better health), have generally comparable psychometric properties and are used to derive two aggregate summary measures: the PCS and the MCS. The SF-12 is a suitable measure where information on the SF-36 PCS and MCS scores is required. Compared to the SF-36, the SF-12 has been less well-studied, thus findings related to the SF-36 may not be transferable to the SF-12. Refer to: Busija L, Pausenberger E, Haines TP, Haymes S, Buchbinder R, Osborne RH. Adult measures of general health and health-related quality of life: Medical Outcomes Study Short Form 36-Item (SF-36) and Short Form 12-Item (SF-12) Health Surveys, Nottingham HealthProfile (NHP), Sickness Impact Profile (SIP), Medical Outcomes Study Short Form 6D (SF-6D), Health Utilities Index Mark 3 (HUI3), Quality of Well-Being Scale (QWB), and assessment of Quality of Life (AQoL). *Arthritis Care Res (Hoboken)*. 2011 Nov;63 Suppl 11:S383–412.
4. Flachenecker P, Henze T, Zettl UK. Spasticity in patients with multiple sclerosis – clinical characteristics, treatment and quality of life. *Acta Neurol Scand*, 2014; 129: 154–162.
5. Thavaneswaran S, Grifi M, Stauble F, Wade AG, Crawford G. Capturing the SCI Patient Voice: A Survey of the burden of spasticity in SCI patients in the UK. The 4th ISCOS and ASIA Joint Scientific Meeting (Montreal, 2015) and Medtronic data on file.
6. Jenkinson C, Stewart-Brown S, Petersen S, Paice C. assessment of the SF-36 version 2 in the United Kingdom. *J Epidemiol Community Health*. 1999 Jan;53(1):46–50.
7. van der Meer P, Post MW, van Leeuwen CM, van Kuppevelt HJ, Smit CA, van Asbeck FW. Impact of health problems secondary to SCI one and five years after first inpatient rehabilitation. *Spinal Cord*. 2017 Jan;55(1):98–104. doi: 10.1038/sc.2016.103. Epub 2016 Jul 19.
8. Welmer et al. Spasticity and Its association with functioning and health-related quality of life 18 months after stroke. *Cerebrovasc Dis* 2006;21:247–53.
9. Gillard PJ, Sucharew H, Kleindorfer D, et al. The negative impact of spasticity on the health-related quality of life of stroke survivors: a longitudinal cohort study. *Health Qual Life Outcomes*. 2015;13:159. doi:10.1186/s12955-015-0340-3.
10. Denno MS, Gillard PJ, Graham GD, et al. Anxiety and depression associated with caregiver burden in caregivers of stroke survivors with spasticity. *Arch Phys Med Rehabil*. 2013;94(9):1731–6. doi:10.1016/j.apmr.2013.03.014.
11. Zettl UK, Henze T, Essner U, Flachenecker P. Burden of disease in multiple sclerosis patients with spasticity in Germany: mobility improvement study (Move I). *Eur J Health Econ*, 2014; 15:953–966.
12. A. Szende et al. (eds.), Self-Reported Population Health: An International Perspective based on EQ-5D. 2014. DOI 10.1007/978-94-007-7596-1_3. Page 30, Table 3.6 EQ-5D index population norms (country-specific TTO value sets) – Germany value of 0.938.


HIGHER HEALTHCARE COSTS ASSOCIATED WITH PATIENTS WHO DEVELOP SPASTICITY

INCREASING
SEVERITY,
INCREASES
COSTS



A Swedish study has shown that in the 12-months **post-stroke**, direct healthcare costs for patients who experience spasticity are **4 times higher** than for those who do not develop spasticity¹.

Mean Direct Healthcare Costs of patients with post-stroke spasticity was PPP \$84,195 vs. PPP \$21,842 for non-spastic stroke patients¹.




*For this study, Swedish currency (SEK) was converted into Purchasing Power Parities, which are rates of currency conversion that eliminate the difference in price levels between countries². US dollar (PPP\$) at the exchange rate of 9.34 SEK for 1 PPP\$ (2003 value)³.

¥ Other direct costs include primary care and medication.

HOSPITALISATIONS,
HOME HELP AND
RESIDENTIAL CARE
WERE THE MAIN
COST DRIVERS
IN THE
POST-STROKE
POPULATION, WITH
HIGHER COSTS
OBSERVED FOR
SPASTIC PATIENTS
COMPARED TO
NON-SPASTIC
PATIENTS¹

HIGHER HEALTHCARE COSTS ARE ASSOCIATED WITH AN INCREASING DEGREE OF SPASTICITY

Data from a German burden-of-disease study on patients with multiple sclerosis (MS)-related spasticity demonstrates **increasing healthcare costs accompany increasing spasticity severity** (N=379)⁴.



⁴2011 costs from Germany

REFERENCES


1. Lundström E, Smits A, Borg J, Terént A. Four-fold increase in direct costs of stroke survivors with spasticity compared with stroke survivors without spasticity: the first year after the event. *Stroke*, 2010;41(2):319-24.
2. Vachris MA, Thomas J International price comparisons based on purchasing power parity. *Monthly Labor Review Online*. 1999;122:3-12.
3. OECD purchasing power parities for GCO, 2003. Available at: http://stats.oecd.org/wbos/Index.aspx?datasetcodeSNA_TABLE4. Accessed July 12, 2009.
4. Zettl UK, Henze T, Essner U, Flachenecker P. Burden of disease in multiple sclerosis patients with spasticity in Germany: mobility improvement study (Move I). *Eur J Health Econ*. 2014; 15:953-966.



TREATMENT OPTIONS

DISABLING SPASTICITY IS CHALLENGING TO CONTROL EFFECTIVELY AND TYPICALLY REQUIRES A COMBINATION OF TREATMENTS

CLICK ON THE ICONS FOR MORE INFORMATION ON EACH TREATMENT OPTION



Treatment options are not mutually exclusive.
They may be used in conjunction with each other.

ORAL MEDICATIONS




ORAL BACLOFEN IS ONE OF THE MOST COMMONLY USED ANTISPASMODIC DRUGS IN THE TREATMENT OF SPASTICITY

For some patients, the amount of drug that penetrates the blood-brain barrier may not provide adequate relief for their disabling spasticity.

Increasing the dose may often result in unacceptable central nervous system (CNS) side effects such as sedation, dizziness, and muscle weakness¹⁻⁴.

Systemic delivery through the circulatory system.




INTRATHECAL BACLOFEN (ITB) THERAPY

AN ALTERNATIVE ROUTE OF DELIVERY

ITB THERAPY DELIVERS THE DRUG DIRECTLY TO THE INTRATHECAL SPACE SURROUNDING THE SPINAL CORD

This allows effective drug concentrations to reach the cerebrospinal fluid, while the plasma remains relatively unaffected, compared with oral administration⁵.

Targeted delivery directly to the cerebral spinal fluid.



Medtronic
Further, Together

REHABILITATION THERAPY



PHYSICAL THERAPY to maximize motor skills, train muscles, improve mobility and promote flexibility.

OCCUPATIONAL THERAPY for upper limb function and fine motor skills to improve self-care and comfort.

SPEECH THERAPY to improve speech or use of alternative means of communication.

INJECTABLE TREATMENTS



PHENOL AND ALCOHOL are neurolytic agents which induce chemodenervation.

BOTULINUM TOXIN (BoNT) agents, administered by intramuscular injection, are available for the treatment of focal spasticity.

Injections must be repeated regularly as their effectiveness wears off over time.

SURGICAL PROCEDURES



General or focal, non-reversible procedures targeting:

NERVES (e.g., selected dorsal rhizotomy, neurotomy).

MUSCLES AND TENDONS (e.g., lengthening transfer).

BONES (eg. osteotomies, fusions).

REFERENCES

1. Hattab JR. Review of European clinical trials with baclofen. In: Feldman RG, Young RR, Koella WP, eds. Spasticity: Disordered Motor Control. Chicago, IL: Year Book; 1980:71-85.
2. Jagatsinh Y. Intrathecal baclofen: Its effect on symptoms and activities of daily living in severe spasticity due to spinal cord injuries: A pilot study. Indian J. Orthop. 2009;43(1):46-49.
3. Lazorthes Y, Sallerin-Caute B, Verdie JC, Bastide R, Carillo JP. Chronic intrathecal baclofen administration for control of severe spasticity. J. Neurosurg. 1990;72(3):393-402.
4. Montane E, Vallano A, Laporte JR. Oral antispastic drugs in nonprogressive neurologic diseases: a systematic review. Neurology. 2004;63(8):1357-1363.
5. Muller H, Zierski J, Dralle D, Krauss D, Mutschler E. Pharmacokinetics of intrathecal baclofen. Localspinal therapy of spasticity. 1988:223-226.

INTRATHECAL BACLOFEN (ITB) THERAPY

94%

Up to 94% of patients would agree to pump implantation again⁵.

88%

More than 88% of caregivers would recommend ITB infusion to others⁶.

ITB is a viable option in patients who experience intolerable side effects or who fail to respond to the maximum recommended dose of oral baclofen¹.



Delivers a liquid form of medicine directly to fluid around spinal cord



Requires **100 to 1000x** smaller doses than oral baclofen²



Expected to produce fewer or more tolerable side effects when compared to oral baclofen³⁻⁴





REFERENCES

1. Ertzgaard et al. Efficacy and safety of oral baclofen in the management of spasticity: a rationale for intrathecal baclofen. Submitted article.
2. Barnes MP, Johnson GR. Upper motor neurone syndrome and spasticity: Clinical management and neurophysiology: Cambridge University Press, 2008.
3. Penn RD. Intrathecal baclofen for spasticity of spinal origin: seven years of experience. J Neurosurg. 1992;77(2):236-240.
4. Meythaler JM, Guin-Renfroe S, Law C, et al. Continuously infused intrathecal baclofen over 12 months for spastic hypertonia in adolescents adults with cerebral palsy. Arch Phys Med Rehabil 2001;82(2):155-161.
5. Sadiq SA, Wang GC. Long-term intrathecal baclofen therapy in ambulatory patients with spasticity J Neurol (2006) 253: 563–569.
6. Campbell WM, Ferrel A, McLaughlin JF et al. Long-term safety and efficacy of continuous intrathecal baclofen. Dev Med Child Neurol. 2002;44(10):660-665.

PHARMACOLOGY

**ITB ONLY REQUIRES
1/100 OF THE ORAL
DAILY DOSAGE AND
RESULTS IN A 50X
HIGHER SPINAL
CONCENTRATION¹⁻⁴**

Compared to oral baclofen, much lower (1/100) doses of ITB are required to achieve a relatively high (50x) concentration in the cerebrospinal fluid (CSF), thus minimizing side effects.



REFERENCES

1. Knutsson E, Lindblom U, Martensson A. Plasma and Cerebrospinal Fluid Levels of Baclofen (Lioresal®) at Optimal Therapeutic Responses in Spastic Paresis. Journal of the neurological Sciences, 1974, 23: 473-484.

2. Kroin JS, Penn RD. Cerebrospinal fluid pharmacokinetics of lumbar intrathecal baclofen. Lakke JPWF (ed) Parenteral drug therapy is spasticity and Parkinson's disease. Parthenon Pub Group. Park Ridge, NJ. 1991. pp 67-77.

3. Muller H, Zierski J, Dralle D, Krauss D, Mutschler E. Pharmacokinetics of intrathecal baclofen. Local-spinal therapy of spasticity. 1988:223-226.

4. Penn RD. Intrathecal Baclofen for Severe Spasticity. Annals of the New York Academy of Sciences. 1988. Volume 531.

BENEFITTING PATIENTS AND CLINICIANS FOR MORE THAN 30 YEARS

Medtronic's global reach enables our intrathecal drug delivery therapy to be available in many countries across the world.



* Design changes to improve performance and reliability; data on file at Medtronic.
Note: This timeline only shows extracts of the Medtronic TDD product portfolio

PERFORMANCE AND RELIABILITY

The SynchroMed™ II Programmable Infusion System:

- is available with a 20ml and 40ml reservoir size
- precisely moves medication through a peristaltic action
- has constant and flexible infusion modes with bolus dosing
- allows patient to self administer a therapeutic bolus
- makes safe 1.5T and 3.0 full body MRI possible*
- is designed to resume therapy immediately after MRI scan

* Under specific conditions for 1.5T and 3.0T MRI scans. Requires interrogation to confirm pump status.

ITB VALUE SUMMARY **INTRATHECAL BACLOFEN THERAPY**

94.2%

Pump survival after 6.5 years¹

Confidence year after year

±0.3%

Repeatability²

Consistency day after day



19.5mm thickness (20ml pump)
26mm thickness (40ml pump)

Medtronic
Further, Together



PERFORMANCE ASCENDA™

This is why Medtronic designed the Ascenda™ catheter³:

Real data shows only
28% of therapy complications are
system related.


Of the 28%,
75.7% are
catheter related.

Of the 75.7%, **82%** are due to:

- Dislodgement/ Migration
 - Fracture/ Break
 - Kink/ Occlusion

DESIGN ENHANCEMENTS⁴

4 layers vs. 1 (silicone catheters)



Designed to reduce catheter dislodgement through more secure anchoring & better stabilization during placement.



PERFORMANCE ASCENDA™

CLINICAL DATA DEMONSTRATES FAVOURABLE RESULTS FOR ASCENDA™ COMPARED TO OTHER CATHETERS^{5,6}

RESULTS (Motta et al 2016)		OTHER CATHETERS	ASCENDA™
Number of ITB Patients	# of Patients	416	92
Catheter related complications	# of Patients (%)	75 (18%)	0 (0%)

RESULTS (Pucks-Faes et al 2018)		OTHER CATHETERS	ASCENDA™
Number of ITB Patients	# of Patients	73	36
Catheter related complications	# of Patients (%)	19 (26%)	3 (8%)

REFERENCES

1. Medtronic Product Performance Report 2017 http://professional.medtronic.com/ppr/wcm/groups/mdtcom_sg/mdt/neuro/documents/documents/mdt_product_performance_2017.pdf
2. Medtronic data on file November 2015. Engineering bench study that included > 1,000 data points from 23 pumps. Typical infusion rates of 100, 300, 500 and 700 mcL/day were included under the following conditions: 37° C, full, mid- and low-reservoir volumes, 300 meters above sea level. "Repeatability" is defined as the extent to which volumes delivered by the SynchroMed II Infusion System varied from refill to refill.
3. Medtronic Product Performance Report 2016 http://professional.medtronic.com/ppr/wcm/groups/mdtcom_sg/mdt/neuro/documents/documents/mdt_product_performance_2016.pdf. Accessed February 2019.
4. Medtronic Technical Design Summary: Ascenda™ Models 8780 and 8781 intrathecal catheter, Medtronic Europe 2010.
5. Motta, F. and Antonello C.E. Comparison between an Ascenda and a silicone catheter in intrathecal baclofen therapy in pediatric patients; analysis of complications. *Journal of Neurosurgery: Pediatrics*, 18(4), 493-498.
6. Pucks-Faes E, Hitznerberger G, Matzak H, et al. Eleven years' experience with Intrathecal Baclofen – Complications, risk factors. *Brain Behav*. 2018 Mar 30;8(5): e00965.

Home

What is
Spasticity?

The Burden
of Spasticity

Treatment
Options

Intrathecal
Baclofen
Therapy

Pharmacology

Technology
Overview

Therapy
Continuum

Clinical
Evidence

Safety

Patient Benefit

Economic Value

Guidelines &
Recommendations

References

Glossary

AVAILABLE IN MULTIPLE REGIONS ACROSS THE GLOBE

PLEASE CONTACT
A MEDTRONIC
REPRESENTATIVE
TO FIND OUT MORE
ABOUT IMPLANT
AND REFILL CENTRES
WORLDWIDE



ITB VALUE SUMMARY **INTRATHECAL BACLOFEN THERAPY**

Medtronic
Further, Together

A STAGED APPROACH TO SUCCESSFUL THERAPY

Patient assessment is an essential step in selecting ITB Therapy candidates.


A thorough assessment should ensure that all facets of the patient's condition are considered.

The goals of treatment vary for each patient and must be realistic given their current level of function and disability.

Patients, caregivers and physicians should work together to develop and agree upon all objectives that may be possible with ITB Therapy.

THERAPY ALWAYS STARTS WITH A SCREENING TEST (TEST DOSE)

A screening test determines if ITB may work for the patient. It involves the administration by lumbar puncture of an ITB bolus injection. The patient's response is monitored for an 8 hour period following the procedure.



If the medical team, together with the patient, decide ITB is the right treatment option, a Programmable Infusion System (pump and catheter) will be implanted in the patient's lower abdomen for long-term management of their spasticity.

Home

What is
Spasticity?

The Burden
of Spasticity

Treatment
Options

Intrathecal
Baclofen
Therapy

Pharmacology

Technology
Overview

Therapy
Continuum

Clinical
Evidence

Safety

Patient Benefit

Economic Value

Guidelines &
Recommendations

References

Glossary

PATIENT
SELECTION

GOAL
SETTING

SCREENING
TEST

IMPLANT



THERAPY
MAINTENANCE

A 5-STEPS IMPLANTATION TECHNIQUE


1. POSITION PATIENT




2. PLACE CATHETER



3. PREPARE PUMP POCKET



4. TUNNEL CATHETER



5. CONNECT PUMP AND CATHETER



INTERDISCIPLINARY THERAPY MANAGEMENT IS ESSENTIAL TO ENSURE OPTIMAL TREATMENT BENEFITS

TITRATION

Schedules are individually set to benefit the patient, depending on activities and needs.

REPLACEMENT

Every 6-7 years¹, an alarm is activated 90 days before battery life ends.

REFILLS

Every 6 weeks-6 months depending on dosing requirements and reservoir size.

TROUBLE SHOOTING

Patients should be monitored for the occurrence of complications.

Home

What is
Spasticity?

The Burden
of Spasticity

Treatment
Options

**Intrathecal
Baclofen
Therapy**

Pharmacology

Technology
Overview

**Therapy
Continuum**

Clinical
Evidence

Safety

Patient Benefit

Economic Value

Guidelines &
Recommendations

References

Glossary

REFERENCE

1. Medtronic Product Performance Report 2017 http://professional.medtronic.com/ppr/wcm/groups/mdtcom_sg/mdt/neuro/documents/documents/mdt_product_performance_2017.pdf

ITB IS SUPPORTED BY CLINICAL EVIDENCE ACROSS INDICATIONS

ITB can help manage spasticity, decrease spasms and pain, as well as improve function and quality of life.

	SPASTICITY OF CEREBRAL ORIGIN			SPASTICITY OF SPINAL ORIGIN	
	CEREBRAL PALSY	BRAIN INJURY	STROKE	SPINAL CORD INJURY	MULTIPLE SCLEROSIS
NO. OF RANDOMISED PATIENTS	Up to 54	Up to 18	Up to 74	Up to 33	Up to 106
LENGTH OF FOLLOW-UP	Up to 9 years	Up to 9 years	Up to 12 months	Up to 12 years	Up to 21 years
MUSCLE TONE	7 Studies ^{1-3,5-9*}	4 Studies ⁹⁻¹²	5 Studies ^{13-16,18}	5 Studies ^{8,19-22}	8 Studies ^{8,20,21,23-27}
SPASMS	2 Studies ^{7,8}	3 Studies ¹⁰⁻¹²		5 Studies ^{8,19-22}	8 Studies ^{8,20,21,23-27}
PAIN	2 Studies ^{2,4-6*}		1 Study ¹⁷		4 Studies ^{23,25-27}
FUNCTION	4 Studies ^{2,3,5,6,8*}	1 Study ¹⁰	3 Studies ¹³⁻¹⁵	2 Studies ^{8,19}	2 Studies ^{8,25}
QUALITY OF LIFE	1 Study ^{4-6*}		3 Studies ^{13,14,17}		3 Studies ^{21,23,25}

Clinical studies including at least 10 patients and showing a significant improvement of the outcome measures.

*Publications 4,5 and 6 refer to the same clinical trial.

CEREBRAL PALSY: 1 retrospective analysis³, 4 prospective studies^{1,2,8,9}, 1 randomized controlled study⁶ followed by 2 prospective studies on the same patient group^{4,5} and 1 prospective controlled study⁷. In the 2 controlled studies, the outcomes have been compared to the conventional medical management.

BRAIN INJURY: 4 prospective studies^{9,10,11,12}.

STROKE: 3 prospective studies^{13,14,18}, 1 case series¹⁵ and 1 randomized controlled study^{16,17}.

SPINAL CORD INJURY: 4 prospective studies^{8,19-21} and 1 retrospective analysis²².

MULTIPLE SCLEROSIS: 2 retrospective study^{24,26}, 5 prospective studies^{8,20,21,25,27} and 1 randomized control study for the first 13 weeks and prospective study for the follow-up at 26 and 52 weeks²³.



REFERENCES

1. Albright AL, Gilmartin R, Swift D et al. Long-term intrathecal baclofen therapy for severe spasticity of cerebral origin. *J Neurosurg.* 2003;98(2):291-295
2. Ramstad K, Jahnson R, Lofterod B, Skjeldal OH. Continuous intrathecal baclofen therapy in children with cerebral palsy – when does improvement emerge?. *Acta Paediatr.* 2010;99(11):1661–1665
3. Motta F, Antonello CE, Stignani C. Intrathecal baclofen and motor function in cerebral palsy. *Dev Med Child Neurol.* 2011;53(5):443–448
4. Vles GF, Soudant DL, Hoving MA et al. Long-term follow-up on continuous intrathecal baclofen therapy in non-ambulant children with intractable spastic cerebral palsy. *Eur J Paediatr Neurol.* 2013;17(6):639–644
5. Hoving MA, van Raak EP, Spincemaille GH et al; Dutch Study Group on Child Spasticity. Safety and one-year efficacy of intrathecal baclofen therapy in children with intractable spastic cerebral palsy. *Eur J Paediatr Neurol.* 2009;13(3):247–256
6. Hoving MA, van Raak EP, Spincemaille GH, Palmans LJ, Becher JG, Vles JS; Dutch Study Group on Child Spasticity. Efficacy of intrathecal baclofen therapy in children with intractable spastic cerebral palsy: a randomised controlled trial. *Eur J Paediatr Neurol.* 2009;13(3):240–246
7. Morton RE, Gray N, Vloeberghs M. Controlled study of the effects of continuous intrathecal baclofen infusion in non-ambulant children with cerebral palsy. *Dev Med Child Neurol.* 2011;53(8):736–741
8. Guillaume D, Van Havenbergh A, Vloeberghs M, Vidal J, Roeste G. A clinical study of intrathecal baclofen using a programmable pump for intractable spasticity. *Arch Phys Med Rehabil.* 2005;86(11):2165–71
9. Yoon YK, Lee KC, Cho HE et al. Outcomes of intrathecal baclofen therapy in patients with cerebral palsy and acquired brain injury. *H. Medicine (Baltimore).* 2017 Aug;96(34):e7472
10. Rawicki B. Continuous intrathecal baclofen delivered via an implantable pump: long-term follow-up review of 18 patients. *J Neurosurg.* 1999;91(5):733–736
11. Meythaler JM, McCary A, Hadley MN. Intrathecal infusion of baclofen for spasticity caused by acquired brain injury: a preliminary report. *J Neurosurg.* 1997;87(3):415–419.
12. Meythaler JM, Guin-Renfroe S, Grabb P, Hadley MN. long-term continuously infused intrathecal baclofen for spastic-dystonic hypertonia in traumatic brain injury: 1-year experience. *Arch Phys Med Rehabil.* 1999;80(1):13–19
13. Schiess MC, Oh IJ, Stimming EF et al. Prospective 12-month study of intrathecal baclofen therapy for poststroke spastic upper and lower extremity motor control and functional improvement. *Neuromodulation.* 2011;14(1):38–45
14. Ivanhoe CB, Francisco GE, McGuire JR et al. Intrathecal baclofen management of poststroke spastic hypertonia: implications for function and quality of life. *Arch Phys Med Rehabil.* 2006;87(11):1509–1515
15. Francisco GE, Boake C. Spastic hemiplegia after intrathecal baclofen therapy: a preliminary study. *Arch Phys Med Rehabil.* 2003;84(8):1194–1199.
16. Creamer M, Cloud G, Kossmehl P, et al. Intrathecal baclofen therapy versus conventional medical management for severe post-stroke spasticity: Results from a multicentre, randomised, controlled, open-label trial (SISTERS). *J Neurol Neurosurg Psychiatry.* 2018 Jun;89(6):642–650
17. Creamer et al. Effect of intrathecal baclofen on pain and quality of life in poststroke spasticity: A randomized trial (SISTERS). *Stroke.* 2018;49:2129–2137
18. Meythaler JM, Guin-Renfroe S, Brunner RC, Hadley MN. Intrathecal baclofen for spastic hypertonia from stroke. *Stroke.* 2001;32(9):2009–2109
19. Azouvi P, Mane, M, Thiebaut, J B, Denys, P, Remy-Neris, O, Bussel, B. Intrathecal baclofen administration for control of severe spinal spasticity: functional improvement and long-term follow-up. *Arch Phys Med Rehabil.* 1996;77(1):35–9
20. Ordia JI, Fischer, Edward, Adamski, Ellen, Chagnon, Kimberly G, Spatz, Edward L. Continuous intrathecal baclofen infusion by a programmable pump in 131 consecutive patients with severe spasticity of spinal origin. *Neuromodulation.* 2002;5(1):16–24
21. Delhaas EM. Long-term outcomes of continuous intrathecal baclofen infusion for treatment of spasticity: A prospective multicenter follow-up study. *Neuromodulation.* 2008;11(3):227–236
22. Burns AS, Meythaler JM. Intrathecal baclofen in tetraplegia of spinal origin: efficacy for upper extremity hypertonia. *Spinal Cord.* 2001;39(8):413–9.
23. Middel B, Kuipers-Upmeyer H, Bourma J, et al. Effect of intrathecal baclofen delivered by an implanted programmable pump on health related quality of life in patients with severe spasticity. *J Neurol Neurosurg Psychiatry.* 1997;63(2):204–9
24. Khan AA, Birks-Agnew I, Bullock P, Rushton D. Clinical outcome and complications of intrathecal baclofen pump in multiple sclerosis patients: a retrospective study. *NeuroRehabilitation.* 2010;27(2):117–20
25. Natale M, D'Oria S, Nero VV, Squillante E, Gentile M, Rotondo M. Long-term effects of intrathecal baclofen in multiple sclerosis. *Clin Neurol Neurosurg.* 2016;143:121–5
26. Lee BS, Jones J, Lang M et al. Early outcomes after intrathecal baclofen therapy in ambulatory patients with multiple sclerosis. *J Neurosurg.* 2018 Oct;129(4):1056–1062
27. Sammariae Y, Yardley M, Keenan L et al. Intrathecal baclofen for multiple sclerosis related spasticity: A twenty year experience. *Mult Scler Relat Disord.* 2018 Oct 15;27:95–100

ITB HAS A DEMONSTRATED SAFETY PROFILE ACROSS INDICATIONS

Adverse effects of different origins may occur with ITB¹⁻⁸.



DRUG

The most common side effects of ITB are related to the drug itself.

E.g., drowsiness, somnolence, nausea, vomiting, muscle weakness, urinary retention.



DEVICE/PROCEDURE

E.g., infection, catheter migration or disconnection, pump dysfunction, cerebrospinal fluid leak, spinal headache.

Almost half of complications related to ITB occur within the first year post surgical implantation, with 25% occurring within the first month.⁹

Most system-related complications are linked to the catheter¹¹. Use of the [Ascenda™ catheter](#) is associated with a more than 3-fold decrease in the risk of catheter complication^{9,10}.

REFERENCES

1. Hsieh JC, Penn RD. Intrathecal baclofen in the treatment of adult spasticity. *Neurosurg Focus* 2006;21(2).
2. Plassat R, Perrouin Verbe B, Menei P, et al. Treatment of spasticity with intrathecal Baclofen administration: long-term follow-up, review of 40 patients. *Spinal Cord* 2004;42(12):686-693.
3. Borriani L, Bensmail D, Thiebaut JB, et al. Occurrence of adverse events in chronic intrathecal baclofen infusion: a one-year follow-up study of 158 adults. *Arch Phys Med Rehabil* 2014;95(6):1032-1038.
4. Hoving MA, van Raak EP, Spincemaille GH, et al. Safety and one-year efficacy of intrathecal baclofen therapy in children with intractable spastic cerebral palsy. *Eur J Paediatr Neurol* 2009;13(3):247-256.
5. Ivanhoe CB, Francisco GE, McGuire JR, et al. Intrathecal baclofen management of poststroke spastic hypertonia: implications for function and quality of life. *Arch Phys Med Rehabil* 2006;87(11):1509-1511.
6. Natale M, Mirone G, Rotondo M, et al. Intrathecal baclofen therapy for severe spasticity: analysis on a series of 112 consecutive patients and future prospectives. *Clin Neurol Neurosurg* 2012;114(4):321-325.
7. Stempien L, Tsai T. Intrathecal baclofen pump use for spasticity: a clinical survey. *Am J Phys Med Rehabil* 2000;79(6):536-541.
8. Creamer M, Cloud G, Kossmehl P, et al. Intrathecal baclofen therapy versus conventional medical management for severe post-stroke spasticity: Results from a multicentre, randomised, controlled, open-label trial (SISTERS). *J Neurol Neurosurg Psychiatry*. 2018 Jun;89(6):642-650.
9. Pucks-Faes E, Hitztenberger G, Matzak H, et al. Eleven years' experience with Intrathecal Baclofen – Complications, risk factors. *Brain Behav*. 2018 Mar 30;8(5): e00965.
10. Motta, F. and Antonello C.E. Comparison between an Ascenda and a silicone catheter in intrathecal baclofen therapy in pediatric patients: analysis of complications. *Journal of Neurosurgery: Pediatrics*, 2016 Oct, 18(4), 493-498.
11. Medtronic Product Performance Report 2017 http://professional.medtronic.com/ppr/wcm/groups/mdtcom_sg/mdt/neuro/documents/documents/mdtproduct_performance_2017.pdf Accessed February 2019.

[Home](#)[What is Spasticity?](#)[The Burden of Spasticity](#)[Treatment Options](#)[Intrathecal Baclofen Therapy](#)[Pharmacology](#)[Technology Overview](#)[Therapy Continuum](#)[Clinical Evidence](#)[Safety](#)[Patient Benefit](#)[Improving QOL](#)[Improving Function](#)[Improving Pain](#)[Economic Value](#)[Guidelines & Recommendations](#)[References](#)[Glossary](#)

PATIENT BENEFIT

7 ITB STUDIES SHOW A SIGNIFICANT IMPROVEMENT IN QUALITY OF LIFE FOR PATIENTS OR CAREGIVERS

CLICK ON EACH AUTHOR NAME FOR MORE INFORMATION ON THE STUDY

	INDICATION	AUTHOR	OUTCOME MEASURE	DOMAINS
CEREBRAL	Stroke	Creamer 2018	EQ-5D	Utility score
	CP	Vles 2013	CHQ-PF50	Pain/discomfort, mental health, parental impact emotional, parental impact time, physical summary
	Stroke	Schiess 2011	SS-QOL	Family role, mobility, social roles, thinking, upper extremity, function work productivity, personality, self-care
	Stroke	Ivanhoe 2006	SIP	Physical dimension, psychological dimension, home management, recreation/pastimes
SPINAL	MS	Natale 2016	SDS	Overall score
	MS and SCI	Delhaas 2008	EQ-5D	Visual analog scale
	MS and SCI	Middel 1997	SIP	Overall score, sleep and rest, recreation and pastimes, mobility, body care and movement, physical dimension
			HSC	Overall score, physical health

CHQ-PF50: Child Health Questionnaire Parent Form 50

CP: Cerebral Palsy

EQ-5D: EuroQol-5D

HSC: Hopkins Symptoms Checklist

MS: Multiple Sclerosis

SCI: Spinal Cord Injury

SDS: Self Rating Depression Scale

SIP: Sickness Impact profile

SS-QOL: Stroke Specific Quality of Life

ITB VALUE SUMMARY **INTRATHECAL BACLOFEN THERAPY**

Medtronic
Further, Together

CREAMER 2018



**PATIENT BENEFIT
IMPROVING
QUALITY OF LIFE**

Intrathecal Baclofen Therapy

Pharmacology

Technology Overview

Therapy Continuum

Clinical Evidence

Safety

Patient Benefit

Improving QOL

Improving Function

Improving Pain

Economic Value

Guidelines & Recommendations

References


Glossary

In stroke patients, after 6 months of ITB treatment a significant improvement of the EQ-5D utility score was observed when compared to patients treated with conventional medical management (oral medication)¹⁷.

EQ-5D UTILITY SCORE



EQ-5D VISUAL ANALOGUE SCORE




CMM= conventional medical management
Change from baseline, between group difference * $p < .05$

VLES 2013



**PATIENT BENEFIT
IMPROVING
QUALITY OF LIFE**

Assessment of the CHQ-PF50 score showed an improvement in pain, mental health, parental impact emotional, parental impact time and physical summary after long-term follow up of ITB treatment (6-9 years) in children with intractable spastic cerebral palsy⁴.




Changes from baseline: * $p < .05$; ** $p \leq .01$

SCHIESS 2011



**PATIENT BENEFIT
IMPROVING
QUALITY OF LIFE**

In stroke patients, the SS-QOL showed significant improvement at 12 months of ITB treatment, specifically in the domains of family roles, mobility, personality, self-care, social roles, thinking, upper extremity function, and work/productivity¹³.



Change from baseline: * $p \leq .05$; ** $p < .01$; *** $p \leq .001$

IVANHOE 2006



**PATIENT BENEFIT
IMPROVING
QUALITY OF LIFE**

Intrathecal Baclofen Therapy

Pharmacology

Technology Overview

Therapy Continuum

Clinical Evidence

Safety

Patient Benefit

Improving QOL

Improving Function

Improving Pain


Economic Value

Guidelines & Recommendations

References

Glossary

In stroke patients, after 12 months of ITB treatment, SIP scores improved significantly overall in both the physical domain (that includes body and care movement, mobility and ambulation) and the psychosocial domains (that includes emotional behaviour, social interest, alertness behaviour and communication), as well as in other categories such as home management or recreation and pastimes¹⁴.



NATALE 2016



**PATIENT BENEFIT
IMPROVING
QUALITY OF LIFE**

Intrathecal Baclofen Therapy

[Pharmacology](#)[Technology Overview](#)[Therapy Continuum](#)[Clinical Evidence](#)[Safety](#)[Patient Benefit](#)[Improving QOL](#)[Improving Function](#)[Improving Pain](#)[Economic Value](#)[Guidelines & Recommendations](#)[References](#)[Glossary](#)

In multiple sclerosis patients, the SDS score varied significantly from pump implantation (baseline value) to last follow-up (74 months in average)²⁵.




DELHAAS 2008



**PATIENT BENEFIT
IMPROVING
QUALITY OF LIFE**

In a group of spasticity patients (71% of which, were multiple sclerosis and spinal cord injured patients), the average EQ-5D visual analog scale was significantly higher after 6 and 12 months of ITB treatment²¹.



Change from baseline: * $p \leq .05$

MIDDEL 1997



**PATIENT BENEFIT
IMPROVING
QUALITY OF LIFE**

Intrathecal Baclofen Therapy

Pharmacology

Technology
OverviewTherapy
ContinuumClinical
Evidence

Safety

Patient Benefit

Improving QOL

Improving Function

Improving Pain



Economic Value

Guidelines &
Recommendations

References

Glossary

In multiple sclerosis and spinal cord injured patients, SIP scale and HSC score improved overall after 12 months of ITB treatment. In particular, sleep and rest, recreation and pastimes, mobility, body care and movement, physical dimensions and physical health showed a significant decrease²³.



Change from baseline: * $p \leq .05$; ** $p \leq .001$

SIP overall score has been calculated without the ambulation subscores as not applicable in this population group.

[Home](#)[What is Spasticity?](#)[The Burden of Spasticity](#)[Treatment Options](#)[Intrathecal Baclofen Therapy](#)[Pharmacology](#)[Technology Overview](#)[Therapy Continuum](#)[Clinical Evidence](#)[Safety](#)[Patient Benefit](#)[Improving QOL](#)[Improving Function](#)[Improving Pain](#)[Economic Value](#)[Guidelines & Recommendations](#)[References](#)[Glossary](#)

PATIENT BENEFIT

11 ITB STUDIES SHOW A SIGNIFICANT IMPROVEMENT IN PATIENT FUNCTION

CLICK ON EACH AUTHOR NAME FOR MORE INFORMATION ON THE STUDY

	INDICATIONS	AUTHOR	OUTCOME MEASURE	DOMAINS
CEREBRAL	CP	Motta 2011	GMFM	Overall score
	Stroke	Schuess 2011	FIM	Transfer, grooming, eating, bathing, toileting, ambulation, dressing-upper body and gait distance
			UEMAL	Gross and fine activities and amount of scale, quality of movement
			Walking Speed	
	CP	Ramstad 2010	PEDI	Self-care, mobility, social function
			GMFM	Overall score
	CP	Hoving 2009	GMFM	Sitting dimension, goal dimension
	Stroke	Ivanhoe 2006	FIM	Overall score, motor , transfer, self-care subscores
	CP	Guillaume 2005	FIM, WeeFIM	Overall score
SPINAL	Stroke	Francisco 2003	Modified FIM Walking Speed	Locomotion-walking, stair, sit to stand, stand to sit
	CP and TBI	Rawicki 1999	FIM, Snow Hygiene Scale	Transfer score, hygiene (ability to clean and self-catheterize)
	MS	Natale 2016	BI	Overall score
	MS and SCI	Guillaume 2005	FIM, WeeFIM	Overall score
	MS and SCI	Azouvi 1996	FIM	Overall score

BI: Barthel Index

CP: Cerebral Palsy

FIM: Functional Independence Measure

GMFM: Gross Motor Function Mobility

MS: Multiple Sclerosis

PEDI: Pediatric Evaluation of Disability Inventory

SCI: Spinal Cord Injury

TBI: Trauma Brain Injury

UEMAL: Upper Extremity Manual Activity Log

WeeFIM: Functional Independence Measure for children

ITB VALUE SUMMARY **INTRATHECAL BACLOFEN THERAPY**

Medtronic
Further, Together

MOTTA 2011




IMPROVEMENT IN PATIENT FUNCTION

Intrathecal Baclofen Therapy

[Pharmacology](#)[Technology Overview](#)[Therapy Continuum](#)[Clinical Evidence](#)[Safety](#)[Patient Benefit](#)[Improving QOL](#)[Improving Function](#)[Improving Pain](#)[Economic Value](#)[Guidelines & Recommendations](#)[References](#)[Glossary](#)

ITB increased the total median GMFM score in children with cerebral palsy after 1 year of treatment³.



Change from baseline: * $p < .001$

SCHIESS 2011



IMPROVEMENT IN PATIENT FUNCTION

Intrathecal Baclofen Therapy

Pharmacology

Technology Overview

Therapy Continuum

Clinical Evidence

Safety

Patient Benefit

Improving QOL

Improving Function

Improving Pain


Economic Value

Guidelines & Recommendations

References

Glossary

After 12 months of ITB treatment, stroke patients demonstrated improvement in all FIM scores (except bed mobility and lower body dressing), in walking velocity speed and in amount of use (AOU) and quality of movement (QOM) of the spastic upper extremities measured with the UEMAL¹³.



Change from baseline: * $p < .05$; ** $p < .01$

SCHIESS 2011



IMPROVEMENT IN PATIENT FUNCTION

Intrathecal Baclofen Therapy

Pharmacology

Technology
OverviewTherapy
ContinuumClinical
Evidence

Safety

Patient Benefit

Improving QOL

Improving Function



Improving Pain

Economic Value

Guidelines &
Recommendations

References

Glossary



Change from baseline: * $p < .05$; ** $p < .01$; *** $p < .001$

RAMSTAD 2010



IMPROVEMENT IN PATIENT FUNCTION

Intrathecal Baclofen Therapy

Pharmacology

Technology Overview

Therapy Continuum

Clinical Evidence

Safety

Patient Benefit

Improving QOL

Improving Function

Improving Pain



Economic Value

Guidelines & Recommendations

References

Glossary

In children with cerebral palsy, ITB treatment improved functional skills as demonstrated by the PEDi individual scores and GMFM-66 total score after 6 and 18 months².



Change from baseline: * p < .05; ** p < .01

HOVING 2009



IMPROVEMENT IN PATIENT FUNCTION

Intrathecal Baclofen Therapy

Pharmacology

Technology Overview

Therapy Continuum

Clinical Evidence

Safety

Patient Benefit

Improving QOL

Improving Function

Improving Pain


Economic Value

Guidelines & Recommendations

References

Glossary

After 12 months of ITB treatment, children with cerebral palsy significantly improved GMFM sitting and goal dimension scores with 3.3 and 4.0 points, respectively⁵.



Change from baseline: * $p < .05$; ** $p < .01$

IVANHOE 2006



IMPROVEMENT IN PATIENT FUNCTION

Intrathecal Baclofen Therapy

Pharmacology

Technology Overview

Therapy Continuum

Clinical Evidence

Safety

Patient Benefit

Improving QOL

Improving Function

Improving Pain


Economic Value

Guidelines & Recommendations

References

Glossary

In stroke patients, after 3 and 12 months of ITB treatment, FIM scores improved significantly overall and particularly in motor, transfer, self-care subscores¹⁴.




Change from baseline, overall (3 and 12 months) p value: * $p \leq .05$; ** $p < .01$

GUILLAUME 2005



IMPROVEMENT IN PATIENT FUNCTION

In patients with intractable spasticity of cerebral and spinal origin, the WeeFIM (used for patients ≤ 7 years of age) and the FIM (for patients over 7 years of age) overall score increased significantly after 12 months of ITB treatment⁸.



Change from baseline: * $p < .001$

FRANCISCO 2003




IMPROVEMENT IN PATIENT FUNCTION

Intrathecal Baclofen Therapy

[Pharmacology](#)[Technology Overview](#)[Therapy Continuum](#)[Clinical Evidence](#)[Safety](#)[Patient Benefit](#)[Improving QOL](#)[Improving Function](#)[Improving Pain](#)[Economic Value](#)[Guidelines & Recommendations](#)[References](#)[Glossary](#)

In stroke patients, ITB significantly improved walking speed and functional mobility rating, measured with a modified FIM score, at an average follow-up time of 8.9 months¹⁵.




Change from baseline: * $p < .05$; ** $p < .01$

RAWICKI 1999



IMPROVEMENT IN PATIENT FUNCTION

In patients suffering from spasticity of cerebral origin (cerebral palsy, cerebrovascular disease, trauma brain injury) ITB improved hygiene score (Snow) and transfer score (FIM) after long-term treatment (from 12 months to 9 years)¹⁰.




Change from baseline: * $p < .05$; ** $p < .001$

NATALE 2016



IMPROVEMENT IN PATIENT FUNCTION

In multiple sclerosis patients, the BI score varied significantly from pump implantation (baseline value) to last follow-up (74 months on average)²⁵.




Change from baseline: * $p < .005$

GUILLAUME 2005



IMPROVEMENT IN PATIENT FUNCTION

In patients with intractable spasticity of cerebral and spinal origin, the WeeFIM (used for patients ≤ 7 years of age) and the FIM (used for patients over 7 years of age) overall score increased significantly after 12 months of ITB treatment⁸.




Change from baseline: * $p < .001$

AZOUVI 1996



IMPROVEMENT IN PATIENT FUNCTION

In patients with severe and disabling spinal spasticity, ITB treatment improved FIM overall score after 6 months. Most dramatic improvements were observed in patients exhibiting a thoracic or low cervical lesion. Functional improvement for the most severely disabled patients, who were nearly totally dependent, was still significant¹⁹.



Change from before ITB: * $p < .05$; ** $p < .01$; *** $p < .001$

[Home](#)[What is Spasticity?](#)[The Burden of Spasticity](#)[Treatment Options](#)[Intrathecal Baclofen Therapy](#)[Pharmacology](#)[Technology Overview](#)[Therapy Continuum](#)[Clinical Evidence](#)[Safety](#)[Patient Benefit](#)[Improving QOL](#)[Improving Function](#)[Improving Pain](#)[Economic Value](#)[Guidelines & Recommendations](#)[References](#)[Glossary](#)

PATIENT BENEFIT

7 ITB STUDIES SHOW A SIGNIFICANT IMPROVEMENT IN PAIN

CLICK ON EACH AUTHOR NAME FOR MORE INFORMATION ON THE STUDY

	INDICATION	AUTHOR	OUTCOME MEASURE
CEREBRAL	Stroke	Creamer 2018	NPRS (0-10 scale, actual, least, worst pain)
	CP	Vles 2013	VAS (satisfaction for decrease in pain)
	CP	Ramstad 2010	Pain frequency during last 4 weeks, pain severity (0-4 scale)
SPINAL	MS	Lee 2018	NPRS (0-10 scale)
	MS	Sammaraiiee 2018	VAS (0-10 scale)
	MS	Natale 2016	VAS (0-10 scale)
	MS	Middel 1997	VAS (0-10 scale)

CP: Cerebral Palsy

MS: Multiple Sclerosis

NPRS: Numerical Pain Rating Scale

VAS: Visual Analogue Scale

CREAMER 2018



PATIENT BENEFIT REDUCING PAIN

Intrathecal Baclofen Therapy

Pharmacology

Technology
OverviewTherapy
ContinuumClinical
Evidence

Safety

Patient Benefit

Improving QOL

Improving Function

Improving Pain


Economic Value

Guidelines &
Recommendations

References

Glossary

In stroke patients, after 6 months of ITB treatment, a significant improvement in NPRS for actual and least pain was observed compared to patients treated with conventional medical management (oral medication)¹⁷.



VLES 2013




PATIENT BENEFIT REDUCING PAIN

Intrathecal Baclofen Therapy

[Pharmacology](#)[Technology Overview](#)[Therapy Continuum](#)[Clinical Evidence](#)[Safety](#)[Patient Benefit](#)[Improving QOL](#)[Improving Function](#)[Improving Pain](#)[Economic Value](#)[Guidelines & Recommendations](#)[References](#)[Glossary](#)

Satisfaction for improvement in pain significantly increased after 6 months of ITB treatment and was stable over time in children with spastic cerebral palsy. Higher VAS score represents less pain⁴.




Change from baseline * $p < .01$

RAMSTAD 2010



PATIENT BENEFIT REDUCING PAIN

Pain frequency (during the last 4 weeks) and pain severity were statistically reduced after 6 and 18 months of ITB treatment in a pediatric population with cerebral palsy².



Change from baseline * $p \leq 0.05$; ** $p \leq 0.01$

LEE 2018




PATIENT BENEFIT REDUCING PAIN

Intrathecal Baclofen Therapy

[Pharmacology](#)[Technology Overview](#)[Therapy Continuum](#)[Clinical Evidence](#)[Safety](#)[Patient Benefit](#)[Improving QOL](#)[Improving Function](#)[Improving Pain](#)[Economic Value](#)[Guidelines & Recommendations](#)[References](#)[Glossary](#)

In ambulatory multiple sclerosis patients, pain score decreased significantly after 6 and 12 months of ITB treatment²⁶.




Change from baseline * $p < .05$, ** $p < .005$

SAMMARAIIEE 2018



PATIENT BENEFIT REDUCING PAIN

In multiple sclerosis patients, the VAS score for pain, stiffness and discomfort decreased significantly after ITB implant and efficacy was sustained over time²⁷.



Change from baseline were significant over time. Pain $p < .005$, Stiffness $p < .0001$, Discomfort $p < .05$

NATALE 2016




PATIENT BENEFIT REDUCING PAIN

Intrathecal Baclofen Therapy

[Pharmacology](#)[Technology Overview](#)[Therapy Continuum](#)[Clinical Evidence](#)[Safety](#)[Patient Benefit](#)[Improving QOL](#)[Improving Function](#)[Improving Pain](#)[Economic Value](#)[Guidelines & Recommendations](#)[References](#)[Glossary](#)

In multiple sclerosis patients, the pain score significantly decreased from pump implantation (baseline value) to last follow-up (74 months on average)²⁵.



Change from baseline * $p < 0.005$


[Home](#)[What is Spasticity?](#)[The Burden of Spasticity](#)[Treatment Options](#)[Intrathecal Baclofen Therapy](#)[Pharmacology](#)[Technology Overview](#)[Therapy Continuum](#)[Clinical Evidence](#)[Safety](#)[Patient Benefit](#)[Improving QOL](#)[Improving Function](#)[Improving Pain](#)[Economic Value](#)[Guidelines & Recommendations](#)[References](#)[Glossary](#)

MIDDEL 1997



**PATIENT
BENEFIT
REDUCING
PAIN**

In multiple sclerosis and spinal cord injured patients, pain score significantly improved after 12 months of ITB treatment²³.



Change from baseline * $p < .01$

REFERENCES

1. Albright AL, Gilmartin R, Swift D et al. Long-term intrathecal baclofen therapy for severe spasticity of cerebral origin. *J Neurosurg.* 2003;98(2):291-295.
2. Ramstad K, Jahnsen R, Lofterd B, Skjeldal OH. Continuous intrathecal baclofen therapy in children with cerebral palsy – when does improvement emerge? *Acta Paediatr.* 2010;99(11):1661-1665.
3. Motta F, Antonello CE, Stignani C. Intrathecal baclofen and motor function in cerebral palsy. *Dev Med Child Neurol.* 2011;53(5):443-448.
4. Vles GF, Soudant DL, Hoving MA et al. Long-term follow-up on continuous intrathecal baclofen therapy in non-ambulant children with intractable spastic cerebral palsy. *Eur J Paediatr Neurol.* 2013;17(6):639-644.
5. Hoving MA, van Raak EP, Spincemaille GH et al; Dutch Study Group on Child Spasticity. Safety and one-year efficacy of intrathecal baclofen therapy in children with intractable spastic cerebral palsy. *Eur J Paediatr Neurol.* 2009;13(3):247-256.
6. Hoving MA, van Raak EP, Spincemaille GH, Palmans LJ, Becher JG, Vles JS; Dutch Study Group on Child Spasticity. Efficacy of intrathecal baclofen therapy in children with intractable spastic cerebral palsy: a randomised controlled trial. *Eur J Paediatr Neurol.* 2009;13(3):240-246.
7. Morton RE, Gray N, Vloeberghs M. Controlled study of the effects of continuous intrathecal baclofen infusion in non-ambulant children with cerebral palsy. *Dev Med Child Neurol.* 2011;53(8):736-741.
8. Guillaume D, Van Havenbergh A, Vloeberghs M, Vidal J, Roeste G. A clinical study of intrathecal baclofen using a programmable pump for intractable spasticity. *Arch Phys Med Rehabil.* 2005;86(11):2165-71.
9. Yoon YK, Lee KC, Cho HE et al. Outcomes of intrathecal baclofen therapy in patients with cerebral palsy and acquired brain injury. *H. Medicine (Baltimore).* 2017 Aug;96(34):e7472.
10. Rawicki B. Continuous intrathecal baclofen delivered via an implantable pump: long-term follow-up review of 18 patients. *J Neurosurg.* 1999;91(5):733-736.
11. Meythaler JM, McCary A, Hadley MN. Intrathecal infusion of baclofen for spasticity caused by acquired brain injury: a preliminary report. *J Neurosurg.* 1997;87(3):415-419.
12. Meythaler JM, Guin-Renfroe S, Grabb P, Hadley MN. long-term continuously infused intrathecal baclofen for spastic-dystonic hypertonia in traumatic brain injury: 1-year experience. *Arch Phys Med Rehabil.* 1999;80(1):13-19.
13. Schiess MC, Oh JJ, Stimming EF et al. Prospective 12-month study of intrathecal baclofen therapy for poststroke spastic upper and lower extremity motor control and functional improvement. *Neuromodulation.* 2011;14(1):38-45.
14. Ivanhoe CB, Francisco GE, McGuire JR et al. Intrathecal baclofen management of poststroke spastic hypertonia: implications for function and quality of life. *Arch Phys Med Rehabil.* 2006;87(11):1509-1515.
15. Francisco GE, Boake C. Spastic hemiplegia after intrathecal baclofen therapy: a preliminary study. *Arch Phys Med Rehabil.* 2003;84(8):1194-1199.
16. Creamer M, Cloud G, Kossmehl P, et al. Intrathecal baclofen therapy versus conventional medical management for severe post-stroke spasticity: Results from a multicentre, randomised, controlled, open-label trial (SISTERS). *J Neurol Neurosurg Psychiatry.* 2018 Jun;89(6):642-650.
17. Creamer et al. Effect of intrathecal baclofen on pain and quality of life in poststroke spasticity: A randomized trial (SISTERS). *Stroke.* 2018;49:2129-2137.
18. Meythaler JM, Guin-Renfroe S, Brunner RC, Hadley MN. Intrathecal baclofen for spastic hypertonia from stroke. *Stroke.* 2001;32(9):2009-2109.
19. Azouvi P, Mane, M, Thiebaut, J B, Denys, P, Remy-Neris, O, Bussel, B. Intrathecal baclofen administration for control of severe spinal spasticity: functional improvement and long-term follow-up. *Arch Phys Med Rehabil.* 1996;77(1):35-9.
20. Ordia JI, Fischer, Edward, Adamski, Ellen, Chagnon, Kimberly G, Spatz, Edward L. Continuous intrathecal baclofen infusion by a programmable pump in 131 consecutive patients with severe spasticity of spinal origin. *Neuromodulation.* 2002;5(1):16-24.
21. Delhaas EM. Long-term outcomes of continuous intrathecal baclofen infusion for treatment of spasticity: A prospective multicenter follow-up study. *Neuromodulation.* 2008;11(3):227-236.
22. Burns AS, Meythaler JM. Intrathecal baclofen in tetraplegia of spinal origin: efficacy for upper extremity hypertonia. *Spinal Cord.* 2001;39(8):413-9.
23. Middel B, Kuipers-Upmeyer H, Bouma J, et al. Effect of intrathecal baclofen delivered by an implanted programmable pump on health related quality of life in patients with severe spasticity. *J Neurol Neurosurg Psychiatry.* 1997;63(2):204-9.
24. Khan AA, Birks-Agnew I, Bullock P, Rushton D. Clinical outcome and complications of intrathecal baclofen pump in multiple sclerosis patients: a retrospective study. *NeuroRehabilitation.* 2010;27(2):117-20.
25. Natale M, D'Oria S, Nero VV, Squillante E, Gentile M, Rotondo M. Long-term effects of intrathecal baclofen in multiple sclerosis. *Clin Neurol Neurosurg.* 2016;143:121-5.
26. Lee BS, Jones J, Lang M et al. Early outcomes after intrathecal baclofen therapy in ambulatory patients with multiple sclerosis. *J Neurosurg.* 2018 Oct;129(4):1056-1062.
27. Sammarraiee Y, Yardley M, Keenan L et al. Intrathecal baclofen for multiple sclerosis related spasticity: A twenty year experience. *Mult Scler Relat Disord.* 2018 Oct 15;27:95-100.

ECONOMIC VALUE

COST-EFFECTIVENESS OF ITB



INTRODUCTION OF ITB DEMONSTRATES COST-SAVINGS OVER TIME, POST IMPLANTATION



A few studies have shown reduction in hospitalization length of stay vs. conventional medical management (CMM)¹⁻³.



Lifetime analysis cost modeling indicates that ITB could afford cost savings of up to \$8,009 USD/patient/year, compared with conventional therapy⁴.

Cost savings from this analysis were derived primarily from reduced hospital admissions and physician visits⁴.


COST-EFFECTIVENESS

ITB IMPROVES
PATIENT OUTCOMES
WHILE MAINTAINING
REASONABLE
HEALTHCARE COSTS,
IN THE LONG TERM.



To date, there are **5 published cost-effectiveness studies** from Europe, Japan and the U.S. which provide cost-effectiveness ratios (ICERs) for ITB vs. comparative treatments that, in most cases, are below, or close to, the willingness-to-pay (WTP) thresholds across different country systems⁵⁻⁹.

These studies demonstrate that ITB provides incremental health benefits at an acceptable incremental cost vs. comparative treatments, when used to treat disabling spasticity in appropriately-selected patients, across time horizons of 1-5 years and patient's lifetime⁵⁻⁹.




[CLICK TO ENLARGE](#)

COST-EFFECTIVENESS OF ITB

Incremental cost-effectiveness ratios for ITB from cost-effective analyses with CMM as comparator

(Original values inflated to 2017 then currency converted to 2017 Euros)



*Range of incremental cost-effectiveness thresholds generally accepted in Europe. Only the UK provides published explicit ICER thresholds.¹⁰

REFERENCES

1. Ordia JI, Fischer E, Adamski E, Spatz EL. Chronic intrathecal delivery of baclofen by a programmable pump for the treatment of severe spasticity. *J Neurosurg*. 1996;85(3):452-457.
2. Becker WJ, Harris CJ, Long ML, Ablett DP, Klein GM, DeForge DA. Long-term intrathecal baclofen therapy in patients with intractable spasticity. *Can J Neurol Sci*. 1995;22(3):208-217.
3. Nance P, Schryvers O, Schmidt B, Dubo H, Loveridge B, Fewer D. Intrathecal baclofen therapy for adults with spinal spasticity: therapeutic efficacy and effect on hospital admissions. *Can J Neurol Sci*. 1995;22(1):22-29.
4. Saulino M, Guillemette S, Leier J, Hinnenthal J. Medical cost impact of intrathecal baclofen therapy for severe spasticity. *Neuromodulation*. 2015;18(2):141-149.
5. de Lissovoy G, Matza LS, Green H, Werner M, Edgar T. Cost-effectiveness of intrathecal baclofen therapy for the treatment of severe spasticity associated with cerebral palsy. *J Child Neurol*. 2007;22(1):49-59. Costs reported in this study were adjusted to base year 2003. Costs in the study were reported in base year 2003. The reported ICER in (USD) inflated to November 2017 rate using the Bureau of Labor Statistics, US Dept. of Labor. CPI Inflation Calculator. Source: http://www.bls.gov/data/inflation_calculator.htm (accessed January 21, 2019). This value was then converted to 2017 Euros (30 November 2017) using Reuters online currency calculator (Exchange rate: USD:Euro/1:0.840) Source: <http://www.reuters.com/finance/currencies> (accessed on January 21, 2019).
6. Hattori N, Hirayama T, Katayama Y. Cost-effectiveness analysis of intrathecal baclofen therapy in Japan. *Neurol Med Chir*. 2012;52(7):482-487. Costs reported in the study were adjusted to base year 2009. The reported ICER of 1,554,428 Japanese Yen/QALY was inflated to 2016 rates (latest available) based on Inflation Rates (all items) year on year from 2009 to 2016 to provide 2016 rates in YEN, using annual inflation data available from OECD country statistics: Japan. Source: http://www.oecd-ilibrary.org/economics/country-statistical-profile-japan-2016-2_csp-jpn-table-2016-2-en (September 21, 2016). This value was then converted to 2017 EUROS (20 November 2018) using Reuters online currency calculator (Exchange rate Euro:Yen/1:0.008) Source: <http://www.reuters.com/finance/currencies> (accessed 21 January 2019).
7. Sampson FC, Hayward A, Evans G, Morton R, Collett B. Functional benefits and cost/benefit analysis of continuous intrathecal baclofen infusion for the management of severe spasticity. *J Neurosurg*. 2002;96(6):1052-1057. Costs in the study were reported in base year 2002. The reported ICERs in (£) were inflated to 2017 rates using the Bank of England Inflation Calculator - Source: <https://www.bankofengland.co.uk/monetary-policy/inflation/inflation-calculator> (accessed on 21/01/2019). This value was then converted to 2017 Euros using Reuters online currency calculator (Exchange rate £:Euro/1:1.136 on 30/11/2017) Source: <http://www.reuters.com/finance/currencies> (accessed on 21/01/2019).
8. Hoving MA, Evers SM, Ament AJ, van Raak EP, Vles JS. Intrathecal baclofen therapy in children with intractable spastic cerebral palsy: a cost-effectiveness analysis. *Dev Med Child Neurol*. 2008;50(6):450-455. Cost in the study were reported in base year 2003. The reported ICER in (€) was inflated using Inflation Rates (all items) year on year from 2003 to 2017 to provide 2017 rates in Euros, using annual inflation data available from OECD country statistics: Netherlands. Source: http://www.oecd-ilibrary.org/economics/country-statistical-profile-netherlands-2011_csp-nld-table-2011-1-enjsessionid=bbjom423ftk95.x-oecd-live-03 (September 21, 2016).
9. Vidal, J, Slof, J, Serrano, D, Marqués, T, Kumru, H & Benito-Penalva, J (2016): Expert Review of Pharmacoeconomics & Outcomes Research, Published online: 03 May 2016, DOI: 10.1080/14737167.2016.118024. Costs in the study were reported in base year 2016. The reported ICERs in (€) were inflated to November 2017 rates using the inflation calculator from the Spanish National Institute of Statistics - Source: <http://www.ine.es/calcula/> (accessed on 21/01/2019).
10. McCabe C, Claxton K, Culyer AJ. The NICE cost-effectiveness threshold: what it is and what that means. *Pharmacoeconomics*. 2008;26(9):733-44. Exchange rate of UK Pounds to Euros of 1:1.136 as of 30/11/2017, source: <http://www.reuters.com/finance/currencies> (accessed on 15/01/2019).

HEALTH TECHNOLOGY ASSESSMENT AND POLICY GUIDANCE

CANADA

HTA/GOVERNMENT BODY

Ontario Ministry of Health and Long-Term Care (Medical Advisory Secretariat)
HTA, 2005: Intrathecal baclofen pump for spasticity¹.

RECOMMENDATIONS

Level 2 evidence supports effectiveness of ITB infusion for short-term reduction of severe spasticity in patients who are unresponsive/cannot tolerate oral baclofen.

Level 3 evidence supports effectiveness of ITB for long-term reduction of severe spasticity in patients who are unresponsive/cannot tolerate oral baclofen.

SPAIN

HTA/GOVERNMENT BODY

Agencia de Evaluación de Tecnologías Sanitarias de Andalucía (AESTA) (2000): Implantable infusion pumps: clinical uses and evidence of effectiveness and safety in pain and spasticity – systematic review³.

RECOMMENDATIONS

ITB can be indicated for selected cases of severe spasticity of spinal origin.

UK

HTA/GOVERNMENT BODY

NHS Clinical Commissioning
Policy: Intrathecal baclofen (ITB)², April 2013

RECOMMENDATIONS

ITB for adults or paediatrics with chronic, severe, diffuse spasticity and/or dystonia of spinal or cerebral origin which renders them a full-time wheelchair user or bed-bound.

The policy supports the use of ITB for the groups for which it is the most cost-effective, where other options are exhausted, and where patient and carer evidence shows a real likelihood of success.

POLAND

HTA/GOVERNMENT BODY

Agency for Health Technology Assessment for provision of healthcare (2011)⁴.

RECOMMENDATIONS

ITB for the treatment of spasticity refractory to pharmacological treatment.

CLINICAL GUIDELINES FOR INTRATHECAL BACLOFEN

INDEPENDENT EUROPEAN GUIDELINES RECOMMEND ITB FOR DISABLING NON-FOCAL SPASTICITY

CLINICAL GUIDELINE & COUNTRY	SYMPTOMS/POPULATION	RECOMMENDATIONS FOR ITB
<p>National Institute for Health & Care Excellence (NICE), United Kingdom. Cerebral palsy in adults. NICE Guideline 119 (2019)⁵</p> <p>Link: https://www.nice.org.uk/guidance/NG119</p>	Cerebral palsy in adults	Consider referral to a tone/spasticity management service offering continuous pump-administered ITB if they still have difficulties with spasticity, despite enteral muscle relaxant drug treatment or botulinum toxin type A treatment.
<p>Federatie Medisch Specialisten: Netherlands Society of Rehabilitation Medicine (NSRM), The Netherlands. Guideline for the treatment of cerebral and/or spinal spasticity in adults (2017)⁶</p> <p>Link: https://revalidatiegeneeskunde.nl/article/richtlijn-behandeling-van-cerebrale-enof-spinale-spasticiteit-bij-volwassenen-gepubliceerd</p>	Spasticity in adults	<p>ITB is efficient in decreasing spasticity compared to placebo in patients with cerebral and/or spinal spasticity with both observable impact on upper and lower extremities.</p> <p>Evidence Grade: High</p>

CLINICAL GUIDELINES FOR INTRATHECAL BACLOFEN

INDEPENDENT EUROPEAN GUIDELINES RECOMMEND ITB FOR DISABLING NON-FOCAL SPASTICITY

CLINICAL GUIDELINE & COUNTRY	SYMPTOMS/POPULATION	RECOMMENDATIONS FOR ITB
<p>The Finnish Medical Society Duodecim, Finland national guidelines for stroke (2016)⁷</p> <p>Link: http://www.kaypahoito.fi/web/kh/suosituksset/suositus;jsessionid=6958396585ECA8BA65BE630C1293EF0C?id=nix00624</p>	Spasticity related to cerebrovascular disease; stroke	<p>ITB is efficient in decreasing spasticity compared to placebo in patients with cerebral and/or spinal spasticity with both observable impact on upper and lower extremities.</p> <p>Evidence Grade: High</p>
<p>German Neurology Society (Deutsche Gesellschaft für Neurologie) S1 Guideline: Therapy of spastic syndromes (2012)⁸</p> <p>(currently under review)</p> <p>Link: https://www.dgn.org/leitlinien/2431-II-89-2012-therapie-des-spastischen-syndroms</p>	Spastic syndromes	In patients with severe, generalized spasticity, spastic tetra- or paraplegia that cannot be adequately treated with physiotherapy, oral antispastic medication, botulinum toxin type A, or ITB should be considered.
<p>National Institute for Health & Care Excellence (NICE), United Kingdom. Spasticity in under 19s: management. Clinical Guideline 145 (2012)⁹</p> <p>(Last updated November 2016)</p> <p>Link: https://www.nice.org.uk/guidance/cg145</p>	Spasticity in people aged 0-19 years old	Consider treatment with continuous pump-administered ITB if, despite the use of non-invasive treatments, spasticity is causing difficulties with pain, muscle spasms, posture, function or self-care (or ease of care by parents or carers).

CLINICAL GUIDELINES FOR INTRATHECAL BACLOFEN

INDEPENDENT EUROPEAN GUIDELINES RECOMMEND ITB FOR DISABLING NON-FOCAL SPASTICITY

CLINICAL GUIDELINE & COUNTRY	SYMPTOMS/POPULATION	RECOMMENDATIONS FOR ITB
<p>Socialstyrelsen, Sweden National guidelines for the treatment of multiple sclerosis (MS) and parkinson's disease: support for governance and management (2016)¹⁰</p> <p>Link: http://www.socialstyrelsen.se/nationellariktlinjermsochparkinsonssjukdom</p>	<p>Multiple sclerosis in adults</p>	<p>ITB pump therapy should be offered for the treatment of severe spasticity associated with MS in the upper or lower limbs (Priority 4*).</p> <p>*Scale of Priority 1-10, where Priority 1 should receive the greatest share of resource allocation, and Priority 10 the least.</p>
<p>British Pain Society, United Kingdom. Intrathecal drug delivery for the management of pain and spasticity in adults; recommendations for best clinical practice (2015)¹¹</p> <p>Link: https://www.britishpainsociety.org/static/uploads/resources/files/itdd_2015_pro_v3.pdf</p>	<p>Pain and spasticity in adults</p>	<p>The working group believes* that ITB is well-established in the management of both cerebral and spinal spasticity and is a cost-effective method of baclofen delivery for spasticity.</p> <p>*These recommendations are primarily evidence-based but where necessary comprise the consensus opinion of the working groups.</p>

REFERENCES

1. Medical Advisory Secretariat. Intrathecal baclofen pump for spasticity: an evidence-based analysis. Toronto: Medical Advisory Secretariat (MAS). Volume 5(7). 2005.
2. National Health Service (NHS) Commissioning Board. NHS Clinical Commissioning Policy: Intrathecal baclofen (ITB). April 2013.
3. Agencia de Evaluación de Tecnologías Sanitarias de Andalucía (AESTA) (2000): Bombas De Infusión Implantables. Implantable infusion pumps: clinical uses and evidence of effectiveness and safety in pain and spasticity – systematic review. Seville: Andalusian Agency for Health Technology Assessment (AETSA) 2000. Accessed at: http://www.aetsa.org/download/publicaciones/antiguas/AETSA_P_2000_Bombas_intratec.pdf.
4. Agency for Health Technology Assessment, Poland. Recommendation No 66/2011 of the President of the Agency for Health Technology Assessment on the eligibility of the provision of health care "baclofen pump implantation in the treatment of spasticity refractory to pharmacological treatment" as a guaranteed provision. Accessed at www.aotm.gov.pl. September, 2011.
5. National Institute for Health and Clinical Excellence (NICE), United Kingdom. Cerebral palsy in adults. NICE Guideline 119 (NG 119). Published: 15 January 2019. Accessed at: <https://www.nice.org.uk/guidance/NG119>.
6. Federatie Medisch Specialisten: Netherlands Society of Rehabilitation Medicine (NSRM), The Netherlands. Guideline for the treatment of cerebral and/or spinal spasticity in adults (2017). Accessed at: <https://revalidatiegeneeskunde.nl/article/richtlijn-behandeling-van-cerebrale-en-of-spinale-spasticiteit-bij-volwassenen-gepubliceerd>.
7. National guidelines for stroke – treatment of spasticity in Aivoverenvuotohäiriö (AVH) patients (2016), The Finnish Medical Society Duodecim, Finland. Accessed at: <http://www.kaypahoito.fi/web/kh/suositus/suositus.jsessionid=6958396585ECA8BA65BE630C1293EFOC?id=nix00624>.
8. German Neurology Society (Deutsche Gesellschaft für Neurologie) (2012) S1 Guideline: Therapy of spastic syndromes8 (currently under review). Accessed at: <https://www.dgn.org/leitlinien/2431-II-89-2012-therapie-des-spastischen-syndroms>.
9. National Institute for Health & Care Excellence (NICE), United Kingdom. Spasticity in under 19s: management. Clinical Guideline 145 (2012). Last updated November 2016. Accessed at: <https://www.nice.org.uk/guidance/cg145>.
10. Socialstyrelsen, Sweden (2016) National guidelines for the treatment of multiple sclerosis (MS) and parkinson's disease: support for governance and management (2016). Accessed at: <http://www.socialstyrelsen.se/nationellariklinjermsochparkinsonssjukdom>.
11. British Pain Society, United Kingdom. Intrathecal drug delivery for the management of pain and spasticity in adults; recommendations for best clinical practice (2015) Accessed at: https://www.britishpainsociety.org/static/uploads/resources/files/itdd_2015_pro_v3.pdf.

GLOSSARY

Acquired Brain Injury (ABI) or Brain Injury (BI): Brain injury is an inclusive term for damage to the brain that occurs after birth and is not related to a hereditary, congenital, or degenerative disease. Traumatic brain injury, hypoxic and anoxic brain injury, and stroke are forms of acquired brain injuries.

Ashworth Scale / Modified Ashworth Scale: The Ashworth Scale tests resistance to passive movement around a joint with varying degrees of velocity. Scores range from 0–4, with 0 indicating no resistance and 4 indicating rigidity. The Modified Ashworth Scale is similar, but adds a +1 score used to indicate resistance throughout less than half of the range of movement.

Baclofen: Drug that inhibits the release of excitatory neurotransmitters, inhibiting spastic response to stretch reflex.

Barthel Index (BI): The Barthel Index is an ordinal scale used to measure the ability of an individual with a neuromuscular or musculoskeletal disorder to care for him/herself and perform activities of daily living. It assesses 10 activities of daily living and mobility activities rated by the amount of assistance required to complete: feeding, bathing, grooming, dressing, bowel control, bladder control, toileting, chair transfer, ambulation, and stair climbing. See also: Functional Independence Measure.

Bolus: Discrete amount of medication, drug or other compound in order to raise its concentration in blood to an effective level.

Case series: Coherent and consecutive set of cases of a disease (or similar problem) which derive from the practice of one or more health care professionals or health care setting. The cases were defined by similar diagnoses/disease state or having undergone the same intervention/procedure.

Central nervous system (CNS): The central nervous system is the part of the nervous system consisting of the brain and spinal cord.

Cerebral palsy (CP): Cerebral palsy is a non-progressive neurological disorder due to brain injury or malformation early in life. It is a group of permanent early childhood movement disorders. Motor symptoms may include spasticity.

Cerebrospinal fluid (CSF): The cerebrospinal fluid is a clear, colorless body fluid found in the brain and spine. It occupies the subarachnoid space (between the arachnoid mater and the pia mater) and the ventricular system around and inside the brain and spinal cord.

Child Health Questionnaire™ Parent Form (CHQ-PF50): The CHQ-PF50 is a standardized quality of life questionnaire designed for children of 5–18 years of age. It measures 14 individual health concepts and can be aggregated into 2 summary component scores for physical functioning and psychosocial health on a 0–100 scale.

Controlled study: A study with 2 or more groups, with different intervention/exposure, and outcomes were measured in all groups and compared between groups.

Conventional Medical Management (CMM): A nonsurgical spasticity treatment protocol that may include physical and occupational therapy, orthotics, mobility aids, oral medications, and chemodenervation. CMM does not include ITB therapy, neurosurgery, or orthopedic surgery.

Cost-effectiveness Analysis (CEA): A form of economic analysis used to support decision making about a treatment choice versus an alternative. CEA compares the relative costs and outcomes (effects) of two courses of action, and its results are expressed as a ratio (ICER).

EuroQol-5D (EQ-5D): The EQ-5D is a non-disease-specific patient self-assessment questionnaire that measures health-related quality of life.

Functional Independence Measure (FIM): The FIM is an 18-item tool used in an inpatient rehabilitation setting to measure the severity of disability and need for assistance in carrying out activities of daily living. Tasks are rated on a 7-point ordinal scale from total assistance/complete dependence to complete independence. A higher score indicates higher level of function. The FIM's 18 items are categorized into 6 motor and cognitive subscales.

Gross Motor Function Measure (GMFM): The GMFM is a standardized clinical tool that evaluates change over time in five dimensions of motor function ability in children aged 5 months to 16 years with CP: lying and rolling, sitting, crawling and kneeling, standing and walking, and running and jumping. It uses a 4-point score for each item, summed to calculate raw and percent scores for each dimension and overall score.

Hopkins Symptoms Checklist (HSCL): The HSCL is a self-report symptom inventory which measures symptoms of anxiety and depression. It is comprised of 58 items which are representative of the symptom configurations commonly observed among outpatients. It is scored on five underlying symptom dimensions—sommatization, obsessive-compulsive, interpersonal sensitivity, anxiety and depression—which have been identified in repeated factor analyses.

Incremental cost-effectiveness ratio (ICER): The key result of a CEA, expressed as the additional cost per QALY gained. The incremental cost per QALY is interpreted as the additional cost for 1 additional year of life in full health.

Intrathecal (subarachnoid) space: The space containing cerebrospinal fluid (CSF) between the arachnoid mater and pia mater of the spinal cord.

Intrathecal baclofen (ITB): Delivery of baclofen directly into the intrathecal space by lumbar puncture or an indwelling intrathecal catheter connected to a drug delivery pump.

Lumbar Puncture: Procedure in which a needle is inserted into the spinal canal.

Multiple sclerosis (MS): Multiple sclerosis is an immune-mediated progressive disease that destroys the myelin membrane of nerves in the brain and spinal cord. Demyelination impedes neuronal communication, which may result in motor symptoms including spasticity.

Patient Health Questionnaire (PHQ-9): The PHQ-9 is a multipurpose instrument for screening, diagnosing, monitoring and measuring the severity of depression.

Pediatric Evaluation of Disability Inventory (PEDI): The PEDI is an assessment tool for children that measures capability and performance of functional activities in 3 content areas: self-care, mobility, and social function. Scores range from 0–100, with higher scores indicating lesser disability.

Peristaltic Pump: Mechanical pump in which pressure is provided by the movement of a constriction along a tube, similar to biological peristalsis.



GLOSSARY

Pharmacology: Branch of medicine concerned with the uses, effects, and modes of action of drugs.

Prevalence: Percentage of a population that is affected with a particular disease at a given time.

Prospective study: The subjects were enrolled in the study after the study had been designed and approved, and the subjects were followed for some period of time.

Quality Adjusted Life Year (QALY): A QALY is a generic measure of effectiveness that encompasses both quality and quantity of life (i.e., survival), providing a consistent and common measure that healthcare funders can use to inform funding decisions.

Randomized Controlled Trial (RCT): An experimental comparative study in which participants are randomly assigned to intervention or control and followed over time to measure differences in outcomes. Concealment of the allocation of the randomized assignment is important. Known and unknown prognostic factors (confounders) should be evenly distributed between the groups.

Retrospective study: The study started after the data/information had been collected.

Rivermead Mobility Index (RMI): 15 items questionnaire for quantifying functional abilities such as gait, balance, and transfers after stroke.

Self Rating Depression Scale (SDS): The Self-Rating Depression Scale is a 20-item self-report questionnaire that is widely used as a screening tool, covering affective, psychological and somatic symptoms associated with depression. Each item is scored on a Likert scale ranging from 1 to 4. A total score is derived by summing the individual item scores, and ranges from 20 to 80. Most people with depression score between 50 and 69, while a score of 70 and above indicates severe depression.

Short Form Survey (SF-36 and SF-12): 36 or 12-question measure of health perception in eight dimensions: physical functioning, role limitation because of physical health, social functioning, vitality or energy, bodily pain, mental health, role limitation because of emotional problems, and general health. Provides Physical Component Scores (PCS) and Mental Component Scores (MCS-36) ranging from 0-100, with higher scores representing better self-reported health.

Sickness Impact Profile (SIP): The SIP is a yes/no patient self-assessment that measures psychosocial and physical quality of life in 12 categories (sleep and rest, eating, work, home management, recreation and pastimes, social interaction, alertness behavior, emotional behavior, communication, ambulation, mobility, and body care and movement). A lower score indicates improvement. It is designed to objectively assess outcomes from health care services. There are 2 versions of the SIP; one with 136 items and one with 68 items.

Snow Hygiene Score: A measurement of ability to clean and self-catheterize first described by Snow et al (1990) within a case series of 9 MS patients receiving botulinum toxin for spasticity. This score was used by Rawicki (1999) to measure hygiene in 18 patients with cerebral origin spasticity receiving ITB therapy.

Spasm: Sudden involuntary muscular contraction or convulsive movement.

Spasticity: Spasticity is an abnormal increase in muscle tone caused by injury of upper motor neuron pathways regulating muscles. Spasticity may be a result of multiple sclerosis, cerebral palsy, stroke, brain injury, or spinal cord injury.

Spinal Cord Injury (SCI): Injury to the spinal cord resulting in a change, either temporary or permanent, in the cord's normal motor, sensory, or autonomic function. SCI symptoms vary widely, and can include pain, paralysis, or spastic hypertonia.

Stretch Reflex: Muscle contraction in response to stretching within the muscle. It is a monosynaptic reflex which provides automatic regulation of skeletal muscle length.

Stroke: Also known as cerebrovascular accident (CVA), or cerebrovascular insult (CVI). Types of stroke include ischemic, resulting from impeded blood flow to a region of the brain, and hemorrhagic, due to bleeding in the brain. Symptoms of stroke may include an inability to move or feel on one side of the body, difficulty understanding or speaking, and spastic hypertonia in the long term.

Stroke-Specific Quality of Life scale (SS-QOL): Patient-centered outcome measure assesses health-related quality of life (HRQOL), specific to patients with stroke.

Traumatic Brain Injury (TBI): Traumatic brain injury, a form of acquired brain injury (ABI), is damage to the brain caused by an external mechanical force, such as a motor vehicle accident, explosive blast injury, or penetrating injury.

Upper Extremity Manual Activity Log (UEMAL): Structured interview intended to examine how much and how well the subject use of their affected arms. Participants are asked standardized questions about the amount of use of their more-affected arm and the quality of their movement during the functional activities indicated uses their more-affected arm.

Upper Motor Neuron Lesion: An upper motor neuron lesion (also known as pyramidal insufficiency) is a lesion of the neural pathway above the anterior horn cell of the spinal cord or motor nuclei of the cranial nerves. This is in contrast to a lower motor neuron lesion, which affects nerve fibers traveling from the anterior horn of the spinal cord or the cranial motor nuclei to the relevant muscle(s).

Upper Motor Neuron Syndrome (UMNS): Upper motor neuron syndrome (UMNS) is the motor control changes that can occur in skeletal muscle after an upper motor neuron lesion.

WeeFIM: The WeeFIM or Pediatric Functional Independence Measure is a validated pediatric outcomes tool that measures disability in children with developmental disorders. It measures how much assistance a child needs to perform activities of daily living. It can be used to track functional improvement and goal attainment.

Willingness-To-Pay Threshold: Willingness-to-pay threshold is a threshold above which treatments are no longer considered cost-effective. An ICER is meaningful with respect to this threshold, which is approximately £20,000–£30,000 in the United Kingdom (Shirowa 2010, Health Economics 19(4):422–37). The probability of payers not paying for a therapy increases significantly with increases in the ICER. (Shirowa T, Sung YK, Fukuda T et al. International survey on willingness-to-pay (WTP) for one additional QALY gained: What is the threshold of cost effectiveness? Health Econ. 2010;19(4):422–437).

Home

What is
Spasticity?

The Burden
of Spasticity

Treatment
Options

Intrathecal
Baclofen
Therapy

Brief Statement

See the device manual for detailed information regarding the instructions for use, the implant procedure, indications, contraindications, warnings, precautions, and potential adverse events. For further information please contact your local Medtronic representative and/or consult the Medtronic website at www.medtronic.com.



www.medtronic.com/manuals

Consult instructions for use on this website. Manuals can be viewed using a current version of major internet browser. For best results, use Adobe Acrobat® Reader in the browser.

Medtronic

Medtronic International Trading Sarl
Route du Molliau 31
Case postale
1131 Tolochenaz
Switzerland
Tel: +41 (0) 21 802 70 00
Fax: +41 (0) 21 802 79 00

medtronic.eu

Medtronic Limited
Building 9
Croxley Park
Hatters Lane
Watford, Hertfordshire
WD18 8WW
Tel: +44 (0) 1923 212213
Fax: +44 (0) 1923 241004

medtronic.co.uk

Medtronic Ireland Limited
Block 3090-3094
Lake Drive
Citywest Business Campus
Dublin
DN24 XN47
Tel: +353 (0) 1 511 1400
Fax: +353 (0) 1 807 7220

medtronic.ie

UC201912633 EE © Medtronic 2019.
All Rights Reserved.

Not for distribution in France

