

Medtronic


SPHERE Per-AF

U.S. IDE

study results

Source: Anter E, Mansour M, Nair DG, et al. Dual-energy lattice-tip ablation system for persistent atrial fibrillation: a randomized trial. *Nat Med.* 2024.

Randomized 1:1 Persistent AF clinical trial



Sphere-9 catheter
with Affera mapping
and ablation system
n = 212

THERMOCOOL SMARTTOUCH® SF
with CARTO®* mapping system
n = 208

Trial results at a glance

Primary safety
event rate[†]

1.4%

Primary
effectiveness[‡]

73.8%

Procedural characteristic
superiority over STSF

25
min less

Superior skin-to-skin
procedural time

27
min less

Superior time between
first and last application

29
min less

Superior energy
application time

[†] For a full list of safety events, review the SPHERE Per-AF manuscript.

[‡] Primary effectiveness endpoint definition: The primary effectiveness endpoint was acute procedure failure, repeat ablation at any time, or after three months: recurrence of AF/AFL/AT, cardioversion, or new/re-initiated/increased AAD usage.

The Sphere-9™ catheter was **noninferior in safety** and provided **quality of life improvements** for patients

1.4%

primary safety event
rate vs. 1% with the
STSF catheter†



0 PV stenosis



0 Phrenic nerve
paralysis



0 Major
vascular access
complications



0 Atrio-esophageal
fistula

3/212 Hospitalizations due to cardiovascular or pulmonary AE‡

0/212 Death	0/212 Thromboembolism
0/212 Myocardial infarction	0/212 Heart block
0/212 Transient ischemic attack (TIA)	0/212 Gastroparesis
0/212 Stroke/Cerebrovascular accident (CVA)	0/212 Severe pericarditis



0 Cardiac
tamponade

Patients treated with the Sphere-9™ catheter experienced improvements to quality of life

**SF-12v2 health survey scores improved by 3.2 points
(mental) and 4.7 points (physical)**
from baseline to 12 months.

AFEQT questionnaire scores improved by 22.3 points
from baseline to 12 months.



† Investigational arm – 212 patients. Control Arm – 208 patients.

‡ The 3 hospitalization events within 1 week post-procedure and were for chronic obstructive pulmonary disease exacerbation, pulmonary edema, and hemoptysis.


The Sphere-9 catheter delivered **higher effectiveness** outcomes for persistent AF patients

73.8%

Primary effectiveness at 12 months vs. 65.8% in the control arm

- Demonstrated noninferior effectiveness
- Higher treatment success compared to the control arm at 3, 6, 9, and 12 months
- 100% acute procedural success

Patients treated with Sphere-9 catheter observed less recurrence of AF, AT, and AFL throughout the entire follow-up period




No. at risk		Days from index ablation procedure				
		0	90	180	270	360
Investigational	212	211	181	166	84	
Control	208	202	159	148	82	


† The **Kaplan-Meier analysis** included early study exits via censoring which accounts for the slightly different effectiveness results.

Primary effectiveness testing

Difference between Sphere-9 catheter and STSF.




The Sphere-9 catheter with all-in-one design delivered **procedural characteristic superiority**



Dual energy wide area focal catheter

- PF lesions are delivered in four-second intervals
- RF lesions are delivered in five-second intervals



Additional procedural characteristic results[†]


Parameter	Sphere-9 (n = 212)	STSF (n = 208)
Fluoroscopy time [†]	4.9 minutes	6.3 minutes
Time between first and last application for pulmonary vein isolation	25.9 minutes	53.6 minutes
Total fluid delivered by ablation catheters	482.0 mL	727.1 mL
Procedures with single transseptal accesses	95.3%	62.0%
Number of mapping and/or ablation catheters used in LA	1 catheter: 97.2% 2 catheters: 2.8%	1 catheter: 0% 2 catheters: 98.5% 3 catheters: 1.5%

[†] For a full list of procedural characteristic results, review the SPHERE Per-AF manuscript.

SPHERE Per-AF trial design and study population

Randomized, controlled clinical trial comparing Sphere-9 catheter with the Affera mapping and ablation system to the THERMOCOOL SMARTTOUCH[®] SF catheter with CARTO[®] mapping system for the treatment of patients with persistent atrial fibrillation refractory or intolerant to drugs.

Randomized 1:1



Treatment strategy

PVI plus additional lesions at physician discretion in both arms

Trial investigator[†] experience with Sphere-9 catheter

6

Average number of patients treated per operator in the investigational arm

87.5%

Physicians using Affera system for the first time

[†] All operators had significant experience with the control arm device (typically > 1,000 patients treated).

By the numbers



3 countries

U.S., Israel, Czechia



23 sites

activated



> 40

operators



420 patients

treated in
primary analysis



< 12 months

to complete enrollment

Brief statement

See the device manual for information regarding the instructions for use, indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative and/or consult the Medtronic website at www.medtronic.eu.

For applicable products, consult instructions for use on manuals.medtronic.com. Manuals can be viewed using a current version of any major internet browser. For best results, use Adobe Acrobat® Reader with the browser.

Medtronic

Europe

Medtronic International Trading Sàrl.
Route du Molliau 31
Case postale
CH-1131 Tolochenaz
medtronic.eu
Tel: +41 0 21 802 70 00
Fax: +41 0 21 802 79 00

medtronic.eu

United Kingdom/Ireland

Medtronic Limited
Building 9
Croxley Green Business Park
Hatters Lane
Watford
Herts WD18 8WW
medtronic.co.uk
Tel: +44 0 1923 212213
Fax: +44 0 1923 241004

UC202411222 EE

CE0123

2024-sphere-9-trial-results-brochure-en-gb-emea-13016877
©2024 Medtronic. All rights reserved. **Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.