


# Clinical evidence summary

Three studies evaluated feasibility, safety, and time savings with the FlexCath Cross transseptal solution:


	Rizzi, et al. <sup>1</sup>	Yap, et al. <sup>2</sup>	Yong Ji, et al. <sup>3</sup>
Overview	<ul style="list-style-type: none"> <li>Nonrandomized study</li> <li>Single center</li> <li>40 patients</li> <li>AcQCross™ system vs. mechanical crossing with Abbott SL0™ sheath and dilator + BRK-1™ needle</li> </ul>	<ul style="list-style-type: none"> <li>Retrospective study</li> <li>Single center</li> <li>40 patients</li> <li>AcQCross system vs. RF crossing with Abbott SL1™ sheath and dilator + Bayliss NRG™ needle</li> </ul>	<ul style="list-style-type: none"> <li>Retrospective analysis</li> <li>Single, high-volume ablation center</li> <li>200 patients</li> <li>AcQCross system vs. mechanical crossing with Abbott SL0 sheath and dilator + BRK needle</li> </ul>
Results	<ul style="list-style-type: none"> <li>Mean transseptal puncture (TSP) time was <math>3.6 \pm 0.7</math> min in the AcQCross system group vs. <math>7.1 \pm 0.6</math> min in the mechanical needle group, <math>P &lt; 0.01</math></li> <li>49.9% relative reduction in the time for TSP with the AcQCross system vs. traditional mechanical needle</li> <li>No adverse events or peri- and postprocedural complications were noted</li> </ul>	<ul style="list-style-type: none"> <li>Mean time for cryoballoon delivery to the left atrium was <math>15.5 \pm 6.8</math> min in the AcQCross system group vs. <math>21.5 \pm 7.4</math> min in NRG RF needle group, <math>P = 0.01</math></li> <li>27.9% relative reduction in time for cryoballoon delivery with the AcQCross system vs. NRG RF needle</li> <li>Procedure time was <math>49.7 \pm 9.0</math> min in the AcQCross system group vs. <math>59.6 \pm 8.1</math> min in the NRG RF needle group, <math>P &lt; 0.001</math></li> <li>16.6% relative reduction in procedure time with the AcQCross system vs. NRG RF needle</li> <li>One temporary phrenic nerve injury in the AcQCross system group that resolved during the procedure</li> </ul>	<ul style="list-style-type: none"> <li>Mean TSP time was <math>9.2 \pm 9.5</math> min for the AcQCross system group vs. <math>16.3 \pm 23.9</math> min for the mechanical needle group</li> <li>43.6% relative reduction in the time for TSP with the AcQCross system vs. traditional mechanical needle</li> <li>No major safety issues were reported</li> </ul>

The FlexCath Cross system is feasible<sup>2</sup> and safe<sup>1-3</sup> and demonstrated improvement in TSP<sup>1,3</sup> and procedure<sup>2</sup> time.


- ✓ No major safety issues were reported<sup>1-3</sup>
- ✓ FlexCath Cross system is more efficient and predictable in time to TS crossing<sup>3</sup>
- ✓ FlexCath Cross system improves efficiency in cryoballoon procedures by reducing the number of exchanges<sup>2</sup>



Relative reduction in **TSP time** with FlexCath Cross system vs. BRK-1<sup>1</sup>



Relative reduction in **TSP time** with FlexCath Cross system vs. BRK<sup>3</sup>



Relative reduction in **procedure time** with FlexCath Cross system vs. NRG<sup>2</sup>

# Bench testing to evaluate coring with RF

## Objective:

To evaluate the risk of coring with FlexCath Cross when RF energy is used

## Experimental setup:

- Swine heart attached to fixture
- Sample submerged in water bath
- FlexCath Cross needle actuated to its maximum exposure while simultaneously pushing the sharp tip into the fossa ovalis of the swine heart and “activating” the ES generator
  - Used “cut mode” at 20 W for approximately two seconds (two-time safety margin used)
- Visually inspected the FlexCath Cross and in the water bath for loose swine heart material debris
- Distal end of the FlexCath Cross was then placed in a glass beaker filled with distilled water and the guidewire was deployed through the distal end to capture any particulate into the beaker; the beaker was visually inspected

## Results:

No loose biological particulate or coring observed (n = 59)

- Visual inspection of the distal tip of the dilator and exposed needle after crossing under 10x magnification resulted in no observations of particulate matter
- Visual inspection after each crossing of the exposed needle before and after pushing a guidewire out of the distal end resulted in no observations of coring
- Each test unit was used to cross the fossa ovalis a total of five (5) times. Visual inspection of the distal tip (n = 29) after each crossing under 10x magnification resulted in no physical damage observed as a result of continued RF energy use as compared to the pre-test condition.



## References

- <sup>1</sup> Rizzi S, Pannone L, Monaco C, et al. First experience with a transseptal puncture using a novel transseptal crossing device with integrated dilator and needle. *J Interv Card Electrophysiol*. December 2022;65(3):731-737.
- <sup>2</sup> Yap SC, Bhagwandien RE, Szili-Torok T. Use of a novel integrated dilator-needle system in cryoballoon procedures: a zero-exchange approach. *J Interv Card Electrophysiol*. November 2022;65(2):527-534.
- <sup>3</sup> Yong Ji S, et al. Use of the Novel AcQCross Transseptal System Improves Procedural Efficiency Compared to Conventional Transseptal Puncture Technique. Presented at AF Symposium 2022.

## FlexCath Cross™ Transseptal Solution

### Indications

Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events, and directions for use. FlexCath Cross™ devices are FDA cleared. Not all FlexCath Cross devices are CE marked.

710 Medtronic Parkway  
Minneapolis, MN 55432-5604  
USA

Toll-free: 800.328.2518  
Tel: +1.763.514.4000

[medtronic.com](https://www.medtronic.com)

©2023 Medtronic. All rights reserved. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic.  
™ Third-party brands are trademarks of their respective owners.  
All other brands are trademarks of a Medtronic company.

UC202400169 EN  
09/2023

**Medtronic**