


# Medtronic

## Get started with the MyCareLink Heart™ mobile app

A personalized  
remote monitoring  
solution that helps  
keep you connected  
and may increase  
your peace of mind




# Why should I choose app-based remote monitoring?

Today, millions of people implanted with heart devices are remotely monitored.

Remote monitoring may

- Give you peace of mind and increase quality of life
- Reduce hospital visits and save you time
- Help you feel assured as you stay connected to your care team

Flexibility



With connectivity information, transmission confirmation and reminders, the MyCareLink Heart™ app helps you to stay connected to your care team.



The app allows you to record symptoms. These symptoms are sent together with heart rhythm data to your clinic.

and connectivity



You can find information about living with your device, as well as resources for additional education.

# How does it work?

## Step 1. Send

Your heart device monitors your rhythm and automatically sends data to the clinic via an app. The clinic is notified when your device detects an irregular rhythm and when you record symptoms.





## Step 2. Review

Your clinic reviews the data on a secure website.



Please visit [medtronic.com/security](https://www.medtronic.com/security) for up-to-date security information.



// Now I feel more comfortable knowing that my device is being monitored! //




// This is so easy to understand. The app shows you the way, step by step. //

Images are not actual patients. Statements reflect individual experiences. Results may vary.

// It's a mobile app,  
so wherever I go in the  
world, it goes with me.  
I feel safe, secure and  
I've got freedom.  
This is brilliant! //



# App overview




## Home page

### Connectivity status

The home page displays the connectivity status between the app, your heart device and your clinic.

\* If prescribed by physician.





## Dashboard page

### ① **RECORD NEW SYMPTOMS**

Records symptoms and heart rhythm data at the time of a symptomatic episode.\*

### ② **Messages**

Shows important messages about your transmission status and using your app.

### ③ **Heart Device**

Displays implant date, ICM device name, model number and serial number.

### ④ **Symptom History**

Provides a log of your recorded symptoms.

### ⑤ **Education**

Provides answers to frequently asked questions about living with a heart device.

# How to record new symptoms


## Step 1.

Select **RECORD NEW SYMPTOMS**




## Step 2.

Record **how you feel**



### Step 3.



Record **your activity level**



9:41


← Back


RECORD SYMPTOMS


 If you have a medical emergency, call your local emergency number. 


What was your activity level?

Select one.

☐ Not at all active 

☐ A little active 

☐ Active 


☒ Very active 

NEXT

CANCEL

## Step 4.

Describe **what you were doing**




The record specific activity feature has a 250 character count limit per entry.



If you have a medical emergency, call your local emergency number.


## Step 5.

Remain in **app range**



## Step 6.

App confirms that **your symptoms have been recorded**




All clinical and patient data are fictitious and are used for demonstration purposes only.

# App install and setup

This process takes  
approximately 15 minutes.

## Step 1. Prepare

- Your smartphone or tablet with **Bluetooth®** turned **on**
- Your Apple App Store® or Google Play® store **ID** and **password**
- An email address you regularly use
- Your **heart device serial number**, which can be found on your medical device ID card



Place heart device  
serial number  
sticker here

## Step 2. Download and install


Scan the QR code or visit [MCLHeart.com/Download](https://MCLHeart.com/Download). If you don't see the **GET** button (on Apple devices) or **Install** button (on Android devices), your device is not compatible. Please contact your physician for other monitoring options.



### Step 3. Setup

Follow the instructions in the app to connect with your heart device. Here are a few **important actions** to keep in mind:

- Create a password as shown on the screen. You will be asked to enter it again in Step 4.
- The **pairing process** may take up to **6 minutes. Please don't interfere** when the app initiates contact with your heart device for the first time. Keep your **Bluetooth enabled at all times.**
- Click **Allow** to receive alerts and important communication from the app. This step can also be done in your smartphone/tablet setting.





## Step 4. Verify email

This last step will complete the setup process. Follow the **VERIFY EMAIL** link as shown, fill in the email address and password from Step 3, then tap **Submit**. The app may take 1 minute to refresh. Now, you are all set. You will be able to access the features of the app as soon as your clinic adds you to their system.

From: <[do-not-reply@crhfcloud.com](mailto:do-not-reply@crhfcloud.com)>

Date: Tue, Dec 21, 2021 at 2:46 PM

Subject: Verify Your Email Address For Your MyCareLink Heart™ Account

To: <[clk.testing.01@gmail.com](mailto:clk.testing.01@gmail.com)>

Sent to: [monica.smith@mail.com](mailto:monica.smith@mail.com)

Hello,

You are almost done with your Medtronic account setup—just verify your email address using the button below.

**VERIFY EMAIL**

# Things to remember

## To ensure you are connected to your clinic:

- Confirm that Bluetooth is turned on
- Stay connected to the internet via Wi-Fi or cellular data
- Keep your smartphone or tablet near you for a minimum of 3 hours throughout the day
- Keep your app open and running in the background

**You can quickly access educational resources in the Education section of the app.**





Do you have more questions  
about your heart device or patient  
monitoring solution?  
Call us, we are here to help!

**BeConnected**

**00800-266-632-82\***

Monday - Friday  
08.00 - 16.00\*\*

Alternative number from  
mobile phones or if travelling abroad  
+44 (0) 1923 202543\*\*\*

\*Free-of-charge number from a landline and some  
mobiles (please check with mobile service providers)

\*\* voicemail service available outside opening hours

\*\*\* number charged at standard call rate (may vary  
between phone providers)

# Frequently asked questions

## Multiple devices

### **Can I use the MyCareLink Heart™ mobile app on multiple mobile devices?**

Yes. Before using a second mobile device, you will need to have paired your first mobile device. You can then use your username and password to log in on your second mobile device. Caregivers can have the app on their phones, too!

## Bluetooth

### **Will my heart device interact with other Bluetooth devices? For example, will it interact with the Bluetooth in my car?**

Your heart device has Bluetooth Low Energy, while cars typically use Bluetooth Classic. Cars will not attempt to establish a connection with your heart device and your heart device will not be able to detect Bluetooth Classic transmissions.

## Battery

### **Will the MyCareLink Heart™ mobile app drain my phone battery?**

If you normally have Bluetooth turned on, you shouldn't experience any change in the battery life of your mobile device. If you do not normally keep Bluetooth turned on, you may notice a decrease in the battery life of your mobile device while using the app between charges.

## App's data consumption

### **How much cellular data does the MyCareLink Heart™ mobile app use each month?**

The app's data consumption (15 MB/month) is equivalent to about two minutes of web surfing per day.

## Security

### **If my mobile device is stolen, what will someone see in the MyCareLink Heart™ mobile app?**

Just like other apps, if your mobile device has Wi-Fi or cellular connection, the data on the app will be available for a viewer to see. It is recommended that you protect your mobile device with a PIN code, fingerprint or other security features.

## Updates

### **Mobile device/operating system compatibility**

Mobile device and operating system (OS) requirements to support the MyCareLink Heart™ mobile app will change over time to ensure reliable monitoring. You may need to update or replace your mobile device or OS to use the app to transfer data between your heart device and the Medtronic CareLink™ network.

Please visit [MCLHeart.com](http://MCLHeart.com) for a list of compatible devices.

Information contained herein is not medical advice and should not be used as an alternative to speaking with your doctor. Discuss indications, contraindications, warnings, precautions, potential adverse events and any further information with your healthcare professional.

# Medtronic

Europe  
Medtronic International Trading Sàrl.  
Route du Molliau 31  
Case postale  
CH-1131 Tolochenaz  
[www.medtronic.eu](http://www.medtronic.eu)  
Tel: +41 (0)21 802 70 00  
Fax: +41 (0)21 802 79 00

United Kingdom/Ireland  
Medtronic Limited  
Building 9  
Croxley Park  
Hatters Lane  
Watford  
Herts WD18 8WW  
[www.medtronic.co.uk](http://www.medtronic.co.uk)  
Tel: +44 (0)1923 212213  
Fax: +44 (0)1923 241004

UC202006845a EE © Medtronic 2022.  
All rights reserved.