

Medtronic

Not all porcine tissue valves are designed the same.

Hemodynamic performance is impacted.



The Mosaic™ bioprosthesis was designed to allow for maximum blood flow. The creation of the valve – from the porcine leaflets and valve geometry (specifically stent height), to how the tissue is mounted in the frame – affects how optimally the blood flows in diastole, and how well the leaflets coapt in systole.

The Mosaic mitral valve demonstrates low gradients and large EOAs.



Clinical study hemodynamics at one year from respective clinical trials^{1,2}

Epic™/Epic Plus mitral ¹					
Size	25	27	29	31	33
Mean gradient (mm) (n =)	NA	6.1 ± 2.9 (30)	5.5 ± 1.7 (41)	4.8 ± 1.4 (26)	4.1 ± 1.6 (24)
Mean EOA (cm ²) (n =)	NA	1.4 ± 0.7 (16)	1.5 ± 0.5 (26)	1.6 ± 0.3 (15)	1.5 ± 0.3 (22)

Mosaic mitral ²					
Mean gradient (mm) (n =)	5.7 ± 1.7 (41)	4.6 ± 2.1 (98)	4.4 ± 1.8 (123)	3.7 ± 1.4 (50)	3.4 ± 1.8 (8)
Mean EOA (cm ²) (n =)	1.6 ± 0.4 (35)	1.7 ± 0.5 (90)	1.8 ± 0.5 (114)	1.7 ± 0.5 (43)	1.9 ± 0.5 (6)

The charts are not intended to be a comparison of the two devices, as there is no head-to-head clinical study, but rather are intended to illustrate the clinical results of the two similar trials. Multiple factors contribute to clinical study outcomes and need to be considered in making any assessments across different studies.

Epic™ tissue mounting technique reduces the space available for flow.³

The tissue mounting technique used for the Epic valve places part of each cusp's aortic wall (Figure 1) into the stent lumen – important space needed for flow and reducing the potential “true” inside diameter (ID). (Figure 2)




Figure 1
Retained excess
aortic wall



Figure 2
Retained excess
aortic wall
protruding into
flow area

Impact of stent height on hemodynamics

The Evin, et al. study⁴ was conducted to establish in vitro nominal reference EOA values and transvalvular gradients independent of manufacturers.

A key finding was that EOAs increased and gradients decreased with increasing bioprosthesis height. EOAs and gradients were also impacted by prosthesis diameter, stent diameter, and heart rate.

Figure 1






Figure 2



¹ Instructions for Use: Stented Porcine Tissue Valves Epic™/Epic™ Supra. St. Jude Medical, Inc. Available at: https://www.accessdata.fda.gov/cdrh_docs/pdf4/P040021S004c.pdf. Accessed April 20, 2023.

² Mosaic Porcine Bioprosthesis. Instructions for Use. Medtronic, Inc. 2018. 1220016001 Rev. 1C.

³ Data on file.

⁴ Evin M, Magne J, Grieve S, Rieu R, Pibarot P. Characterization of Effective Orifice Areas of Mitral Prosthetic Heart Valves: An In-Vitro Study. *J Heart Valve Dis*. November 2017;26(6):677-687.

This material should not be considered the exclusive source of information, it does not replace or supersede information contained in the device manual(s).

Please note that the intended use of a product may vary depending on geographical approvals.

See the device manual(s) for detailed information regarding the intended use, the (implant) procedure, indications, contraindications, warnings, precautions, and potential adverse events.

For a MRI compatible device(s), consult the MRI information in the device manual(s) before performing a MRI.

If a device is eligible for eIFU usage, instructions for use can be found at Medtronic's website manuals.medtronic.com.

Manuals can be viewed using a current version of any major internet browser. For best results, use Adobe Acrobat® Reader with the browser.

Medtronic products placed on European markets bear the CE mark and the UKCA mark (if applicable).

For any further information, contact your local Medtronic representative and/or consult Medtronic's websites.

Medtronic

Europe
Medtronic International Trading Sàrl.
Route du Molliou 31
Case postale
CH-1131 Tolochenaz
Tel: +41 (0)21 802 70 00
Fax: +41 (0)21 802 79 00

[medtronic.eu](https://www.medtronic.eu)

2025-mosaic-mitral-lifetime-management-series-hemodynamic-brochure--
emea-20004655 ©2025 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. TM* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.