


# Medtronic

## Esoflip™ balloon dilation catheter

# Standard protocol quick reference guide



### Esoflip™\* ES-310

1. Placed trans-orally during sedated upper endoscopy
2. Connect a pressure monitor with a range of at least 0 to 5ATM (507kPa) to the Balloon Inflation Monitor Port. See statement from product labeling. Do not exceed the maximum inflation pressure of 1.5 ATM (152 kPa)  
*8 sensors 0.5 cm spaced*
3. Positioned with Stricture waist mid-balloon
4. Guide wire port in tip (0.035")  
*Catheter markings referenced from center of balloon*  
*Initial fill 5 ml to measure Stricture diameter pre-dilation*
5. Inflate to desired dilation diameter

### Precautions

- Balloon Rated Burst Pressure (RBP): **3 ATM (304 kPa)**
- Maximum inflate volume: **16 ml**  
Do not exceed the maximum inflation pressure of 1.5 ATM (152 kPa)

### Esoflip™\* ES-320

1. Placed trans-orally during sedated upper endoscopy
2. Connect a pressure monitor with a range of at least 0 to 5ATM (507kPa) to the Balloon Inflation Monitor Port. See statement from product labeling. Do not exceed the maximum inflation pressure of 2ATM (203kPa).  
*16 sensors 0.5 cm spaced*
3. Positioned with Stricture waist mid-balloon
4. Guide wire port in tip (0.035")  
*Catheter markings referenced from center of balloon*  
*Initial fill 20 ml to measure Stricture diameter pre-dilation*
5. Inflate to desired dilation diameter

### Precautions

- Balloon Rated Burst Pressure (RBP): **4 ATM (405 kPa)**
- Maximum inflate volume: **42 ml**  
Do not exceed the maximum inflation pressure of 2ATM (203kPa).

### Esoflip™\* ES-330

1. Placed trans-orally during sedated upper endoscopy
2. Connect a pressure monitor with a range of at least 0 to 5ATM (507kPa) to the Balloon Inflation Monitor Port. See statement from product labeling. Do not exceed the maximum inflation pressure of 1.5 ATM (152 kPa).  
*14 sensors 0.5 cm spaced*
3. Positioned with EGJ sphincter waist mid-balloon
4. Guide wire port in tip (0.035")  
*Catheter markings referenced from center of balloon*  
*Initial fill 30 ml to measure EGJ sphincter diameter pre-dilation*
5. Inflate to desired dilation diameter

### Precautions

- Balloon Rated Burst Pressure (RBP): **3 ATM (304 kPa)**
- Maximum inflate volume: **75 ml**  
Do not exceed the maximum inflation pressure of 1.5 ATM (152 kPa).

### Indications for use

The Esoflip™ ES-310 and the ES-320 balloon catheter are indicated for use to dilate esophageal strictures due to esophageal surgery, primary gastro-esophageal reflux and radiation therapy. The Esoflip™ ES-330 balloon dilation catheter is indicated for use in a clinical setting for dilating the gastroesophageal junction of a patient with Achalasia. **Note:** The Esoflip™ catheter is to be used only with the Endoflip™ system.

### Contraindications

The Esoflip™ catheter is contraindicated where endoscopy is contraindicated. Do not use the Esoflip™ catheter on patients with actively-bleeding varices in the esophagus or with active esophageal perforation. The Esoflip™ ES-310 catheter is not suitable for diameter measurements and dilation of strictures smaller than 6 mm or greater than 10 mm. The Esoflip™ ES-320 catheter is not suitable for diameter measurements smaller than 8 mm or greater than 20 mm. The Esoflip™ ES-330 catheter is not suitable for diameter measurements smaller than 8 mm or greater than 30 mm.

### Potential complications

Potential complications include allergic reaction, anaphylaxis, bleeding, cardio-respiratory complications, dental trauma, infection, pain, perforation, pulmonary aspiration, vasovagal response.

### References

1. Esoflip™ ES-310 Catheter Instructions for Use, PT00140542 Rev A 2021-02
2. Esoflip™ ES-320 Catheter Instructions for Use, PT00140523 Rev A 2021-02
3. Esoflip™ ES-330 Catheter Instructions for Use, PT00140543 Rev A 2021-02

Important: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2025 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.  
esoflip-standard-protocol-qrg-18934273

[medtronic.com/covidien/uk](https://www.medtronic.com/covidien/uk)

\*If the catheter pre-use check procedure detects a problem, a message notifies you that the catheter may be faulty and should be replaced.