

# Medtronic


Engineering the extraordinary

## HEALTH SYSTEM & PAYER VALUE


SynchroMed™ II Drug Infusion System

# Targeted drug delivery for refractory cancer pain

Together, let's help patients  
take control of their  
cancer-related pain







# Cancer pain is a significant problem - And it's growing

Pain persists throughout the cancer lifecycle, reported in:

- 55% of patients undergoing active treatment;
- 39% of patients with cancer in remission;
- 66% of patients with advanced, metastatic, or terminal disease<sup>1</sup>.

**Refractory pain is suffered by people who have inadequate pain control, despite optimized use of systematic analgesics<sup>3</sup>.**

Cancer pain can be episodic with no pain in between or with lesser background pain.

**With better treatments, cancer patients are surviving longer - however chronic pain can persist in patients many years from diagnosis<sup>4-6</sup>.**

Refractory cancer pain

Prevalence

Breakthrough pain


Burden of cancer pain

Current treatment management

Intrathecal drug delivery

References

About  
**25%**  
of cancer patients  
die in pain<sup>2</sup>




Breakthrough cancer pain

[Click for more](#)


**Cancer pain has a high incidence in the cancer population<sup>9</sup>**

Prevalence of pain by cancer type (% of patients)



[Click to enlarge](#)

Medtronic



Refractory cancer pain

Prevalence

Breakthrough pain


Burden of cancer pain

Current treatment management


Intrathecal drug delivery

# Cancer pain has a high incidence in the cancer population<sup>9</sup>

Prevalence of pain by cancer type (% of patients)



References




# Breakthrough cancer pain

Breakthrough cancer pain (BTcP) is a complex, heterogenous symptom described as a transitory flare of pain in the setting of chronic pain<sup>10,11</sup>.

In a European cohort of 1,000 cancer pain patients, a daily median of three BTcP episodes was reported, with the majority (61.8%) describing their BTcP as 'severe'<sup>12</sup>.

**BTcP Intensity: European cohort of 1,000 patients**




Study Demographics: 51% male, median age of 62 years. Tumor location: gastrointestinal (26.4 %), lung (17.2%), urological (16%), breast (12.5 %), other body locations (23.9%), unknown (1.6%), not stated (2.4%)<sup>12</sup>.

A 2018 study from 17 centers\* across 16 provinces throughout Spain examined the impact of BTcP and found<sup>13</sup>:

- Patients experienced a mean of 3.1 episodes of BTcP/day (mean duration of 30.6 min<sup>13</sup>).
- 20% described their BTcP intensity as "unbearable"<sup>13</sup>.
- Almost two-thirds indicated that their BTcP was "unpredictable" in its occurrence<sup>13</sup>.


**Intensity of BTP in cancer patients (N=152)**



Study demographics: 65.8% male with mean age of 66.8 years. Tumor location: gastrointestinal tract (23%), lung (22.4%), breast (9.2%), prostate gland (5.3%), other body locations (40.1%)<sup>13</sup>.

\*8 pain units; 8 palliative care units;  
1 oncology department<sup>13</sup>





## Quality of life impact

### Cancer patients perceptions of their pain:

67%

describe it  
as distressing<sup>9</sup>

32%

describe their pain  
as so bad they wish  
they would die<sup>9</sup>

36%

describe their pain  
as an unbearable aspect  
of their disease<sup>9</sup>

Pain is often the most  
tangible sign of disease  
that cancer patients and  
their families perceive<sup>7</sup>

Refractory  
cancer pain

Burden of  
cancer pain

Patient impact


Treatment barriers

Cost impact

Current  
treatment  
management

Intrathecal  
drug delivery


References



# Cancer pain has a detrimental impact on patients & their families

Cancer pain contributes to impaired performance of daily activities, disability, negatively affects mental health, and reduces health-related quality of life<sup>9</sup>.

The 10 most common symptoms experienced by colorectal cancer (CRC) survivors with pain are almost twice as prevalent in comparison to those without pain who report these symptoms<sup>14</sup>.



Refractory cancer pain

Burden of cancer pain

Patient impact

Treatment barriers


Cost impact

Current treatment management

Intrathecal drug delivery

References





# Barriers to seeking cancer pain analgesia

Barriers to seeking adequate cancer pain relief include:

- **Reluctance** to report pain, due to the belief that healthcare professionals must be focused on the cancer itself<sup>15</sup>;
- **Beliefs** that pain is an inherent part of cancer and fears that higher pain levels, increased analgesia needs, and the symbolic significance of opioids indicate disease progression, deterioration and approaching death<sup>16</sup>;
- **Fears** of side effects and addiction associated with perceptions of opioids use<sup>15</sup>.

“... I don't want to let them see I'm in pain, so I have to hide it.”

Patient-reported extract from a study of narrative interviews with severe chronic cancer pain patients, from a single centre in Italy<sup>23</sup>



Family members of end-of-life cancer patients report anxiety about giving correct and timely dosages of analgesia, and concerns about keeping the patient comfortable without causing harm<sup>16</sup>.

Refractory cancer pain

Burden of cancer pain

Patient impact

Treatment barriers


Cost impact

Current treatment management

Intrathecal drug delivery

References





# The cost of unmanaged cancer pain

## High dose opioids

- Cancer pain patients prescribed conventional medical management (CMM) with strong opioids have an increased prevalence of opioid-induced constipation (OIC) or bowel dysfunction, despite the use of laxatives<sup>17,18</sup>.
- 2017 data from a Danish national registry found total healthcare costs were 25% higher for cancer patients with OIC vs. without OIC\*<sup>19</sup>.

\* adjusted for age, gender, opioid usage, marital status, (p<0.001)

## Ability to work

- 78% of cancer pain patients reported inability to work, compared to 40% of cancer patients without pain<sup>14</sup>.
- A Norwegian study found increased pain to be significantly associated with long-term cervical cancer survivors holding disability pensions, compared to those holding paid work<sup>20</sup>.

Refractory  
cancer pain

Burden of  
cancer pain

Patient impact

Treatment barriers


Cost impact

Current  
treatment  
management

Intrathecal  
drug delivery

References






# The cost of unmanaged cancer pain

## Breakthrough cancer pain

- The monthly cost per breakthrough cancer (BTcP) patient in Spain was estimated at €2,941 (2017 values)<sup>13</sup>.
- Assuming the annual distribution of cost is proportional, this translates to a cost of €35,000/year/patient for the management of BTcP alone.

## Monthly cost drivers for a BTcP patient




\* primary care, specialist physician, emergency unit, day hospital, visit by home hospitalization, radiotherapy session;

§ nonhealthcare resources, psychotherapy, physiotherapy, caregiver (paid and unpaid);

¥ patient and caregiver leave<sup>13</sup>.

Multicentre study (17 centers) - 8 pain units; 8 palliative care units; 1 oncology department - across 16 provinces throughout Spain (2017)<sup>13</sup>






# Management of cancer pain is frequently suboptimal<sup>21,22</sup>

## With undertreatment for those with moderate to severe pain<sup>9,24-27</sup>

- Since 1996, The World Health Organization's (WHO) 3-step analgesic ladder (Figure 3, left) has been the stalwart of cancer pain management<sup>28</sup>.
- However, in real-life clinical practice, the likelihood of clinical success if the WHO steps are followed is low (12%, 74% and 14% at Steps 1, 2 and 3, respectively)<sup>29</sup>
- For patients at stage 3 (with moderate-severe pain), 25% do not respond to recommended treatment, or are poor responders to pain relief with transdermal and oral opioids<sup>30</sup>.



A. WHO 3-step analgesic ladder reported in 1996. Adapted from<sup>28</sup>.

Proposed  
4<sup>th</sup> step

Click  
for more

With the rise in cancer survival rates<sup>4-6</sup> the treatment of cancer pain is shifting from short-term analgesia to long-term pain management<sup>8</sup>. Yet, there is little evidence of the safety or long-term efficacy of opioid therapy in cancer patients with chronic pain<sup>33</sup>.

Refractory  
cancer pain

Burden of  
cancer pain

Current  
treatment  
management

WHO 3-step  
ladder


Impact of high  
dose opioids

Inequality  
of access

Intrathecal  
drug delivery

References






# Management of cancer pain is frequently suboptimal<sup>21,22</sup>

## With undertreatment for those with moderate to severe pain<sup>9,24-27</sup>

- As the management of cancer pain is now considered to be more complex, the WHO ladder is under review for adults<sup>34</sup>, and has been abandoned for the management of children<sup>35</sup>.
- A fourth step, including interventional techniques, has been proposed (Figure B, right)<sup>31</sup>.




A. WHO 3-step analgesic ladder reported in 1996. Adapted from<sup>28</sup>.

B. The proposed 4<sup>th</sup> step, as reported by Miguel et al, 2000<sup>31</sup>, Adapted from<sup>31</sup>.

In the last decade, intrathecal drug delivery has been recognized as a treatment option for cancer patients in whom the expected clinical outcome of systemic treatments is not obtained, or who present intolerance to low doses of these treatments<sup>8,9</sup>.






# Management of cancer pain is frequently suboptimal<sup>12,13</sup>

## High dose opioids: patient impact

Evidence supports pain as a key prognostic variable associated with a shorter survival in cancer patients<sup>36-39</sup>. Furthermore, high level opioid use is significantly associated with shorter overall survival<sup>40</sup>.

High dose opioids prescribed for cancer pain have well-established adverse events<sup>41</sup>:

- Increased prevalence of opioid-induced constipation or bowel dysfunction (OIC, OIBD) despite the use of laxatives - which can negatively impact on cancer patients' already reduced quality of life<sup>17,42</sup>.
- OIC may be indicative of poor opioid tolerability, which could lead patients lowering their opioid dose, resulting in inadequate pain control<sup>42</sup>.



Refractory cancer pain

Burden of cancer pain

Current treatment management

WHO 3-step ladder


Impact of high dose opioids

Inequality of access

Intrathecal drug delivery

References






# Inequalities of access to optimized cancer pain treatments

A 2018 survey of European healthcare professionals and cancer patients suffering from breakthrough cancer pain (BTcP) found significant time lags of (58 and 13 weeks in Germany and Spain, respectively) between diagnosis of BTcP and initiation of specific treatment.<sup>44</sup>

## Time with cancer pain (weeks)



Refractory cancer pain

Burden of cancer pain

Current treatment management


WHO 3-step ladder

Impact of high dose opioids

Inequalities of access


Intrathecal drug delivery

References



# Inequalities of access to optimized cancer pain treatments

Multidisciplinary management and collaboration between oncologists and other specialists remains uncommon.<sup>44</sup>



\* Nationwide survey of medical oncologists in Spain, carried out through a self-administered, written questionnaire (conducted in two waves over Sept-Oct 2015 and Dec 2015-Jan 2016. Total of 73 and 82 oncologists participated in the first and second wave, respectively. The study participants were from 15 of the 17 Autonomous Communities of Spain. All respondents (155 in total) fulfilled the questionnaires completely.

Refractory cancer pain

Burden of cancer pain

Current treatment management

WHO 3-step ladder


Impact of high dose opioids

Inequalities of access

Intrathecal drug delivery

References





# Inequalities of access to optimized cancer pain treatments

Access to effective palliative pain management for cancer patients is a European priority<sup>41</sup>.

**In the 2018 Resolution for the Provision of Palliative Care in Europe<sup>45</sup>**, The Council of Europe (the continent's leading human rights organization) calls on member states to:

"ensure access to pain treatment and management as a crucial component of palliative care.... including address educational and attitudinal barriers by raising awareness of appropriate and effective pain management, including opioid-based treatments, among health-care professionals and the general public"<sup>45</sup>

- Integration of palliative medicine in oncology improves patient outcomes and decreases healthcare costs<sup>46</sup>.
- Early involvement of palliative medicine after the cancer diagnosis is supported by national guidelines <sup>46</sup>.



**Multi-disciplinary collaboration between Oncology, Pain Management, Nursing & Palliative Care are essential to achieving effective and timely pain relief for cancer pain patients and survivors<sup>43,47</sup>.**

Refractory cancer pain

Burden of cancer pain

Current treatment management

WHO 3-step ladder


Impact of high dose opioids

Inequalities of access

Intrathecal drug delivery

References





# The benefits of intrathecal targeted drug delivery for cancer pain patients


- A lower morphine equivalent dose is needed via intrathecal application<sup>48</sup>.
- Targeted drug delivery (TDD) via the intrathecal route offers significant benefits to cancer pain patients who are not optimally managed on oral or transdermal analgesia. The smaller doses may reduce the systemic side effects<sup>49,50</sup>.

Route of Administration	Conversion Ratio
Oral	300
Intravenous	100
Epidural	10
Intrathecal	1

**Morphine equivalent dose across different routes of morphine administration<sup>48</sup>**


## Intrathecal

Targeted drug delivery (TDD) directly to the cerebral fluid



## Oral

Systemic delivery through the circulatory system



Refractory cancer pain

Burden of cancer pain

Current treatment management

Intrathecal drug delivery

Clinical evidence

Economic evidence

Therapy guidelines


Pump benefits

Patient selection

SynchroMed™ II drug infusion system

References





# The procedure: clinical considerations

## Trial implant and refill management<sup>62</sup>

Refractory  
cancer pain

Burden of  
cancer pain

Current  
treatment  
management

Intrathecal  
drug delivery

Clinical  
evidence

Economic  
evidence

Therapy  
guidelines

Pump benefits

Patient  
selection

Synchromed™ II  
drug infusion  
system


References

### Trial

- Simple lumbar puncture
- Medications delivered via needle or catheter
- Evaluate therapy efficacy


### Implant


- Time: 1-2 hours
- May be performed under general/local MAC
- Two incisions: Catheter **1** and Pump **2**
- Relative size of the pump



### Refill

- 22 gauge needle
- No sedation required






# Clinical evidence supports TDD for chronic cancer pain


## International, multicenter RCT assessing TDD + conventional medical management

The safety and efficacy of targeted drug delivery (TDD) for the treatment of cancer-related pain have been demonstrated in randomized controlled clinical trials<sup>49,50</sup>.

### More patients with increased pain control and reduced toxicity<sup>49</sup>.



### % of patients with a reduction in both pain & toxicity



Study Details: International, multicenter RCT (Evidence Level II), assessing TDD + Conventional Medical Management (CMM) for treating chronic cancer pain vs. CMM. A total of 202 patients (Randomization: 101 in TDD arm; 99 in CMM arm) were enrolled at 21 centers (16 in the United States, 4 in Europe, and 1 in Australia). Clinical success was defined as >20% reduction in VAS scores, or equal scores with > 20% reduction in toxicity<sup>49,50</sup>.

Refractory cancer pain

Burden of cancer pain

Current treatment management

Intrathecal drug delivery

Clinical evidence

Economic evidence

Therapy guidelines

Pump benefits


Patient selection

SynchroMed™ II drug infusion system

References








# Clinical evidence supports TDD for chronic cancer pain

## International, multicenter RCT assessing TDD + conventional medical management

More patients with less fatigue and reduces consciousness<sup>49</sup>.

Reduction in individual toxicities from baseline to 4-week follow-up.



Study Details: International, multicenter RCT (Evidence Level II), assessing TDD + Conventional Medical Management (CMM) for treating chronic cancer pain vs. CMM. A total of 202 patients (Randomization: 101 in TDD arm; 99 in CMM arm) were enrolled at 21 centers (16 in the United States, 4 in Europe, and 1 in Australia). Clinical success was defined as >20% reduction in VAS scores, or equal scores with > 20% reduction in toxicity<sup>49,50</sup>.

Refractory cancer pain

Burden of cancer pain

Current treatment management

Intrathecal drug delivery

Clinical evidence

Economic evidence

Therapy guidelines


Pump benefits

Patient selection

SynchroMed™ II drug infusion system

References





# Clinical evidence supports TDD for chronic cancer pain

## International, multicenter RCT assessing TDD + conventional medical management

The safety and efficacy of targeted drug delivery (TDD) for the treatment of cancer-related pain have been demonstrated in randomized controlled clinical trials<sup>49,50</sup>.

### Safety outcomes

- The 194 serious adverse events (SAEs) reported were evenly distributed between the two study groups; 49% were reported in the CMM group and 51% in the TDD+CMM group<sup>49</sup>.
- The frequency of two common adverse events associated with opioid use, fatigue and depressed level of consciousness, was significantly lower in the TDD+CMM group compared to the CMM group from baseline to 4 weeks ( $p < 0.05$ )<sup>49</sup>.

Study Details: International, multicenter RCT (Evidence Level II), assessing TDD + Conventional Medical Management (CMM) for treating chronic cancer pain vs. CMM. A total of 202 patients (Randomization: 101 in TDD arm; 99 in CMM arm) were enrolled at 21 centers (16 in the United States, 4 in Europe, and 1 in Australia). Clinical success was defined as >20% reduction in VAS scores, or equal scores with > 20% reduction in toxicity<sup>49,50</sup>.

Refractory  
cancer pain

Burden of  
cancer pain

Current  
treatment  
management

Intrathecal  
drug delivery

Clinical  
evidence

Economic  
evidence

Therapy  
guidelines


Pump benefits

Patient  
selection

SynchroMed™ II  
drug infusion  
system

References





# Clinical evidence supports TDD for chronic cancer pain

## Systematic Review & Meta-Analysis: Management of intrathecal drug delivery (Perruchoud et al. 2021)<sup>32</sup>

### Objectives

To examine the efficacy of managing cancer-related pain with IDD with external pump or implanted infusion systems. Secondary objectives included the effects of IDD on systemic opioid use (oral morphine equivalent [OME]) and infection rates.

### Results

Twenty-nine studies, 17 for calculating pain levels and 13 for weighted mean morphine dose, were identified.

Pain levels significantly decreased postintervention from baseline. Mean differences (on a 0 to 10 scale) were:

- -4.34 ( $p < 0.001$ ) at 4 to 5 weeks (short-term)
- -4.34 ( $p < 0.001$ ) at 6 to 12 weeks (mid-term)
- -3.32 ( $p < 0.001$ ) at >6 months (long-term)

Mean systemic opioid use (OME) was reduced by 308.24 (SE = 22.72) mg/d. (54% reduction of OME).

Mean infection rates were ~3% for intrathecal pumps, which is similar to that found in other populations treated with IDD implanted systems.

Survival - The pooled weighted mean survival time was 130.68 (SE = 1.50) days. This may have been possible due to reduced toxicity and a general improvement in the patients' QoL post-IDD.

### Study Design

A systematic literature search according to the Preferred Reporting Items for Systematic Reviews and Meta-Analyses (PRISMA) guidelines.

### Limitations

Meta-analysis could not be performed based on control group values because of the heterogeneous designs of the selected studies.

A substantial intertrial heterogeneity was found, as anticipated from the variability in study designs.

Complications, adverse events, and drug side effects were not consistently collected and/or analyzed across selected studies.

### Key Insights

This 2021 published meta-analyses showed a statistically significant and sustained decrease in cancer pain with IDD, compared with baseline.

Systemic opioid consumption was reduced on average by >50% after IDD. Infection rates were comparable with other indications.

Refractory  
cancer pain

Burden of  
cancer pain

Current  
treatment  
management

Intrathecal  
drug delivery

Clinical  
evidence

Economic  
evidence


Therapy  
guidelines

Pump benefits

Patient  
selection

Synchromed™ II  
drug infusion  
system

References



# Clinical evidence supports TDD for chronic cancer pain

## Effectiveness and Safety of Intrathecal Drug Delivery Systems for the Management of Cancer Pain: A Systematic Review and Meta-Analysis (Duarte et al. 2022)<sup>51</sup>

### Objectives

To evaluate the effectiveness and safety of IDD and SCS for cancer pain.

### Results

Twenty-two studies were included on either IDD or SCS for cancer pain. Eight studies were included in the meta-analysis of pain intensity. Pain levels significantly decreased postintervention from baseline.

Mean differences (on a 0 to 10 scale) were:

- $-3.31$  ( $p < 0.001$ ) at the latest posttreatment follow-up time
- $-3.53$  ( $p < 0.001$ ) up to one month after treatment.

Studies with either an IDDS or an SCS device showed similar results.

Improvements were also observed in survival, HRQoL or functional outcomes, and use of systemic opioids.

Postdural puncture headache was the most reported complication, whereas urinary retention, nausea, and vomiting were commonly reported side effects.

### Study Design

A systematic review method according to the general principles outlined in the Centre for Reviews and Dissemination (CRD) guidance for conducting reviews in health care and reported in accordance with the Preferred Reporting Items for Systematic Reviews and Meta-Analyses (PRISMA).

### Limitations

- Substantial unexplained heterogeneity was present within the meta-analyses, meaning that the magnitude of pain intensity reduction with IDDS therapy is uncertain.
- Limited evidence was identified evaluating the effectiveness of SCS for patients with cancer pain.

### Key Insights

This 2022 published meta-analyses showed a statistically and clinically significant reduction in pain intensity up to one month and the latest posttreatment follow-up, compared with baseline.

The evidence also suggests that IDD results in improvements in other outcomes, such as survival, HRQoL or functional outcomes, and use of systemic opioids.

Refractory cancer pain

Burden of cancer pain

Current treatment management

Intrathecal drug delivery

Clinical evidence

Economic evidence

Therapy guidelines


Pump benefits

Patient selection

SynchroMed™ II drug infusion system

References





# Budget impact of ziconotide through IDD for cancer pain

## Ziconotide for the Management of Cancer Pain: A Budget Impact Analysis (Lambe et al. 2022)<sup>52</sup>

### Objectives

To conduct a budget impact analysis from the NHS perspective over a five-year time horizon that reflects a patient population receiving IDD ziconotide or morphine therapy for cancer pain in England.

### Study Design

Markov-like analytic modelling to estimate the budget impact of ziconotide monotherapy vs morphine monotherapy through intrathecal drug delivery (IDD) for the management of cancer pain. The perspective adopted was that of the UK National Health Service, with a five-year time horizon.


Patients requiring an ITDD were expected to incur the device cost, device implantation cost, refill procedure and the drug costs. The time horizon was five years.

### Results

The total costs of ziconotide monotherapy and morphine monotherapy for the first five years are shown in Figure 1. The estimated five-year cumulative budget impact of treatment with ziconotide monotherapy for the five-year time horizon was £2,355,675 whereas that of morphine monotherapy was £913,804.

The results of this study suggest that the use of ziconotide is associated with higher costs to the health care system. However, the additional costs in any of the first five years are below the resource impact significance level of £1 million for medical technologies in England.

### Budget impact of Morphine vs Ziconotide monotherapy



### Key Insights

This budget impact analysis suggest that although intrathecal ziconotide is associated with higher costs to the health care system in England, the incremental costs are not significant.

Routine commissioning of ziconotide would provide an alternative for a population with limited IDD treatment options.

Refractory cancer pain

Burden of cancer pain

Current treatment management

Intrathecal drug delivery

Clinical evidence

Economic evidence


Therapy guidelines

Pump benefits

Patient selection

SynchroMed™ II drug infusion system

References




# Clinical evidence supports TDD for chronic cancer pain

## Real world clinical effectiveness

A 2019 publication based on data from a prospective, long-term multicenter registry (Medtronic's Product Surveillance Registry) report on the **largest cohort of TDD cancer pain patients to date** (1,403 patients included).


### Significant Reduction of Pain at 6- and 12-months post-implant<sup>53</sup>

Paired pain scores




### Improvement in Health-related Quality of Life at 6 months, and maintained at 12 months<sup>53</sup>

Paired EQ-5D utility scores



Results support the benefits of TDD as an effective therapeutic option with a positive benefit-risk ratio in the treatment of chronic cancer pain<sup>53</sup>.





# Clinical evidence supports TDD for chronic cancer pain

## Real world safety

- Infection rates following TDD system implants in cancer pain patients are low, ranging from 0.9%\* to 3.2%\*\*


Medtronic actively tracks and annually publishes the performance reports for the Synchromed II™ drug infusion system.<sup>66</sup>

Find these available at:  
<https://professional.medtronic.com/ppr>



\* Retrospective chart review of 217 TDD implants from a single center. Surgical Site infection risk after TDD placement for cancer pain is low, despite frequent concurrent antineoplastic therapy and leukopenia in the perioperative period.<sup>54</sup>

\*\* Prospective, long term, multicenter registry on 1,403 patients suffering from cancer-related pain implanted with TDD. 3.2% experienced infection requiring surgical intervention (e.g., explant, replacement, revision, or debridement).<sup>53</sup>




# Economic value of TDD for cancer pain


## Compared to conventional medical management

Recent published evidence has shown that compared to CMM alone, TDD+CMM significantly reduced healthcare utilization as early as 2 months and through 12 months ,in the form of lower inpatient visits, inpatient days and emergency department visits<sup>55</sup>.


### Inpatient visits



### Emergency room visits



### Total Inpatient days (across all admissions)



● TDD + CMM ● CMM

- From a separate study, in selected patients prescribed high-cost conventional opioid regimens, TDD management of cancer-related pain resulted in cost savings as early as 7.6 months compared with conventional opioid therapy, from a US payer perspective<sup>56</sup>.

Refractory  
cancer pain

Burden of  
cancer pain

Current  
treatment  
management

Intrathecal  
drug delivery

Clinical  
evidence

**Economic  
evidence**

Therapy  
guidelines


Pump benefits

Patient  
selection

SynchroMed™ II  
drug infusion  
system

References





Refractory cancer pain

Burden of cancer pain

Current treatment management

Intrathecal drug delivery

Clinical evidence

Economic evidence

Therapy guidelines

Pump benefits

Patient selection


Synchromed™ II drug infusion system

References

# Recommendations:

## Intrathecal drug delivery for cancer-related pain

Guideline and society (year)	Recommendations for intrathecal drug delivery
<b>European Pain Federation (EFIC) (2018)<sup>58</sup></b>	<p>Patients should be referred to specialist advice and treatment if pain is not improving within a short time or if they are experiencing intolerable side effects of analgesia [GRADE 1C]</p> <p>This includes access to advanced pain management techniques such as intrathecal pumps</p>
<b>European Society for Medical Oncology (ESMO) (2018)<sup>21</sup></b>	<p>Level II B - Intraspinal techniques delivered and monitored by a skilled team should be included as part of the cancer pain management strategy</p>
<b>Polyanalgesic Consensus Conference (PACC) convened by the International Neuromodulation Society (2017)<sup>59</sup></b>	<p>Localized pain can be adequately covered with intrathecal therapy (II B Strong)</p> <p>Diffuse pain can be adequately treated with intrathecal therapy (III C Moderate)</p> <p>Global pain can be adequately treated with intrathecal therapy (III D Moderate)</p>
<b>European Association for Palliative Care (EAPC) (2012)<sup>60</sup></b>	<p>Spinal (epidural or intrathecal) administration of opioid analgesics in combination with local anaesthetics or clonidine should be considered for patients in whom analgesia is inadequate or who have intolerable adverse effects despite the optimal use of oral and parenteral opioids and non-opioid agents. (Weak Recommendation)</p>



Refractory cancer pain

Burden of cancer pain

Current treatment management

Intrathecal drug delivery

Clinical evidence

Economic evidence

Therapy guidelines

Pump benefits

Patient selection

SynchroMed™ II drug infusion system


References

# Recommendations:

## Intrathecal drug delivery for cancer-related pain

Guideline and society (year)	Recommendations for intrathecal drug delivery
<b>Federatie Medisch Specialisten - The Netherlands: Pijn bij patiënten met kanker (2019)</b> <b>(Pain in Patients with Cancer, 2019)<sup>57</sup></b>	<p>Neuraxial (intrathecal or epidural) administration of opioids in cancer pain</p> <p>Recommendation: In cancer patients in the palliative phase: Consider intrathecal opioid administration, if necessary in combination with a local anesthetic and / or clonidine, for the treatment of pain in cancer patients when oral, transdermal or parenteral treatment with opioids has insufficient analgesic effect and / or is associated with severe side effects (1D).</p> <p>If neuraxial administration of opioids is expected to be necessary for longer than a few weeks, it is preferable to opt for administration via a Fully Implantable Delivery System (FIDS). This reduces the risk of catheter dislocation and infection.</p>
<b>British Pain Society (2015)<sup>61</sup></b> <b>Intrathecal drug delivery for the management of pain and spasticity in adults; recommendations for best clinical practice</b>	<p>The working group believes that there is reasonable evidence supporting the use of ITDD in pain in patients with cancer where this is not controlled by systemic analgesia or where systemic analgesia causes intolerable side effects.</p> <p>ITDD has also been found to be a cost-effective alternative to systemic, intravenous or external infusion devices for cancer patients who require pain management for 3 months or more.</p>






# Gaps in access to intrathecal drug delivery for cancer pain

Despite demonstrated clinical evidence and medical society recommendations, there is a substantial gap between the need and provision of intrathecal drug delivery for patients with refractory cancer pain

**The recent 2019 UK collaborative Framework for Provision of Pain Services for Adults Across the UK with Cancer or Life-limiting Disease<sup>47</sup> has stated the need for its guidance is timely because of:**

“the importance to the public of pain associated with cancer and life-limiting disease, the evidence of under-treatment or poor access to care, the need to show evidence of better pain management in CQC inspections from 2016, and the need to meet new commissioning requirements for managing complex interventions (e.g., intrathecal pumps)”<sup>47</sup>



<sup>§</sup> Based on NHS England estimate that 5-15% of cancer patients have refractory pain and require advanced techniques which may include: chordotomy, spinal injections, nerve root ablation, ITDD, radiofrequency ablation, chemical destruction of nerve, cryotherapy (7)

\*2014-2019 includes ITDD for management of spasticity; 2014 and 2015 includes ITDD for noncancer pain.

<sup>†</sup> Patients with an ICD10 C code up to 180 days prior to receiving an ITD

<sup>‡</sup> Up to January 2020.

HES, Hospital Episode Statistics; ITDD, intrathecal drug delivery device

Refractory cancer pain

Burden of cancer pain

Current treatment management

Intrathecal drug delivery

Clinical evidence

Economic evidence

Therapy guidelines

Pump benefits

Patient selection

Synchromed™ II drug infusion system

References



Refractory  
cancer pain

Burden of  
cancer pain

Current  
treatment  
management

Intrathecal  
drug delivery

Clinical  
evidence

Economic  
evidence

Therapy  
guidelines

Pump benefits

Patient selection

SynchroMed™ II  
drug infusion  
system

References

# Considerations for appropriate patients

## Patient selection recommendations

1. **Medium to long-term life expectancy**  
(typically  $\geq 3$  month)<sup>49,65</sup>
2. **Visual analogue score (VAS)  $\geq 5$  despite 200 mg/day**  
of oral morphine or analgesic equivalent<sup>49,50</sup>
3. **Consider those on lower dose analgesics** if opioid  
side effects are not well tolerate<sup>49,50,59</sup>
4. **In your assessment, consider the patient perspective  
on mobility**  
(i.e. preferred place of treatment - home or  
specialized care service)


Try the therapy  
with your patients

Test with  
injection




or

Test with external pump



Medtronic





# Benefitting patients and clinicians for more than 30 years

Refractory cancer pain

Burden of cancer pain

Current treatment management

Intrathecal drug delivery

Clinical evidence

Economic evidence

Therapy guidelines


Pump benefits

Patient selection

**Synchromed™ II drug infusion system**

References


Medtronic's global reach enables our intrathecal drug delivery therapy to be available in many countries across the world.



\* Design changes to address causes of motor stalls; Product Performance Report 2021, v1.0 03Mar2022. The implementation of these changes does not imply equivalent percent reductions of motor stall.

Note: This timeline only shows extracts of the Medtronic TDD product portfolio.






# Performance and reliability

The SynchroMed™ II Programmable Infusion System<sup>62</sup>:

- available with a 20ml and 40ml reservoir size
- precisely moves medication through a peristaltic action
- constant and flexible infusion modes with bolus dosing
- allows patient to self administer a therapeutic bolus with myPTM™
- safe 1.5T and 3.0 full body MRI possible<sup>\*64</sup>
- designed to resume therapy shortly after MRI scan
- replacement after 6-7 years<sup>66</sup>
- an alarm is activated 90 days before battery life ends




19.6mm thickness (20ml pump)  
26.1mm thickness (40ml pump)




Confidence year after year

## Design Enhancements Ascenda™

### Anchor dispenser tool with bi-wing anchor



### Catheter design




### Guide wire interlock handle



\* Under specific conditions for 1.5T and 3.0T MRI scans.  
Requires interrogation to confirm pump status.






# MRI performance

Medtronic's Product Surveillance Registry, reporting on the **largest cohort of TDD cancer pain patients to date** (1,403 patients included) tracked MRI activity amongst patients implanted with the Synchromed II™ programmable infusions system for cancer pain, and found:

- 3 MRIs were reported in 51 patients<sup>53</sup>
- All MRI-induced motor stalls recovered as expected within 24 hours<sup>53</sup>
- There were no reports of post-MRI drug withdrawal or sequelae<sup>53</sup>

## Synchromed II™ offers a simple protocol to manage MRI in cancer pain patients



Refractory cancer pain

Burden of cancer pain

Current treatment management

Intrathecal drug delivery

Clinical evidence

Economic evidence

Therapy guidelines

Pump benefits

Patient selection


Synchromed™ II drug infusion system

References




# Benefits for breakthrough cancer pain

- Onset of breakthrough cancer pain (BTcP) is sudden and reaches a maximum intensity within 1 second to 30 min (median  $\approx$  3 min)<sup>67</sup>
- Patients can average 4-6 episodes/day, with a median duration of 15-30 min/episode<sup>67</sup>
- Episodes may, or may not, be associated with a precipitating factor and therefore **may not be predictable**. Precipitating factors can be identified in 55-80% of all episodes<sup>67</sup>



**More patients able to better control breakthrough pain and 3x faster with myPTM™ compared to conventional medication.**<sup>63</sup>



- Breakthrough pain reduction was **46.8% before** implant and **65.2% after** implant and with myPTM™<sup>63</sup>
- Patients with severe pain were **84.2% before** implant and **35.2% after** implant and myPTM™<sup>63</sup>
- With myPTM™ patients can manage breakthrough pain **quicker compared to systemic opioids** (10 vs 30 min on average)<sup>63</sup>

Refractory cancer pain

Burden of cancer pain

Current treatment management

Intrathecal drug delivery

Clinical evidence

Economic evidence

Therapy guidelines


Pump benefits

Patient selection

Synchromed™ II drug infusion system

References






# Easy programming of myPTM™ with clinician programmer


## Compare with ease

Easily compare current and pending updates while programming simple continuous mode or myPTM™.



## Managing Pain. Managing options.

Targeted drug delivery (TDD) with myPTM™ allows for patient control of breakthrough pain – within programmed parameters – ensuring 100% compliance. When pain is breakthrough, your prescription is right at hand.



myPTM™  
Personal Therapy Manager

Refractory cancer pain

Burden of cancer pain

Current treatment management

Intrathecal drug delivery

Clinical evidence

Economic evidence

Therapy guidelines


Pump benefits

Patient selection

Synchromed™ II drug infusion system

References





# References

Refractory  
cancer pain

Burden of  
cancer pain


Current  
treatment  
management

Intrathecal  
drug delivery

1. Van Den Beuken-Van Everdingen M H J, Hochstenbach L M J, Joosten E A J, Tjan-Heijnen V C G and Janssen D J A 2016 Update on Prevalence of Pain in Patients with Cancer: Systematic Review and Meta-Analysis J. Pain Symptom Manage. 51 1070-1090.e9
2. Stearns L, Boortz-Marx R, Du Pen S, Friehs G, Gordon M, Halyard M, Herbst L K J 2005 Intrathecal Drug Delivery for the Management of Cancer Pain: A Multidisciplinary Consensus of Best Clinical Practices J Support Oncol , 3 399-408
3. O'Brien T and Kane C M 2014 Pain services and palliative medicine - an integrated approach to pain management in the cancer patient Br. J. Pain 8 163-71
4. Malvezzi M, Carioli G, Bertuccio P, Boffetta P, Levi F, Vecchia C La and Negri E 2018 European cancer mortality predictions for the year 2018 with focus on colorectal cancer Ann. Oncol. 29 1016-22
5. Malvezzi M, Carioli G, Bertuccio P, Boffetta P, Levi F, La Vecchia C and Negri E 2019 European cancer mortality predictions for the year 2019 with focus on breast cancer Ann. Oncol. 30 781-7
6. Malvezzi M, Carioli G, Bertuccio P, Boffetta P, Levi F, La Vecchia C and Negri E 2017 European cancer mortality predictions for the year 2017, with focus on lung cancer Ann. Oncol. Off. J. Eur. Soc. Med. Oncol. 28 1117-23
7. Ovayolu Ö, Ovayolu N, Aytaç S, Serçe S, Sevinc A. Pain in cancer patients: pain assessment by patients and family caregivers and problems experienced by caregivers. Supportive Care in Cancer 2015;23(7): 1857-64.
8. Bruel BM, Burton AW. Intrathecal Therapy for Cancer-Related Pain. Pain Med. 2016;17(12):2404-21.
9. Breivik H, Cherny N, Collett B, de Conno F, Filbet M, Foubert A J, Cohen R and Dow L 2009 Cancer-related pain: A pan-European survey of prevalence, treatment, and patient attitudes Ann. Oncol. 20 1420-33
10. R K Portenoy; N A Hargen 1990 Breakthrough Pain: Definition, Prevalence and Characteristics Pain 41 271-81
11. Deandrea S, Corli O, Consonni D, Villani W, Greco M T and Apolone G 2014 Prevalence of breakthrough cancer pain: A systematic review and a pooled analysis of published literature J. Pain Symptom Manage. 47 57-76
12. Davies A, Buchanan A, Zeppetella G, Porta-Sales J, Likar R, Weismayr W, Slama O, Korhonen T, Filbet M, Poulain P, Mystakidou K, Ardavanis A, O'Brien T, Wilkinson P, Caraceni A, Zucco F, Zuurmond W, Andersen S, Damkier A, Vejlgård T, Nauck F, Radbruch L, Sjolund K F and Stenberg M 2013 Breakthrough cancer pain: An observational study of 1000 european oncology patients J. Pain Symptom Manage. 46 619-28
13. Pérez-Hernández C, Jiménez-López A J, Sanz-Yagüe A, Mar-Medina J, Larrañaga I and Soler-López B 2018 Observational Study Evaluating the Economic Impact of Breakthrough Pain in Cancer Patients in Clinical Practice in Spain: The IMDI Study Pain Ther. 7 227-40
14. Drury A, Payne S and Brady A M 2017 The cost of survival: an exploration of colorectal cancer survivors' experiences of pain Acta Oncol. (Madr). 56 205-11
15. Maltoni M 2008 Opioids, pain, and fear Ann. Oncol. 19 5-7
16. Payne S, Turner M, Seamark D, Thomas C, Brearley SG, Wang X, Blake S M C 2015 Managing end of life medications at home - accounts of bereaved family carers: a qualitative interview study BMJ Support Palliat Care 5 181-8
17. Abramowitz L, Béziaud N, Labreze L, Giardina V, Caussé C, Chuberre B, Allaert F A and Perrot S 2013 Prevalence and impact of constipation and bowel dysfunction induced by strong opioids: A cross-sectional survey of 520 patients with cancer pain: DYONISOS study J. Med. Econ. 16 1423-33
18. Annemans L 2011 Pharmacoeconomic impact of adverse events of long-term opioid treatment for the management of persistent pain Clin. Drug Investig. 31 73-86
19. Søndergaard J, Christensen H N, Ibsen R, Jarbøl D E and Kjellberg J 2017 Healthcare resource use and costs of opioid-induced constipation among non-cancer and cancer patients on opioid therapy: A nationwide register-based cohort study in Denmark Scand. J. Pain 15 83-90
20. Dahl A A, Bentzen A G, Fosså S D, Hess S L, Steen R and Kiserud C E 2020 Long-term cervical cancer survivors on disability pension : a subgroup in need of attention from health care providers 5816
21. Fallon M, Giusti R, Aielli F, Hoskin P, Rolke R, Sharma M and Ripamonti C I 2018 Management of cancer pain in adult patients: ESMO Clinical Practice Guidelines Ann. Oncol. 29 iv166-91
22. Kwon J H 2014 Overcoming barriers in cancer pain management J. Clin. Oncol. 32 1727-33
23. Torresan M Ma, Garrino L, Borraccino A, Macchi G, De Luca A, Dimonte V. Adherence to treatment in patient with severe cancer pain: A qualitative enquiry through illness narratives. Eur J Oncol Nurs 2015 Aug;19(4):397-404
24. Reis-Pina P, Lawlor P G and Barbosa A 2017 Adequacy of cancer-related pain management and predictors of undertreatment at referral to a pain clinic J. Pain Res. 10 2097-107
25. Deandrea S, Montanari M, Moja L and Apolone G 2008 Prevalence of undertreatment in cancer pain. A review of published literature Ann. Oncol. 19 1985-91
26. Meuser T, Pietruck C, Radbruch L, Stute P, Lehmann K A and Grond S 2001 Symptoms during cancer pain treatment following WHO-guidelines: A longitudinal follow-up study of symptom prevalence, severity and etiology Pain 93 247-57
27. Greco M T, Roberto A, Corli O, Deandrea S, Bandieri E, Cavuto S and Apolone G 2014 Quality of cancer pain management: An update of a systematic review of undertreatment of patients with cancer J. Clin. Oncol. 32 4149-54
28. World Health Organisation 1996 Cancer pain relief, 2nd Edn vol 34 (Geneva)
29. Orhan M E, Bilgin F, Ergin A, Dere K and Güzeldemir M E 2008 Pain treatment practice according to the WHO analgesic ladder in cancer patients: Eight years experience of a single center Agri 20 37-43
30. Corli O, Floriani I, Roberto A, Montanari M, Galli F, Greco M T, Caraceni A, Kaasa S, Dragani T A, Azzarello G, Luzzani M, Cavanna L, Bandieri E, Gamucci T, Lipari G, Di Gregorio R, Valenti D, Reale C, Pavesi L, Iorno V, Crispino C, Pacchioni M and Apolone G 2016 Are strong opioids equally effective and safe in the treatment of chronic cancer pain? A multicenter randomized phase IV "real life" trial on the variability of response to opioids Ann. Oncol. Off. J. Eur. Soc. Med. Oncol. 27 1107-15
31. Miguel R 2000 Interventional treatment of cancer pain: The fourth step in the world health organization analgesic ladder? Cancer Control 7 149-56
32. Perruchoud C., Dupoirion D., Papi B., Calabrese A., Brogan S.E. 2021. Management of Cancer-Related Pain With Intrathecal Drug Delivery: A Systematic Review and Meta-Analysis of Clinical Studies. Neuromodulation 2022 Jan 21;S1094-7159
33. Bennett M, Paice J A and Wallace M 2017 Pain and Opioids in Cancer Care: Benefits, Risks, and Alternatives Am. Soc. Clin. Oncol. Educ. B. 37 705-13
34. Yang J, Bauer B A, Wahner-Roedler D L, Chon T Y and Xiao L 2020 The modified WHO analgesic ladder: Is it appropriate for chronic non-cancer pain? J. Pain Res. 13 411-7
35. World Health Organisation 2012 WHO guidelines on the pharmacological treatment of persisting pain in children with medical illness
36. Mantyh P W 2006 Cancer pain and its impact on diagnosis, survival and quality of life Nat. Rev. Neurosci. 7 797-809
37. Armstrong A J, Garrett-Mayer E, Yang Y C O, Carducci M A, Tannock I, De Wit R and Eisenberger M 2007 Prostate-specific antigen and pain surrogacy analysis in metastatic hormone-refractory prostate cancer J. Clin. Oncol. 25 3965-70

References





# References

Refractory  
cancer pain


Burden of  
cancer pain

Current  
treatment  
management

Intrathecal  
drug delivery

38. Herndon J E, Fleishman S, Kornblith A B, Kosty M, Green M R and Holland J 1999 Is quality of life predictive of the survival of patients with advanced nonsmall cell lung carcinoma? Cancer 85 333-4
39. Quinten C, Coens C, Mauer M, Comte S, Sprangers M A, Cleeland C, Osoba D, Bjordal K and Bottomley A 2009 Baseline quality of life as a prognostic indicator of survival: a meta-analysis of individual patient data from EORTC clinical trials Lancet Oncol. 10 865-71
40. Zylla D, Steele G, Shapiro A, Richter S and Gupta P 2018 Impact of opioid use on health care utilization and survival in patients with newly diagnosed stage IV malignancies Support. Care Cancer 26 2259-66
41. Wanderer J 2016 Anesthesia & Analgesia Infographic Anesth. Analg. 123 2016
42. Fine P G, Chen Y, Wittbrodt E and Datto C 2019 Impact of opioid-induced constipation on healthcare resource utilization and costs for cancer pain patients receiving continuous opioid therapy 687-96
43. García-Mata J, Álamo C, de Castro J, Contreras J, Gálvez R, Jara C, Llombart A, Pérez C, Sánchez P, Traseira S and Cruz J J 2018 A survey of perceptions, attitudes, knowledge and practices of medical oncologists about cancer pain management in Spain Clin. Transl. Oncol. 20 1061-71
44. O'Hagan P and Mercadante S 2018 Breakthrough cancer pain: The importance of the right treatment at the right time Eur. J. Pain (United Kingdom) 22 1362-74
45. Council of Europe-Parliamentary Assembly 2018 The provision of palliative care in Europe 2249 1-14
46. Finn L, Green A R and Malhotra S 2017 Oncology and palliative medicine: Providing comprehensive care for patients with cancer Ochsner J. 17 393-7
47. Faculty of Clinical Oncology, The Royal College of Radiologists U K 2019 Framework for Provision of Pain Services for Adults Across the UK with Cancer or Life-limiting Disease
48. Bhatia G, Lau M E and Gulur P 2013 Intrathecal Drug Delivery (ITDD) systems for cancer pain F1000Research 2 96
49. Smith T J, Staats P S, Deer T, Stearns L J, Rauck R L, Boortz-Marx R L, Buchser E, Català E, Bryce D A, Coyne P J and Pool G E 2002 Randomized clinical trial of an implantable drug delivery system compared with comprehensive medical management for refractory cancer pain: Impact on pain, drug-related toxicity, and survival J. Clin. Oncol. 20 4040-9
50. Smith T J, Coyne P J, Staats P S, Deer T, Stearns L J, Rauck R L, Boortz-Marx R L, Buchser E, Català E, Bryce D A, Cousins M, Pool G E, Wallace M, Yaksh T, Magnuson S, Leung A, Ahadian F, Bullock S, McBeth M, Hoyer A, Miguel R, Weitzner M A, Balducci L, Follett K, Hitchon P, Minore W S, Weiss H, Jaworowicz J, Croy S, Davis S, Grossman S, Grieb M, Carson L, Mitchell S, Stuckey M, Charapata S, McCracken M, Sorensen R, Calava J, Dunn J, Kosek P, Weinstein S, Banerjee T K, Caraway D, Kim C, Serafini M, McNeil K, Frame J, Orlandini G, Siragusa A, Aliaga F, Lopez-Pousa A, Pericay C, Perey L, Muriel C, Boyle F, Molloy A, Brooker C and Walker S E 2005 An implantable drug delivery system (IDDS) for refractory cancer pain provides sustained pain control, less drug-related toxicity, and possibly better survival compared with comprehensive medical management (CMM) Ann. Oncol. 16 825-33
51. Duarte R., Copley S., Nevitt S., Maden M., Al-Ali A.M., Dupoirion D., Eldabe S. 2022. Effectiveness and Safety of Intrathecal Drug Delivery Systems for the Management of Cancer Pain: A Systematic Review and Meta-Analysis. Neuromodulation 2022 Apr 11;S1094-7159
52. Lambe T., Duarte R., Eldabe R., Copley S., Kansal A., Black S., Dupoirion D., Eldabe S. 2022. Ziconotide for the Management of Cancer Pain: A Budget Impact Analysis. Neuromodulation 2022 Oct 3;S1094-7159(22)
53. Stearns L M, Abd-Elseyed A, Perruchoud C, Spencer R, Hammond K, Stromberg K and Weaver T 2020 Intrathecal Drug Delivery Systems for Cancer Pain: An Analysis of a Prospective, Multicenter Product Surveillance Registry Anesth. Analg. 130 289-97
54. Sindt J E, Larsen S D, Dalley A P, Collier W H and Brogan S E 2020 The Rate of Infectious Complications After Intrathecal Drug Delivery System Implant for Cancer-Related Pain Is Low Despite Frequent Concurrent Anticancer Treatment or Leukopenia Anesth. Analg. XXX 1
55. Stearns L J, Narang S, Albright R E, Hammond K, Xia Y, Richter H B, Paramanandam G K, Haagenson K K and Doth A H 2019 Assessment of Health Care Utilization and Cost of Targeted Drug Delivery and Conventional Medical Management vs Conventional Medical Management Alone for Patients With Cancer-Related Pain JAMA Netw. open 2 e191549
56. Brogan S E, Winter N B, Abiodun A and Safarpour R 2013 A Cost Utilization Analysis of Intrathecal Therapy for Refractory Cancer Pain: Identifying Factors Associated with Cost Benefit Pain Med. (United States) 14 478-86
57. Nederlandse Vereniging voor Anesthesiologie 2019 Pijn bij Patiënten met Kanker
58. Bennett M I, Eisenberg E, Ahmedzai S H, Bhaskar A, O'Brien T, Mercadante S, Krčevski Škvarč N, Vissers K, Wirz S, Wells C and Morlion B 2019 Standards for the management of cancer-related pain across Europe—A position paper from the EFIC Task Force on Cancer Pain Eur. J. Pain (United Kingdom) 23 660-8
59. Deer T R, Pope J E, Hayek S M, Bux A, Buchser E, Eldabe S, De Andrés J A, Erdek M, Patin D, Grider J S, Doleys D M, Jacobs M S, Yaksh T L, Poree L, Wallace M S, Prager J, Rauck R, DeLeon O, Diwan S, Falowski S M, Gazelka H M, Kim P, Leong M, Levy R M, McDowell G, McRoberts P, Naidu R, Narouze S, Perruchoud C, Rosen S M, Rosenberg W S, Saulino M, Staats P, Stearns L J, Willis D, Krames E, Huntoon M and Mekhail N 2017 The Polyanalgesic Consensus Conference (PACC): Recommendations on Intrathecal Drug Infusion Systems Best Practices and Guidelines Neuromodulation 20 96-132
60. Caraceni A, Hanks G, Kaasa S, Bennett M I, Brunelli C, Cherny N, Dale O, De Conno F, Fallon M, Hanna M, Haugen D F, Juhl G, King S, Klepstad P, Laugsand E A, Maltoni M, Mercadante S, Nabal M, Pigni A, Radbruch L, Reid C, Sjogren P, Stone P C, Tassinari D and Zeppetella G 2012 Use of opioid analgesics in the treatment of cancer pain: Evidence-based recommendations from the EAPC Lancet Oncol. 13 58-68
61. NHS England 2015 Clinical Commissioning Policy: Intrathecal Pumps for Treatment of Severe Cancer Pain Clin. Comm. Policy Intrathecal Pumps Treat. Sev. Cancer Pain
62. Medtronic Synchromed II Product Manual. www.medtronic.com/manuals
63. Brogan S E, Winter N B and Okifuji A 2015 Prospective Observational Study of Patient-Controlled Intrathecal Analgesia: Impact on Cancer-Associated Symptoms, Breakthrough Pain Control, and Patient Satisfaction Reg. Anesth. Pain Med. 40 369-75
64. Lecouvet F E, Talbot J N, Messiou C, Bourguet P, Liu Y and De Souza N M 2014 Monitoring the response of bone metastases to treatment with Magnetic Resonance Imaging and nuclear medicine techniques: A review and position statement by the European Organisation for Research and Treatment of Cancer imaging group Eur. J. Cancer 50 2519-31
65. Burton A W, Rajagopal A, Shah H N, Mendoza T, Cleeland C, Hassenbusch S J and Arens J F 2004 Epidural and intrathecal analgesia is effective in treating refractory cancer pain Pain Med. 5 239-47
66. Medtronic Product Performance Report 2018
67. Mishra, Seema, Bhatnagar, Sushma, Chaudhary, Prakash and S P S R 2009 Breakthrough Cancer Pain: Review of Prevalence, Characteristics and Management Indian J. Palliat. Care 15 14-8

References



Refractory cancer pain

Burden of cancer pain

Current treatment management

Intrathecal drug delivery

**Brief Statement**

See the device manual for detailed information regarding the instructions for use, the implant procedure, indications, contraindications, warnings, precautions, and potential adverse events. For further information please contact your local Medtronic representative and/or consult the Medtronic website at [www.medtronic.eu](http://www.medtronic.eu).

For applicable products, consult instructions for use on [www.medtronic.com/manuals](http://www.medtronic.com/manuals). Manuals can be viewed using a current version of major internet browser. For best results, use Adobe Acrobat® Reader in the browser.

When ITB is mentioned, we are considering Intrathecal baclofen (an antispasmodic) administered by an intrathecal drug delivery pump therapy. Medtronic provides only the intrathecal drug delivery pump and the catheter; the baclofen is provided by an external company.

**Medtronic**

**Europe**  
Medtronic International Trading Sàrl.  
Route du Molliau 31  
Case postale  
CH-1131 Tolochenaz  
[www.medtronic.eu](http://www.medtronic.eu)  
Tel: +41 (0)21 802 70 00  
Fax: +41 (0)21 802 79 00

**United Kingdom/Ireland**  
Medtronic Limited  
Building 9  
Croxley Park  
Hatters Lane  
Watford  
Herts WD18 8WW  
[www.medtronic.co.uk](http://www.medtronic.co.uk)  
Tel: +44 (0)1923 212213  
Fax: +44 (0)1923 241004

UC202312121EE © Medtronic 2023.  
All rights reserved.

**medtronic.eu**