

## Pediatric cardiac critical care

### Cerebral oxygen saturation monitoring with the INVOS™ system

Ronald A. Bronicki, MD and Paul A. Checchia, MD

#### Introduction

Congenital heart disease (CHD) is the most common birth defect, with an incidence of moderate to severe forms of 6 in 1000 live births.<sup>1</sup> Since the first surgical repair of congenital heart disease in 1953, the overall mortality in pediatric cardiac surgery has progressively decreased to <3%.<sup>2</sup> Despite this, surviving patients continue to display a relatively high incidence of neurological, developmental, and psychiatric disabilities.

During the early postoperative period following surgical repair of congenital heart defects, the vulnerable neonatal brain is subject to unstable hemodynamics with reduced systemic oxygen delivery ( $\text{DO}_2$ ) at a time when oxygen demand ( $\text{VO}_2$ ) is elevated.<sup>3-5</sup> Studies to date suggest that the first 12-24 hours after surgery represents a critical period for managing the risk of ischemic neurological injury. Reductions in cerebral tissue oxygenation during this time are associated with adverse neurodevelopmental outcomes and interventions designed to optimize hemodynamics and improve  $\text{DO}_2$  may thus reduce neurological injury.<sup>5-9</sup>

In a cohort of 13 patients, Hoffman et al. demonstrated that decreased systemic oxygen delivery in the early postoperative period was significantly associated with hypoxic-ischemic brain injury and neurodevelopmental abnormalities while these adverse outcomes were not

associated with either deep hypothermic circulatory arrest or cardiopulmonary bypass (CPB) time.<sup>9</sup> Similarly, in 22 neonates treated for the hypoplastic left heart syndrome (HLHS), 73% of patients had new or worsened ischemic lesions upon repeat MRI analysis and there was a significant association between prolonged reductions in postoperative cerebral oxygenation (<45% for >180 cumulative minutes) and the development of MRI lesions.<sup>6</sup> Additionally, long-term follow-up studies have indicated that reduced perioperative cerebral oxygenation is associated with impaired neurodevelopmental outcomes and brain abnormalities on MRI 4-6 years following surgery.<sup>7,8</sup> Together, these data indicate that perioperative reductions in cerebral oxygenation adversely impact neurodevelopmental outcomes and effectively managing the postoperative period may be critical for reducing the development or progression of neurological injury.

#### Current clinical practice and challenges

With these data in mind, simply targeting mortality following congenital heart surgery is no longer an acceptable goal and current treatment guidelines now focus on improving neurodevelopmental outcomes. The early identification of perioperative factors that result in neurological injury allow for the rapid implementation of corrective strategies and improved outcomes. A key component to critical care monitoring of the pediatric cardiac surgery patient is the rapid and accurate assessment of cardiovascular function, cardiac output, and tissue oxygenation. Unfortunately, estimations of these variables based on routine clinical data such as the physical examination and standard hemodynamic parameters, such as heart rate, blood pressure, and central venous pressure, are often unreliable.<sup>10-12</sup> In children, there is a poor correlation between estimations of cardiac output and systemic vascular resistance based on peripheral pulses, capillary refill, and peripheral/core body temperatures and actual measurements of these parameters.<sup>13,14</sup> Blood pressure is also a poor correlate of cardiac output, as systemic vascular resistance rises to compensate for a decrease in perfusion.

#### Key Learning Points

- Cerebral hypoperfusion is a common complication of complex thoracic aortic operations that adversely affects clinical outcome and the cost of care.
- Cerebral oxygenation and perfusion can be monitored continuously and noninvasively during thoracic aortic operations using the INVOS™ cerebral/somatic oximeter.
- Monitoring with the INVOS™ system may aid in the early detection of cerebral hypoperfusion or malperfusion.

## The INVOS™ monitor

Near-infrared spectroscopy (NIRS) provides a non-invasive tool for continuous monitoring of regional tissue oxygen saturation ( $rSO_2$ ).<sup>15</sup> NIRS analysis is based on the fact that hemoglobin has characteristic infrared absorption spectra that shift when bound to oxygen. This spectral shift permits the concentration and ratio of oxygenated to deoxygenated hemoglobin to be measured.

The INVOS™ cerebral/somatic oximeter employs disposable sensors with an integrated near infrared light source and photodetector that can be applied to each side of the forehead for monitoring blood in the brain or to other sites for monitoring tissue beneath the sensor. The INVOS™ oximeter uses 2 wavelengths (730 and 810 nm) of near-infrared light, which are emitted from a light emitting diode, and detected at 2 separate distances (3 and 4 cm) so that the signal arising from superficial tissues can be separated from the signal arising from the cerebral cortex (or somatic tissue of interest). More specifically, for cerebral oximetry, the 3 cm detector is designed to measure the signal generated by light passing through shallow structures such as skin, skull, and soft tissue while the 4 cm detector is designed to measure the signal generated by light passing through both these shallow structures and the frontal cerebral cortex. With this design, approximately 85% of the signal is derived from cortical tissue.<sup>16</sup>

In contrast to pulse oximetry, NIRS monitoring evaluates the non-pulsatile signal, reflecting the oxygen saturation of the microcirculation. Because the majority of blood in the microcirculation is venous (upwards of 75-80%), the oxygen saturation is used as a surrogate for venous oxygen saturations. The INVOS™ oximeter calculates the cerebral oxygen saturation using an algorithm that corrects the signal for extracranial noise and assumes a ratio of 25% arterial to 75% venous blood in the sampled tissue. The pediatric/neonatal INVOS™ system is designed for patients weighing <40 kg and differs from the adult system in that it incorporates a specific algorithm – taking into account the thinner skull and extracranial tissue of infants and children.<sup>17</sup>

## Monitoring oxygen transport balance


Overall, NIRS analysis of  $rSO_2$  reflects the balance between oxygen supply and demand and the ratio of oxyhemoglobin to total hemoglobin reflects the surplus of oxygen remaining after tissue extraction (see Figure 1 for summary). In brief, the Fick equation dictates that  $VO_2$  is determined by the product of cardiac output (CO) and the arterial venous oxygen content difference ( $CaO_2 - CvO_2$ ) or the volume of oxygen extracted per minute.<sup>11</sup> The Fick equation may be simplified to:  $SaO_2 - SvO_2 / SaO_2 = VO_2 / DO_2$ , where  $SaO_2$  is the arterial oxygen saturation and  $SvO_2$  is the venous oxygen saturation. Thus, by determining the  $O_2$  extraction ratio ( $SaO_2 - SvO_2 / SaO_2$ ) the relationship between  $DO_2$  and  $VO_2$  can be assessed.

## The Fick Principle

$$VO_2 = CO \times CaO_2 - CvO_2$$

$$SaO_2 - SvO_2 / SaO_2 = VO_2 / DO_2 = \text{Oxygen Transport Balance}$$

$$\text{Oxygen Extraction Ratio (OER)} = SaO_2 - SvO_2 / SaO_2$$



**FIGURE 1.** Physiological relationship between oxygen transport balance and oxygen extraction ratio. When oxygen content is decreased, either due to hypoxemia or anemia, or when  $O_2$  demand is increased, global (increase in CO) as well as regional (vasodilation) compensatory circulatory mechanisms maintain oxygen delivery and a relatively stable regional and global oxygen extraction ratio. It is not until these compensatory mechanisms are exhausted that oxygen delivery falls and the oxygen extraction ratio begins to increase. As the oxygen extraction ratios rise above 50-60%, tissue hypoxia ensues and lactate begins to accumulate.

When oxygen content is decreased, either due to hypoxemia or anemia, or when  $O_2$  demand is increased, global (increase in CO) as well as regional (vasodilation) compensatory circulatory mechanisms maintain  $DO_2$  and a relatively stable regional and global oxygen extraction ratio. It is not until these compensatory mechanisms are exhausted that  $DO_2$  falls and the oxygen extraction ratio begins to increase (Figure 1) in order to maintain tissue oxygenation. In another words, it is when  $DO_2$  becomes limited that the only compensatory mechanism available to maintain tissue oxygenation is an increase in oxygen extraction. The normal global and cerebral oxygen extraction ratios are 25-30%. As the oxygen extraction ratio rises above 50-60%, tissue hypoxia ensues and lactate begins to accumulate (Figure 1).

## Clinical value of perioperative oxygen saturation monitoring with INVOS™ technology

Over the past ~15 years, a number of studies have evaluated  $rSO_2$  monitoring as a surrogate measure of either systemic or cerebral oxygenation in children with critical heart disease.<sup>18-23</sup> Overall, these studies have indicated that  $rSO_2$  correlates well with jugular venous ( $SjO_2$ ) and central venous oxygen saturations ( $SvO_2$ ) (see Table 1).

**Table 1. Summary of Key Studies Comparing  $rSO_2$  with  $SjO_2$  and/or  $SvO_2$**

Daubeney et al. (1996)	<ul style="list-style-type: none"> <li>• Compared cerebral <math>rSO_2</math> and <math>SjO_2</math> in 40 children with diverse forms of congenital heart disease in the cardiac catheter laboratory or during cardiac surgery.</li> <li>• Across patients, the correlation between <math>rSO_2</math> and <math>SjO_2</math> was 0.69 (<math>p &lt; 0.0001</math>).</li> <li>• For surgery patients, correlations ranged from 0.78 to 0.96 with a median value of 0.91.</li> <li>• <math>rSO_2</math> tended to overestimate <math>SjO_2</math> at low <math>SjO_2</math> values and tended to underestimate <math>SjO_2</math> at high <math>SjO_2</math> values.</li> </ul>
Tortoriello et al. (2005)	<ul style="list-style-type: none"> <li>• Compared cerebral <math>rSO_2</math> with <math>SvO_2</math> in 20 children following surgical repair of congenital heart disease.</li> <li>• <math>rSO_2</math> was well correlated with <math>SvO_2</math> across patients (<math>r = 0.67</math> and <math>r^2 = 0.45</math>; <math>p &lt; 0.0001</math> for both).</li> <li>• The overall bias for <math>rSO_2</math> as a predictor of <math>SvO_2</math> was 3.3% with a precision of 16.6%; <math>rSO_2</math> was a more accurate predictor of <math>SvO_2</math> in biventricular repair patients (bias = -0.3%, precision = 11.8%).</li> </ul>
Li et al. (2006)	<ul style="list-style-type: none"> <li>• Evaluated cerebral and somatic <math>rSO_2</math>, <math>VO_2</math>, and systemic hemodynamic and oxygen transport variables in 11 children undergoing the Norwood procedure.</li> <li>• Cerebral <math>rSO_2</math> was correlated with <math>SvO_2</math> (<math>r = 0.43</math>), <math>SaO_2</math> (<math>r = 0.61</math>), <math>PaO_2</math> (<math>r = 0.99</math>), and MAP (<math>r = 0.52</math>; <math>p &lt; 0.0001</math> for all); somatic <math>rSO_2</math> showed similar but weaker correlations.</li> </ul>
Nagdyman et al. (2008)	<ul style="list-style-type: none"> <li>• Compared cerebral <math>rSO_2</math> with jugular bulb oxygen saturation (<math>SjO_2</math>) and central <math>SvO_2</math> in 31 children with congenital heart defects.</li> <li>• Cerebral <math>rSO_2</math> was well correlated with <math>SjO_2</math> (<math>r = 0.83</math>, <math>p &lt; 0.0001</math>) and <math>SvO_2</math> (<math>r = 0.93</math>, <math>p &lt; 0.0001</math>).</li> <li>• The mean bias was -5.2% for <math>rSO_2</math> and <math>SjO_2</math> and 5.6% for <math>rSO_2</math> and <math>SvO_2</math>.</li> </ul>
Ricci et al. (2010)	<ul style="list-style-type: none"> <li>• Retrospectively compared <math>SvO_2</math> and <math>rSO_2</math> in 100 newborn patients between 24 hours preoperative and 72 hours postoperative for cardiac surgery.</li> <li>• The correlation between <math>rSO_2</math> and <math>SvO_2</math> was 0.37 (<math>p &lt; 0.001</math>); when analyzed by time (pre/postoperative), the correlations were 0.24 (<math>p = 0.05</math>) for preoperative, 0.34 (<math>p &lt; 0.01</math>) for early postoperative, and 0.50 (<math>p &lt; 0.001</math>) for late postoperative.</li> <li>• During the preoperative period, <math>rSO_2</math> tended to overestimate <math>SvO_2</math> (bias = 1.5); during the early postoperative period, <math>rSO_2</math> tended to underestimate <math>SvO_2</math> (bias = -1.2).</li> <li>• The correlation between average <math>rSO_2</math> values at consecutive time points and average <math>SvO_2</math> at the same time points was 0.8 (<math>p &lt; 0.005</math>).</li> </ul>
Moreno et al. (2013)	<ul style="list-style-type: none"> <li>• Compared <math>SvO_2</math> and <math>rSO_2</math> (cerebral and flank) in 23 newborns and infants younger than 45 days undergoing heart surgery.</li> <li>• The mean difference (<math>\pm</math>SD) between <math>SvO_2</math> and cerebral <math>rSO_2</math> was <math>10.45 \pm 14.12\%</math>.</li> <li>• The mean difference (<math>\pm</math>SD) between <math>SvO_2</math> and somatic (flank) <math>rSO_2</math> was <math>7.16 \pm 16.67\%</math>.</li> <li>• The correlation between <math>SvO_2</math> and cerebral <math>rSO_2</math> was 0.48 (<math>p &lt; 0.001</math>); for <math>SvO_2</math> and somatic <math>rSO_2</math>, the correlation was 0.35 (<math>p &lt; 0.001</math>).</li> </ul>

**Table 2. Summary of Key Cerebral INVOS™ Studies**

Dent et al. (2005)	<ul style="list-style-type: none"> <li>• Twenty-two neonates with HLHS were studied to compare preoperative and postoperative brain MRI findings.</li> <li>• Continuous cerebral oxygen saturation monitoring with INVOS™ technology was employed with a critical rSO<sub>2</sub> value of 45% used to indicate cerebral desaturation.</li> <li>• 73% of patients had new or worsened ischemic lesions upon repeat MRI analysis.</li> <li>• 70% of patients had prolonged (&gt;180 minutes) cerebral desaturation (rSO<sub>2</sub> &lt;45%) during the postoperative period.</li> <li>• There was a significant association (p=0.029) between prolonged postoperative cerebral low rSO<sub>2</sub> (&gt;180 cumulative minutes &lt;45%) and the development of MRI lesions, with a sensitivity and specificity of 82% and 75%, respectively.</li> <li>• The positive predictive value of rSO<sub>2</sub> for new MRI lesions was 90% and the negative predictive value was 60%.</li> </ul>
Li et al. (2008)	<ul style="list-style-type: none"> <li>• Sixteen neonates undergoing the Norwood procedure were studied to examine the effects of systemic hemodynamics and O<sub>2</sub> transport variables on rSO<sub>2</sub> during the first 72 hours postoperatively.</li> <li>• rSO<sub>2</sub> was 66% at baseline, 88% at the end of cooling, 90% during SACP, 64% during rewarming, 54% after separation from CPB, and 51% on arrival in ICU.</li> <li>• By 72 hours postoperative, rSO<sub>2</sub> was 56%, which was still significantly below baseline (p=0.02).</li> <li>• Overall ICU rSO<sub>2</sub> values ranged from 27% to 79% with 28% of readings lower than 48%.</li> <li>• rSO<sub>2</sub> was positively correlated with SBP, DBP, MAP, Qs, Qp, and DO<sub>2</sub>; rSO<sub>2</sub> was negatively correlated with ERO<sub>2</sub> and SVR.</li> <li>• Overall, critical rSO<sub>2</sub> values were seen frequently in these patients and changes in rSO<sub>2</sub> were significantly related to changes in systemic hemodynamics and oxygen transport variables.</li> <li>• Hemodynamic interventions to modify systemic O<sub>2</sub> transport may help reduce the risk of cerebral ischemia and improve neurodevelopmental outcomes.</li> </ul>
Horvath et al. (2010)	<ul style="list-style-type: none"> <li>• Thirty-six postoperative children who had undergone delayed sternal closure after cardiac surgery were studied to determine if rSO<sub>2</sub> correlated with indicators of hemodynamic compromise.</li> <li>• Cerebral and somatic (renal) rSO<sub>2</sub> was reduced after delayed sternal closure suggesting mild and transient hemodynamic compromise after sternal closure.</li> </ul>
Bronicki et al. (2013)	<ul style="list-style-type: none"> <li>• Twenty-three patients were monitored for 30 minutes before and after extubation following repair of tetralogy of Fallot to determine the effects of extubation and loading of respiratory muscles on the distribution of cardiac output.</li> <li>• With extubation, central venous oxygen saturation increased from 65% to 70% (p=0.003) and cerebral rSO<sub>2</sub> increased from 67% to 72% (p=0.0001).</li> <li>• Before extubation, 4 patients were below or near the critical cerebral rSO<sub>2</sub> value of 50%.</li> <li>• rSO<sub>2</sub> data suggest that the brain was in or approaching an oxygen supply-dependent state with a subgroup of patients near or below critical rSO<sub>2</sub> values.</li> </ul>
Hansen et al. (2013)	<ul style="list-style-type: none"> <li>• Thirty-two patients with congenital heart disease who underwent superior cavopulmonary anastomosis were studied to assess the clinical value of NIRS monitoring.</li> <li>• Cerebral rSO<sub>2</sub> was 44% at the end of surgery and reached a minimal value of 40% at 2 hours postoperative.</li> <li>• Cerebral rSO<sub>2</sub> gradually increased during the postoperative period to 57% by the end of the 48 hours postoperative monitoring period and increased following extubation (49% before, 53% after).</li> <li>• Cerebral rSO<sub>2</sub> was moderately correlated with SvO<sub>2</sub> and SaO<sub>2</sub> (r=0.686 and r=0.547, respectively; p&lt;0.001 for each).</li> <li>• Postoperative complications were observed in 7 patients and the mean cerebral rSO<sub>2</sub> values during the first 4 hours postoperative were lower in patients with complications compared to those without (29% vs. 45%; p&lt;0.001).</li> <li>• Reduced postoperative cerebral rSO<sub>2</sub> may be predictive of postoperative complications.</li> </ul>

**Table 2. Summary of Key Cerebral INVOS™ Studies (cont'd.)**

Hoffman et al. (2013)	<ul style="list-style-type: none"> <li>• Twenty-one patients who had undergone stage 1 palliation of HLHS were studied to determine if early postoperative rSO<sub>2</sub> was related to later neurodevelopmental performance.</li> <li>• Patients with low to abnormal visual-motor integration scores had significantly reduced postoperative rSO<sub>2</sub> (64% vs. 68%; p&lt;0.05).</li> <li>• The time at rSO<sub>2</sub> &lt;45% and &lt;55% was related to low visual-motor integration and neurodevelopmental index scores.</li> <li>• Perioperative cerebral rSO<sub>2</sub> monitoring can detect hypoxic-ischemic conditions associated with neurological injury and reduced neurodevelopmental performance.</li> <li>• Efforts to avoid cerebral hypoxia are likely to improve outcomes in this high-risk population.</li> </ul>
-----------------------	---

### Cerebral oxygen saturation monitoring with INVOS™ technology

Oxygen saturation monitoring with INVOS™ technology may aid the clinician in the detection of inadequate cerebral perfusion. Inadequate cerebral perfusion may lead to neurological injury and adverse outcomes. A number of studies have examined the clinical value of cerebral oxygen saturation monitoring with INVOS™ technology in the postoperative congenital heart patient.<sup>5,6,8,24-27</sup>

In general, these studies have suggested that these patients are highly susceptible to reduced cerebral oxygenation and that reductions in cerebral rSO<sub>2</sub> are associated with increased incidence of ischemic lesions, neurological injury, and adverse neurodevelopmental outcome. Postoperative manipulations designed to improve cerebral perfusion and thus rSO<sub>2</sub> may help reduce the likelihood of poor clinical outcome. Key studies describing the use of perioperative cerebral oxygen saturation monitoring with INVOS™ technology in the congenital heart patient are shown in Table 2.

### Cerebral oxygen saturation monitoring with INVOS™ technology as an indicator of global / systemic tissue oxygenation

Cerebral NIRS oximetry may be used as a surrogate for mixed or central venous oxygen saturations and in doing so used to assess the global oxygen supply and demand relationship. However, the correlation between cerebral NIRS oximetry and central or mixed venous oxygen saturations may vary. Changes in PaCO<sub>2</sub> / pH “uncouple” cerebral blood flow from metabolism, altering cerebral oxygenation irrespective of changes in systemic hemodynamics. In a low cardiac output state blood flow is redistributed to maintain perfusion of the most vital organs such as the brain. Thus, cerebral rSO<sub>2</sub> may underestimate the degree of global disturbances in the oxygen supply and demand relationship based on an assessment of central or mixed oxygen saturations.

### Somatic oxygen saturation monitoring with INVOS™ technology

In addition to cerebral monitoring, rSO<sub>2</sub> monitoring of peripheral tissue may provide insight into the global balance between oxygen saturation supply and demand. The rationale behind this strategy lies in the fact that as CO falls and becomes limited, blood flow is redistributed away from less vital organs (egs., renal and mesentery) to maintain perfusion of the most vital organs. As a result, renal and mesenteric oxygen saturations fall earlier and to a greater extent than cerebral oxygen saturations. A number of studies have used somatic monitoring (abdominal or splanchnic and renal) to demonstrate peripheral changes in tissue oxygenation following congenital heart surgery.<sup>26-31</sup> In general, these studies have indicated that perioperative monitoring of somatic oxygen saturation with INVOS™ technology in the congenital heart patient may allow for an even earlier appreciation of a falling CO. Key studies describing the use of INVOS™ technology in this setting are shown in Table 3.

### Incorporating INVOS™ into clinical practice

Existing clinical experience and published studies support the benefits of incorporating oxygen saturation monitoring with INVOS™ technology into the routine perioperative management of pediatric cardiac critical care patients. The INVOS™ monitor provides an easy to use, noninvasive, real-time measurement of cerebral/somatic oxygen saturation. Overall, strategies to improve dropping INVOS™ values focus on improving oxygen balance by increasing perfusion pressure and/or arterial oxygen content and reducing oxygen demand (see Table 4).<sup>32</sup>

### Limitations of NIRS oximetry

While NIRS oximetry appears to provide a valuable tool for real-time analysis of tissue oxygenation, there remain several unanswered questions regarding the overall clinical utility of rSO<sub>2</sub> monitoring. In a recent study comparing



**Table 3. Summary of Key Somatic INVOS™ Studies**

Kaufman et al. (2008)	<ul style="list-style-type: none"> <li>• Twenty neonates and infants were studied within 48 hours of cardiac surgery to determine the correlation between abdominal NIRS oximetry and gastric tonometry.</li> <li>• Evaluation of anterior abdominal and renal rSO<sub>2</sub> revealed strong correlation between abdominal rSO<sub>2</sub> and intramucosal gastric pH (<math>r=0.79</math>; <math>p&lt;0.0001</math>) as well as between abdominal rSO<sub>2</sub> and SvO<sub>2</sub> (<math>r=0.89</math>; <math>p&lt;0.0001</math>).</li> <li>• There was a strong negative correlation between abdominal rSO<sub>2</sub> and serum lactate (<math>r=-0.77</math>; <math>p&lt;0.0001</math>).</li> </ul>
Chakravarti et al. (2009)	<ul style="list-style-type: none"> <li>• Twenty-three children were studied to determine the relationship between rSO<sub>2</sub> and blood lactate levels following surgery.</li> <li>• Measures included cerebral, splanchnic, renal, and muscle rSO<sub>2</sub> monitoring for 24 hours postoperatively with blood lactate levels measured at 0, 2, 4, 6, and 24 hours.</li> <li>• Cerebral rSO<sub>2</sub> showed the strongest negative correlation with blood lactate (<math>r=-0.74</math>; <math>p&lt;0.0001</math>) followed by splanchnic, renal, and muscle (<math>r=-0.61</math>, <math>-0.57</math>, and <math>-0.48</math>, respectively [all <math>p&lt;0.0001</math>]).</li> <li>• The correlation between rSO<sub>2</sub> and blood lactate was further improved by averaging cerebral and renal values (<math>r=-0.82</math>; <math>p&lt;0.0001</math>).</li> <li>• Averaged cerebral and renal rSO<sub>2</sub> <math>\leq 65\%</math> predicted a lactate level <math>\geq 3.0</math> mmol/L with a sensitivity of 95% and specificity of 83% (<math>p=0.0001</math>).</li> <li>• Averaged cerebral and renal rSO<sub>2</sub> predicts hyperlactatemia in acyanotic children after congenital heart surgery.</li> <li>• Noninvasive, continuous oxygen saturation monitoring may facilitate the identification of global hypoperfusion due to low cardiac output.</li> </ul>
Owens et al. (2011)	<ul style="list-style-type: none"> <li>• Forty infants undergoing biventricular repair were studied to determine if persistent low renal rSO<sub>2</sub> was associated with acute kidney injury.</li> <li>• Continuous renal rSO<sub>2</sub> monitoring was employed for 48 hours postoperatively.</li> <li>• Subjects with low renal oximetry (<math>&lt;50\%</math> for <math>&gt;2</math> hours) had significantly higher postoperative peak creatinine levels (0.8 vs. 0.5; <math>p=0.003</math>) by 48 hours along with higher incidence of acute kidney injury (50% vs. 3%; <math>p=0.003</math>).</li> <li>• Patients with reduced renal oximetry required more ventilator days and greater vasoactive support and had higher lactate levels.</li> <li>• Prolonged reductions in renal oxygen saturation appeared to correlate with renal dysfunction, decreased systemic oxygen delivery, and overall postoperative course.</li> </ul>
Bronicki et al. (2013)	<ul style="list-style-type: none"> <li>• Twenty-three patients were monitored for 30 minutes before and after extubation following repair of tetralogy of Fallot to determine the effects of extubation and loading of respiratory muscles on the distribution of cardiac output.</li> <li>• Mesenteric rSO<sub>2</sub> fell from 74% to 72% (<math>p=0.04</math>); there were no changes in renal rSO<sub>2</sub>.</li> <li>• Following extubation, the presumed increase in respiratory pump perfusion, as well as the concurrent increase in cerebral perfusion, came at the expense of mesenteric perfusion.</li> </ul>
Hansen et al. (2013)	<ul style="list-style-type: none"> <li>• Thirty-two patients with congenital heart disease who underwent superior cavopulmonary anastomosis were studied to assess the clinical value of NIRS monitoring.</li> <li>• Somatic SO<sub>2</sub> (sSO<sub>2</sub>) decreased from 77% during the early postoperative course to 68% at the end of the 48-hour monitoring period.</li> <li>• Postoperative complications were observed in 7 patients and the mean sSO<sub>2</sub> values during the first 4 hours postoperative were lower in patients with complications compared to those without (70% vs. 80%; <math>p=0.004</math>).</li> </ul>
Oyaizu et al. (2013)	<ul style="list-style-type: none"> <li>• Thirty-four patients with congenital heart disease were studied to evaluate the change in liver oxygen saturation following pulmonary artery banding (PAB).</li> <li>• PAB was considered effective (<math>n=26</math>) or ineffective (<math>n=8</math>) based on the need for secondary PAB correction following the first surgery.</li> <li>• Regional liver tissue oxygenation was significantly increased compared to baseline in patients with effective PAB (72% vs. 60%).</li> <li>• Changes in liver rSO<sub>2</sub> reflected the effectiveness of PAB; monitoring liver rSO<sub>2</sub> may be a good indicator to decide the tightness of the PAB.</li> </ul>

NIRS oximetry devices in healthy adult volunteers, Bickler et al. showed that while the 5 NIRS monitors studied were all responsive to hypoxia-induced desaturation there was significant subject-to-subject variability and measurement bias across devices. These authors concluded that while all marketed NIRS oximeters assume a constant ratio of arterial to venous blood in the sampled region (though different devices assume different ratios); this ratio might in fact be relatively dynamic in the face of changing oxygen balance. In fact, knowledge of this dynamic physiologic environment prevents any tissue oximeter from being “absolute” in the mixed clinical environment. Given this variability, it is difficult to establish a critical cutoff for  $rSO_2$ , either in terms of absolute percent saturation or reduction from baseline. There are data to suggest, however, that the  $rSO_2$  threshold for neurological injury may be approximately 45-50% (i.e., an oxygen extraction ratio of 50-55%). In fact, brain lactate levels begin to rise and functional disturbances on EEG begin to be seen as cerebral  $rSO_2$  falls below 50%.<sup>33,34</sup> Based on these data, interventions designed to maintain  $rSO_2$  above this critical point may lessen the incidence of ischemic neurological injury.

### Conclusions

Despite advances in techniques and improved mortality, surgical repair of congenital heart defects remains associated with significant risk of poor neurodevelopmental outcomes. Studies to date suggest that one of the most important factors in determining the clinical outcome of these complex surgical procedures is the incidence and severity of cerebral hypoperfusion during the early postoperative period. The INVOS™ oximeter provides noninvasive and continuous monitoring of cerebral/somatic oxygen supply and demand and thus aids in the clinician’s real-time evaluation of the adequacy of cerebral/somatic tissue perfusion. Use of oxygen saturation monitoring with INVOS™ technology during this critical period allows for the rapid identification of inadequate cerebral perfusion that might otherwise go undetected. Interventions based on  $rSO_2$  analysis have been shown to improve cerebral  $rSO_2$  values and may lead to improved clinical outcomes in this vulnerable patient population. Similarly, monitoring of peripheral tissue oxygen saturation with INVOS™ technology may help identify impaired systemic oxygen delivery, which can lead to organ dysfunction and contribute to adverse postoperative outcomes.

**Table 4. Summary of  $rSO_2$  and Oxygen Transport Balance**

$rSO_2$			
<ul style="list-style-type: none"> <li>Increases with rise in oxygen delivery or fall in oxygen demand.</li> <li>Decreases when oxygen delivery falls or there is an uncompensated rise in oxygen demand.</li> </ul>			
Oxygen delivery/supply		Oxygen demand/consumption	
Influenced by:		Increased by:	Decreased by:
<ul style="list-style-type: none"> <li>Oxygen content               <ul style="list-style-type: none"> <li>– Hemoglobin concentration</li> <li>– Hemoglobin saturation</li> </ul> </li> </ul>	<ul style="list-style-type: none"> <li>Cardiac output               <ul style="list-style-type: none"> <li>– Optimize HR</li> <li>– Optimize preload</li> <li>– Intropy</li> <li>– Reduce afterload</li> </ul> </li> </ul>	<ul style="list-style-type: none"> <li>Fever, shivering</li> <li>Seizures</li> <li>Wakefulness, anxiety, and pain</li> </ul>	<ul style="list-style-type: none"> <li>Normo- or hypothermia</li> <li>Sedation and muscle relaxants</li> <li>Mechanical ventilation</li> </ul>

## REFERENCES

- Hoffman JI, Kaplan S. The incidence of congenital heart disease. *J Am Coll Cardiol.* 2002;39:1890-1900.
- Bronicki RA, Checchia PA. Advances in pediatric cardiac intensive care. *Pediatr Clin North Am.* 2013;60:655-667.
- Li J, Zhang G, Holtby HM, et al. Inclusion of oxygen consumption improves the accuracy of arterial and venous oxygen saturation interpretation after the Norwood procedure. *J Thorac Cardiovasc Surg.* 2006;131:1099-107.
- Li J, Zhang G, McCrindle BW, et al. Profiles of hemodynamics and oxygen transport derived by using continuous measured oxygen consumption after the Norwood procedure. *J Thorac Cardiovasc Surg.* 2007;133:441-448.
- Li J, Zhang G, Holtby H, et al. The influence of systemic hemodynamics and oxygen transport on cerebral oxygen saturation in neonates after the Norwood procedure. *J Thorac Cardiovasc Surg.* 2008;135:83-90,90 e81-82.
- Dent CL, Spaeth JP, Jones BV, et al. Brain magnetic resonance imaging abnormalities after the Norwood procedure using regional cerebral perfusion. *J Thorac Cardiovasc Surg.* 2005; 130(6): 1523-30.
- Kussman BD, Wypij D, Laussen PC, et al. Relationship of intraoperative cerebral oxygen saturation to neurodevelopmental outcome and brain magnetic resonance imaging at 1 year of age in infants undergoing biventricular repair. *Circulation.* 2010;122:245-254.
- Hoffman GM, Brosig CL, Mussatto KA, et al. Perioperative cerebral oxygen saturation in neonates with hypoplastic left heart syndrome and childhood neurodevelopmental outcome. *J Thorac Cardiovasc Surg.* 2013;146:1153-64.
- Hoffman GM, Mussatto KA, Brosig CL, et al. Systemic venous oxygen saturation after the Norwood procedure and childhood neurodevelopmental outcome. *J Thorac Cardiovasc Surg.* 2005;130:1094-1100.
- Connors AF, McCaffree DR, Gray BA. Evaluation of right-heart catheterization in the critically ill patient without acute myocardial infarction. *N Engl J Med.* 1983;308:263-267.
- Bronicki RA. Venous oximetry and the assessment of oxygen transport balance. *Pediatr Crit Care Med.* 2011;12:S21-26.
- Tibby SM, Hatherill M, Marsh MJ, Murdoch IA. Clinicians' abilities to estimate cardiac index in ventilated children and infants. *Arch Dis Child.* 1997;77:516-518.
- Tibby SM, Hatherill M, Murdoch IA. Capillary refill and core-peripheral temperature gap as indicators of haemodynamic status in paediatric intensive care patients. *Arch Dis Child.* 1999;80:163-166.
- Lobos AT, Lee S, Menon K. Capillary refill time and cardiac output in children undergoing cardiac catheterization. *Pediatr Crit Care Med.* 2012;13:136-140.
- Owen-Reece H, Smith M, Elwell CE, Goldstone JC. Near infrared spectroscopy. *Br J Anaesth.* 1999;82:418-426.
- Murkin JM, Arango M. Near-infrared spectroscopy as an index of brain and tissue oxygenation. *Br J Anaesth.* 2009;103 Suppl 1 i3-13.
- Dullenkopf A, Frey B, Baenziger O, et al. Measurement of cerebral oxygenation state in anaesthetized children using the INVOS 5100 cerebral oximeter. *Paediatr Anaesth.* 2003;13:384-391.
- Daubeney PE, Pilkington SN, Janke E, et al. Cerebral oxygenation measured by near-infrared spectroscopy: comparison with jugular bulb oximetry. *Ann Thorac Surg.* 1996;61:930-934.
- Tortoriello TA, Stayer SA, Mott AR, et al. A noninvasive estimation of mixed venous oxygen saturation using near-infrared spectroscopy by cerebral oximetry in pediatric cardiac surgery patients. *Paediatr Anaesth.* 2005;15:495-503.
- Li J, Van Arsdell GS, Zhang G, et al. Assessment of the relationship between cerebral and splanchnic oxygen saturations measured by near-infrared spectroscopy and direct measurements of systemic haemodynamic variables and oxygen transport after the Norwood procedure. *Heart.* 2006;92:1678-1685.
- Nagdyman N, Ewert P, Peters B, et al. Comparison of different near-infrared spectroscopic cerebral oxygenation indices with central venous and jugular venous oxygenation saturation in children. *Paediatr Anaesth.* 2008;18:160-166.
- Ricci Z, Garisto C, Favia I, et al. Cerebral NIRS as a marker of superior vena cava oxygen saturation in neonates with congenital heart disease. *Paediatr Anaesth.* 2010;20:1040-1045.
- Moreno GE, Pilan ML, Manara C, et al. Regional venous oxygen saturation versus mixed venous saturation after paediatric cardiac surgery. *Acta Anaesthesiol Scand.* 2013;57:373-379.
- Phelps HM, Mahle WT, Kim D, et al. Postoperative cerebral oxygenation in hypoplastic left heart syndrome after the Norwood procedure. *Ann Thorac Surg.* 2009;87:1490-1494.
- Horvath R, Shore S, Schultz SE, et al. Cerebral and somatic oxygen saturation decrease after delayed sternal closure in children after cardiac surgery. *J Thorac Cardiovasc Surg.* 2010;139:894-900.
- Bronicki RA, Checchia PA, Anas NG, et al. Cerebral and somatic oxygen saturations after repair of tetralogy of Fallot: effects of extubation on regional blood flow. *Ann Thorac Surg.* 2013;95:682-686.
- Hansen JH, Schlagen J, Armbrust S, et al. Monitoring of regional tissue oxygenation with near-infrared spectroscopy during the early postoperative course after superior cavopulmonary anastomosis. *Eur J Cardiothorac Surg.* 2013;43:e37-43.
- Kaufman J, Almodovar MC, Zuk J, et al. Correlation of abdominal site near-infrared spectroscopy with gastric tonometry in infants following surgery for congenital heart disease. *Pediatr Crit Care Med.* 2008;9:62-68.
- Chakravarti SB, Mitnacht AJ, Katz JC, et al. Multisite near-infrared spectroscopy predicts elevated blood lactate level in children after cardiac surgery. *J Cardiothorac Vasc Anesth.* 2009;23:663-667.
- Owens GE, King K, Gurney JG, Charpie JR. Low renal oximetry correlates with acute kidney injury after infant cardiac surgery. *Pediatr Cardiol.* 2011;32:183-188.
- Oyaizu T, Hatori E, Katori N, et al. Change in regional liver oxygen saturation may reflect adequacy of pulmonary artery banding. *Am Soc Anesth Annual Meeting.* 2013: A4187 (Abstract).
- Hoffman GM, Ghanayem NS, Tweddell JS. Noninvasive assessment of cardiac output. *Semin Thorac Cardiovasc Surg Pediatr Card Surg Annu.* 2015:12-21.
- Levy WJ, Levin S, Chance B. Near-infrared measurement of cerebral oxygenation: Correlation with electroencephalographic ischemia during ventricular fibrillation. *Anesthesiology.* 1995;83:738-746.
- Kurth CD, Levy WJ, McCann J. Near-infrared spectroscopy cerebral oxygen saturation thresholds for hypoxia-ischemia in piglets. *J Cereb Blood Flow Metab.* 2002;22:335-341.

The INVOS™ monitoring system should not be used as the sole basis for diagnosis or therapy and is intended only as an adjunct in patient assessment. Reliance on the INVOS™ system alone for detecting cerebral desaturation events is not recommended.

# Medtronic

© 2022 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. EMEA-PM-2200051-invos-pediatric-cardiac-critical-care-white-paper-en-we-7463037

[medtronic.co.uk/pmri](https://www.medtronic.co.uk/pmri)