



Medtronic

Quick reference guide

# Monitoring premature neonates

INVOS™ regional oximeter





Premature neonates may have immature organs with suboptimal functioning – putting them at risk for poor oxygen exchange and perfusion – which may be missed by routine vital sign monitoring that does not capture oxygen supply, demand, or content at the organ level.

Information from INVOS™ near-infrared spectroscopy (NIRS) monitoring of the cerebral, peri-renal, and splanchnic sites – individually, or in combination – can integrate with routine vital signs and labs to provide a more comprehensive picture of organ-specific oxygen extraction and utilization to inform hemodynamic management, ventilation, and resuscitation.


It is intended for use as an adjunct trend monitor of regional hemoglobin oxygen saturation of blood in tissue beneath the sensor in any individual, including neonates weighing less than 2.5 kg.

# Blood pressure and pulse oximetry do not always reflect perfusion and oxygenation at the tissue level

|            | <div>Hemodynamic instability</div> <div>  </div>   | <div>Suboptimal respiratory support and resuscitation</div> <div>  </div>  |
|------------|---|---|
| Causes     | <ul style="list-style-type: none"> <li>• Poor transition</li> <li>• Low cardiac output</li> <li>• Impaired autoregulation</li> <li>• Anemia</li> <li>• Shock, sepsis</li> </ul> | <ul style="list-style-type: none"> <li>• Impaired lung function</li> <li>• Poor gas exchange capabilities</li> <li>• Shunting</li> </ul>  |
| Management | <ul style="list-style-type: none"> <li>• Volume expansion</li> <li>• Vasopressors / dilators</li> <li>• Inotropes</li> <li>• Blood transfusion</li> </ul>                       | <ul style="list-style-type: none"> <li>• Supplemental oxygen</li> <li>• Surfactant administration</li> <li>• Bag mask support in L&amp;D</li> <li>• Non-invasive ventilation</li> <li>• Invasive ventilation</li> </ul> |
| Challenges | <p>Maintaining an arbitrary gestational age-based blood pressure target <b>does not guarantee optimal organ perfusion</b></p>   | <p>Pulse oximetry is an indirect measure of oxygen delivery to the organs; <b>too much or too little can cause long term damage</b></p>   |

## The INVOS™ monitoring system provides continuous noninvasive measurement of organ-specific oxygenation

Measures the percentage of oxyhemoglobin or the venous reserve capacity following tissue oxygen extraction in a specific region/tissue under the sensor.



There is both a shallow pathway (shorter distance photodetector) and a deep pathway (farther distance photodetector). The short pathway is subtracted from the deep pathway so that the tissue of interest at 2.5 cm below the probe is interrogated, avoiding contamination from the skin, bone, and dura when monitoring  $rSO_2$  of the brain.<sup>1</sup>



Provides noninvasive cerebral/somatic monitoring that measures regional, real-time tissue oxygen extraction and utilization



Provides oxygen saturation from vascular beds to assess organs individually, or in combination to track brain/body perfusion shifts



Can help clinicians identify ischemic threats in the cerebral and peripheral circulatory systems so they can intervene earlier<sup>2</sup>

# Regional oxygen saturation (rSO<sub>2</sub>) reflects post-extraction oxygen balance of the tissue

Allows for more personalized, patient-centered approaches – helping delineate the underlying pathophysiology more effectively when added to clinical assessments.<sup>2</sup>

## Regional oxygen saturation (rSO<sub>2</sub>)


- Venous-weighted :
  - ~ 25% arterial contribution
  - ~ 75% venous (post-extraction) contribution
- The rSO<sub>2</sub> value reflects the balance between oxygen supply, demand, and content in the tissue – how much oxygen is left after the tissue extracts a portion of what is available<sup>1</sup>

## Factors affecting rSO<sub>2</sub> values

- Blood pressure (supply)
- Cardiac output (supply)
- Shunting (supply, content)
- Carbon dioxide (supply)
- Transfusion (content)
- Brain activity (demand)
- Temperature (demand)



# Standard of care vital signs and labs are incomplete markers of tissue perfusion and oxygenation



Courtesy of Scott Duncan, M.D.







The  $rSO_2$  measurement may be an **early and more sensitive warning** to a perfusion change

Routine clinical monitors do not detect changes in systemic vascular resistance which result in redistribution of cardiac output

Monitoring site-specific perfusion often provides an earlier warning of developing pathology and deterioration than systemic measures or lab tests<sup>2</sup>


# Monitoring rSO<sub>2</sub> at single sites or in combination can reveal changes in cardiac output and impending shock

|  | What do rSO <sub>2</sub> values and trends look like at this site?  | What can monitoring rSO <sub>2</sub> tell you?   |
|--|---|--|
| <b>Cerebral</b><br>     | <ul style="list-style-type: none"> <li>• High flow</li> <li>• High oxygen extraction</li> <li>• Stable values</li> <li>• Typical rSO<sub>2</sub> = 60 to 80%<sup>4</sup></li> </ul>                               | Monitors for <b>changes in cerebral perfusion, oxygen content, and oxygen extraction</b>                   |
| <b>Peri-renal</b><br>  | <ul style="list-style-type: none"> <li>• Lower flow</li> <li>• Lower oxygen extraction</li> <li>• Variable values</li> <li>• Typical rSO<sub>2</sub> = 5 to 15 points higher than cerebral<sup>4</sup></li> </ul> | Reveals changes in organ-area perfusion related to <b>low cardiac output or compensated shock states</b>   |
| <b>Splanchnic</b><br> |   |  |
| <b>Multi-site</b><br> | Somatic rSO <sub>2</sub> trending lower than cerebral rSO <sub>2</sub> may indicate abnormal pathology  | Combining cerebral and somatic monitoring can indicate <b>changes in cardiac distribution</b> <sup>1</sup> |

To avoid pressure sores, do not apply external pressure (e.g., headbands, wraps, or tape) to the sensor. Medtronic recommends using a new sensor every 24 hours or if adhesive does not seal the sensor to the skin adequately.

# rSO<sub>2</sub> changes can alert clinicians to investigate changes in cardiovascular status, oxygen transport, respiratory status, and blood glucose<sup>3</sup>

Although the NIRS monitor may indicate an abnormality, **the clinician should understand all the possible underlying causes** of rSO<sub>2</sub> changes



## Abbreviations:

rSO<sub>2</sub>, regional tissue oxygen saturation of hemoglobin

CO, cardiac output

SVC, superior vena cava

P<sub>aw</sub>, mean airway pressure

MV, minute ventilation

PDA, patent ductus arteriosus

[Hb], blood hemoglobin concentration

RBC, red blood cells

SaO<sub>2</sub>, arterial hemoglobin saturation

FiO<sub>2</sub>, inspired oxygen fraction




PCO<sub>2</sub>, partial pressure of carbon dioxide



## Published case reports to illustrate use in a variety of clinical circumstances

There is a learning curve when using NIRS initially.

As you begin selecting babies to study, **start by choosing a neonate that has an identified and understood physiological process**, which allows for confirmation by NIRS. This helps clinicians develop an understanding and comfort level with this technology.

|        |  |  |
|--------|--|--|
| Case 1 | <b>Hemodynamic management</b><br>   | <b>NIRS as a guide to manage hemodynamic instability and shock</b><br><br>25-week gestational age neonate with early onset sepsis                          |
| Case 2 | <b>Ventilation management</b><br> | <b>NIRS as a signal of decreased perfusion due to excessive mean airway pressure</b><br><br>28-week gestational age neonate on HFOV                        |
| Case 3 | <b>Resuscitation</b><br>          | <b>NIRS as a guide to cardiopulmonary resuscitation after delivery</b><br><br>24-week gestational age male neonate requiring immediate support after birth |


# Hemodynamic management

## Case 1: NIRS as a guide to manage hemodynamic instability and shock<sup>5</sup>



Simultaneous near-infrared spectroscopy (NIRS) and amplitude-integrated electroencephalography (aEEG): Dual use of brain monitoring techniques improves our understanding of physiology


Variane GFT, Chock VY, Netto A, Pietrobon RFR, Van Meurs KP.  
(BRAZIL, USA)



|           |   |
|-----------|---|
| Patient   | <ul style="list-style-type: none"><li>• 25-week gestational age neonate weighing 610 g</li><li>• Showed <b>early indicators of hemodynamic instability</b></li><li>• Diagnosed with early onset sepsis and intubated to manage respiratory distress syndrome</li></ul>  |
| Situation | <ul style="list-style-type: none"><li>• <b>3 days post-birth: experienced severe hypotension and septic shock</b></li><li>• Peri-renal rSO<sub>2</sub> decreased 30% from baseline</li><li>• 40 minutes later, cerebral rSO<sub>2</sub> decreased from 75% to 15%</li><li>• <b>After intervening with fluids and inotropes, both peri-renal and cerebral rSO<sub>2</sub> improved</b></li></ul> |
| NIRS      | <ul style="list-style-type: none"><li>• Although mean arterial blood pressure remained steady at 32mmHg, peri-renal rSO<sub>2</sub> began to drop steadily indicating impending instability</li></ul>   |

# Ventilation management

## Case 2: NIRS as a signal of decreased perfusion due to excessive mean airway pressure ( $\bar{P}_{aw}$ )



|           |  |
|-----------|--|
| Patient   | <ul style="list-style-type: none"><li>• 28-week gestational age neonate admitted to NICU with <b>respiratory distress syndrome</b> and to rule out sepsis</li></ul>  |
| Situation | <ul style="list-style-type: none"><li>• High-frequency oscillatory ventilation (HFOV) support was initiated with a mean airway pressure (<math>\bar{P}_{aw}</math>) of 16 cmH<sub>2</sub>O and subsequently increased to 17 cm H<sub>2</sub>O</li><li>• Steady decline in peri-renal rSO<sub>2</sub> over the first hour and a half</li><li>• <b>Decrease in <math>\bar{P}_{aw}</math> from 17cmH<sub>2</sub>O to 15cmH<sub>2</sub>O lead to an immediate increase in peri-renal rSO<sub>2</sub></b></li></ul> |
| NIRS      | <ul style="list-style-type: none"><li>• NIRS indicated the peri-renal tissue was not adequately perfused</li><li>• <b>With that information, they determined that the increased <math>\bar{P}_{aw}</math> was putting too much pressure on the heart decreasing oxygen delivery to the peri-renal tissue and clinicians were able to decrease ventilation support</b></li></ul>  |


# Resuscitation

## Case 3: NIRS as a guide to cardiopulmonary resuscitation after delivery<sup>6</sup>



Respiratory function and near infrared spectroscopy recording during cardiopulmonary resuscitation in an extremely preterm newborn

Li ES, Cheung PY, Pichler G, Aziz K, Schmölzer GM. (Canada, Austria)




|           |   |
|-----------|---|
| Patient   | <ul style="list-style-type: none"><li>• 24-week gestational age male neonate weighing 650 g</li><li>• Apgar scores were 3, 2, and 7 at 1, 5, and 10 minutes post-birth</li><li>• <b>Neonate was “floppy”, cyanotic, and apneic</b> after delayed cord clamping (DCC)</li></ul>  |
| Situation | <ul style="list-style-type: none"><li>• Positive pressure mask ventilation did not improve the neonate’s condition</li><li>• <b>CPR was initiated and neonate was successfully intubated after two attempts between chest compressions (CC)</b></li><li>• Return of spontaneous circulation (ROSC) was observed after the second intubation attempt</li></ul>   |
| NIRS      | <ul style="list-style-type: none"><li>• <b>Despite an adequate arterial oxygen saturation of ~ 95% to 100%, cerebral rSO<sub>2</sub> remained very low at ~15% until ROSC</b></li><li>• Cerebral rSO<sub>2</sub> may provide a more complete picture of whether resuscitation efforts are sufficient for restoring cerebral perfusion and oxygenation</li></ul> |

# Summary

All three of these cases are examples of  $rSO_2$  as an early alert of changes in oxygen supply and demand which could be indicative of possible impending deterioration or the impact of interventions.

Monitoring of cerebral, peri-renal, and splanchnic  $rSO_2$  – individually, or in combination – may be used in the NICU to inform hemodynamic management, ventilation, and resuscitation.





The INVOS™ monitoring system should not be used as the sole basis for diagnosis or therapy and is intended only as an adjunct in patient assessment. Reliance on the INVOS™ system alone for detecting cerebral desaturation events is not recommended.

INVOS™ System and device discussed in this guide is referring to the INVOS™ 7100 system. The sensor shown in this guide may also be used with the INVOS™ 5100C system.

## References

1. Marin T, Moore J. Understanding Near-Infrared Spectroscopy *J. Adv Neonatal Care* 2011 Dec;11(6):382-8
2. Mintzer JP, Moore JE. Regional tissue oxygenation monitoring in the neonatal intensive care unit: evidence for clinical strategies and future directions. *Pediatr Res*. 2019;86(3):296-304.
3. Pellicer A, Greisen G, Benders M, et al. The SafeBoosC phase II randomised clinical trial: a treatment guideline for targeted near-infrared-derived cerebral tissue oxygenation versus standard treatment in extremely preterm infants. *Neonatology*. 2013;104(3):171-178.
4. McNeill S, Gatenby JC, McElroy S, Engelhardt B. Normal cerebral, renal and abdominal regional oxygen saturations using near-infrared spectroscopy in preterm infants. *J Perinatol*. 2011;31(1):51-57. doi:10.1038/jp.2010.71
5. Variane GFT, Chock VY, Netto A, Pietrobon RFR, Van Meurs KP. Simultaneous Near-Infrared Spectroscopy (NIRS) and Amplitude-Integrated Electroencephalography (aEEG): Dual Use of Brain Monitoring Techniques Improves Our Understanding of Physiology. *Front Pediatr*. 2020;7:560.
6. Li ES, Cheung PY, Pichler G, Aziz K, Schmölzer GM. Respiratory function and near infrared spectroscopy recording during cardiopulmonary resuscitation in an extremely preterm newborn. *Neonatology*. 2014;105(3):200-204.

This resource is intended for educational purposes only. It is not intended to provide comprehensive or patient-specific clinical practice recommendations for rSO<sub>2</sub> monitoring technology. The clinical choices discussed in this text may or may not be consistent with your own patient requirements, your clinical practice approaches, or guidelines for practice that are endorsed by your institution or practice group. It is the responsibility of each clinician to make his/her own determination regarding clinical practice decisions that are in the best interest of patients. Readers are advised to review the current product information, including the Indications for use currently provided by the manufacturer. Neither the publisher, authors, nor Covidien LP, a Medtronic company, assumes any responsibility for any injury and or damage to persons or property resulting from information provided in this text.

# Medtronic

Important: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2022 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. EMEA-PM-2100039 -invo- premature-neonates-quick-reference-guide-en-we-7455372

[medtronic.co.uk/pmri](https://www.medtronic.co.uk/pmri)