

Medtronic

Symplicity Spyral™
renal denervation system






Setting the
standard
in renal denervation

The Symplicity™ blood pressure procedure delivers
significant, safe, and sustained BP reductions.¹⁻⁴

When typical hypertension treatments aren't enough

Lifestyle changes and medications have defined hypertension (HTN) treatment, but they are not always enough to help patients achieve control.



How renal denervation with Symplicity Spyral works:

- The kidneys modulate the sympathetic tone via the renal nerves to control blood pressure
- Symplicity Spyral supplies precisely controlled and targeted radiofrequency energy to the renal nerves⁶
- The procedure safely disrupts the overactive sympathetic signaling between the kidneys and brain, to reduce blood pressure⁶

Turn to an approach recommended by societies

2024 European Society of Cardiology (ESC) Guidelines and 2023 European Society of Hypertension (ESH) Guidelines now recommend renal denervation as a safe and effective treatment option in uncontrolled and resistant hypertension^{7,8}.

European Guidelines: Patient selection criteria^{7,8}

Uncontrolled BP



Patients who have **uncontrolled hypertension** with eGFR >40 ml/min/1.73m² on fewer than 3 antihypertensive drugs with increased CV Risk

Resistant Hypertension



Patients with resistant hypertension with eGFR >40 ml/min/1.73m² who have uncontrolled BP despite a 3-drug combination*

* Including a thiazide or thiazide-like diuretic

Additional recommendations^{7,8}



Selection of patients should be done in a **shared decision-making process** taking into account **patient preference**



Renal denervation should be performed in experienced RDN centers after a **multidisciplinary assessment****

** A multidisciplinary team is defined as two or more healthcare professionals working together. This can include, but is not limited to, a proceduralist and a non-invasive physician.



Based on definitions, ESC classification IIb and ESH classification I vary slightly on shared decision making and experienced centers.

Setting the standard in renal denervation

Most used
RDN
system⁹

>30,000

patients treated globally with
the Symplicity RDN system⁹


The Symplicity™ blood pressure procedure is a simple, minimally-invasive approach to achieve BP reductions, independent of patient adherence. Precisely controlled and targeted **radiofrequency (RF) energy** is delivered to the renal nerves⁶ with no permanent implant left behind.

Symplicity Spyral™
Multi-Electrode Renal
Denervation Catheter

- 4 F catheter,[†] compatible with 6 F guide catheter, 0.014" guidewire
- Four independently controlled electrodes
- Shaped atraumatic tip with radiopaque marker

Simple and versatile

- One catheter size fits vessels 3-8mm¹⁰
- Easy-to-use, **plug-and-play** design
- **Nonocclusive design** allows for continuous blood flow to naturally protect the vessel wall⁶



Safe and precise

Safe

- Radiofrequency energy preferentially heats fat tissue where renal nerves are located and **avoids non-target structures**¹¹
- Unique, real-time and **responsive algorithm** automatically adjusts power by monitoring temperature and impedance for safe energy distribution⁶

Precise ablation

- Multi-electrode, helical design covers four quadrants simultaneously for a circumferential ablation^{6,10}
- Symplicity Spyral™ allows for denervation throughout the renal anatomy, including the distal branches,¹⁰ where late arriving nerves are accessible and total nerve density is highest^{12,13}

Symplicity G3™ Renal Denervation RF Generator




Significant, safe, and sustained blood pressure reductions¹⁻⁴

Supported by clinical evidence **unmatched** in scope,
quantity and quality.

>5,000 patients enrolled in the
global clinical program^{‡18}

Significant

Consistent blood pressure reductions in 4 sham-controlled trials,
in both the presence and absence of medications^{1,2,14,15}



9
mmHg

OSBP reduction in
patients off and on
medications^{1,2,14,15}

Safe

Excellent safety profile^{1,2}


Pooled data indicated low incidence of procedural and clinical adverse events²

0.4%

major adverse event rate at composite endpoint, including no new incidence of renal artery stenosis (>70%) at 6 months²

Stable kidney function

Minimal impact to kidney function at 3-year follow-up¹⁶



Sustained

Confidence in durable, real-world outcomes with Symplicity™ system, demonstrating lasting blood pressure reductions at 3 years in over 1,300 real-world patients.¹⁷

17
mmHg

reduction in OSBP at 3 years. n=1337¹⁷



RDN016 Symplicity Spyral multi-electrode renal denervation catheter: single-use RDN catheter for use only with the Symplicity G3 renal denervation RF generator

RDNG3A/RDN017 (as approved in your country) Symplicity G3 generator: reusable RF generator, power cable, remote control, and DVI-D output for use with the Symplicity Spyral catheter

RDN019 Mobile cart for Symplicity G3 generator (optional accessory)

Guidewire: 0.014" non-hydrophilic guidewire with a supportive shaft and a floppy tip

Guide catheter: 6 F guide catheter (90 cm max length)

Recommended models:

SB6RDND1K - Sherpa NX™ Balanced 6Fr 55-cm guide catheter with RDND1 curve

LA6RDCK - Launcher™ 6Fr 55-cm guide catheter with RDC curve

LA6IMAK - Launcher™ 6Fr 55-cm guide catheter with IMA curve

SB6IMAK - Sherpa NX™ Balanced 6Fr 55-cm guide catheter with IMA curve

Dispersive Electrode for RDNG3A: Compatible with ValleyLab REM Polyhesive Adult Patient Return Electrode (Model **E7507, E7507-DB**)

Dispersive Electrode for RDN017: A disposable dispersive electrode with a molded connector specified for adult use that complies with IEC 60601-2-2 is required

¹Catheter dimension of 0.052" is average diameter determined during design verification. Upper bound allowable is 0.061".

²Study follow-up is ongoing. Data does not represent follow-up for all patients.

³Not all subjects analyzed have completed follow-up through 3 years. All available data included at the time of this analysis.

¹ Böhm M, Kario K, Kandzari DE, et al. Efficacy of catheter-based renal denervation in the absence of antihypertensive medications (SPYRAL HTN-OFF MED Pivotal): a multicentre, randomised, sham-controlled trial. *Lancet*. May 2, 2020;395(10234):1444-1451.

² Kandzari DE, Townsend RR, Kario K, et al. Safety and Efficacy of Renal Denervation in Patients Taking Antihypertensive Medications. *J Am Coll Cardiol*. 2023;82(19):1809-1823

³ Mahfoud F, Mancia G, Schmieder RE, et al. Outcomes Following Radiofrequency Renal Denervation According to Antihypertensive Medications: Subgroup Analysis of the Global SYMPPLICITY Registry DEFINE. *Hypertension*. 2023;80(8):1759-1770

⁴ Mahfoud F, Kandzari DE, Kario K, et al. Long-term efficacy and safety of renal denervation in the presence of antihypertensive drugs (SPYRAL HTN-ON MED): a randomised, sham-controlled trial. *Lancet*. April 9, 2022;399(10333):1401-1410.

⁵ NCD Risk Factor Collaboration (NCD-RisC). Worldwide trends in hypertension prevalence and progress in treatment and control from 1990 to 2019: a pooled analysis of 1201 population-representative studies with 104 million participants. *Lancet*. September 11, 2021;398(10304):957-980.

⁶ Coates P, Tuney S, Trudel J, Hettrick DA. Time, temperature, power, and impedance considerations for radiofrequency catheter renal denervation. *Cardiovasc Revasc Med*. September 2022;42:171-177.

⁷ Mancia G, Kreutz R, Brunstrom M, et al. Authors/Task Force Members: 2023 ESH Guidelines for the management of arterial hypertension The Task Force for the management of arterial hypertension of the European Society of Hypertension Endorsed by the European Renal Association (ERA) and the International Society of Hypertension (ISH). *J Hypertens*. June 21, 2023.

⁸ McEvoy et al. European Heart Journal 2024 Oct 7;45(38):3912-4018. doi: 10.1093/eurheartj/ehae178

⁹ Medtronic Data on File, RDN Catheter Historic Data, Feb 2025. Data includes both Symplicity Flex and Symplicity Spyral.

¹⁰ Medtronic Symplicity Spyral™ Instructions for Use.

¹¹ Sato Y, Sharp A, Mahfoud F, et al. Translational value of preclinical models for renal denervation: a histological comparison of human versus porcine renal nerve anatomy. *EuroIntervention*. February 6, 2023;18(13):e1120-e1128.

¹² Garcia-Touchard A, Maranillo E, Mompeo B, Sañudo JR. Microdissection of the human renal nervous system : Implications for performing renal denervation procedures. *Hypertension*. October 2020;76(4):1240-1246.

¹³ Struthoff H, Launder L, Hohl M, et al. Histological examination of renal nerve distribution, density, and function in humans. *EuroIntervention*. September 18, 2023;19(7):612-620.

¹⁴ Townsend RR, Mahfoud F, Kandzari DE, et al. Catheter-based renal denervation in patients with uncontrolled hypertension in the absence of antihypertensive medications (SPYRAL HTN-OFF MED): a randomised, sham-controlled, proof-of-concept trial. *Lancet*. November 11, 2017;390(10108):2160-2170.

¹⁵ Kandzari DE, Böhm M, Mahfoud F, et al. Effect of renal denervation on blood pressure in the presence of antihypertensive drugs: 6-month efficacy and safety results from the SPYRAL HTN-ON MED proof-of-concept randomised trial. *Lancet*. June 9, 2018;391(10137):2346-2355.

¹⁶ Symplicity Spyral Renal Denervation System. Sponsor Executive Summary. US FDA Circulatory Systems Devices Panel. Meeting date, August 23, 2023. Pg 121; Figure 58.

¹⁷ Mahfoud F, et al. Outcomes following radiofrequency renal denervation according to antihypertensive medications: subgroup analysis of the Global SYMPPLICITY Registry DEFINE. EuroPCR 2023. Data includes both Symplicity Flex and Symplicity Spyral.

¹⁸ Medtronic data on file, Patients in clinical program, Feb 2025. Data includes both Symplicity Flex and Symplicity Spyral.

This material should not be considered the exclusive source of information, it does not replace or supersede information contained in the device manual(s).

Please note that the intended use of a product may vary depending on geographical approvals.

See the device manual(s) for detailed information regarding the intended use, the implant procedure, indications, contraindications, warnings, precautions, and potential adverse events.

For a MRI compatible device(s), consult the MRI information in the device manual(s) before performing a MRI.

If a device is eligible for eIFU usage, instructions for use can be found at Medtronic's website manuals.medtronic.com.

Manuals can be viewed using a current version of any major internet browser. For best results, use Adobe Acrobat® Reader with the browser.

Medtronic products placed on European markets bear the CE mark and the UKCA mark (if applicable).

For any further information, contact your local Medtronic representative and/or consult Medtronic's websites.

For distribution only in markets where the Symplicity Spyral™ multi-electrode renal denervation catheter and Symplicity G3™ renal denervation RF generator have been approved. Not for distribution in the USA, Japan or France.

Europe

Medtronic International Trading Sàrl.
Route du Molliau 31
Case postale
CH-1131 Tolochenaz
Switzerland
www.medtronic.eu
Tel: +41 (0)21 802 70 00
Fax: +41 (0)21 802 79 00

United Kingdom/Ireland

Medtronic Limited
Building 9
Croxley Park
Hatters Lane
Watford
Herts WD18 8WW
www.medtronic.co.uk
Tel: +44 (0)1923 212213
Fax: +44 (0)1923 241004

medtronic.eu/RDN

17021189-en-gb-emea.01 @2025 Medtronic. All rights reserved. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. TMThird party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

Learn more on renal denervation

Medtronic