

The Sonicision™ Ultrasonic Dissection Devices

**The Science of Performance.
The Ease of a Full Portfolio.**

Product Overview



The industry's first cordless ultrasonic portfolio

Specifications

Ultrasonic Dissector

- Single Patient Use
- 5 mm diameter, 13 cm, 26 cm, 39 cm & 48 cm lengths
- 14.5 mm active blade
- Dual-mode energy button (minimum and maximum)

Reusable Generator

- Reusable for up to 100 sterilization cycles
- Sterilizable using low temperature hydrogen peroxide gas, such as STERRAD and STERRIS V-PRO™*
- Integrated LED indicator provides device feedback

Reusable Battery

- Reusable for up to 100 sterilization cycles
- Advanced lithium polymer battery
- Sterilizable using low temperature hydrogen peroxide gas, such as STERRAD and STERRIS V-PRO™*

Order Information

Catalog #	Quantity Per Package	Product Description
SCD13	6 each	13 cm Cordless Ultrasonic Dissector (1 torque wrench included with each dissector)
SCD26	6 each	26 cm Cordless Ultrasonic Dissector (1 torque wrench included with each dissector)
SCD391	1 each	39 cm Cordless Ultrasonic Dissector (Select countries only; 1 torque wrench included)
SCD396	6 each	39 cm Cordless Ultrasonic Dissector (1 torque wrench included with each dissector)
SCD48	6 each	48 cm Cordless Ultrasonic Dissector (1 torque wrench included with each dissector)
SCG	1 each	Reusable Generator
SCB	1 each	Reusable Battery
CBC	1 each	Battery Charger
SCST	1 each	Sterilization Tray
SCST1	1 each	Sterilization Tray (For Outside of the U.S. Customers)

Change the energy in your OR. Contact your Covidien Sales Representative at 1-800-255-8522 or visit Covidien.com/Sonicision today!