



DESIGNED TO SIMPLIFY **COMPLEX** **INSERTIONS**

Tips for using Dover™ specialty catheters


Your patients have specialized needs. They need specialized products. That's why we developed the Dover™ specialty catheters.

We offer catheters:

- With coude, urethral coude, and council tips
- Designed for transurethral resection of the prostate (TURP)
- With temperature sensor

Our TURP catheter offers:

- An extra strong balloon that can be inflated up to 90 cc
- Four or six large staggered drain eyes
- Embedded nylon spirals to provide stability to catheter wall



Our catheter with coude tip offers:

- Curved tip designed to help navigate the compromised urethral track



Our 400 series catheter with temperature sensor:


- Is MRI compatible
- Is available in a closed system tray



Our coude tip catheter features a funnel bump, indicating the direction of the curved tip.



Our catheter with council tip allows for placement over a guidewire or stylet.



Medtronic
Further, Together

SPECIALIZED PRODUCTS FOR UNIQUE NEEDS

Your patient might need a specialty catheter because of:

- Acute urinary retention
- Postoperative procedures
- Negotiation of prostatic urethra
- Continuous irrigation
- Continuous temperature monitoring

Choose the best way to fit your patient's need

	Size	Capacity		Material				2-way Foley	3-way Foley	Hydrogel coating	Radiopaque
		5 cc	30 cc	Silver	Silicone	Latex	Red latex				
Coudé	12–24 Fr	■	■	■	■	■		■		■	■
Urethral coudé	14–16 Fr						■			■	
Council	16–24 Fr	■			■			■			■
Temperature sensing	14–18 Fr	■		■	■			■		■	
TURP	22–24 Fr	Up to 90 cc				■		■	■		

For complete insertion instructions, consult the product IFU.

© 2017 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.
05/2017 – H9622 – [WF# 1497472]

15 Hampshire Street
Mansfield, MA 02048

800.962.9888
508.261.8000

medtronic.com/covidien

Medtronic