
















Compatibility of patient cables with temporary leads and catheters

Compatible patient cables	Lead or catheter used		
<p>5492AL 5492VL</p> 	 <p>Myocardial wires</p>	 <p>Unprotected (exposed) connectors</p>	
<p>5487L</p> 	 <p>Myocardial wires</p>	 <p>Stimulation probe with IS1 connector</p>	
<p>5832S</p> 	 <p>Myocardial wires</p>	 <p>Unprotected (exposed) connectors</p>	 <p>Stimulation probe with IS1 connector</p>
<p>5833SL</p> 	 <p>Myocardial wires</p>	 <p>Unprotected (exposed) connectors</p>	 <p>Stimulation probe with IS1 connector</p>
<p>5846AL 5846VL</p> 	 <p>Protected (shrouded) connectors</p>		

General Indication statements:

For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the Instructions for Use. Consult instructions for use on this website <http://manuals.medtronic.com/manuals/main/region>. Manuals can be viewed using a current version of any major Internet browser. For best results, use Adobe Acrobat Reader® with the browser.

Medtronic

Medtronic Canada

99 Hereford Street
Brampton, Ontario, L6Y 0R3
Toll-free: 800.268.5346
Tel: 905.460.3800

medtronic.ca

©2024 Medtronic. All rights reserved. Medtronic, Medtronic logo and Engineering the extraordinary are trademarks of Medtronic.
™*Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
UC202313699EC CA-CTL-0289-E Rev. 02/2024