






SPECIMEN RETRIEVAL BAGS

A comprehensive
portfolio from
Medtronic

Bag Dimensions and Volume	Product Name	Order Code	Bag Material	Trocar Size
 <p>Volume 275 ml</p> <p>2.44" 6.53"</p>	ReliaCatch™ device 10 mm	CATCH10	Ripstop nylon	11 mm or 12 mm
 <p>Volume 220 ml</p> <p>2.28" 5.9"</p>	Endo Catch™ Gold device 10 mm	173050G	Polyurethane	11 mm or 12 mm
 <p>Volume 1500 ml</p> <p>3.83" 9.3"</p>	Endo Catch™ II device 15 mm	173049	Polyurethane	15 mm
 <p>Volume 150 ml</p> <p>3" 6"</p>	Endo Bag™ device	25030	Polyurethane	10 mm, 11 mm, or 12 mm
 <p>Volume 750 ml</p> <p>5" 8"</p>	Endo Bag™ device	25040	Polyurethane	10 mm, 11 mm, or 12 mm

For more information, call **800-772-8772**
or visit us at **medtronic.com/covidien**