


Medtronic

TruClear™ hysteroscope set duo

Integrate hysteroscopy into your office.



Perform see-and-treat hysteroscopic procedures with our expanded portfolio of TruClear™ hysteroscopes.

The latest to join the TruClear™ system family of products, the TruClear™ hysteroscope set duo provides the option of a diagnostic or operative sheath for the office setting.

It's the complete hysteroscopic offering you need to expand diagnostic and therapeutic procedures for your patients in the comfort and convenience of your office.

Diagnostic benefits

- Easy insertion into the cervical canal and uterine cavity with slanted tip design¹
- Angled view you need to see the uterine cavity with a 30° scope
- Evaluate the uterine cavity with a small 4-mm diagnostic sheath
- Little-to-no dilation with diagnostic sheath^{1,2}

Operative benefits


- Maintains clear operative field with built-in continuous flow^{2,3}
- Access to the uterine cavity with long scope length
- Can be used with your preferred hand instrument†

Practice benefits


- Perform see-and-treat hysteroscopic procedures with the flexibility of using either a diagnostic and operative sheath

†Hand instruments must be 5Fr diameter and 18 inches in length

TruClear™ operative sheath duo




TruClear™ diagnostic sheath duo



Ordering information

Product number	Description	Unit of measure
8300057945	TruClear™ Hysteroscope Duo	Each
8300057946	TruClear™ Diagnostic Sheath Duo	Each
8300057947	TruClear™ Operative Sheath Duo	Each
8300055246	Rubber Replacement Seal	Each
8300059605	TruClear™ Stopcock Replacement Set, Diagnostic (for use with 8300057946)	Each
8300059594	TruClear™ Threaded Ring, Duo (for use with 8300057946)	Each
8300059604	TruClear™ O-Ring, Duo (for use with 8300057946)	Pack of 5
8300054627	Polishing Paste (for use with 8300057945)	Each
8300055288	Small Stopcock Replacement Set (for use with 8300057947)	Each



Discover more ways to streamline your in-office procedures with the TruClear™ hysteroscopic tissue removal system.

For more information, call your Medtronic sales representative or visit [medtronic.com/covidien/en-gb/clinical-solutions/office-hysteroscopy](https://www.medtronic.com/covidien/en-gb/clinical-solutions/office-hysteroscopy)

1. Internal test report RE00440188 Rev Duo Insertion Claim Test Report, 26 January 2023
2. Internal test report RE00192361 Dynamite Hysteroscope Evaluation Testing Report, 29 March 2019
3. Internal test report RE00192361-1 Appendix A - OEM-Private Label Screening Test Results, 29 March 2019

© 2023 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. EMEA-UW-2300003-truclear-duo-info-sheet-en-we-8673516

[medtronic.com/covidien/uk](https://www.medtronic.com/covidien/uk)

Medtronic