

US Returned Product Analysis Summary

Model Number	Family	US Market Release	Estimated US Implants	Estimated US Active	Implant Damage	Electrical Malfunction	Other
3830	SelectSecure	Aug-05	12,200	10,600	21	5	1
4003, 4003M	CapSure	Jul-86	38,000	5,300	24	60	2
4004, 4004M	CapSure	Feb-89	72,600	5,800	56	688	19
4011	Target Tip	Nov-82	58,400	6,200	29	152	5
4012	Target Tip	Jul-83	93,700	6,500	50	827	34
4023	CapSure SP	Aug-91	41,200	9,600	47	21	6
4024	CapSure SP	Oct-91	222,100	56,100	264	135	34
4033	CapSure Z	Not US released	NA	NA	2	0	0
4057, 4057M	Screw-in	Aug-88	10,100	1,800	39	6	4
4058, 4058M	Screw-in	Jan-89	101,900	20,800	388	261	23
4067	CapSureFix	Jan-97	1,000	300	3	1	1
4068	CapSureFix	Mar-96	124,800	43,900	406	111	11
4073	CapSure Sense	Jun-02	600	400	1	0	0
4074	CapSure Sense	Jun-02	63,000	45,200	13	8	1
4076	CapSureFix Novus	Feb-04	199,200	165,500	96	13	8
4081	Target Tip	Jul-89	3,900	800	4	5	0
4092	CapSure SP Novus	Sep-98	151,000	83,100	39	19	5
4503, 4503M	CapSure	Jul-86	8,000	1,500	2	12	0
4504, 4504M	CapSure	Mar-90	15,400	1,700	5	172	4
4512	Target Tip	Jul-83	10,300	1,200	4	85	8
4523	CapSure SP	Aug-91	11,200	3,200	5	2	1
4524	CapSure SP	Oct-91	101,700	32,300	47	31	8
4533	CapSure Z	Not US released	NA	NA	0	0	0
4557, 4557M	Screw-in	Aug-88	19,700	4,400	53	14	4
4558M	Screw-in	Nov-94	20,000	5,400	111	12	1
4568	CapSureFix	Jan-97	69,700	31,700	198	18	4
4574	CapSure Sense	Jun-02	39,900	30,300	7	2	0
4592	CapSure SP Novus	Oct-98	74,400	43,500	13	5	0
5023, 5023M	CapSure SP	Nov-88	9,900	2,300	15	7	0
5024, 5024M	CapSure SP	Mar-90	201,500	56,200	723	115	29
5026	CapSure	Feb-88	7,400	1,100	60	7	1
5033	CapSure Z	Feb-96	2,400	700	6	1	3
5034	CapSure Z	Feb-96	56,300	17,000	85	31	11
5054	CapSure Z Novus	Jun-98	85,400	44,900	43	16	6
5068	CapSureFix	Jan-97	103,100	40,600	455	75	15
5072	SureFix	Jun-98	8,900	4,500	28	5	1
5076	CapSureFix Novus	Aug-00	1,025,600	723,800	949	276	84
5092	CapSure SP Novus	Jun-98	111,800	61,600	48	28	11
5524, 5524M	CapSure SP	Mar-90	60,600	21,000	67	25	7
5534	CapSure Z	Feb-96	26,300	9,300	29	8	5
5554	CapSure Z Novus	Jun-98	54,800	31,500	8	12	4
5568	CapSureFix	Jan-97	70,400	45,000	264	15	12

continued

US Returned Product Analysis Summary continued

Model Number	Family	US Market Release	Estimated US Implants	Estimated US Active	Implant Damage	Electrical Malfunction	Other
5592	CapSure SP Novus	Jun-98	27,800	18,100	6	3	0
5594	CapSure SP Novus	Jun-01	10,908	1,700	0	4	0
6940	CapSureFix	Oct-98	25,500	10,100	114	21	3
6957	Spectraflex	Jul-79	29,100	2,000	85	39	25
6957J	Spectraflex	Sep-80	30,000	2,100	74	30	30
6961	Tenax	Jan-78	44,700	2,100	103	27	0
6962	Tenax	Jan-78	70,600	3,200	170	84	0