



FACT SHEET

Contact:  
Deanne McLaughlin  
818-576-4325

## Paradigm™ 511 Insulin Infusion Pump

### What is the Paradigm<sup>®</sup> pump?

The Paradigm insulin infusion pump is Medtronic MiniMed's newest product. This small, battery-operated device (about the size of a pager) delivers insulin to patients with diabetes and makes pump therapy even easier. The Paradigm pump features a simplified menu system, making it easy to teach, learn and use. With an emphasis on ease-of-use, the pump is expected to appeal to a broader range of diabetes patients and to physicians prescribing pump therapy.

### Key features:

<u>Easy to use</u>	The Paradigm insulin infusion pump introduces <i>E-Z Path Programming<sup>®</sup></i> , a simplified and intuitive menu system.
<u>AAA battery</u>	The Paradigm pump uses a standard AAA alkaline battery, which is affordable and widely available.
<u>Small</u>	The Paradigm is the smallest, full-featured pump available and is easy to wear.
<u>Safe</u>	The Paradigm pump conducts more than 6 million safety checks each day. Programming block options prevent tampering by small children.
<u>Water Use</u>	Protection against splashes or a temporary dunk in water.
<u>Remote control</u>	Diabetes patients can use a small, hand-held remote programmer to discreetly program insulin delivery without accessing the pump.
<u>Multiple Bolus &amp; Basal Options</u>	The Paradigm pump contains various insulin delivery options to meet patients' needs, including Normal, Square™ and Dual Wave™ bolus capabilities, and multiple basal patterns and a temporary basal rate. In addition, the pump uses BioPulse™ delivery for precise insulin control.
<u>Color options</u>	Clear, blue and charcoal

### **Obtaining a Medtronic MiniMed pump**

An insulin pump requires a prescription from a physician. Physicians and diabetes patients can contact Medtronic MiniMed to obtain information about a pump by calling 1-866-MiniMed, (818) 362-5958 outside the United States., or on our website at [www.minimed.com](http://www.minimed.com). Medtronic MiniMed is dedicated to customer service, providing pump initiation training, and educational support for physicians and diabetes patients.

### **Insurance coverage**

Medtronic MiniMed has an extensive payor contract base, and the Paradigm pump is covered by most insurance plans. Medtronic MiniMed provides insurance verification, billing and reimbursement assistance to its customers, making the adoption of pump therapy as easy as possible.

**For a brief statement about the use of Insulin Pump Therapy, see the enclosed insert.**