



CARDIOSIGHT™ SERVICE


Powered by Medtronic CareLink® Network



1 Patient checks in at the clinic receptionist desk and proceeds to the waiting room until called by a clinic staff member.

2 Clinic staff member performs a CardioSight reading to obtain information from the patient's Medtronic device.

3 Clinic staff member escorts the patient to the exam room.



4 A staff member retrieves the faxed report and places it in the patient record.

5 A clinician reviews the cardiac trends to help assess the patient's status.

6 The clinician meets with the patient to discuss status and management.

Direct Access. Unique Insight.

Brief Statement

Intended Use

The Model 2020A CardioSight™ Reader is intended for use in a clinical setting. The CardioSight Reader is indicated for use to interrogate compatible Medtronic implantable devices to collect patient and device data and send the information to the clinician. The Model 2020A CardioSight Reader cannot be used to change therapy.

Contraindications

There are no contraindications for the Model 2020A CardioSight Reader.

Warnings and Precautions

The CardioSight Reader must only be used for interrogating compatible Medtronic implantable devices. Do not use a cellular phone while the antenna is positioned over the implanted device. The CardioSight Reader is designed for use in the continental United States, Alaska, and Hawaii.

See the appropriate technical manual for detailed information regarding instructions for use, indications, contraindications, warnings and precautions, and potential adverse events. For further information, please call Medtronic at 1-800-328-2518 and/or consult Medtronic's website at www.medtronic.com.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

World Headquarters

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432-5604
USA
Tel: (763) 514-4000
Fax: (763) 514-4879
www.medtronic.com

Medtronic CardioSight Service
Tel: 1 (800) 929-4043
(24-hour technical support for
physicians and medical professionals)
www.medtronic.com/cardiosight

Medtronic USA, Inc.
Toll-free: 1 (800) 328-2518

UC200600477 EN
© Medtronic, Inc. 2005
All Rights Reserved
Printed in USA
August 2005

